

OCT 26 2001



K01335.5

AMERICAN MEDICAL SYSTEMS

## 510(k) SUMMARY

**Submitter's Name:** American Medical Systems, Inc.

**Address:** 10700 Bren Road West  
Minnetonka, MN 55343

**Tel:** 952-930-6541

**Fax:** 952-930-6496

**Contact Person:** Ginger Sackett Glaser

**Date of Summary Preparation:** October 5, 2001

**Device Common Name:** Surgical Mesh, Sling, Urethral Sling

**Device Trade Name:** SPARC™ Sling System

**Device Classification Name:** Surgical Mesh, polymeric

**Predicate Device:** SPARC™ Sling System – K011251

### Device Description

The SPARC™ Sling System is a sterile, single use procedure kit consisting of:

- Two stainless steel, curved, 22-cm long, needle passers (also called insertion tools).
- One piece of AMS Polypropylene sling mesh with attached dilating connectors. The AMS Polypropylene sling mesh is constructed of polypropylene monofilament that is precut to 1.0cm width x 50cm length. A fixed blue polypropylene anchoring suture runs through the middle of the sling mesh. Two plastic sheaths that overlap in the center of the sling mesh, cover the sling mesh and protect it during placement.

Dilating connectors are attached to either end of the plastic sheaths. The dilating connectors are used to attach to the vaginal ends of the SPARC™ needle passers during the procedure to facilitate sling placement.

- Two blue colored plastic cystoscopy aids are included in the kit in order to facilitate cystoscopic viewing of the bladder. The use of these cystoscopy aids is optional.



### **Indications for Use**

The indication for use for the SPARC™ Sling System is not changing. It continues to read as follows:

“The SPARC™ Sling System is intended for the placement of a pubourethral sling for the treatment of female stress urinary incontinence (SUI) resulting from urethral hypermobility and/or intrinsic sphincter deficiency.”

### **Comparison to Predicate Device**

The fundamental scientific technology of the device will not change with the proposed alternative configuration of the device. A polypropylene suture identification loop with an attached silicone tab has been added for easier location of the portion of the tensioning suture that should be snipped after sling tensioning is complete.

### **Supporting Information**

A risk analysis of the proposed modification and bench test data reported in this 510(k) application substantiate equivalence to the predicate and did not raise any new questions of safety or effectiveness.

### **Conclusion**

The proposed modification is equivalent to the predicate with respect to intended use, technological characteristics, and performance.



OCT 26 2001

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Ms. Ginger Sackett Glaser  
Senior Regulatory Affairs Specialist  
American Medical Systems, Inc.  
10700 Bren Road West  
Minnetonka, Minnesota 55343

Re: K013355

Trade/Device Name: SPARC™ Sling System  
Regulation Number: 21 CFR 878.3300  
Regulation Name: Surgical Mesh  
Regulatory Class: Class II  
Product Code: FTL  
Dated: October 9, 2001  
Received: October 10, 2001

Dear Ms. Glaser:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 21 CFR Part 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,



Celia M. Witten, Ph.D., M.D.

Director

Division of General, Restorative  
and Neurological Devices

Office of Device Evaluation

Center for Devices and  
Radiological Health



Enclosure

**INDICATIONS FOR USE ENCLOSURE**

**510(k) Number:**

K013355

**Device Name:**

SPARC™ Sling System

**Indications for Use:**

The SPARC™ Sling System is intended for the placement of a pubourethral sling for the treatment of female stress urinary incontinence (SUI) resulting from urethral hypermobility and/or intrinsic sphincter deficiency.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE  
IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)  
(Division Sign-off)  
Division of General and Restorative Devices

510(k) Number

K013355

Prescription Use \_\_\_\_\_  
(Per 21 CFR801.109)

OR

Over the Counter Use \_\_\_\_\_



(Division Sign-Off)

Division of General, Restorative  
and Neurological Devices

510(k) Number

K013355