

TRANSCRIPT OF PROCEEDINGS

8158 '00 FEB -8 P2:38

DEPARTMENT OF HEALTH AND HUMAN SERVICES
FOOD AND DRUG ADMINISTRATION
CENTER FOR DEVICES AND RADIOLOGICAL HEALTH

OBSTETRICS AND GYNECOLOGY DEVICES PANEL

62nd MEETING

This transcript has not
been edited and FDA
makes no representation
regarding its accuracy

Pages 1 thru 293

Gaithersburg, Maryland
January 24, 2000

MILLER REPORTING COMPANY, INC.

507 C Street, N.E.
Washington, D.C. 20002
(202) 546-6666

sgg

TCB

1

DEPARTMENT OF HEALTH AND HUMAN SERVICES
FOOD AND DRUG ADMINISTRATION
CENTER FOR DEVICES AND RADIOLOGICAL HEALTH

OBSTETRICS AND GYNECOLOGY DEVICES PANEL

62nd MEETING

Monday, January 24, 2000

9:00 a.m.

Gaithersburg Holiday Inn
Gaithersburg, Maryland

MILLER REPORTING COMPANY, INC.
507 C Street, N.E.
Washington, D.C. 20002
(202) 546-6666

P A R T I C I P A N T S

Jorge Blanco, M.D., Panel Chair
Elisa Harvey, D.V.M., Ph.D., Executive Secretary

Panel Voting Members:

Donald Chatman, M.D.
Subir Roy, M.D.
Nancy Sharts-Hopko, Ph.D.
Machelle Allen, M.D.
Ralph D'Agostino, Ph.D.
Mike Diamond, M.D.
Gary Eglinton, M.D.
Jay Iams, M.D.
Michael Neuman, Ph.D., M.D.
Mary Jo O'Sullivan, M.D.
Robert Wolfson, Ph.D., M.D.

Industry Representative

Gary Jarvis

Consumer Representative

Diony Young

C O N T E N T S

Meeting Comes to Order, Jorge Blanco, M.D.	4
Introductory Comments, Colin Pollard, Chief, Obstetrics and Gynecology Devices Branch	10
Regulatory Issues, Elisa Harvey, D.V.M., Ph.D.	14
Open Public Hearing:	
Susan M. Ramin, M.D., FACOG, ACOG	17
Ms. Susan Meikle, NIH	25
Sponsor Presentation	
Nellcor N-400 Fetal Oxygen Saturation Monitoring System:	
Introduction and Proposed Indication for Use, Ms. Donna Page	27
Technology and Critical Threshold, David B. Swedlow, M.D.	30
Pivotal Study Design and Results, Thomas J. Garite, M.D.	36
Cesarean Section for Dystocia and Clinical Utility, Richard P. Porreco, M.D.	49
Impact on Clinical Practice, Frank Boehm, M.D.	55
Nursing Perspective, Nancy Townsend, RN, MSN	61
Patient Perspective, Ms. Lucy Woods	67
Restatement of Proposed Indication for Use, Ms. Donna Page	69
FDA Presentations:	
Preclinical Aspects, Ms. Kathryn S. Daws-Kopp	90
Clinical Aspects, Diane Mitchell, M.D.	95
Panel Deliberations	129
Open Public Hearing:	
Michael Ross	228

P R O C E E D I N G S

1
2 CHAIRMAN BLANCO: Let's go ahead and get started.
3 There are a couple of panel members who are on the way in,
4 so we will go through some of the preliminaries so that we
5 can stay on time. This is a lot of information and a lot of
6 things that we are going to be doing today, so we need to
7 try to make sure that we get going. We want to make sure
8 that we stay on time, and it is very important that we give
9 all the appropriate time to the PMA that we are evaluating
10 today.

11 I would like to go ahead and formally call the
12 meeting to order. I want to remind everyone that there is a
13 sign-in sheet by the door. If you would please sign in, let
14 us know who you are and who was here.

15 Before we have a strict agenda, and there is a
16 time period for comments from the audience. We know of two
17 organizations, people, that want to speak before the panel.
18 If you would like to speak before the panel, that is the
19 time to do it. If you feel like you would like to make a
20 comment during the panel deliberations, you must be
21 recognized. We do not accept outbursts from the audience,
22 and any time that you are coming forward to speak, you need
23 to come to the mike and identify yourself. You need to
24 identify whether you have any conflicts of interest. That
25 means if any organization, company, etc. has funded any or

1 part of your trip, research, or any other possible conflict
2 of interest, you really need to let us know for the record
3 what that relationship is.

4 At this time, I would like to just go around
5 through the panel and have everyone introduce themselves to
6 the audience, and if we can go ahead and start from this
7 left side.

8 MR. JARVIS: Gary Jarvis, the industry
9 representative.

10 MS. YOUNG: I am Diony Young, the consumer
11 representative from Genesco, New York.

12 DR. ROY: Subir Roy, from the University of
13 Southern California.

14 DR. SHARTS-HOPKO: Nancy Sharts-Hopko, from
15 Villanova University.

16 DR. DIAMOND: Michael Diamond, Professor of
17 Obstetrics and Gynecology, Wayne State University in
18 Detroit.

19 DR. IAMS: Jay Iams, obstetrician from Ohio State
20 University.

21 CHAIRMAN BLANCO: I am Jorge "George" Blanco, the
22 University of Florida at Pensacola.

23 DR. HARVEY: I am Elisa Harvey, from the Center
24 for Devices. I am the Executive Secretary for the
25 Obstetrics and Gynecology Devices Panel.

1 DR. EGLINTON: Gary Eglinton, New York Hospital,
2 Queens.

3 DR. O'SULLIVAN: Mary Jo O'Sullivan, University of
4 Miami, OB-GYN.

5 DR. WOLFSON: Robert Wolfson, Colorado Springs,
6 OB-GYN/Perinatology.

7 MS. ALLEN: Machele Allen, OB-GYN, NYU, Bellevue.

8 DR. D'AGOSTINO: Ralph D'Agostino, Boston
9 University, biostatistician.

10 DR. CHATMAN: Donald Chatman, obstetrician-
11 gynecologist, Northwestern University.

12 DR. NEUMAN: Michael Neuman, from the Joint
13 Program in Biomedical Engineering of the University of
14 Tennessee and the University of Memphis.

15 DR. SCHULTZ: I am Dan Schultz. I am the Acting
16 Director of the Division of Reproductive, Abdominal, and
17 Radiological Devices, Office of Device Evaluation, Center
18 for Devices, FDA.

19 CHAIRMAN BLANCO: All right. Thank you very much.
20 I also would like to let everyone know that Dr. Dan Schultz,
21 the Acting Division Director, is the FDA press contact, and
22 if anyone would like some press information, he is the
23 person to get in touch with. Now Dr. Harvey is going to do
24 some more of the preliminaries.

25 DR. HARVEY: I would like to start by reading a

1 couple of documents. One is the appointment to temporary
2 voting status for today:

3 Pursuant to the authority granted under the
4 Medical Devices Advisory Committee Charter, dated October
5 27, 1990 and amended April 20, 1995, I appoint the following
6 people as voting members of the Obstetrics and Gynecology
7 Devices Panel for the duration of this panel meeting on
8 January 24, 2000: Dr. Machelles Allen; Dr. Ralph D'Agostino;
9 Dr. Michael Diamond; Dr. Gary Eglinton; Dr. Jay Iams; Dr.
10 Michael Neuman; Dr. Mary Jo O'Sullivan; and Dr. Robert
11 Wolfson. For the record, these people are special
12 government employees and are consultants to this panel.
13 They have undergone the customary conflict of interest
14 review and they have reviewed the material to be considered
15 at this meeting.

16 And it is signed by Dr. David Feigal, the Director
17 of the Center for Devices and Radiological Health.

18 The second document I would like to put into the
19 record is the conflict of interest statement for this
20 meeting:

21 The following announcement addresses conflict-of-
22 interest issues associated with this meeting and is made a
23 part of the record to preclude even the appearance of an
24 impropriety. To determine if any conflict existed, the
25 agency reviewed the submitted agenda and all financial

1 interests reported by the committee participants. The
2 conflict-of-interest statutes prohibit special government
3 employees from participating in matters that could affect
4 their or their employer's financial interest.

5 However, the agency has determined that
6 participation of certain members and consultants, the need
7 for whose services outweighs the potential conflict of
8 interest involved is in the best interest of the government.
9 A waiver has been granted for Dr. Donald Chatman for his
10 interest in firms that could potentially be affected by the
11 panel's deliberations. The waiver allows him to participate
12 fully in all matters before the panel today.

13 Copies of this waiver may be obtained from the
14 agency's Freedom of Information Office, Room 12A-15 of the
15 Parklawn Building.

16 We would like to note for the record that the
17 agency took into consideration certain matters regarding
18 Drs. Michael Neuman and Robert Wolfson. These individuals
19 reported past or current interests in firms at issue, but in
20 matters not related to the topics for today's session.
21 Therefore, the agency has determined that they may
22 participate fully in the deliberations.

23 In the event that the discussions involve any
24 other products or firms not already on the agenda for which
25 an FDA participant has a financial interest, the participant

1 should excuse him or herself from such involvement and the
2 exclusion will be noted for the record.

3 With respect to all other participants, we ask, in
4 the interest of fairness, that all persons making statements
5 or presentations disclose any current or previous financial
6 involvement with any firm whose products they may wish to
7 comment upon.

8 The other things I would like to point out are
9 that there is information on getting transcripts and videos
10 of today's meeting at the table at the back of the room.

11 Anyone who has any comments to make to the panel,
12 if you could provide a hard copy with your remarks, that
13 would be helpful. Mr. Mike Kuchinsky, at the podium, will
14 take those from you.

15 The last thing I would like to point out is, for
16 the panel's sake, what the panel folder contents are so that
17 they can follow through today's proceedings. You should
18 have a copy of the agenda, the discussion questions, and the
19 panel roster, as well as the presentations for today. I
20 will be giving some information on regulatory definitions.
21 You have the sponsor's presentation in your folder. You
22 have presentations from FDA, by Kathy Daws-Kopp and Diane
23 Mitchell. You have some information from a representative
24 of ACOG, and you have a previous statement that was made to
25 the OB-GYN Devices Panel in 1996 by Dr. Larry Gilstrap.

1 Thank you.

2 CHAIRMAN BLANCO: Thank you, Dr. Harvey. All
3 right. It's my pleasure now to begin the meeting by
4 introducing Mr. Colin Pollard, Chief of the Obstetrics and
5 Gynecologic Devices Branch, who will give us some
6 information.

7 **Introductory Comments**

8 MR. POLLARD: Thank you, Dr. Blanco. Good
9 morning, ladies and gentlemen of the panel, distinguished
10 audience.

11 We have brought you together today to consider the
12 premarket approval application, or PMA, submitted by the
13 Nellcor Perinatal Business of Mallinckrodt, for its
14 intrapartum fetal oxygen saturation monitoring system, the
15 Nellcor N-400.

16 The sponsor has proposed that this monitor be
17 indicated for women who are in labor with term pregnancies
18 when the strip-chart tracing from conventional intrapartum
19 fetal monitoring is non-reassuring. Immediately following
20 my opening remarks, Dr. Harvey, Executive Secretary of your
21 panel, will go over the basic ground rules of your panel
22 deliberations on this PMA, especially with respect to what
23 constitutes valid scientific evidence, safety, and
24 effectiveness.

25 I don't think I need to tell you that we consider

1 this a very important PMA. With nearly four million births
2 in the U.S. each year, the currently proposed indication for
3 this new sensor may impact more than a quarter of that
4 number.

5 Some of you will recall that we convened this
6 panel three-and-a-half years ago, in July of 1996, to
7 consider this technology and others like it in a general
8 way. We, at FDA, wanted to develop a guidance document that
9 would help manufacturers and clinical researchers put
10 together cogent clinical development plans for products like
11 this.

12 At that meeting, Nellcor shared its plan with the
13 panel for the pivotal clinical study that would support its
14 future PMA. Besides the FDA-invited guest speakers, several
15 other manufacturers and researchers also addressed the
16 panel.

17 You should know that we have put issuance of this
18 draft guidance document aside for the moment so that we can
19 digest how this first PMA goes and the panel input on it.

20 From a regulatory viewpoint, our meeting today is
21 a natural progression from our 1996 meeting, because we are
22 now going to look at the data from that clinical study to
23 see whether it supports approval of that PMA. I bring to
24 your attention that two other PMAs will hinge on the outcome
25 of this PMA before you today. These two secondary PMAs are

1 from GE\Marquette and Agilant Technologies, who will
2 integrate the Nellcor fetal pulse oximetry technology into
3 their currently marketed fetal monitors, the Corometrics and
4 Hewlett-Packard monitors, respectively.

5 Although there certainly is important information
6 and data for FDA to review in these two secondary PMAs to
7 ensure that the technology is integrated properly, we do not
8 plan to bring either of them before the panel and they are
9 not the subject of today's agenda.

10 We have tried our best to bring together a truly
11 top-notch panel, with experts in instrumentation and
12 biostatistics, not to mention extensive representation by
13 perinatologists. We are fortunate to have here today many
14 of the same panel participants from that 1996 meeting,
15 including our panel chair, Dr. Blanco, as well as Dr.
16 Eglinton, Dr. Neuman, Dr. Allen, and Dr. Diamond. We
17 believe that this will add some regulatory continuity to our
18 review process.

19 Because this product potentially represents a big
20 step for intrapartum clinical management in the U.S., we
21 also strengthened the perinatology expertise on the panel by
22 adding maternal fetal medicine specialists, Dr. Iams, Dr.
23 O'Sullivan, and Dr. Wolfson.

24 You should also know that, using one of our newer
25 PMA approaches, Mallinckrodt/Nellcor submitted this PMA in a

1 shell/module configuration. As our review team will explain
2 to you later this morning, this allowed us to review and
3 close out a number of modules of preclinical information.

4 So, as you well know, the PMA at this point is
5 based primarily on a couple of key clinical studies,
6 including a pivotal randomized control trial that looks at
7 the effect of the new monitor on intervention.

8 My only point here is that there is a lot of data
9 here and the study results are complex. The panel is a
10 little larger than usual, and you have to work your way
11 completely through the agenda in the limited amount of time
12 we have today. We have impressed upon your panel chair, Dr.
13 Blanco, the importance of due process and the need to arrive
14 at a panel recommendation at the end of the day that
15 conforms to one of the three formats that Dr. Harvey will
16 explain to you in a minute.

17 To achieve that end, I only ask that we all stay
18 focused, keep our remarks succinct, and respect each other's
19 views. If we succeed with those three things, Dr. Blanco's
20 job will be a lot easier, and we have a very good chance for
21 a successful outcome for this meeting.

22 I don't want to take up any more of your time, so
23 my comments are concluded. Thank you, Dr. Blanco.

24 CHAIRMAN BLANCO: Thank you, Colin. We will try
25 to rise to the occasion and ensure that we fulfill your

1 expectations. I think Dr. Harvey is now going to speak to
2 us on some regulatory issues.

3 **Regulatory Issues**

4 DR. HARVEY: Mike should be putting this up on the
5 slide, if he can, but you have the full handout in your
6 folder as well. I want to provide for the panel and the
7 audience reminders of the regulatory definitions that we are
8 obliged to adhere to today.

9 The first is valid scientific evidence. Valid
10 scientific evidence is evidence from well-controlled
11 investigations, partially-controlled studies, studies and
12 objective trials without matched controls, well-documented
13 case histories conducted by qualified experts, and reports
14 of significant human experience with the marketed device
15 from which it can fairly and responsibly be concluded by
16 qualified experts that there is reasonable assurance of the
17 safety and effectiveness of a device under its conditions of
18 use.

19 The definition of safety: There is reasonable
20 assurance that a device is safe when it can be determined,
21 based upon valid scientific evidence, that the probable
22 benefits to health from the use of the device for its
23 intended uses and conditions of use, when accompanied by
24 adequate directions and warnings against unsafe use,
25 outweigh any probable risks.

1 The definition of effectiveness is that there is
2 reasonable assurance that a device is effective when it can
3 be determined, based upon valid scientific evidence, that in
4 a significant portion of the target population, the use of
5 the device for its intended uses and conditions of use, when
6 accompanied by adequate directions for use and warnings
7 against unsafe use, will provide clinically significant
8 results.

9 I also want to point out that your PMA review
10 should be independent of cost, any previous regulatory
11 difficulties, clinical data submitted in any other PMAs, or
12 the medical-legal climate and its effect on the standard of
13 care.

14 At the end of the day, you will be voting on this
15 premarket approval application, and your options will be one
16 of the three: The first is you will vote for recommendation
17 of approval with no conditions attached to the approval. A
18 second option will be approvable subject to specified
19 conditions, and these are such as resolution of very clearly
20 identified deficiencies cited either by you, the panel, or
21 FDA staff. Examples could include resolutions of questions
22 concerning some of the data or changes in the draft
23 labeling.

24 You may conclude that post-approval requirements
25 should be imposed as a condition of approval. These

1 conditions may include a continuing evaluation of the device
2 and submission of periodic reports. If you believe that
3 such requirements are necessary, your recommendation must
4 address the following points: The reason or purpose of the
5 requirement; the number of patients to be evaluated; and the
6 reports required to be submitted.

7 Your third voting option will be not approvable,
8 and if you vote in that way, you must have one of the
9 following reasons for recommending not approvable: Either
10 safety -- and that is that the data do not provide
11 reasonable assurance that the device is safe under the
12 conditions of use prescribed, recommended, or suggested in
13 the proposed labeling; or effectiveness -- reasonable
14 assurance has not been given that the device is effective
15 under the conditions of use in the labeling; and third,
16 labeling -- based on a fair evaluation of all the material
17 facts in your discussions, you believe the proposed labeling
18 to be false or misleading.

19 Thank you, Dr. Blanco.

20 CHAIRMAN BLANCO: Thank you, Dr. Harvey. I just
21 would like to add to that, having been at some of these a
22 few times before, for the new panel members especially, that
23 after your vote, you are also asked to spend a few minutes
24 justifying why you voted the way you did, or at least trying
25 to explain for the record the way that you voted.

1 We are moving along very nicely, and it is now
2 time for the public comments. At this time, I have two
3 individuals who have requested time for public comment. I
4 again would like to remind these individuals to please note
5 any conflict of interest. I also would like to remind them
6 to keep their notes to the allotted five minutes for each
7 individual.

8 The first individual that I have is Dr. Susan
9 Ramin, who I believe is representing the American College of
10 Obstetricians and Gynecologists. Is that correct?

11 **Open Public Hearing**

12 MS. RAMIN: Good morning. My name is Dr. Susan
13 Ramin from the University of Texas, Houston Medical Center,
14 and I do not have any conflict of interest.

15 Let's begin with basically just a brief background
16 in fetal pulse oximetry. It is a new technology that,
17 hopefully, has the potential to aid in evaluating the fetus.
18 It specifically measures -- or at least in the simplest
19 terms, it measures the level of oxygen in the fetus. And
20 this is important because lack of oxygen during labor can
21 result in neurologic damage and cerebral palsy in a newborn.

22 Now, currently, our ability to monitor or to
23 evaluate or assess fetal well-being includes monitoring the
24 fetal heart rate, either by auscultation or electronically,
25 and also the use of fetal scalp blood pH.

1 When we look at electronic fetal heart rate
2 monitoring, which was developed in the late 1960s, it was
3 hoped that it would decrease the incidence of cerebral palsy
4 and neonatal mortality but, unfortunately, as we have
5 discovered over the last three decades, this has not done
6 so. The reason for this is because the electronic fetal
7 heart rate monitor is sensitive but it is not very specific.

8 In other words, a normal fetal heart rate pattern
9 does predict a good neonatal outcome, however, an abnormal
10 pattern is a poor predictor of fetal acidosis, with a 50
11 percent predictive value. More importantly, most newborns
12 who have an abnormal pattern will be normal.

13 Electronic fetal heart rate monitoring, however,
14 has been associated with an increase in the cesarian
15 delivery rate, and thus the reason for the development of
16 fetal pulse oximetry, in order to try to help determine
17 which fetus is actually compromised and whether or not
18 intervention needs to be done.

19 I would like to state a quote by Benson and
20 colleagues, back in 1968, regarding fetal heart rate
21 auscultation: Naivete and wishful thinking inspired our
22 hope for a simple rule-of-thumb estimate of fetal distress.
23 Obviously, the problem is much too complex for such an early
24 appraisal. And I think this holds true for electronic fetal
25 heart rate monitoring, as well.

1 Moving on with fetal scalp blood pH as a means of
2 assessing fetal well-being, this, too, is a cumbersome
3 technique. It is invasive. It does require multiple
4 determinations, and it has been abandoned by many
5 clinicians. And thus, pulse oximeters have been utilized
6 recently, with past studies and ongoing research.

7 Now, there are two different types of sensors.
8 There is the transmission sensor and the reflectance sensor,
9 and both measure the amount of oxyhemoglobin absorbed and
10 not absorbed.

11 When we look at the transmission pulse oximetry,
12 this is utilized with both the adults and also children.
13 This technique utilizes a light-emitting diode that's placed
14 directly across from the photo detector, and it is now
15 currently used routinely in anesthesia, in critical care,
16 and in newborn nurseries. And, I think there is little
17 question that this has been a significant impact on
18 decreasing morbidity and mortality in these settings.

19 Now, this is just a picture from the article by
20 Lien (phonetic) and Tom Garite in Contemporary OB-GYN,
21 illustrating the transmission sensor, where the device is
22 basically put over a patient's finger. It has been shown to
23 be accurate and has proven efficacy and safety.

24 Now, for obvious reasons, one cannot utilize
25 transmission pulse oximetry, and therefore reflectance pulse

1 oximetry has been developed, where both the light emitters
2 and the photo detectors are on the same surface and they
3 measure the amount of light that is reflected back. Again,
4 from the Contemporary OB-GYN article, this illustrates where
5 both are on the same surface.

6 And, this is again a picture just illustrating
7 where placement of the fetal pulse oximeter, the sensor, is
8 placed through the cervix and fits up against the fetal
9 cheek. This device has been tested in several studies and
10 appears to provide an adequate signal in at least 50 to 60
11 percent of the time.

12 So, in simple terms, the fetal pulse oximetry
13 measures not only the arterial oxygen saturation, but it
14 also measure fetal heart rate and peripheral perfusion.

15 I would like to go over the background or at least
16 the animal studies, especially using the sheep model that
17 has looked at the fetal pulse oximeter. What the animal
18 studies have shown is that the arterial oxygen saturation
19 does seem to correlate with the oxygen that is measured
20 directly in the blood. And, more importantly, the cut-off
21 value of 30 percent, which we are going to discuss again in
22 a few minutes, seems to be less than 30 percent. In other
23 words, in these fetuses, aerobic metabolism is maintained
24 until the arterial oxygen saturation falls below 30 percent,
25 and then metabolic acidosis begins.

1 When we look at the human studies, the first one
2 reported by Dildy and colleagues in 1994, looking at 160
3 women who had a normal, spontaneous vaginal delivery, there
4 does appear to be a wide range of arterial oxygen values,
5 and, as we would expect, there is a decrease in the arterial
6 oxygen saturation during labor.

7 In the first stage of labor, the mean arterial
8 oxygen saturation was 59 percent, and during the second
9 stage of labor, the mean is 53 percent. And if we look at
10 two standard deviations below the mean, the arterial oxygen
11 saturation is 33 percent. Moreover, over 95 percent of all
12 the arterial oxygen saturation values were greater than 30
13 percent.

14 In another study, by Steelbach and Gobel in 1995,
15 of 122 women in labor, they found that the duration of
16 decreased arterial oxygenation is also important in
17 predicting newborn outcomes. In fact, what they reported
18 was that when the arterial oxygen saturation fell below 30
19 percent for more than 10 minutes, the umbilical artery pH
20 was less than 7.20 in more than half of the cases.

21 So this leads us to the clinical efficacy data,
22 and the big question is does it impact upon detecting the
23 compromised fetus and does it decrease the cesarian delivery
24 rate? Recently, Steve Bloom and Ken Levino in Dallas
25 reported their results of using the fetal pulse oximetry for

1 intrapartum outcome in the Obstetrics and Gynecology
2 Journal, in 1999.

3 They utilized this device in 129 fetuses from
4 uncomplicated pregnancies that were 36 weeks gestation or
5 greater, and found that 53 percent had at least one or more
6 episodes of an arterial oxygen saturation of less than 30
7 percent. They found no difference, however, between the
8 high oxygen saturation and the low oxygen saturation group
9 as far as the rate of cesarian delivery, 13 percent versus
10 nine percent, and there was no difference in the umbilical
11 artery pH of less than 7.20, 10 percent versus 9 percent.

12 However, they looked at the duration of an
13 arterial oxygen saturation of less than 30 percent for
14 greater than or equal to 2 minutes, and they found that it
15 was associated with an increase in fetal compromise, i.e.,
16 cesarian delivery for non-reassuring fetal heart rate
17 tracing and umbilical artery pH of less than 7.2, admission
18 to the special care nursery and 5-minute Apgar scores of 3
19 or less.

20 Now, when you look at Nellcor's data, they had 472
21 women in their baseline phase and then 1190 women in their
22 pilot study and randomized clinical trial. And, it is
23 important to note that in enrollment two-thirds of the women
24 had at least one or more risk factors for subsequent
25 cesarian delivery, and one-third of the fetuses had a risk

1 factor.

2 They found that there was no difference in the
3 overall cesarian delivery rate, however, the rate for
4 cesarian delivery for a non-reassuring fetal heart rate
5 tracing was decreased by 50 percent, from 10 percent to 5
6 percent. They also found that there was a decrease by 43
7 percent for non-reassuring fetal heart rate tracing and
8 dystocia.

9 So, the company's data does appear to show that
10 there is a decrease in the cesarian delivery rate for a non-
11 reassuring fetal heart rate tracing.

12 Then, we have to ask the question about the safety
13 data. Although it is difficult to ascertain safety from the
14 available literature, there are no reports of significant
15 morbidity to either the mother or the fetus.

16 This technology does not appear to be associated
17 with an increase in infection, as reported by Johnson and
18 colleagues in 1994. Looking at 112 women with a fetal pulse
19 oximeter compared to matched controls, the infectious
20 morbidity rate was the same.

21 When we look at the company's data, the device,
22 again, does not appear to be associated with significant
23 morbidity, and adverse events appear to be minor.
24 Importantly, there were no maternal deaths and, as reported
25 by Nellcor, maternal adverse effects include fever, mucus

1 membrane disorder, urinary retention, endometrial disorder,
2 postpartum hemorrhage and anemia.

3 As far as fetal adverse effects are concerned,
4 again, they were considered mild but included ecchymosis,
5 accidental injury, jaundice, perinatal disorder and dyspnea.

6 There were 10 adverse events that related to the
7 device, primarily accidental injury, and there were 6
8 neonatal deaths. Again, the investigators thought that
9 there was no causal relationship to these neonatal deaths
10 and the use of the device. Four of these babies had a
11 cardiac congenital anomaly.

12 According to Nellcor, significantly more mothers
13 in the fetal heart rate and the fetal pulse oximetry group
14 in the randomized clinical trial had no adverse events, 70
15 percent versus 66 percent, and there were no adverse
16 maternal events that were considered by the investigators to
17 have a causal relationship with this device.

18 I think the questions and areas of needed research
19 include peer review of the entire set of data, and then we
20 need to ask the question of whether or not this pulse
21 oximetry device has an impact on neonatal morbidity and
22 mortality, and importantly, doesn't have an impact on long-
23 term neurological outcome.

24 Another question we must ask is does the reduction
25 in the cesarian delivery rate for non-reassuring fetal heart

1 tracing -- is that due to other factors, other than the
2 device itself? What is the cost of the technology, and is
3 it only beneficial for high risk patients?

4 The current recommendation would that if the data
5 holds up to peer review, the reduction in cesarian delivery
6 for non-reassuring fetal heart tracing would appear to
7 justify its approval and use at least in the high risk
8 population because it would have a significant impact on the
9 practice of obstetrics, especially considering the high
10 cesarian delivery rate in this country. It might also
11 potentially have an impact on the high rate of litigation
12 cases for babies that have neurologic dysfunction.

13 Thank you very much.

14 CHAIRMAN BLANCO: Dr. Ramin, are you representing
15 ACOG?

16 DR. RAMIN: I am representing ACOG.

17 CHAIRMAN BLANCO: Thank you.

18 DR. RAMIN: That is correct.

19 CHAIRMAN BLANCO: The next speaker that I have on
20 the list is Susan Meikle -- I apologize if I mispronounced
21 your name -- representing the National Institutes of Health.

22 MS. MEIKLE: Good morning. My name is Susan
23 Meikle. I am currently the Acting Program Officer for the
24 Maternal-Fetal Medicine Unit Network for NICHD. I have no
25 conflicts of interest.

1 The Network is a group of thirteen academic
2 institutions, two of which are represented by panel members,
3 who perform protocols. The main objective of the Network
4 was originally to study the etiology of prematurity. This
5 has since evolved over almost fifteen years of evaluation
6 and studies and funding from NICHD to include the areas of
7 low birth weight, medical complications such as asthma, and
8 interventions during labor and delivery.

9 Some of the studies that are currently undergoing
10 are using magnesium to reduce cerebral palsy; multiple dose
11 steroids versus single dose steroids is something we will do
12 in the future; and we do have a cesarian section registry to
13 look at current rates of cesarian section and those
14 outcomes.

15 These studies are chosen by the steering committee
16 which is composed of principal investigators, and at our
17 next steering committee meeting the results of this meeting
18 will be presented to the PIs, and I assume that there will
19 be some discussion about interest in that work and looking
20 into some of the questions that were presented by ACOG.

21 Thank you.

22 CHAIRMAN BLANCO: Thank you very much. There have
23 been a couple of questions for Dr. Ramin. Would you mind
24 going back to the podium? Did you participate in the study
25 at all?

1 DR. RAMIN: I did not.

2 CHAIRMAN BLANCO: Do you have any connection?

3 DR. RAMIN: I have no connection.

4 CHAIRMAN BLANCO: Okay, thank you. I would also
5 like to remind the panel that the data being presented today
6 -- the sponsor and the FDA will be the responsible folks to
7 present the data that we need to assess for the panel
8 deliberations.

9 I do not have anyone else registered to speak.
10 Does anyone from the audience care to make any comments?
11 This is the time to make the public comments.

12 [No response]

13 Then, we will use the time wisely and proceed on
14 with our sponsor presentation.

15 **Sponsor Presentation**

16 **Introduction and Proposed Indication for Use**

17 MS. PAGE: Good morning. My name is Donna Page.
18 I am employed by Mallinckrodt as a Senior Regulatory Affairs
19 Manager for their Perinatal Division, located in Pleasanton,
20 California.

21 It is had been a long time coming and we are
22 really pleased to be here today, to present the N-400 Fetal
23 Oxygen Saturation Monitoring System for the panel's
24 consideration.

25 We have provided the panel with an outline of our

1 program. We have six key speakers following me. They will
2 each introduce themselves fully at the beginning of their
3 presentation. We also have a lot of information to share
4 with you today so we respectfully request that you hold all
5 questions until the speakers have finished their
6 presentations. To facilitate the question and answer
7 period, we have provided the panel with a list of our
8 attendees and their areas of expertise.

9 For the record, the N-400 is a pulse oximetry
10 system. It is designed to be used during labor and delivery
11 to continuously monitor the oxygen saturation of the fetus.
12 It is to be used as an adjunct to standard fetal heart rate
13 monitoring.

14 The system consists of a monitor, a patient module
15 and a sterile sensor. We do have a system here. We are
16 hoping that you will all take the opportunity to take a
17 closer look at it.

18 For the past four years, the N-400 has been the
19 subject of a clinical investigation under IDE G95106. The
20 results of this clinical investigation will be the focus of
21 our presentation today.

22 While under investigation in the U.S., the N-400
23 has been commercially available in the international
24 marketplace. It was introduced to Europe in 1996 with the
25 CE Mark. More recently, it has obtained TGA approval in

1 Australia, CSA approval in Canada, and is the subject of a
2 pending Shonin application in Japan. It has been well
3 received by the obstetrical community. It has also been
4 thoroughly researched and is the subject of over 300
5 publications. Since we introduced the product, we have
6 shipped over 35,000 sensors, and we believe that it has been
7 used in a similar number of labors.

8 Based on the results of Mallinckrodt's clinical
9 investigation, and supported by its history of use in other
10 markets, Mallinckrodt is requesting marketing approval for
11 the N-400 under PMA P990053, with labeling that we believe
12 will be supported by the data and information presented here
13 today.

14 The next slide shows the components of our
15 indications for use statement. The statement is equivalent
16 to the statement that was in your packet, however, it has
17 been reformatted to more clearly identify the essential
18 elements. The essential elements are the Nellcor N-400
19 system is intended for use as an adjunct to fetal heart rate
20 monitoring. It is not intended to replace conventional
21 fetal heart rate monitoring during labor.

22 The patient population for which the N-400 system
23 is intended consists of term infants in active labor, with
24 ruptured membranes, with a non-reassuring fetal heart
25 pattern. The purpose of the N-400 system is to improve the

1 physician's ability to assess the fetal status. The N-400
2 directly measures the fetal oxygen saturation. This permits
3 the safe continuation of labor during periods of non-
4 reassuring fetal heart and reassuring FSpO2, reducing the
5 rate of C-sections performed for the indication of non-
6 reassuring fetal status without causing injury to the mother
7 or fetus.

8 Finally, the addition of the N-400 to conventional
9 fetal heart rate monitoring improves the sensitivity and
10 specificity for matching the delivery indication to
11 immediate neonatal condition.

12 This concludes my introductory remarks. I would
13 now like to turn the podium over to Dr. David Swedlow. Dr.
14 Swedlow will be providing you with more information on the
15 N-400 technology and discussing the rationale for the 30
16 percent critical threshold.

17 **Technology and Critical Threshold**

18 DR. SWEDLOW: Good morning. My name is David
19 Swedlow. I am a pediatrician and an anesthesiologist, and I
20 specialize in critical care medicine, until about 1987, when
21 I left academic medicine to join Nellcor as a Senior Vice
22 President of Medical Affairs and Technology Development. I
23 was there for nine years, working on this project and
24 others, until I retired two years ago, and I am currently a
25 paid technical adviser to the company for this project.

1 The journey for me to today, to fetal pulse
2 oximetry, actually begins 28 years ago, when I was a
3 pediatric house officer at Johns Hopkins. I remember being
4 called in the middle of the night from the neonatal nursery
5 to run down four flights of stairs, run over two buildings,
6 and run up five flights of stairs, which is not easy for me,
7 to attend to the delivery as a pediatrician of a patient
8 being sectioned for fetal distress. I arrived to see
9 anxious, concerned parents and anxious obstetrical house
10 staff, only to be rewarded by a totally normal child.

11 I thought to myself at that time that there must
12 be a better way of determining who is in trouble and who is
13 not in trouble, but at the time there really wasn't any.
14 Six years later, as an anesthesia resident at the University
15 of Pennsylvania, I was again called in the middle of the
16 night but this time to administer anesthesia for a crash C-
17 section for an emergency delivery. Once again, the baby
18 came out fine but at that time there was no technology
19 available to better define and improve fetal assessment.

20 In 1982, ten years after the beginning of this
21 story, I saw for the first time a pulse oximeter. It was a
22 device that wrapped a little band-aid around the finger,
23 shone light through the finger and measured the color of
24 blood and, in so doing, measured the oxygen saturation of
25 the patient. At that point in time, I was absolutely

1 convinced that it would change the practice of medicine, and
2 I thought to myself when I first saw it that this is what we
3 need for the babies, but it has been a long time coming.

4 Pulse oximetry addresses a very serious medical
5 problem, and that is uncertain patient oxygenation. It has
6 been applied in surgery, anesthesia and in the ICU unit. It
7 provides an objective, continuous, and direct measurement of
8 oxygen in the adult, child and neonate and, in so doing, I
9 certainly feel and I imagine most other physicians would
10 feel as well, it has truly transformed the practice of
11 medicine in that population. However, until now
12 obstetricians had no way of measuring fetal oxygenation
13 directly. They were forced to rely on indirect measures,
14 such as the fetal heart rate.

15 We believe that the N-400, fetal oxygen saturation
16 system, brings for the first time an objective, continuous
17 and direct measurement -- and that is the important thing
18 from my point of view, that it is objective and direct
19 measurement of fetal oxygen to obstetrics.

20 What we needed to do, starting back in 1990, or
21 so, was to extend the conventional oximetry technology to
22 the laboring fetus. We had several special fetal issues to
23 deal with. Most importantly, was that the fetal environment
24 was extraordinarily difficult for that technology. We had a
25 wet, unseen patient, with no accessible appendages -- no

1 hands, feet, fingers, ears or toes that we could wrap a
2 band-aid around.

3 So, we had to develop a device that could be
4 inserted through the cervical os and come to lie alongside
5 the fetal face. This is such a sensor. The sensor is
6 inserted during a vaginal exam, is inserted gently and
7 allowed to advance until it comes to lie alongside the fetal
8 face or temple. Most kids don't have beards so that is not
9 a problem.

10 [Laughter]

11 It shines light into the fetal skin. It measures
12 the color of the reflected light coming from the blood
13 cells, and the monitor itself computes the saturation and
14 displays it in real time. It is non-invasive to the fetus.
15 It is, as you will hear later, quite easy to insert and
16 comfortable for the mother. So, for the first time we are
17 able to provide the obstetrician with a direct measurement.

18 The next problem we had was that we had to
19 discover the threshold for clinical reassurance. That is,
20 we had to define or discover the value above which we could,
21 with assurance, say this patient or this child -- this fetus
22 -- is adequately oxygenated, and below which there might be
23 a risk of development of metabolic acidosis due to hypoxia.

24 We went about this in a methodological approach.
25 We went to the literature and found studies from Brian

1 Richardson that indicated that 30 percent was a reasonable
2 target. We did prospective animal studies which defined 30
3 percent as the critical threshold for the development of
4 acidosis. I will explain that in a minute.

5 Then we did human studies, looking at the
6 relationship between what is currently the gold standard for
7 fetal assessment, scalp pH, and the value of saturation.
8 That too indicated 30 percent.

9 This slide is from a prospective study of
10 instrumented, near-term fetal lambs in non-laboring,
11 unanesthetized ewes. The maternal ewe was exposed to graded
12 hypoxia, and also the fetus was exposed to graded ischemia
13 with common iliac artery occlusion.

14 On the vertical axis, on the left, you see the
15 fetal saturation value from a catheter in the fetal lamb.
16 In the vertical axis on the right, you see the value of the
17 base excess from a blood gas drawn from that same catheter.
18 At time zero the mother was made hypoxic; 30 or 40 minutes
19 later the fetal oxygen saturation had fallen to a level
20 below 30 percent, at which the lamb began to accumulate
21 excess acid. As long as the saturation remained below 30
22 percent the acid continued to accumulate. When the
23 saturation was allowed to rise above 30 percent, the lamb
24 recovered the acid base status, and we found a critical
25 threshold at 30 percent. That is, no animal above 30

1 percent accumulated acid; all animals below 30 percent
2 accumulated acid at varying rates.

3 We repeated that -- this wasn't us; this was a
4 perinatal research group. They repeated this same study
5 with graded ischemia and found the same result.

6 We then went to Germany, where we did a multi-
7 center clinical trial, looking at the relationship between
8 scalp pH and saturation. In a series of 46 patients with
9 non-reassuring fetal heart, they drew blood gases and
10 compared the fetal scalp pH, shown on the vertical axis, to
11 the fetal saturation, shown on the Y axis.

12 Using an ROC analysis, we determined that the
13 critical threshold here seemed to be between 30 percent and
14 40 percent. We chose 30 percent for the clinical study.
15 The reasons were quite simple. Animal studies had
16 demonstrated that the critical threshold was, indeed, 30
17 percent, not 40 percent. We wanted clinically to reduce the
18 number of false positives and, thereby, reduce the number of
19 unnecessary interventions. We wanted to, therefore,
20 maximize specificity and that way indicate choosing a
21 threshold value on the left-hand side. In this case
22 maximizing specificity would suggest a value of about 30
23 percent.

24 We also realized that this scalp pH versus
25 saturation is inherently a conservative approach because 7.2

1 as a definition of fetal acidosis from scalp pH is quite
2 conservative. Some people feel that the real threshold for
3 concern ought to be lower and that, too, would suggest a
4 lower threshold.

5 Finally, to be honest, there was a practical
6 issue. We needed a value that had a line printed on the
7 uterine activity charts, which is where we were going to
8 display this thing, so that we could unequivocally in the
9 protocol say you did reach the threshold or you did not
10 reach the threshold. From a clinical use point of view,
11 that turns out actually to be a pretty important issue.

12 For all those reasons, we used 30 percent as the
13 threshold in the randomized controlled trial, and I am going
14 to turn the podium over now to Dr. Tom Garite, who is not
15 only the principal investigator of this study but actually
16 the architect of the study.

17 **Pivotal Study Design and Results**

18 DR. GARITE:

19 Good morning, Dr. Blanco, distinguished members of
20 the panel and distinguished interested members of the
21 audience.

22 My name is Tom Garite. I am an obstetrician and
23 specialist in maternal-fetal medicine. I am a professor and
24 chairman of the Department of Obstetrics and Gynecology at
25 the University of California, Irvine and, as David

1 mentioned, the principal investigator of this study.

2 I have no financial interest in Mallinckrodt, nor
3 have I received any compensation from the company other than
4 the support I received for the research to do the study and
5 the compensation of my travel expenses to attend this
6 meeting.

7 The design of the study that we will present today
8 was as a result, as you might imagine, of extensive
9 negotiation among the investigators and the sponsors to get
10 a group of investigators to agree on a uniform protocol for
11 interpretation and intervention for electronic fetal heart
12 rate monitoring, as you all know, was no small feat. But we
13 did eventually come to a uniform agreement which was
14 presented and agreed upon by the Food and Drug
15 Administration and, as previously mentioned, again presented
16 and discussed a great deal at the advisory panel meeting of
17 this agency in July of 1996.

18 The overall goal for this technology is to improve
19 the accuracy and reliability of intrapartum fetal
20 assessment. Perhaps the ideal study to do would be to
21 design a study wherein we could demonstrate that this device
22 improved fetal outcome. However, we concluded that this
23 goal was unrealistic in that fetal damage due to intrapartum
24 asphyxia is a rare event, and that the numbers required for
25 such a study would be unapproachable.

1 Alternatively, we considered two other approaches.
2 One would be a comparative study between fetal pulse
3 oximetry and scalp pH. We rejected this for two reasons.
4 First, it was just a physiologic surrogate which really did
5 not test how the device performed in actual practice.
6 Second, in reality fetal scalp pH is not something commonly
7 done in the United States.

8 We ultimately chose to evaluate the endpoint of
9 reduction of cesarian section for non-reassuring fetal
10 status. This approach had the advantage of testing how the
11 implementation of this device affected actual clinical
12 behavior, and was a direct reflection of the hoped for
13 improvement in the accuracy of fetal assessment by this new
14 technology.

15 The study, then, was designed to test improvement
16 in fetal assessment by measuring impact on physician
17 behavior and neonatal outcome. Therefore, the study had
18 three goals. First was to test to see if we could use fetal
19 pulse oximetry, together with conventional fetal heart rate
20 monitoring, to reduce the rate of cesarian section for the
21 specific indication of non-reassuring fetal status. Next,
22 was to be sure that continuing labor in the face of a non-
23 reassuring heart rate pattern but a reassuring oxygen
24 saturation was, indeed, safe. Finally, to demonstrate that
25 the device itself was safe for mother and child.

1 The study we conducted was a prospective,
2 randomized, controlled, unblinded interventional clinical
3 trial comparing fetal assessment with electronic fetal heart
4 rate monitoring alone versus electronic fetal heart rate
5 monitoring backed up by fetal oxygen saturation monitoring.
6 In both arms of the trial a study nurse was present from
7 randomization until delivery to ensure protocol compliance
8 and to enhance data collection.

9 To monitor study compliance and evaluate the
10 accuracy of fetal heart rate interpretation and
11 intervention, 100 percent of the cases and all of the fetal
12 heart rate tracings were reviewed by an independent
13 reviewer, Dr. Mike Nageotte, who did not participate as a
14 clinical investigator. For the final analysis we used an
15 intent-to-treat analysis with all patients included.

16 The study was conducted in three phases. The
17 baseline phase was observational, with no pulse oximetry
18 used. The purpose of this was to get an estimate of
19 baseline clinical practice and to screen for sites willing
20 and able to perform the study. As you can see, 472 patients
21 were evaluated.

22 The second was the pilot phase. Here,
23 randomization and data collection was practiced, as was the
24 placement and use of the sensor, and for the sites to
25 familiarize themselves with the use of the clinical

1 management protocol. Each site used a minimum of 15
2 patients and, as you can see, we had a total of 180
3 patients.

4 Finally, the randomized, controlled trial was
5 begun. In this trial we approached a total of 4545
6 patients, of whom 2996 were consented. Ultimately, 1010
7 patients met entry criteria, were enrolled and randomized,
8 502 in the control and 508 in the test arms.

9 This slide shows all the study sites and their
10 principal investigators. As you can see, these sties are
11 geographically dispersed and well mixed between community
12 and university hospitals, and those with and without
13 teaching services -- nine sites in all.

14 Eligible patients included those in active labor,
15 with a cephalic presentation at or below a minus 2 station,
16 with ruptured membranes, and at or beyond 36 weeks of
17 gestation, and also a singleton.

18 Patients were excluded if they were entered in any
19 other intrapartum research study; if they were planning to
20 have an elective cesarian section; had a placenta previa;
21 any other need for immediate delivery, or an infection that
22 precluded internal monitoring.

23 In general, patients were consented on admission
24 but they were only actually enrolled if they developed a
25 specifically defined abnormal heart rate tracing. These

1 fetal heart rate entry criteria included mild to moderate
2 non-reassuring fetal heart rate patterns defined to allow
3 early enrollment of patients who were at risk for developing
4 more severe and concerning fetal heart rate patterns later.

5 In general, patient management was driven by the
6 pattern on the electronic fetal monitor. These patterns
7 were categorized into one of three classes. All of these
8 classes are shown in detail in the clinical report.
9 Patients with a class 3 or ominous pattern were delivered
10 immediately in both groups. Patients with class 2 or non-
11 reassuring pattern, such as persistent, late or non-
12 reassuring variable decelerations. In both groups non-
13 operative corrective measures were applied at the discretion
14 of the managing physician and nurse to try to correct the
15 fetal heart rate pattern. If these measures were
16 unsuccessful and a non-reassuring heart rate persisted, then
17 management differed between the control and test groups. In
18 the control groups the physicians used either accelerations,
19 spontaneous or elicited, or scalp pH to rule out acidosis.
20 If acidosis could not be ruled out and the pattern
21 persisted, the patient was delivered.

22 In the pulse oximetry group the clinician used
23 oxygen saturation. If the oxygen saturation was above 30
24 percent at any time between contractions the physician was
25 reassured and labor was continued. If not, if it remained

1 below 30 percent, then the patient was managed or, if a
2 signal could not be obtained, then the patient was managed
3 as if she had a heart rate monitor alone, i.e.,
4 accelerations or scalp pH was used, and if we could not be
5 reassured with those, the patient was delivered.

6 The remainder of the patients with class 1
7 reassuring heart rate patterns were allowed continuance of
8 labor.

9 Evidence to support that this protocol was clear
10 and reasonable is supported by the fact that only four
11 patients in each group had a significant protocol violation
12 as determined by the independent reviewer.

13 The next four slides provide the pre-randomization
14 clinical characteristics in the control and test groups. In
15 this slide you can see that there no difference in maternal
16 age, in racial distribution, in parity, in source of
17 funding, or in the frequency of previous cesarian section.

18 There were minimal differences in maternal and
19 fetal risk factors between the groups. I apologize for
20 going through these kind of quickly. The labor
21 characteristics were also quite similar in the two groups,
22 with the exception that the frequency of labor induction and
23 its associated pre-induction prostaglandin ripening were
24 more common in the pulse oximeter group.

25 Critical to establishing appropriate randomization

1 for this study is the observation that the frequency of the
2 defined entry fetal heart rate patterns was virtually
3 identical between the two groups. You can see that the
4 specific type of heart rate patterns that allowed enrollment
5 and randomization were virtually identical. Therefore, we
6 concluded that the two groups were well matched and that
7 there was no evidence of meaningful selection bias between
8 the two groups.

9 We also wanted to assure ourselves that there was
10 no evidence of investigator bias or clinical management
11 bias. We addressed the issue of investigator bias in two
12 ways. First, management of two groups was reviewed by the
13 independent reviewer for protocol violation, as I previously
14 mentioned, with four significant violations in each group.
15 Including or excluding these eight patients had no impact on
16 any conclusion in the analysis.

17 Secondly, we examined the pattern and frequency of
18 labor interventions and fetal evaluations between the two
19 groups. Remember that both groups with a non-reassuring
20 pattern were allowed standard labor interventions before
21 proceeding to trying to reassure yourself, and you can see
22 that things like oxygen administration, repositioning, etc.
23 were quite similar between the groups, with the only
24 difference being a slightly increased frequency of need for
25 correction of hypotension in the pulse oximeter group.

1 Similarly, though not shown on this slide, and it
2 is available in your clinical report, the frequency of both
3 the most severely abnormal fetal heart rate pattern and the
4 specific types of pattern, which led to the definition of a
5 non-reassuring or ominous pattern and led to intervention,
6 were also virtually identical between the two groups. So,
7 there was an equal frequency of the worst fetal heart rate
8 pattern in the groups. We concluded, therefore, that there
9 was no evidence of investigator or management bias between
10 the two groups.

11 The primary results of the study are shown on this
12 summary slide. We found that the overall distribution of
13 the route of delivery was similar between the two groups.
14 As you can see, spontaneous vaginal delivery -- similar;
15 assisted vaginal delivery -- similar; and the overall
16 cesarian section rates of 26 and 29 percent were not
17 statistically difficult. However, the indications for
18 cesarian section were substantially altered. The rate of
19 cesarian section performed for non-reassuring fetal status
20 was reduced from 10.2 percent to 4.5 percent, more than a 50
21 percent reduction. This was the outcome for which the study
22 was designed.

23 A logistic regression analysis demonstrated that a
24 strong independent effect of pulse oximetry test group
25 assignment resulted in this decreased risk for non-

1 reassuring fetal status. Unexpectedly and surprisingly, the
2 C-section rate for dystocia increased to offset the
3 reduction for C-section for non-reassuring fetal status.

4 Before discussing this dystocia issue, I would
5 like to first discuss the reasons for the reduction for non-
6 reassuring fetal status and the safety issues for mother and
7 baby.

8 We were able to demonstrate that continuation of
9 labor during periods of non-reassuring fetal heart rate but
10 with reassuring fetal pulse oximetry was safe for the baby.
11 This slide compares the immediate neonatal condition between
12 the two groups, and shows that the reduction of C-section
13 performed for non-reassuring fetal status was achieved
14 without increase in adverse outcome.

15 There were five deaths. Two in each group were
16 for complex congenital cardiac malformations. The third
17 death in the test group is detailed in the clinical report.
18 This baby's demise appeared to result from a delayed
19 appreciation of attention pneumothorax, recognized about one
20 or two hours post-natally, which resulted in severe
21 metabolic acidosis. The baby had a normal 5-minute Apgar
22 and borderline cord pH and went to the newborn nursery.

23 There were no statistically significant
24 differences between the groups in the frequency of low Apgar
25 scores, low cord pH, or base excess in the cord, or in need

1 for resuscitation. However, note that all of the neonates
2 with extremely low cord arterial base excess are in the
3 fetal monitor group, as are two-thirds of the low five-
4 minute Apgar scores.

5 There were no significant differences in maternal
6 outcome between the two groups either. Specifically, there
7 is no difference in intrapartum fever or in postpartum fever
8 or postpartum endometritis.

9 Thus, we conclude that both the use of the fetal
10 oximetry sensor itself, as well as the management protocol
11 tested in the study, is safe for mother and child.

12 This slide is very complex but very important,
13 this and the next slide, and let me try to explain it as
14 best I can. It is important to take time to describe we
15 were able to reduce the cesarian section for non-reassuring
16 fetal status without causing increased injury to the baby.

17 In this study, reduction of cesarian section for
18 non-reassuring fetal status is a behavioral surrogate for
19 improved accuracy of fetal assessment. Another method of
20 evaluating accuracy of fetal assessment is to examine the
21 degree of agreement between the physician's choice for
22 operative intervention for non-reassuring fetal status and
23 the actual immediate newborn condition. If we have provided
24 a method of improved assessment for the clinician, then we
25 should see an improvement in matching between the decision

1 to proceed with or avoid operative intervention and the
2 actual immediate neonatal condition.

3 As previously pointed out, too often we perform an
4 urgent cesarian section for concern over fetal oxygenation
5 only to deliver an extremely vigorous and well oxygenated
6 baby. Thus, the real hope for pulse oximetry is that we
7 will be able to avoid such unnecessary operative
8 intervention without missing babies who would really benefit
9 from immediate delivery.

10 The most direct way to examine this agreement
11 between behavior and immediate neonatal condition is to
12 construct a series of paired 2 X 2 tables, one for the
13 control group and one for the test group, comparing
14 operative delivery or continued expectant management and
15 actual depression versus no depression.

16 As you can see, as you saw on the previous slide,
17 the results are similar for immediate adverse condition.
18 However, this slide shows that numerous neonatal descriptors
19 of immediate condition, and the number of neonates in each
20 group with that condition, and the sensitivity and
21 specificity for fetal heart rate monitoring versus pulse
22 oximetry for that device.

23 In each case shown here, and with nearly any
24 threshold for pH or base excess chosen, the sensitivity and
25 specificity of the decision for operative intervention

1 versus neonatal condition is better in the test group than
2 in the control group. Here you see sensitivity for a low pH
3 compared with fetal heart rate monitoring versus oximetry
4 statistically improved, and specificity for fetal heart rate
5 monitoring versus oximetry statistically improved in the
6 oximetry group. The final column gives the significant
7 value of the difference in sensitivity and specificity given
8 by the Mantel-Haenzel test for homogeneity of odds ratios.

9 I want to point out that although the numerical
10 differences in the specificity may not appear to be as
11 impressive as the statistical significance, it is important
12 to point out that, for example, for pH the 78 percent
13 specificity represents a false-positive rate of 22 percent
14 for fetal monitoring versus a false-positive rate of only 14
15 percent for pulse oximetry. Thus, the numerical difference
16 and the statistical significance are both not only
17 statistically significant but clinical relevant.

18 We, therefore, draw the following conclusions:
19 This is a large, well executed, multi-center randomized,
20 controlled trial, with a high degree of study compliance.
21 The addition of fetal oximetry monitoring to conventional
22 fetal heart rate monitoring improves the accuracy of fetal
23 assessment and permits safe reduction in the number of
24 cesarean sections performed for the specific indication of
25 non-reassuring fetal status without causing injury to mother

1 or child.

2 On a personal note, I just want to state how
3 gratifying it has been to be involved in a large, well-
4 conducted, randomized, controlled trial of a diagnostic
5 device, using an intervention protocol, before the device
6 was introduced into clinical practice.

7 Thank you for the opportunity to present this
8 study. Now we need to examine and understand the reason why
9 the overall number of cesarean sections because of the
10 increase of cesarean sections for dystocia in the study
11 group occurred, and my colleague, Dr. Rich Porreco, will
12 present this study analysis.

13 **Cesarean Section for Dystocia and Clinical Utility**

14 DR. PORRECO: Thank you, Dr. Garite. My name is
15 Rich Porreco. I am a perinatologist in Denver, and I was
16 the principal investigator at Presbyterian St. Luke's
17 Medical Center for this trial. I have no financial
18 relationship to Mallinckrodt, although I am told they will
19 cover my expenses to come here, to Gaithersburg.

20 You have seen this slide before. This is the
21 primary outcome slide that Dr. Garite showed you, showing
22 that the overall cesarean birth rate was identical between
23 the control and test groups, largely as a result of the fact
24 that dystocia was increased among the test patients compared
25 to the control group.

1 This study, as the panel may be reminded, was not
2 designed to investigate the incidence of dystocia and,
3 therefore, I think our retrospective analysis of these
4 observations have to be qualified. However, the observation
5 was compelling enough that I think all the investigators
6 felt that we should carefully review this data
7 retrospectively and try to come to some reasonable
8 conclusions.

9 We considered four possibilities for these
10 observations: One, unbalanced patient risk factors, despite
11 the randomized methodology, might account for it.
12 Mislabeling of dystocia due to bias of the investigators,
13 that is, investigators rooting for the device might mislabel
14 their true indication for cesarean section as dystocia
15 rather than fetal distress. Thirdly, did the sensor use
16 itself somehow manage to slow labor? Finally, was the
17 finding of reassuring fetal oximetry information, permitting
18 the continuation and the natural evolution of labor
19 unmasking some underlying risk of dystocia?

20 Let me try to flesh out some of these
21 considerations. As Dr. Garite has pointed out, the study
22 entry characteristics were well balanced. The increase in
23 inductions and prostaglandin use among study patients washed
24 out with a logistic regression analysis. Also, the same
25 analysis showed that C-section for dystocia was an

1 independent effect of being assigned to the test group.

2 On this slide, we considered secondly whether
3 there was mislabeling of dystocia due to investigator bias.
4 In a retrospective, blinded analysis of partograms defined
5 dystocia, as you see noted here -- arrest of dilatation in
6 the active phase of labor for more than three hours, arrest
7 of descent in the second stage for more than two hours, or
8 failed induction, that is, oxytocin administration for more
9 than twelve hours in the presence of ruptured membranes was
10 looked at.

11 On this slide you see that retrospectively
12 assigned this definition of dystocia, both 90 percent of the
13 control patients and 90 percent of the test patients who
14 were sectioned for dystocia truly had dystocia, and those
15 that did not were represented equally in either group. So,
16 that is evidence of no mislabeling by the investigators.

17 Additionally, there was this concern that the
18 investigators, again rooting for the device, might let
19 fetuses with true distress languish somewhat and section
20 them belatedly for dystocia. Looking at patients in the
21 test group who were sectioned for dystocia, you can see that
22 occurrence of depressed fetuses was uncommonly seen as
23 opposed to those in the test group sections for non-
24 reassuring fetal status. Again, corroborative evidence that
25 there was no mislabeling due to investigator bias.

1 On this slide, the third concept was whether the
2 use of the sensor itself slowed labor. In this Kaplan-Meier
3 analysis of time to delivery, you can see that the curves
4 are superimposable. Indeed, the test group had slightly
5 faster labors and, if there was any slowing by the device,
6 it should occur irrespective of the mode of delivery. These
7 Kaplan-Meier curves show no slowing irrespective of the mode
8 of delivery.

9 Lastly then, the issue of whether improved
10 knowledge of fetal status would not only explain the
11 decreased occurrence of cesarean section for fetal distress
12 for non-reassuring fetal status, but also explain the
13 increased occurrence of cesarean birth for dystocia. You
14 have already seen that patients who were delivered for
15 dystocia actually did have dystocia, and in patients with
16 non-reassuring heart rate traces, reassuring oxygen
17 information permitted the natural evolution of labor,
18 continuation of labor, potentially unmasking underlying risk
19 of dystocia. And, if this is the case, we should see an
20 increased rate of these class 2 fetal heart rate patterns,
21 the ones that require some intervention on behalf of the
22 investigators.

23 This is exactly what we have seen. You can see
24 that in patients who were sectioned for dystocia, in the
25 test group, the class 2 non-reassuring patterns, especially

1 variable decelerations, are largely segregated here, in this
2 group, much more so than the fetal heart rate control group
3 alone. This, in balance, was seen only among patients who
4 were sectioned for dystocia, not among other cesarean birth
5 indications.

6 So in summary, we would suggest to you that our
7 retrospective analysis of dystocia and non-reassuring heart
8 rate traces -- that patients in the control group with non-
9 reassuring heart rate traces are delivered by cesarean
10 section for presumed fetal distress, syphoned off, if you
11 will, and sometimes inappropriately, and delivered by
12 cesarean section for that indication, whereas, in the test
13 group the reassuring information of fetal oximetry permits
14 labor to continue, to evolve naturally, unmasking, if you
15 will, their underlying risk for dystocia. Therefore, the
16 non-reassuring fetal heart rate that formulates our
17 inclusion criteria, especially variable decelerations, may
18 be an inherent marker for dystocia.

19 Now, this is not a new finding. The early studies
20 of fetal heart rate traces, from the '70s by Havercam,
21 showed an increased occurrence of cesarean birth for fetal
22 distress and for dystocia. And, we know that positions of
23 the fetal vertex, especially persistent occipitoposterior
24 positions, are associated also with severe variable
25 deceleration and also associate strongly with cesarean birth

1 for dystocia.

2 So in summary, abnormal fetal heart rate patterns
3 are common, and even with expert interpretation and lots of
4 experience we find some built-in ambiguity. This ambiguity
5 causes a lot of "medicalization" in a birth environment.
6 Nurses and physicians look at the monitor strip, wrinkle
7 their brow, cause a lot of anxiety in the birth mother, they
8 turn her on her side, turn the IV up, put oxygen on, and the
9 whole event becomes "medicalized." I think the addition of
10 oxygen information can objectively and unambiguously tell
11 that family that their fetus is well, and remove that
12 "medicalization" from the birth scene.

13 Finally, intervention for the right indication at
14 the right time is extremely valuable to these families. A
15 calm and confident cesarean birth, when it is done in a
16 program methodically for dystocia, is very difficult from
17 one that may be done urgently, potentially under general
18 anesthesia, with concern by the family about the well-being
19 of their child, and we should not underestimate the value of
20 being able to reassure these families that their fetus is
21 well even if ultimately their underlying risk for dystocia
22 is unmasked three or four hours later down the line.

23 Additional comments about the clinical utility of
24 this device are going to be presented to you by Dr. Frank
25 Boehm, and I will turn it over to Dr. Boehm.

1 **Impact on Clinical Practice**

2 DR. BOEHM: Good morning. My name is Frank Boehm.
3 I am a Professor of Obstetrics and Gynecology at Vanderbilt
4 Medical Center, and Director of Maternal-Fetal Medicine at
5 Vanderbilt Hospital.

6 For the record, I have no financial interest in
7 Mallinckrodt Company. Neither myself nor any of my
8 immediate family own stock in the company. In addition, I
9 have received no financial remuneration for this particular
10 randomized study, however, my nurse's salary was paid for
11 Mallinckrodt during the study period. I will receive travel
12 expenses for this trip to make this presentation to the FDA
13 panel.

14 I believe that the clinical significance of this
15 new technology is in its ability to aid the healthcare
16 provider in caring for laboring patients whose electronic
17 fetal monitor tracing reveals non-reassuring fetal status.
18 It will allow physicians to do the right thing for the right
19 reason at the right time. It will allow us to reassure not
20 only ourselves but our patients as well that while the fetal
21 heart rate monitor may be somewhat confusing or not
22 specific, the Nellcor pulse oximeter sensor indicates that
23 the fetus is not hypoxic. A reassuring fetal oxygen
24 saturation is reassuring not only to the physician but also
25 to the patient. It acts as an adjunct to electronic fetal

1 monitoring which may be somewhat confusing or non-specific,
2 and fetal oxygen saturation of 30 percent or greater
3 indicates that the fetus is not hypoxic.

4 Importantly, this device will help physicians and
5 nurses in their continued endeavor to have women empower
6 themselves in decision-making processes during labor. When
7 physicians are responding only to electronic fetal monitor
8 tracings, they take away some of the options that women have
9 in making decisions. With the Nellcor device patients will
10 be better able to more objectively understand the issues and
11 risks to their fetus.

12 In addition, because of the ability to be more
13 sensitive and specific, the Nellcor device will allow
14 physicians to approach the need for cesarean section in a
15 more intelligent fashion. In the past we, physicians,
16 approached cesarean section in either an urgent manner or an
17 elective manner. This increases some of the risks to
18 patients during the surgical procedure. By being able to
19 reassure ourselves that surgery can be done in a more
20 relaxed and planned manner, a non-emergent procedure can,
21 therefore, reduce not only anxiety but, more importantly,
22 can reduce risks to patients undergoing cesarean section.

23 Doctors and nurses are under tremendous pressures
24 in taking care of the laboring patient. One of these
25 pressures is to properly interpret electronic fetal monitor

1 tracings so as to make appropriate decisions in how to
2 manage labor, as well as to explain to patients and their
3 families the status of the unborn child. However, because
4 of many ambiguities of electronic fetal monitor tracings, we
5 need a mechanism to allow for objective data that will
6 remove these ambiguities in electronic fetal monitor
7 interpretation. The fetal pulse oximeter is a device which
8 will give healthcare providers and patients that needed
9 objective data.

10 In summary, I believe the Nellcor fetal pulse
11 oximeter will have a significant impact on the field of
12 obstetrics. Particularly, the approximately 30 percent of
13 patients with non-reassuring fetal heart rate patterns will
14 have more objective data that will allow clinicians to
15 ascertain whether, in fact, the fetus needs to be delivered
16 quickly or whether the patient can continue in the laboring
17 process. This will have a significant impact on the
18 cesarean section rate for fetal stress or distress, and will
19 lead to more appropriate decision-making by all involved.

20 Now, to highlight these comments I would like to
21 present one case from our institution at Vanderbilt. While
22 each investigator had their own case, this is ours that I
23 would like to share with the group. This is a private
24 patient of mine, a 42-year old, gravida 2, para 1 at 38 and
25 6, 7 weeks. The patient underwent a low transverse cesarean

1 section for failure to progress at approximately 8 cm, with
2 an occiput transverse position. The baby weighed 8 lbs. 9
3 oz. This was a delivery that I attended.

4 In the patient's second delivery the estimated
5 fetal weight at term was 7.5 lbs. and the patient strongly
6 desired a vaginal birth after cesarean section. Prenatally,
7 the patient underwent an amnio for normal carrier type and
8 was treated with Synthroid for hypothyroidism.

9 She was admitted to the hospital for an induction.
10 The patient was 2 cm, 70 percent, minus 2 station, vertex.
11 The fetal heart rate was 140 and reactive. There were no
12 contractions. This was the fetal heart rate upon admission
13 to the hospital. As you can see, the fetal heart rate is
14 stable and reactive with these accelerations.

15 At 9:45, rupture of the membranes was performed.
16 Clear fluid was noted, and the patient was now 3 cm, 50 pc
17 and minus 1 station. At 10:45, an epidural was placed
18 because of painful contractions. The patient now was 3-4
19 cm, more effaced and now engaged. This was the fetal heart
20 rate at approximately this time, 10:45 in the morning, and
21 you can see again the acceleration, moderate variability,
22 and a normal fetal heart rate tracing, with uterine
23 contractions occurring down here.

24 At one o'clock, the patient dropped her blood
25 pressure slightly, had an increase in intravenous fluids and

1 the blood pressure improved, and she was making progress.
2 This was the fetal heart rate tracing at the time, again a
3 reassuring fetal heart rate pattern.

4 At 5:30, however, things changed. The fetal heart
5 rate baseline went to 175 to 180. There were persistent
6 late decelerations according to the nurses. The blood
7 pressure was normal. The patient once again hit that magic
8 8 cm, full effacement, and was now plus 1 and left occiput
9 transverse. A fetal scalp electrode was placed, and because
10 of a mildly elevated temperature, ampicillin was started;
11 oxygen was given and Pitocin was continued, and the patient
12 was explained the process of randomization and the
13 availability of the Nellcor sensor device.

14 This was the tracing at 5:30 when the discussions
15 began and, as you can see, the heart rate is 180. There is
16 minimal variability and there are these decelerations that
17 were interpreted as late decelerations -- certainly, a quite
18 non-reassuring fetal heart rate pattern.

19 At 6:30, the patient consented, was randomized
20 fortunately to the sensor group, and her heart rate showed
21 persistent late decelerations with an advancing tachycardia
22 now, as high as 195 and, interestingly, the pulse oximeter
23 reassured us with a 35 percent reading.

24 This was the initial heart rate tracing, and you
25 can see the heart rate here has not really changed. There

1 is still considerable tachycardia, persistent decelerations
2 after the peak of the contractions and minimal variability.
3 At this particular point, the sensor is not reading
4 appropriately here and the nurse has written down the pulse
5 oximetry of 35 percent and 33 percent. We will get to the
6 appropriate reading in a second.

7 This is immediately following, and you can see
8 this tracing of a patient with a previous cesarean section
9 still at 8 cm -- it was very intriguing to consider a
10 cesarean section but because the pulse oximeter reading was
11 in the 33 percent to 41 percent we continued the labor.

12 At seven o'clock, you can see that the heart rate
13 was still up, minimal beat to beat, persistent lates,
14 however, the pulse oximetry was in the reassuring range and
15 the patient was now 9 cm and plus 1, with a temperature of
16 100.7.

17 At this point, the pulse oximeter is working
18 appropriately to register onto the paper, and you can see
19 that it is above 30 percent during this time. There are
20 these continued decelerations. Again, the pulse oximeter
21 reading is down here. I believe there is one point where it
22 drops for a few minutes below 30 but the rest of the time it
23 is above 30. The patient is allowed to continue her labor.

24 At 7:40, the patient was now complete, plus 2 and
25 pushing, and again the pulse oximeter was in the 32 percent

1 to 40 percent range. This is that particular period. You
2 can see the pulse oximeter is well above 30. The patient is
3 pushing, and this is the fetal heart rate that we are
4 obtaining with the fetal scalp electrode in place; minimal
5 variability and decelerations.

6 Finally, at the end the scalp electrode is removed
7 and delivery occurs. This is the last portion of the
8 tracing prior to taking everything off. The pulse oximeter
9 was in the reassuring above 30 range.

10 The outcome is seen on this slide, at 8:04. No
11 forceps delivery. Apgar of this male child was 9 and 9, and
12 the cord blood gases, artery and vein, you can see are a
13 reassuring 7.20 and base deficit of minus 6.8 in the artery
14 which I think is most reassuring.

15 At this time, I would like to introduce Nancy
16 Townsend, clinical nurse specialist, who was the Vanderbilt
17 Nellcor project research nurse. Thank you very much for
18 your attention.

19 **Nursing Perspective**

20 MS. TOWNSEND: Good morning. My name is Nancy
21 Townsend. I am an advance practice nurse at Vanderbilt
22 Medical Center, in Nashville, Tennessee, and was the study
23 coordinator for the fetal pulse oximetry research project at
24 our institution.

25 For the record, I have no financial interest in

1 the Mallinckrodt Company. Neither I nor any of my family
2 members own stock in the company. During the course of the
3 research project, my salary and benefits package was
4 entirely supported by Nellcor and Mallinckrodt. I am not
5 being paid for this presentation, however, Mallinckrodt will
6 reimburse my travel expenses.

7 Before I begin the formal part of my presentation,
8 I would like to make a note about the case study which Dr.
9 Boehm just presented. In the very beginning part of the
10 case where the sensor was placed, I was the nurse who was
11 present during that case. I handwrote the fetal oxygen
12 saturation readings on the monitor tracing. We had adequate
13 signal quality and that was an appropriate reading, however,
14 it just was not tracing on the monitor because I had not
15 plugged in at the back of the monitor the proper plug to
16 allow it to trace, and I figured that out and plugged it in
17 and that is why it was handwritten in the beginning. Again,
18 we had adequate signal quality and it was an appropriate
19 reading.

20 I have polled many nurses with experience in fetal
21 oximetry, both in my own institution and around the country.
22 The responses have been overwhelmingly positively from a
23 nursing standpoint. There is no greater challenge to the
24 perinatal nurse than coordinating a safe, satisfying birth
25 experience for the child-bearing family. To provide a safe

1 birth experience, the priority in nursing care is to ensure
2 that the mother and fetus are adequately oxygenated. In
3 performing this duty, nurses face a similar dilemma as
4 physicians in that standard electronic fetal heart rate
5 monitoring is only an indirect means of assessing fetal
6 oxygenation status.

7 Furthermore, since nurses have a tremendous impact
8 on patient satisfaction during the labor and delivery
9 experience, nurses must do everything within their power to
10 act as advocates for their patients.

11 Within the context of nursing's key role in the
12 birthing process, my belief is that the fetal oximeter is an
13 important tool for nurses in providing the best possible
14 outcome for families during childbirth in the frequent
15 presence of a non-reassuring fetal heart rate tracing. The
16 entire healthcare team and, thus, the childbearing family is
17 reassured by certain measures of fetal well-being on an
18 electronic fetal heart rate tracing, such as moderate
19 baseline variability or accelerations above the baseline
20 fetal heart rate.

21 Based on scientific evidence, nurses know that
22 these signs are almost always indicative of an adequately
23 oxygenated fetus. However, nurses are also well aware that
24 in the presence of non-reassuring signs, such as fetal
25 tachycardia or late decelerations, the evidence is non-

1 specific. In other words, a non-reassuring tracing may or
2 may not indicate deterioration in fetal oxygenation status.
3 Up until now, the only method of directly assessing fetal
4 oxygenation has been fetal scalp blood sampling.

5 This process is uncomfortable and undignified for
6 the patient, time consuming for the obstetrical staff,
7 costly and, most importantly, invasive to the baby. When
8 fetal status continues to be non-reassuring as assessed by
9 fetal heart rate monitoring, the fetal scalp pH test must
10 often be repeated several times. Furthermore, in the
11 presence of a vaginal infection, laceration of the fetal
12 scalp may potentially expose the fetus to harmful
13 microorganisms.

14 The fetal oxygen saturation monitor is a non-
15 invasive device that gives a clear and direct reading of
16 fetal oxygenation status. The sensor is easy to place at
17 different dilations, with no awkward positioning required
18 for the mother. If the signal is lost, a simple
19 readjustment restores the signal usually without a vaginal
20 examination.

21 Patients are reassured by the presence of a
22 comfortable, non-invasive device which gives clear and
23 continuous information about their unborn child's oxygen
24 status on a monitor that is easy to understand. A concise
25 explanation about the technology and the data allows the

1 nurse to effectively communicate key information about fetal
2 status and, thus, decrease patient anxiety. Objective
3 information from the fetal pulse oximeter is helpful since
4 the nurse can refocus the patient's attention to reassuring
5 fetal oximetry data. Plus, since there is a specific number
6 above which metabolic acidosis may be ruled out, clear-cut
7 reassurance is available to both the patient and the nurse.

8 Information is important to both nurses and
9 patients at every stage of labor. There is a clear link
10 between a patient's feeling of control and with their
11 satisfaction with the birth experience. It has also been
12 suggested that with greater birth experience satisfaction, a
13 woman is better able to mother her child. If the nurse is
14 able to provide clear, concise information from the fetal
15 pulse oximetry system, more information and ultimately more
16 control are given to the patient. Information must be
17 shared with families no matter what the situation. A new
18 technology that is easy to use, non-invasive and provides
19 reliable and objective information is, therefore, a very
20 good thing.

21 Nurses favor fetal pulse oximetry for another
22 major reason. When the fetal heart rate tracing is non-
23 reassuring but fetal oximetry data is reassuring, as if
24 often the case, the nurse is reassured. Thus, the patient
25 is given every opportunity to have a vaginal delivery. This

1 reassurance is documented on the fetal heart rate strip and
2 in the nursing notes, and can be clearly communicated to the
3 physician. The nurse is then free to provide more direct
4 labor support to the laboring mother and her family instead
5 of focusing too much energy on the often ambiguous fetal
6 heart rate tracing. This clear communication of fetal
7 status to other members of the healthcare team, particularly
8 to the physician, is another reason why nurses favor this
9 new technology.

10 Because fetal heart rate tracing interpretation is
11 subjective, it is not unusual for nurses and physicians to
12 disagree when reviewing the same fetal monitor tracing.
13 This disagreement increases the tension and stress among all
14 those involved in the labor-delivery experience. Another
15 advantage of the oximeter is that nurses, with a physician's
16 order and with proper competency validation, may potentially
17 place the sensor independently.

18 In summary, fetal pulse oximetry is easy to use.
19 It is non-invasive to the fetus. It is comfortable for the
20 mother, and it provides valuable, objective information that
21 reassures both the childbearing family, especially the
22 mother, and the obstetrical staff. For all these reasons,
23 nurses like this technology. Most important, it is good for
24 the patients whom we serve.

25 Next, I would like to introduce Ms. Lucy Woods.

Patient Perspective

MS. WOODS: Good morning. My name is Lucy Woods. I have no financial interest in the Mallinckrodt Company. I am not being paid for this presentation, although they are reimbursing me my travel expenses.

By the way, I am the patient in the case that Dr. Boehm discussed earlier. I did have a C-section with my first little boy and, along with having a very healthy big baby boy, I was also very uncomfortable for the first few weeks. I really couldn't be the kind of mom I wanted to be. It was uncomfortable holding him. The baby and I were both uncomfortable when he was trying to nurse. So, it was a few weeks down the road before I really felt that special bond because I think he sensed that I was not myself, and I think the C-section really made a difference in how quick we got comfortable with each other.

I know when I went into labor with my second baby boy, Seth, how much I wanted to deliver him, and I distinctly remember, I was lying there, watching the monitors and I thought that doesn't look right. Well, about the time I felt that my nurse looked at me and I said, "something's not right, is it?" And, she said, "well, I'm going to call Dr. Boehm." Well, he came in and he explained the situation. He told me the possibility that I could have another C-section and immediately I was just, like,

1 terrified. I didn't want another C-section. I wanted to be
2 able to go home to my two-year old, take my baby home, you
3 know, be able to pick them up, love them just like a mom
4 ought to do.

5 Well, Dr. Boehm informed me of a study that
6 Vanderbilt was doing at that time, and it was the pulse
7 oximeter. Well, I said, "please call Nancy in; I would love
8 to talk to her." Well, Nancy, the study nurse, came in and
9 she explained what the pulse oximeter was. She explained
10 how it worked, everything involved with it. Of course, my
11 main question was would there be any harmful effects that
12 the baby could have by attaching the sensor, or anything
13 like that. Well, she reassured me. Dr. Boehm seemed very
14 confident with the instrument and, after talking with both
15 of them, I felt very confident in using the pulse oximeter.

16 There was no discomfort at all when it was placed
17 on the baby's face. A few hours later I delivered a very
18 healthy 7 lb. boy. I remember when Dr. Boehm first showed
19 him to me, I saw tiny little red mark on his cheek. Well,
20 by the time they brought him back to me the red mark was
21 totally gone. You know, I thought, "wow, this is
22 wonderful." Within an hour, I told my husband, I said, "I'm
23 ready to get up." I said, "I feel great." And the nurse
24 that was in my room, she said, "are you sure you can get
25 up?" I said, "yes, I feel wonderful."

1 Well, I got to nurse the baby. He and I bonded
2 immediately. I mean, I could hold him and love on him.
3 When I went home, of course, my two-year old, he missed his
4 mommy; he came running up and I got to hug him and hold him
5 instead of saying, "oh, no, don't touch mommy." And, this
6 just made all the difference in the world. I got to, you
7 know, hug and love my two-year old. I got to nurse my baby,
8 you know, and we were both comfortable doing it. From the
9 very first time he nursed, everything went perfect, which
10 was a definite change from the first time. It was just
11 wonderful being able to go home and be the mom that I wanted
12 to be for both my children, and I just hope -- I think this
13 instrument is just a wonderful thing and if it saves anyone,
14 any woman from having an unnecessary C-section it is well
15 worth it.

16 Thank you.

17 **Restatement of Proposed Indication for Use**

18 MS. PAGE: I would like to recap the results
19 presented here today. We have demonstrated that the use of
20 the N-400 system with conventional fetal heart rate
21 monitoring allows the safe continuation of labor during
22 periods of non-reassuring fetal heart rate and reassuring
23 FSpO2.

24 Use of FSpO2 improves the quality of fetal
25 assessment and results in better matching of delivery

1 indication and immediate neonatal condition.

2 Finally, the improvement in fetal assessment
3 results in a clinically meaningful 50 percent reduction in
4 the rate of C-sections performed for non-reassuring fetal
5 status even though, as we have seen, the overall rate of C-
6 sections remains unchanged.

7 Furthermore, we believe that the results support
8 our indication for use statement as previously presented.
9 The system is to be used as an adjunct to fetal heart rate
10 monitoring. It is to be used in a population of term
11 infants with ruptured membranes who have a non-reassuring
12 fetal heart rate pattern, and the purpose of the N-400
13 system is to improve the physician's ability to assess the
14 fetal status.

15 In summary, we believe that the study we have
16 discussed today does, indeed, constitute a valid
17 scientifically sound study; that the results are clinically
18 significant and, most importantly, we believe that we have
19 demonstrated the safety and effectiveness of the N-400 Fetal
20 Oxygen Saturation Monitoring System for its stated intended
21 use.

22 This concludes our presentation. We thank you for
23 your attention and we do welcome your questions at this
24 time.

25 CHAIRMAN BLANCO: Thank you very much to the

1 company for the presentation. We have a few minutes left
2 over in the schedule and I will see if any of the panel
3 members have any questions that they would like to direct to
4 any of the representatives that presented at this time.
5 Diony?

6 MS. YOUNG: Thank you. I am Diony Young. It has
7 been very interesting listening to all of the information
8 about the study. I want to talk about the real world.
9 First of all, the study assumes that all women are put on an
10 electronic fetal monitor. That is not the case, though the
11 majority are in this country. No mention has been made of
12 osculation, and the use of osculation with the use of this
13 device.

14 My first question would be if a non-reassuring
15 heart rate pattern is found on osculation, what is the
16 sequence of obstetrical interventions? Electronic fetal
17 monitoring, external or internal fetal scalp sampling, then
18 followed by oximetry? That is my first question.

19 My second question is in the study it was
20 mentioned that all the women were attended by a nurse
21 throughout the study. In the real world, with nursing
22 cutbacks, the situation I think is very difficult in most
23 hospitals and few women are lucky enough to have a nurse
24 with them, one-on-one, throughout their labor. So, does
25 this device require one-on-one staff attendants at the

1 mother's bedside throughout her labor? So, could I have
2 answers to those questions, please?

3 DR. GARITE: Yes, ma'am. I think your question is
4 important but I think both questions can be answered simply,
5 and then I will get to the details, with one statement and
6 that is, this device is only intended to be used at present
7 in patients who already have evolved to the point where they
8 have a sufficiently non-reassuring fetal heart rate tracing
9 that the clinician needs further information to assure
10 himself or herself that intervention is not required.
11 Therefore, in the osculation setting we, at our hospital,
12 for a number of years had a free-standing birthing center
13 where our midwives used osculation for the primary means of
14 fetal surveillance. But we, as I believe virtually everyone
15 on this country, in the presence of a non-reassuring
16 osculation backed it up with an electronic monitor. So, you
17 go through the same sequence. If the electronic monitor is
18 non-reassuring then, depending on your mode of practice, you
19 use either accelerations or scalp pH to back that up. If
20 you can't be reassured at that point, you end up with a
21 cesarean section, but this device would then fall into place
22 in that sequence, allowing further more accurate fetal
23 assessment.

24 In the monitor protocol, the research nurse was
25 not the patient nurse. Okay? We weren't trying to create

1 an unreal situation. The research nurse was there to
2 monitor protocol compliance and the patient had a clinical
3 nurse. But, nonetheless, your question is cogent, and the
4 reality is again answered by the same question, you don't
5 have a one-on-one nurse in real life. I wish we did, but we
6 don't, especially in my institution. We can't afford it.
7 But when we do have a patient who has a sufficiently
8 concerning fetal heart rate pattern that you are to the
9 point where you might be intervening for non-reassuring
10 fetal status by either a section or a forceps, at that time
11 most of the time you do have a nurse that is pretty much
12 with the patient.

13 So, that is where we are in the evolution. We are
14 not there at the beginning or in the normal labor, or
15 anything like that. This device is intended to be used when
16 you are at the point where you need further information, and
17 our best estimate is that it is somewhere around 25 percent
18 of patients in this country.

19 MS. YOUNG: Okay. One more quick question, and it
20 relates to maternal position. This device is an adjunct so
21 the woman is actually hooked up to a lot of different
22 things. I sort of have this picture of her flat on her back
23 which, in and of itself, can cause fetal distress. She has
24 a lot of wires coming out of her body. You know, how long
25 is she going to stay there, in that position? I understand

1 that there is a problem if she actually moves. Excessive
2 movement is mentioned in the materials. I don't know quite
3 what excessive movement is but presumably she can't get up
4 and ambulate which also itself may, you know, increase the
5 fetal well-being, if she is allowed to do that because
6 changes in maternal position can definitely have an effect
7 on fetal well-being. So, is she flat on her back for many
8 hours, hooked up to all of these things?

9 DR. GARITE: There is no need to be flat on your
10 back to have this device. Mothers labor on their sides when
11 they have non-reassuring fetal heart rate patterns, and with
12 this device they labor on their sides.

13 Yes, it precludes ambulation. I am not aware of
14 data that suggests that ambulation improves fetal
15 oxygenation. Being on their side or their back does. But,
16 again, I will go back to my statement on my previous two
17 questions, this isn't intended for use on the average
18 laboring patient; it is intended to use when you get to the
19 point where you are so concerned about fetal status that you
20 are now being interventive. In reality, instead of a
21 patient going urgently for cesarean section, the intent is
22 that a woman who can continue laboring at least on her side
23 and, you know, in the spectrum of ambulation, laboring on
24 your side and having a knife on your abdomen, I would choose
25 laboring on my side. But that is just a little editorial

1 comment.

2 CHAIRMAN BLANCO: Let's try to be succinct. We
3 are still in the information stage so other questions?

4 DR. SHARTS-HOPKO: Yes, I understand that this
5 would have been done for control purposes in the study but
6 it is sort of reinforced in all the training materials, is
7 this device intended only when the fetus is in a vertex
8 presentation?

9 DR. GARITE: That is correct.

10 DR. DIAMOND: I had a question about the protocol
11 itself. The patients were consented into the study at what
12 point?

13 DR. GARITE: There were three different phases.
14 In many, and most of the institutions, if possible a
15 brochure was distributed in the prenatal period to describe
16 the study. That was not required. On admission, whenever
17 possible, all patients were consented whether they were
18 eligible for enrollment by the fetal heart rate pattern or
19 not. In other words, if they had a vertex, if they were at
20 term, in labor, they were consented but they were not
21 enrolled or randomized until they met all the eligibility
22 requirements, including abnormal heart rate pattern. Some
23 patients were not enrolled on admission; they were consented
24 at the last point as well.

25 DR. DIAMOND: Do you have the information on the

1 breakdown of patients between the latter two categories that
2 you just described?

3 DR. GARITE: No, I don't have that available.

4 DR. DIAMOND: Okay. The ominous heart rate
5 pattern, the class 3 -- maybe I missed it but I didn't see
6 where that data was -- how many of the patients in each
7 group had that, or how many patients who were consented were
8 then not enrolled because they developed that, precluding
9 time for their randomization to one group or the other? Do
10 you have any of that data?

11 DR. GARITE: Yes, we do.

12 DR. SWEDLOW: This is Dave Swedlow. My
13 recollection is that there were 15 or 16 ominous patterns in
14 both groups, no difference between the two. Could you
15 repeat the first question?

16 DR. GARITE: Actually, the study criteria, in
17 essence, excluded some of those patients. If they developed
18 an ominous pattern before randomization, they weren't
19 eligible for randomization because requiring immediate
20 delivery was an ineligible criterion. I know that could
21 inherently bias things, but you can't randomize someone --

22 DR. DIAMOND: Sure. That is why I was asking. As
23 you indicated, that could potentially bias --

24 DR. GARITE: They couldn't even be enrolled if
25 they developed that pattern as their first pattern.

1 DR. DIAMOND: But they could already have been
2 consented potentially.

3 DR. GARITE: They could have been consented but
4 they would not have been randomized. Those ominous patterns
5 could only be included in patients who were already
6 randomized and already either in the monitor or control
7 group. That could not have been an entry criterion.

8 DR. DIAMOND: Yes, but there could have been some
9 that were consented that developed that, which precluded
10 just what you described.

11 DR. GARITE: Absolutely.

12 DR. DIAMOND: Did that happen once? Did that
13 happen in a hundred patients?

14 DR. GARITE: I have no idea.

15 DR. DIAMOND: You don't have that information?

16 CHAIRMAN BLANCO: Okay. Subir?

17 DR. ROY: I was trying to understand better Dr.
18 Porreco's presentation. In slide 38 we have increased rate
19 of non-reassuring fetal heart rate seen only in test
20 patients sectioned for dystocia, and it is broken down into
21 class 2 variable decelerations, both of which occurred more
22 frequently in the test group. Then, in the next slide, it
23 is concluding that test group reassuring fetal oxygen
24 permits labor to continue unmasking dystocia, therefore,
25 non-reassuring FHR is a marker for dystocia. I mean, did

1 these people have a true identification of fetal distress or
2 not?

3 DR. PORRECO: No. Well, they met inclusion
4 criteria to be randomized and that is what we are
5 suggesting. When we go back and look at our database and
6 observe this retrospectively, the increased occurrence of
7 dystocia. The fourth possibility is that these patients
8 unwittingly were selected not only to have an increased risk
9 of fetal distress or non-reassuring patterns by their
10 inclusion but unwittingly we selected a group, we believe,
11 that had an increased underlying risk of dystocia, and being
12 randomized to the sensor group permits the evolution of
13 their labor and unmask, if you will, that increased
14 predisposition to dystocia. In the absence of the monitor,
15 those patients would have been syphoned off and sectioned
16 inappropriately in some cases for fetal distress had they
17 been randomized to the control group.

18 DR. ROY: So the reason for the cesarean section
19 in the control group would have been for these tracing
20 patterns.

21 DR. PORRECO: Correct.

22 DR. ROY: And for the monitored group for
23 dystocia.

24 DR. PORRECO: The sensor group for dystocia.

25 DR. ROY: But they also had the tracing

1 abnormalities?

2 DR. PORRECO: Yes.

3 DR. ROY: I guess that is what I am getting that.

4 DR. PORRECO: Yes, and their labor, instead of
5 being truncated, if you will, continued because the
6 reassuring improved assessment allowed it to proceed and
7 then ultimately if they, indeed, were in that group of an
8 underlying risk of dystocia, it became manifest two, three
9 or four hours later. That is our interpretation of the
10 occurrence or the observation that these patients had a
11 higher rate of cesarean birth for dystocia.

12 DR. ROY: Are we saying that that is a good thing?

13 DR. PORRECO: Are we saying that is a good thing?

14 DR. ROY: Yes.

15 DR. PORRECO: We are saying that improved fetal
16 assessment is a good thing, doing the right thing at the
17 right time for the right reason. That is, knowing that
18 patients who are intervened upon sometimes in urgent status
19 by cesarean birth is appropriate only when, indeed, the
20 fetus requires that kind of urgent intervention. The fact
21 that they may have an underlying increased risk for
22 dystocia, such that cesarean birth down the line is
23 necessary for that indication, in a different program
24 setting is still appropriate. I mean, that is a good thing.

25 CHAIRMAN BLANCO: Let me go ahead. We need to

1 keep going; the panel has a lot of other questions. Jay?

2 DR. IAMS: I just want to keep going on that same
3 theme here, Dr. Porreco. Sooner or later there would have
4 to be, if I am understanding you correctly, a decline in the
5 number of cesarean sections if this technology were applied
6 to a large population. I don't quite figure out how the C-
7 section rate can stay the same if what you are saying is
8 that your monitor identifies a group of women who may be at
9 increased risk for dystocia, as well as for fetal
10 compromise.

11 DR. PORRECO: Our inclusion criteria --

12 DR. IAMS: The inclusion criteria. Then, once
13 that additional technology is in place, the fetal pulse
14 oximeter should decrease the number of cesareans performed
15 for non-reassuring fetal status, but that you are already
16 operating on a group that has an increased risk of dystocia.
17 We know that if you let labor that is a little dysfunctional
18 go on -- at least, we think we know from the Alabama data
19 and other places, that if you give women a little more time
20 with a slow labor they may reduce their section rate there
21 also. So, somehow or other, I am stumped on the notion that
22 the C-section rate isn't going to change; we are just going
23 to know who needs what. Shouldn't it ultimately go down?
24 If you allow labor to progress safely we should see a
25 reduced number of false diagnoses of dystocia plus fetal

1 intolerance to labor, and we should see a reduced number of
2 sections for fetal compromise, no matter how you define it.
3 So, is it just simply a matter of the N here? We just need
4 to have --

5 DR. PORRECO: Yes --

6 DR. IAMS: -- ten thousand women and then we will
7 see that? Is that what you are saying?

8 DR. PORRECO: Yes, and also the fact that if it
9 was dystocia that we were after, as you suggest, we would
10 need other information -- how much uterine activity were
11 these women subjected to? What were the positions of the
12 vertex in patients who weren't progressing? I think in our
13 ability to look at our database -- since that wasn't what we
14 were after, I think if one were after that, that is the kind
15 of information that would have to be at hand to prove what
16 you have just suggested.

17 CHAIRMAN BLANCO: We are starting to kind of get
18 into a discussion and we really need to get the FDA
19 presentation in before we do that. So, are any other
20 questions of fact that we can either leave with the company
21 to try to get to us by the time we start discussion, and so
22 forth? And, let's try to do that quickly because we are now
23 running behind instead of running ahead. Machele first and
24 then Dr. D'Agostino.

25 DR. ALLEN: I actually have a number of fairly

1 minor questions. I apologize for the number of questions
2 and I hope they are simple to answer.

3 In our reading material, in volume 1 on page 10,
4 talking about the accuracy of the monitor, they say a true
5 fetal pO saturation of 30 percent -- 67 percent of the SpO2
6 readings can be expected to fall between 25 percent and 34
7 percent. I was wondering were you comfortable with the
8 remaining 33 percent that are not accurate?

9 CHAIRMAN BLANCO: I think that is kind of part of
10 the discussion, one of the questions that we are going to be
11 addressing. If there is something that you need to know
12 that they can look up and give us back before we start the
13 discussion --

14 DR. ALLEN: I am just want to go right through
15 them and see if they are appropriate questions or not. I
16 wanted to know of the ten patients where we had difficulty
17 in placing the monitor, any idea what caused the difficulty?

18 DR. GARITE: The most common reason was the
19 vernix. Vernix is the one thing -- a really thick vernix --
20 that interferes with this. Some patients, however, with
21 advanced dilation or low station --

22 DR. ALLEN: You just can't get it in?

23 DR. GARITE: -- also will have difficulty in
24 getting it in, just like an intrauterine pressure catheter.

25 DR. ALLEN: On page 17 there is a referral that

1 there might be an interaction between placing epidurals at
2 less than 5 cm along with the monitor. Is there a real
3 interaction?

4 DR. GARITE: I will give you my best answer.
5 Since the vast majority of patients had epidurals, what
6 happens in the patients who got epidurals at greater than 5
7 cm, which is the minority, is that your N gets so small that
8 you just may not have a substantial size. There is not
9 enough data there to conclude in the other arm that there
10 really is a difference or not.

11 DR. ALLEN: And then with the breaches, is that
12 just because they were excluded so you have no data to
13 support using this with breaches?

14 DR. GARITE: Well, one we didn't include breach as
15 an inclusion criterion and there is really even
16 internationally minimal data on malpresentations. Now, you
17 may have some -- I know there is some.

18 CHAIRMAN BLANCO: Well, let's remember that the
19 indication was to lower cesarean section in vertex with non-
20 reassuring fetal heart rate patterns. So, the breach is
21 really not an issue.

22 DR. ALLEN: Okay. Then, just philosophically, are
23 you comfortable with non-reassuring fetal heart rate
24 tracings being a surrogate marker for eventual dystocia?

25 DR. GARITE: I think we should be very careful how

1 far we go with this conclusion.

2 CHAIRMAN BLANCO: Let's leave that. That is
3 really a discussion question, not a point of fact. Dr.
4 D'Agostino?

5 DR. D'AGOSTINO: Let me ask and maybe they can
6 respond quickly, but it probably would be more appropriate
7 to respond later, after the break or something, the
8 statistical analysis is extremely complex and, being a
9 clinical trial even though open-label and so forth, and
10 unblind once the randomization was done, it still is
11 randomized and I guess one of the analysis techniques that
12 one could have done was to look at overall analysis, like
13 the Fisher exact test or something, very simply chi square
14 or what-have-you, and then, after you have done that,
15 instead of building a big logistic regression model, look to
16 see what happens in subsets.

17 I guess my logic is that I would really worry
18 about subsets in the sense of do you think that you have
19 identified some groups where the procedure doesn't work, the
20 monitor doesn't work? You mentioned site 3, for example,
21 but then there is the obesity and the hypotensive. I think
22 I would like some assurance, and maybe the panel also, is
23 that even though you have identified all these variables
24 that, in fact, you have no reason to think the procedure is
25 going to not work in those subsets. Maybe you could get

1 that for us later.

2 Also, with regard to this dystocia question, let
3 me just throw out an idea that would give some comfort to me
4 and give me some understanding, if you have the monitor
5 added to the monitoring, can you identify individuals that,
6 without the new monitoring, would have gotten a C-section
7 for the non-reassurance and then see what happens to those
8 individuals? You may have said that but it is not clear.
9 It would be very interesting to say where are the dystocia
10 cases coming from.

11 DR. GARITE: I can answer that last question
12 pretty quickly. It is not the exact same patients who had
13 non-reassuring fetal heart rate patterns that ultimately
14 developed dystocia. I mean, you heard from one that didn't.
15 There are a lot of examples of those. It is a population
16 effect; it is not the exact same patient. So, like I said,
17 I think we need to be careful where we go with our ultimate
18 conclusion.

19 DR. D'AGOSTINO: Then lastly, again just to think
20 about it, in slide 27 you have 108 individuals versus 78
21 individuals. It is not completely clear to me where the 108
22 and 78 are coming from.

23 DR. GARITE: I will answer that real quickly.
24 Those are the patients who had either operative vaginal
25 delivery or operative cesarean section for the specific

1 individual of non-reassuring fetal status in each group.

2 DR. D'AGOSTINO: It would be interesting to tie
3 them back to the rates that you give when you are looking at
4 your efficacy analysis because they are not the same --
5 there is an overlap but not exactly the same.

6 DR. GARITE: I think you are confusing C-section
7 and operative vaginal delivery. If you add the two of them
8 together, they should add up. If they don't, we will check
9 it.

10 CHAIRMAN BLANCO: Any other questions of fact over
11 on this side? Dr. O'Sullivan?

12 DR. O'SULLIVAN: Tom, can you tell me how many
13 patients actually had pO2's less than 30 during the
14 monitoring process with the device?

15 DR. GARITE: Persistently less than 30?

16 DR. O'SULLIVAN: That that was what made you
17 decide to do the cesarean section.

18 DR. GARITE: Right, how many actually had that and
19 how many had no signal -- I will have to break it down. So,
20 if you will allow us --

21 CHAIRMAN BLANCO: We will go ahead and maybe after
22 lunch you can give us the answer.

23 DR. GARITE: Yes.

24 DR. CHATMAN: Dr. Garite, the epidural rate in
25 these institutions generally is 95 percent and, if it is

1 not, how that might affect the data, if it does or if it
2 could. I would like to know, secondly, how you arrived at
3 the definition for dystocia because, as you know, there are
4 other criteria that people use for the definition of
5 dystocia.

6 DR. GARITE: We will check. What you really want
7 to know though is, is the epidural rate that high in
8 nulliparous patients, in patients who meet the same criteria
9 who are not study patients because epidurals are clearly
10 high in nullips and mulips so you can't take the hospital's
11 gross rate to compare. So, I will have to get that for you.

12 I can answer the second question though, if you
13 can repeat it again. I am sorry.

14 DR. CHATMAN: Just how you arrived at the
15 definition of dystocia.

16 DR. GARITE: The definition of dystocia was solely
17 based on partogram definitions. We did not prospectively
18 know that dystocia was going to be a problem so we didn't
19 include data on pressure catheters, etc., etc. So, we had
20 to take some arbitrary retrospective definition purely based
21 on labor progress, and it is arbitrary.

22 DR. CHATMAN: Okay.

23 CHAIRMAN BLANCO: Mike?

24 DR. DIAMOND: The first question is a little bit
25 like the one that was asked. As I understood the protocol,

1 for those patients randomized to the test group if the pO2
2 was below 30, that is when they got evaluation for
3 accelerations and scalp pH. How many patients actually had
4 that testing done and, therefore, did not end up with a
5 cesarean section based on the pO2?

6 DR. GARITE: Well, I can tell you the scalp pH.
7 It was 15. It was about 15 in each group.

8 DR. DIAMOND: And did not end up with cesarean
9 section because of that?

10 DR. GARITE: No, not necessarily. That is how
11 many were evaluated. Because, you know, what you have is a
12 dynamic problem. Sometimes the SpO2 falls below 30 percent,
13 you revert to the heart rate, you have accelerations or
14 normal scalp pH, and now your oximetry is normal. Or, the
15 same thing happens with the heart rate, it goes back and
16 forth.

17 CHAIRMAN BLANCO: But he is asking how many had it
18 and what the result was. So, if you can get that data,
19 after lunch we will expect that data.

20 DR. DIAMOND: The point is that the device is
21 being used for that purpose. Was that information helpful?
22 The second question is, based on the pilot data that you
23 collected from other sources, what was the power calculation
24 that you performed? What was the difference you expected
25 between the groups and how did that turn out? How did you

1 make your calculations?

2 DR. GARITE: Well, the calculation was, indeed,
3 based on the pilot study. We determined an overall rate of
4 cesarean section for non-reassuring fetal status, which was
5 --

6 DR. SWEDLOW: We expected a rate of 10-12 percent
7 and we wanted a 50 percent reduction.

8 DR. GARITE: And we calculated 1000 patients, and
9 the numbers were right on the button.

10 CHAIRMAN BLANCO: Dr. Eglinton?

11 DR. EGLINTON: Do you have your statistician here
12 and your database on disc? Can you do another regression
13 analysis here today?

14 [Dr. Garite nods in agreement]

15 CHAIRMAN BLANCO: I think we are going to go ahead
16 and have a break. We are going to have a five-minute break.
17 So, that is all we are going to get so we can get back on
18 track. Thank you.

19 [Brief recess]

20 CHAIRMAN BLANCO: Let's get back and convene the
21 panel again, please. We want to try and keep on time, and
22 we are going to shorten lunch so that we make sure that we
23 have enough time for the discussion agenda. Let's go ahead
24 and get started.

25

FDA Presentations

Preclinical Aspects

MS. DAWS-KOPP: Good morning, ladies and gentlemen, distinguished panel members and guests. I am Kathy Daws-Kopp, the lead reviewer for FDA on this PMA. I have been working with this device since the IDE stage, which started over four years ago, and I am here to give you a brief overview of the review process that we have gone through on the preclinical portions of the PMA.

First, I would like to acknowledge the review team. As you can see, a number of people have been involved in the review of this PMA application. In addition to the basics of clinical and statistical portions that you are most familiar with, our review addressed engineering, materials safety, animal data, optics, manufacturing, human factors, sterilization and patient labeling.

In my presentation, I am going to focus on what we have looked at in the course of our preclinical review. Then I will outline a couple of ongoing issues that we are still working on with the company to resolve.

As Mr. Pollard mentioned in his opening remarks, this PMA was submitted as a modular PMA, where the company is allowed to submit portions of the information required for a PMA ahead of time. In this case, much of the preclinical information was submitted by the company before their clinical data was ready to be submitted.

1 In this slide, I have listed the modules that are
2 "accepted and closed," a term we use with modular
3 submissions to indicate that the issues addressed within the
4 module have been addressed to our satisfaction. Most of the
5 modules listed here were "accepted and closed" prior to the
6 submission of the formal PMA.

7 I will briefly go through the list: General
8 information gave us the basic administrative information
9 required for a PMA.

10 Device characteristics was a detailed device
11 description, how it works, the principles of operation, etc.

12 Pulse Simulator II describes a test tool that the
13 company developed that uses recorded oxygen saturation
14 information to verify and validate device modifications.

15 Software covers the software information,
16 including the software hazard analysis, performance
17 requirements, specifications, design and verification and
18 validation information.

19 Biocompatibility provided information about the
20 safety of the materials used, including test data as
21 necessary.

22 Sterilization provided the sterilization methods,
23 validation, etc.

24 Reports and other information provided a
25 comprehensive literature search and analysis.

1 The PMA was officially submitted in September, and
2 there were still some preclinical review issues outstanding
3 at that time. These were rolled into the PMA. Since
4 September, some of the issues have been resolved. The three
5 items listed here are the only remaining issues for the
6 preclinical review. I will briefly explain each of them.

7 Manufacturing has not been completely resolved
8 because we have not yet conducted our inspection of the
9 manufacturing facilities, a requirement for all PMAs.
10 However, we have reviewed the documentation that was
11 submitted and did not find any problems.

12 The product safety module addressed issues such as
13 optical, thermal and electrical safety. A couple of
14 questions remain in this area.

15 The animal, bench and non-IDE testing covered
16 preclinical and clinical studies other than the pivotal
17 trial discussed here today. This is where the concept of
18 device accuracy is first addressed in the PMA submission.

19 I will now summarize the status of these last two
20 items. In regard to product safety, one of our purposes has
21 been to evaluate the information provided in regard to the
22 thermal and optical safety of the sensor.

23 Because the optical and thermal output of the
24 sensor is so low, we do not expect this to pose a problem.
25 Additionally, very few adverse events that can be

1 specifically linked to the device, such as erythema, were
2 recorded during the pivotal clinical study. However, we are
3 still working with the company to close the gaps on the
4 scientific basis to corroborate the clinical experience.

5 Human factors is intended to look at the usability
6 of the device, that is, to evaluate the ergonomics or how
7 user-friendly the system is. Occasionally, this type of
8 review can turn up troubling juxtaposition of controls or
9 sequences of actions that might lead to unintended results.

10 During a much earlier review of this system, our
11 human factors specialist identified some concerns in regard
12 to ease of use of the device that we communicated to the PMA
13 sponsor. These were in the context of a review of the
14 device in an early draft of the labeling and the safety
15 report. We are now reviewing findings from the pivotal
16 clinical study to see if anything else can be learned about
17 this aspect of device use.

18 We expect that any remaining human factors
19 concerns can be resolved with labeling improvements or
20 changes to the training program for the product. We also
21 typically see improvements in human factors features as
22 products hit the marketplace and the PMA sponsor gets
23 feedback from customers. We would review these changes as
24 PMA supplements down the road.

25 Finally, I would like to speak briefly about

1 device accuracy. The sponsor provided data for the various
2 aspects of accuracy. This was derived from piglet and human
3 neonate and infant studies. As you know, the fetal oxygen
4 saturation value taken from this system can't be compared to
5 a fetal co-oximeter reading because it is not practical to
6 obtain a fetal arterial blood sample. Therefore, reading
7 comparisons could not be conducted in humans at the 20-50
8 percent oxygen saturation levels of interest. Accuracy
9 estimates are, thus, limited in this respect. We accept
10 this as unavoidable.

11 We also continued to work with the PMA sponsor to
12 translate bias and precision data into a clinically
13 meaningful expression that is understandable by the
14 healthcare provider, for example obstetricians, family
15 practice physicians, nurse and nurse midwives, who may not
16 understand these statistical terms.

17 Our clinical reviewer, Dr. Mitchell, whose talk
18 will immediately follow mine, will continue the discussion
19 of accuracy in her presentation.

20 So to conclude, you should know that we will
21 continue to work with the company to resolve these issues in
22 areas of product safety, including human factors and
23 accuracy.

24 I will take questions when the panel Chair would
25 like to entertain them, now or at the end of Dr. Mitchell's

1 talk.

2 CHAIRMAN BLANCO: Let's go ahead and have the
3 questions now, if there are any. There don't appear to be
4 any. We will continue. Thank you.

5 MS. DAWS-KOPP: Okay. Now I would like to
6 introduce Dr. Mitchell, the clinical reviewer on this PMA.

7 **Clinical Aspects**

8 DR. MITCHELL: Good morning. Thank you, Ms. Daws-
9 Kopp, for the review of the preclinical issues that we
10 continue to look at.

11 My job today is to review the clinical study
12 results. The FDA review of the clinical data is ongoing.
13 The purpose of my discussion is to point out the outstanding
14 issues and ask the panel to comment on them within the
15 context of the data as presented by the company.

16 The issues that I will be discussing today include
17 the clinical use characteristics, the study objectives, the
18 unanticipated finding of the cesarean section rates for
19 dystocia. I will do a brief comparison of the baseline
20 versus the randomized, controlled trial, and then I will
21 touch on the labeling.

22 The clinical use characteristics that I am going
23 to discuss today include the 30 percent FSpO2 threshold, and
24 FSpO2 means the saturation as measured by the sensor; the
25 bias and precision; the management matrix; and the

1 registration time.

2 Dr. Swedlow did a nice review of how the 30
3 percent threshold was determined. I will just go back over
4 that quickly. First, there was a literature search of
5 animal investigations, followed by the prospective
6 demonstration of the 30 percent threshold in a fetal sheep
7 model, and then the prospective study in human fetuses.

8 This prospective human study was a measurement of
9 the relationship between fetal oxygen saturation and fetal
10 scalp pH, with acidosis being defined as a pH of less than
11 7.2.

12 There were 46 subjects, laboring patients, and to
13 paired data points. Some subjects had more than one data
14 point. As a result of this study, they found that the
15 sensitivity of the device -- I have a laser pointer but I am
16 afraid to use it because of the way the room is structured,
17 so bear with me a little bit. They found that the
18 sensitivity of the device was 81 percent so that in the 16
19 data points where the pH was less than 7.2, in 13 of those
20 pairs the oximeter reading also was less than 30 percent.
21 The specificity was 100 percent, which meant that of the 34
22 pHs that were greater than 7.2, all 34 oximeter readings
23 were greater than 30 percent.

24 Next, the company went on to look at the bias and
25 precision of the device. SaO2 is blood saturation. And, in

1 the first study listed up there, in an animal piglet model
2 the observed average bias between individual readings of the
3 SaO2 and the FSpO2, in the range of 15-40 percent, was minus
4 0.6 percent, with a standard deviation of 4.8 percent.

5 They did the same study in sick infants and
6 children and found that the average bias was in the same
7 direction, minus 1.9 percent, with a standard deviation of
8 5.4 percent.

9 Then, the precision of the device was examined,
10 and the precision was measured as the standard deviation
11 between two devices measuring oxygenation at the same time
12 from the same fetus, and the precision was 4.7 percent.

13 The sponsor concludes with this summarizing
14 statement: The average difference between FSpO2 from the N-
15 400 and blood SaO2, in the range of 15 and 40 percent is 0.6
16 percent, and the typical variation between these readings is
17 4.7 percent. Thus, for example, at a true fetal SaO2 of 30
18 percent, two-thirds of the FSpO2 readings can be expected to
19 fall between 25 percent and 34 percent.

20 This also means that at a fetal SaO2 of 30
21 percent, 95 percent of the FSpO2 readings can be expected to
22 fall between 20 and 40 percent.

23 Next, I would like to briefly review for you the
24 management matrix for the algorithm for evaluating the fetal
25 heart rate tracings that was used during the study and that

1 is also in the labeling.

2 The company defined the FSpO2 as non-reassuring
3 when the FSpO2 remains below 30 percent between
4 contractions, or no value is available despite sensor
5 adjustment. And, FSpO2 was defined as reassuring when it
6 returns to a value of greater than or equal to 30 percent
7 between contractions.

8 This is a histogram showing the distribution of
9 registration times. Each bar represents 5 percent. So, if
10 you look at the pink bar at the far left, 5 of the patients
11 had the sensor give a reading during 5 percent of their
12 labor. If you look at the light yellow bar on the far right
13 of the tracing, 15 of the patients had the sensor produce a
14 reading 100 percent of the time during their labor. The
15 median registration time then was 67 percent. As you can
16 see, it is not a normal distribution.

17 So, this meant that 50 percent of the patients had
18 a reading present two-thirds of the time during their labor.
19 Conversely, 50 percent of the patients had a reading present
20 less than two-thirds of the time during their labor.

21 So in summary, the clinical use characteristics I
22 have highlighted include the 30 percent FSpO2 threshold;
23 precision and bias; management matrix; and registration
24 time. That will be part of a discussion question.

25 In the next part of my talk I will review the

1 study objectives, as well as the results from the study
2 objectives. There were three objectives to the study. The
3 primary objective was for effectiveness, and it was to
4 reduce the cesarean sections for non-reassuring fetal
5 status, as already stated.

6 Then, there were two safety objectives, that labor
7 would be safe to continue if the oxygen saturation was above
8 30 percent, and that use of the sensor was safe for the
9 mother and the fetus.

10 Before I describe the results as they relate to
11 these three specific objectives, I would just like to take a
12 look at the overall cesarean section rate. As you can see,
13 these are rates so it is reported as a percentage between
14 the test and control group. The test group is on the left
15 and reads as FHR plus FSpO2. The control group is on the
16 right. And, the cesarean section rates were essentially the
17 same, at 26 percent and 29 percent.

18 This is a histogram. It doesn't report rates; it
19 reports the actual number of patients on the vertical axis
20 for the test and control groups, and this is for cesarean
21 sections for non-reassuring fetal status. As we can see,
22 there is a difference. There were 23 cesarean sections for
23 NRFS in the test group and 51 in the control group.

24 The sponsor then spent some time discussing the
25 sensitivity and specificity of the device. The comment I