

TRANSCRIPT OF PROCEEDINGS

1152 60 SEP 20 AM 53

DEPARTMENT OF HEALTH AND HUMAN SERVICES

PUBLIC HEALTH SERVICE

FOOD AND DRUG ADMINISTRATION

CIRCULATORY SYSTEM DEVICES PANEL

This transcript has not
been edited and FDA
makes no representation
regarding its accuracy

Pages 1 thru 243

Gaithersburg, Maryland
September 11, 2000

MILLER REPORTING COMPANY, INC.

735 8th Street, S.E.

Washington, D.C. 20003

(202) 546-6666

mc

TCB

1

DEPARTMENT OF HEALTH AND HUMAN SERVICES

PUBLIC HEALTH SERVICE

FOOD AND DRUG ADMINISTRATION

CIRCULATORY SYSTEM DEVICES PANEL

10:10 a.m.

Monday, September 11, 2000

Gaithersburg Marriott Washingtonian
Washingtonian Boulevard
Gaithersburg, Maryland

P A R T I C I P A N T SPanel Participants

Cynthia M. Tracy, M.D., Acting Chairperson
Megan Moynahan, Executive Secretary

Voting Members

Michael D. Crittenden, M.D.
Julie Freischlag, M.D.

Consultants (with Temporary Voting Status)

Kent R. Bailey, Ph.D.
Melvin L. Griem, M.D.
Geoffrey Ibbott, M.D.
Mitchell Krucoff, M.D.
Kenneth E. Najarian, M.D.
Tony W. Simmons, M.D.
J. Frank Wilson, M.D.

Industry Representative

Gary Jarvis

Guest

Robert L. Ayers, Ph.D.

FDA Participants

James E. Dillard III
Donna-Bea Tillman, Ph.D.
Bram Zuckerman, M.D.
Chris M. Sloan
Kimberly B. Peters
Henry T. Heaton II

C O N T E N T S

<u>AGENDA ITEM</u>	<u>PAGE</u>
Call to Order - Cynthia Tracy, M.D., Acting Chairperson	4
Conflict of Interest Statement - Megan Moynahan, Executive Secretary	4
Sponsor Presentation: Novoste Corporation, PMA P000018, Beta-cath Intravascular Brachytherapy System	
- Andrew Green, Director, Regulatory Affairs, Novoste Corporation	8
- Burton Speiser, M.D., Director, Radiation Oncology, St. Luke's Regional Medical Center	9, 29
- Jeffrey J. Popma, M.D., Principal Investigator, START Trial, Brigham and Women's Hospital, Harvard Medical School	15, 41
- Richard E. Kuntz, M.D., Chief, Clinical Biometrics Division, Brigham and Women's Hospital, Harvard Medical School	32
FDA Presentation	
- Kimberly Peters, Leader Reviewer, FDA	43
- Bram Zuckerman, M.D., Clinical Reviewer, FDA	49
Open Committee Discussion	56
Lunch	87
Open Committee Discussion (Continued)	88
Open Public Hearing	--

P R O C E E D I N G S

1
2 ACTING CHAIRPERSON TRACY: Good morning. While
3 we're waiting for our panel to assemble, I'd like to call to
4 order this meeting of the Circulatory System Devices Panel,
5 and our Executive Secretary will read the conflict of
6 interest statement.

7 MS. MOYNAHAN: The following announcement
8 addresses conflict of interest issues associated with this
9 meeting and is made part of the record to preclude even the
10 appearance of an impropriety. The agency reviewed the
11 submitted agenda for this meeting and all financial
12 interests reported by the committee participants to
13 determine if any conflict exist.

14 The conflict of interest statutes prohibit special
15 government employees from participating in matters that
16 could affect their or their employer's financial interest.
17 However, the agency has determined that the participation of
18 certain members and consultants, the need for whose services
19 outweighs a potential conflict of interest involved, is in
20 the best interest of the government.

21 Therefore, a waiver has been granted for Dr.
22 Mitchell Krucoff for his interest in a firm that could
23 potentially be affected by the panel's recommendations. A
24 copy of this waiver may be obtained from the agency's
25 Freedom of Information Office, Room 12A15 of the Parklawn

1 Building.

2 We would like to note for the record that the
3 agency also took into consideration other matters regarding
4 Drs. Krucoff, Cynthia Tracy, Julie Freischlag, Frank Wilson,
5 and Kenneth Najarian. These panelists reported interests in
6 firms at issue, but in matters that are not related to
7 today's agenda. The agency has determined, therefore, that
8 they may participate fully in all discussions. In the event
9 that the discussions involve any other products or firms not
10 already on the agenda for which an Food, Drug and Cosmetic
11 Act participant has a financial interest, the participant
12 should excuse him- or herself from such involvement, and the
13 exclusion will be noted for the record.

14 With respect to all participants, we ask in the
15 interest of fairness that all persons making statements or
16 presentations disclose any current or previous financial
17 involvement with any firm whose products they may wish to
18 comment upon.

19 ACTING CHAIRPERSON TRACY: All right. At this
20 time I'd like to ask the panel members if they could briefly
21 introduce themselves.

22 DR. BAILEY: I'm Kent Bailey, Biostatistics at
23 Mayo Clinic.

24 DR. CRITTENDEN: Michael Crittenden, cardiac
25 surgeon, West Roxbury VA, Harvard Medical School.

1 DR. SIMMONS: Tony Simmons, Cardiology, Wake
2 Forest University.

3 DR. IBBOTT: I'm Geoff Ibbott, medical physicist
4 at the University of Kentucky in Lexington.

5 MS. MOYNAHAN: I'm Megan Moynahan, Executive
6 Secretary of the Circulatory System Devices Panel.

7 ACTING CHAIRPERSON TRACY: I'm Cynthia Tracy. I'm
8 from Georgetown Hospital, cardiology.

9 DR. FREISCHLAG: I'm Julie Freischlag, a vascular
10 surgeon from UCLA Medical Center.

11 DR. KRUCOFF: I'm Mitch Krucoff. I'm a
12 cardiologist at Duke University Medical Center and the
13 director of device clinical trials at the Duke Clinical
14 Research Institute.

15 DR. WILSON: Frank Wilson. I'm a radiation
16 oncologist, Medical College of Wisconsin, Milwaukee.

17 DR. NAJARIAN: Ken Najarian, interventional
18 radiologist, University of Vermont.

19 DR. GRIEM: Mel Griem, University of Chicago,
20 radiologist and radiation biologist.

21 MR. DILLARD: Jim Dillard. I'm the Director of
22 the Division of Cardiovascular and Respirator Devices at the
23 Food and Drug Administration.

24 MS. MOYNAHAN: I'd like to briefly mention that
25 Robert Dacy, our consumer representative, won't be

1 participating today. He was hospitalized late last week.
2 And we attempted to find a replacement for him but were
3 unsuccessful. So we'll be proceeding today without Robert
4 Dacy.

5 I'd like to also read the appointment to temporary
6 voting status for today:

7 Pursuant to the authority granted under the
8 Medical Devices Advisory Committee Charter, dated October
9 27, 1990, as amended April 18, 1999, I appoint the following
10 people as voting members of the Circulatory System Devices
11 Panel for this meeting on September 11, 2000: Cynthia
12 Tracy, Tony Simmons, Kent Bailey, Kenneth Najarian, Frank
13 Wilson, Mitchell Krucoff, Melvin Griem, Geoffrey Ibbott. In
14 addition, I appoint Dr. Cynthia Tracy to act as temporary
15 Chair for the duration of this meeting. For the record,
16 these people are special government employees and are
17 consultants to the panel under the Medical Devices Advisory
18 Committee. They have undergone the customary conflict of
19 interest review and have reviewed the material to be
20 considered at this meeting.

21 It's signed by David W. Feigald, Director of the
22 Center for Devices and Radiological Health.

23 ACTING CHAIRPERSON TRACY: Okay. We'll move on to
24 the open public hearing. Are there any parties present who
25 would like to make a presentation at this time?

1 [No response.]

2 ACTING CHAIRPERSON TRACY: If not, then we'll
3 begin with the sponsor's presentation, and I'd like to
4 remind the speakers to introduce yourselves and to state any
5 conflict of interest.

6 MR. GREEN: Good morning. My name is Andrew
7 Green, and I'm the Director of Regulatory Affairs for
8 Novoste Corporation. Madam Chairman, panel members,
9 representatives of the FDA, we are pleased today to present
10 to you a review of the data that supports the safety and
11 effectiveness of the Beta-cath system in the treatment of
12 in-stent restenosis.

13 Today's presentation will include: this overview;
14 a device and procedure summary by Dr. Burton Speiser, a
15 radiation oncologist and investigator in the START trial, as
16 well as the principal investigator in the START 4020 trial;
17 a review of clinical results by Dr. Jeffrey Popma, principal
18 investigator for the START trial; a device performance and
19 training review, again, by Dr. Speiser, investigator in the
20 START trial; a discussion of specific topics by Dr. Kuntz,
21 Director of the Cardiovascular Data Analysis Center, who did
22 the analysis for the START trial; and some concluding
23 remarks by Dr. Popma, the principal investigator.

24 The Novoste PMA P000018 was submitted on April 17,
25 2000, and requests approval for the 30 millimeter Beta-cath

1 system, specifically developed for intravascular
2 brachytherapy in the cath lab, in the treatment of in-stent
3 restenosis of native coronary arteries 2.7 to 4 millimeters
4 in diameter. This was a randomized, multi-center, placebo,
5 triple-masked trial for in-stent restenosis for the Beta-
6 cath system.

7 We believe that the data that will be reviewed
8 today by Drs. Speiser, Popma, and Kuntz will show that the
9 Beta-cath system has demonstrated effectiveness in
10 significant reductions in all clinical and angiographic
11 outcomes, demonstrated safety and significance reductions in
12 major adverse cardiac events without increased risk of
13 thrombosis to patients, and demonstrated ease of use, short
14 treatment times, minimal exposure to patients and staff,
15 allowing the clinicians to stay in close contact with the
16 patients through the duration of the treatment.

17 At this time I'd like to have Dr. Speiser come up
18 and present the device and procedure summary. He is, again,
19 a radiation oncologist and an investigator in the START
20 trial.

21 DR. SPEISER: Thank you. Novoste has reimbursed
22 normal travel expenses and as well as paid an honorarium for
23 my meeting attendance.

24 I'd like to first discuss the use of Strontium-90.
25 Overall, radiation in any form has a fairly long history for

1 proliferative diseases. External radiation has been used to
2 prevent keloid formation and heterotrophic bone formation.
3 Brachytherapy, specifically Strontium-90, has approximately
4 a 50-year history for the treatment of pterygia, a
5 proliferative disorder of the eyes.

6 The Strontium-90 has been used, as I mentioned,
7 for benign proliferative conditions in the past. The
8 primary mechanism of in-stent restenosis is a proliferative
9 problem in intimal hyperplasia, and the Strontium-90 has an
10 excellent therapeutic ratio, that is, a dose-to-target which
11 is much higher than the dose-to-non-target tissue.

12 Now, while I'm talking about Strontium-90, in
13 effect, Strontium-90 decays to Yttrium-90, and in the decay
14 of Yttrium-90 to Zirconium, we have an energetic beta
15 particle of 2.27 MeV, and that in effect is the issue(?)
16 that we're using. However, for simplicity, I will be
17 referring only to Strontium-90.

18 Strontium-90 has some very advantageous features.
19 It has a dose rate which is quite high, providing very short
20 treatment times, in the range of 3 to 5 minutes, a long
21 half-life of 28.8 years, which eliminates the problem
22 associated with frequent source replacement. Its dose
23 penetration is limited, which matches the dose profile
24 needed in the coronary arteries. It then also means that
25 there's minimal exposure to non-target tissues, as defined

1 as tissues greater than 1 cm from the source axis, and the
2 physician is able to stay with the patient during the entire
3 procedure.

4 This is a graph showing the ISO dose curve with--
5 iridium-192 is the isotope that I most frequently use for
6 brachytherapy and which I have the greatest experience, and
7 in red is the Strontium-90. The bottom is in centimeters,
8 not in millimeters. And what it shows in the shaded area is
9 that in the range of the arteries that we'll be treating,
10 both isotopes provide excellent delivery. The added value
11 for the Strontium-90 is that there's a very rapid dose fall-
12 off, providing extra safety.

13 Now, the exposure on Strontium-90 to the patient
14 is approximately three-tenths of a millirem per procedure,
15 and that's contrasted to an average dose that they received
16 during the fluoroscopy during the procedure of 350 millirem.

17 The radiation oncologist and the interventional
18 radiologist or interventional cardiologist receives
19 approximately two-tenths a millirem and 4 to 16 millirem for
20 the fluoroscopy. In addition, the radiation oncologist
21 receives another 4 millirem per procedure hand dose by
22 handling the device both pre- and post-procedure.

23 Now, the radiation exposure for the patient is
24 quite low. The Strontium-90 component is less than one-
25 tenth of 1 percent of the total dose received during the

1 procedure. The radiation oncologist and interventionalist
2 receive an extremely low dose of the yearly maximum
3 allowable, and the cath lab personnel receive a very tiny
4 dose of their total yearly maximum allowable.

5 I'd like to go into the Beta-cath system and the
6 procedure as it's used.

7 Now, the system is an integrated system that
8 contains a source train of Strontium/Yttrium-90 within a
9 transport device which is then mated to a specific delivery
10 catheter, the Beta-cath catheter, and is complemented with
11 different system accessories.

12 Now, the features of the system is that it's a
13 completely closed system, which allows for controlled
14 delivery and return of the source train such that the
15 sources never make contact with the patient's blood or
16 tissues. In addition, the device has safety interlocks to
17 prevent the sources from inadvertently being discharged from
18 the device unless everything is hooked up completing the
19 closed system. Once again, the treatment time is very
20 short, 3 to 5 minutes, and that allows the physician to
21 remain with the patient during the entire procedure.

22 Now, the State of Georgia has performed a safety
23 evaluation of the Beta-cath system and has issued a sealed
24 source and device registration certificate August 4th of
25 this year, and this certificate has been included in the

1 Nuclear Regulatory Commission sealed source and device
2 registry.

3 Now, the team consists of radiation oncologists,
4 the interventionalist, either the radiologist or
5 cardiologist, a medical physicist, and the complementary
6 cath lab staff.

7 Now, while the interventionalist is completing the
8 angioplasty, the radiation oncologist, with the aid of the
9 physicist, will prepare the Beta-cath system. And this
10 includes over here putting the device in a sterile bag and
11 attaching the syringe--this is a hydraulic system--then next
12 attaching the delivery catheter, and you prime the system to
13 ensure that there's sufficient fluid in the system before
14 the start of the procedure.

15 Next is the prescription of the dose and treatment
16 time based on individual assessment of the reference vessel
17 diameter.

18 Now, the dose is prescribed at 2 millimeters from
19 the center of the source axis, and this is based on
20 individual assessment of the reference vessel diameter. So
21 for an RVD that's equal to or greater than 2.7 millimeters
22 or less than or equal to 3.3, the dose delivered was 18.4
23 Gray and for the larger vessel it was 23 Gray.

24 Next, the delivery catheter is placed across the
25 injury site, and this shows the marker bands here and here

1 that on the delivery catheter to aid the interventionalist
2 to place it across the appropriate area. And this is done
3 with the aid of fluoroscopy.

4 And next is the delivery of radiation, and this is
5 an animation showing the sources coming out hydraulically.
6 Now they're here, there are markers here and here to further
7 assess the placement of the sources. And the source train
8 remains for approximately 3 to 5 minutes, and the placement
9 is checked periodically with fluoroscopy, so you can see the
10 marker here and the marker band on the catheter just outside
11 of that area to ensure that the source train is at the
12 proper position.

13 In addition, what I'd like to do is, if you can
14 see this arrow over here, there's a light that indicates the
15 amount of pressure. As long as that light or any of the two
16 lights above it are lit, there's proper pressure to maintain
17 the source train at the right position without any source
18 drift. Following the completion, the sources are
19 hydraulically removed back into the transport device.

20 At that point the radiation oncologist removes the
21 system, which is both the transport device and catheter, and
22 the interventionalist completes the procedure.

23 I'd like to introduce Dr. Jeffrey Popma, who is
24 the principal investigator for the START trial, who will
25 discuss the clinical results of the trial.

1 DR. POPMA: Dr. Tracy, panel members--next slide,
2 please?--I have no equity interest in Novoste. I will
3 receive travel expenses but no honorarium for today's
4 meeting.

5 Next slide, please?

6 Again, I apologize to Dr. Krucoff, but what I want
7 to do for the non-interventionalist group is to put a little
8 bit of a perspective on where we stand with in-stent
9 restenosis.

10 We estimate this year that over 725,000
11 procedures, coronary interventional procedures will be
12 completed in the United States, and in these 725,000
13 procedures, 80 percent of the patients will receive one or
14 more stents.

15 Now, stents have been very useful for us in the
16 cath lab to prevent restenosis, but, nevertheless, clinical
17 restenosis still occurs in 10 to 20 percent of patients.
18 And what that means is that over 100,000 patients will
19 develop recurrent symptoms due to in-stent restenosis in the
20 United States this year.

21 We have several existing treatment options. Most
22 commonly we perform balloon angioplasty, repeat dilatation
23 within the segment. We tried for a while using a stent
24 within a stent, and we were disappointed that that did not
25 prevent recurrence for in-stent restenosis. We've tried

1 rotational atherectomy, directional atherectomy, ex-(?) or
2 angioplasty, but to date we've had no randomized trials that
3 have demonstrated that this lowers the frequency for
4 recurrence. And, of course, oftentimes, the patient is left
5 with the option of bypass surgery or in many cases a repeat
6 coronary artery bypass operation because the first one has
7 had some limitations.

8 Next slide, please?

9 Now, we have also learned that, depending upon the
10 pattern of restenosis, the recurrence rate will vary. In
11 very focal, discrete lesions, we know from data from Roxanna
12 Mayron (ph) at the Washington Hospital Center that the
13 recurrence rate is 20 percent. However, restenosis often
14 recurs in a much more proliferative pattern, and when the
15 lesion recurs with length more than 10 millimeters, when
16 it's proliferative, or when it's totally occluded, the
17 recurrence rates after treatment of in-stent restenosis
18 range from 35 to 83 percent. It's still a problem.

19 Now, we have some good randomized trial data that
20 has been done. This is from the ARTIST trial that was
21 reported at the ESC last year, and this was a randomized
22 trial of patients coming into these European investigators'
23 clinical practices where patients either received balloon
24 angioplasty, conventional way of treating patients, or a
25 debulking device, rotational atherectomy. And what was

1 found in this I think was relatively important. The
2 restenosis rates ranged between 51 percent to 64 percent, a
3 slight increase in the restenosis rate with rotational
4 atherectomy, but target vessel revascularizations occurred
5 between a third and a half of patients. So it's still a
6 clinical problem, and we have not fixed this with aggressive
7 debulking therapies.

8 So the purpose of the START trial was to assess
9 the safety and efficacy of intracoronary beta radiation
10 using a Strontium-90 source train following successful
11 coronary interventions in patients with in-stent restenosis.

12 I should also emphasize that this is a unique
13 trial design for device trials in that it is a large-scale
14 trial, which was prospectively constructed, including 50
15 centers, was triple-masked so that the patients, the
16 investigators, as well as the core laboratories did not know
17 which treatment strategy the patient received, and this
18 randomized trial included 476 patients with successfully
19 treated in-stent restenosis. Two hundred and forty-four of
20 these patients were randomized to treatment with Strontium-
21 90, and 232 of these patients were randomized to treatment
22 with placebo.

23 We'll talk a lot about endpoints, and I'm going to
24 give you some definitions for these endpoints in just a
25 moment. But the primary efficacy endpoint of this trial was

1 8-month target vessel failure. The secondary efficacy
2 endpoint was 8-month angiographic restenosis, the occurrence
3 of in-stent--the in-stent minimal lumen diameter, and the
4 degree of the late lumen loss within the vessel.

5 The safety endpoint was the 8-month major adverse
6 cardiac event rate, and the occurrence of a new aneurysm
7 formation.

8 Let's talk about these definitions, and they will
9 be presented as the primary and secondary endpoints. So
10 I'll ask your patience to move through these.

11 First of all, the primary endpoint will be target
12 vessel failure, and that's defined as target vessel
13 revascularization, the clinical need for a repeat
14 revascularization procedure, myocardial infarction or death
15 that could not be clearly attributed to a vessel other than
16 the target vessel. Major adverse cardiac events were
17 defined as death, Q wave and non-Q wave myocardial
18 infarction, emergency CABG, and target vessel
19 revascularization.

20 For this study and all the studies that we do that
21 evaluate restenosis, we used target vessel revascularization
22 as an endpoint, which was defined as any clinically driven
23 repeat percutaneous intervention of the target vessel or
24 bypass surgery of the target vessel. What this means for
25 the patient is the following: They've had a successful

1 treatment. They go home from the hospital. Somewhere
2 between 3 and 6, maybe 8 months later, they develop
3 recurrent symptoms. They have an exercise test that shows
4 that there's a problem with the distribution in the area
5 that was treated. They come back in and they have an
6 angiogram, and the angiogram shows that there's been a
7 significant renarrowing within the area of either the vessel
8 for target vessel revascularization or of the lesion itself
9 when we use the endpoint of target lesion revascularization.

10 Again, clinically driven, exercise test recurrent
11 symptoms associated with the appropriate anatomy that shows
12 that this is, in fact, clinical restenosis.

13 This trial was supported by the following
14 individuals: Rick Kuntz ran the data coordinating center;
15 Alexander Lansky ran the core laboratory; Peter Fitzgerald
16 served as the IVUS intravascular core laboratory director;
17 Peter Zimentbaum did the EKG core laboratory.

18 I think it's very important to emphasize that this
19 study was watched over very carefully with Tom Ryan as the
20 head of the Data and Safety Monitoring Committee. Tom Ryan
21 is a professor of medicine at Boston University and has
22 chaired up some very important Data and Safety Monitoring
23 Committees. I think it's fair to say that he really--with
24 this study and with the Beta-cath trial--has done a very
25 good job, I think, overseeing this very important committee.

1 The Clinical Events Committee was chaired by Dave
2 Cohen at the Beth Israel Hospital in Boston.

3 Now, let's put this talk about which patients are
4 going to be included in the study. This study included
5 single lesion, single intervention, where there was a
6 greater than 50 percent narrowing within the previously
7 placed stent. The target vessel diameter needed to be
8 between 2.7 and 4 millimeters in diameter. And we'll talk
9 about this in a bit, but the target lesion length was
10 treatable with a 20-millimeter balloon, in which case we
11 used a 30-millimeter source train, or treatable with a 30-
12 millimeter balloon, in which case we used a 40-millimeter
13 source train, in order to have adequate margins outside the
14 injured area.

15 The major exclusion criteria were multi-vessel
16 coronary intervention, a target lesion residual stenosis of
17 greater than 30 percent. The patients needed to have a
18 successful procedure before they were entered into this
19 study. Other exclusion criteria included unprotected left
20 main disease and prior chest radiotherapy.

21 Dr. Speiser has reviewed the dose prescription.
22 It was done visually, and it was 18.4 Gray in reference
23 vessel diameters between 2.4 and 3.3, and 23 Gray in
24 reference vessel diameters between 3.3 and 4 millimeters.

25 We'll talk about antiplatelet therapy at the

1 initiation of this trial. We felt and left the adjunct
2 antiplatelet therapy to the discretion of the physician.

3 Now, I also want to stop and say this is adjunct
4 antiplatelet therapy we're talking about. Aspirin was
5 standard therapy in all patients after their intervention.
6 The adjunct antiplatelet therapy that we'll be discussing
7 are drugs like ticlopidine(ph) or clopidogril(ph), and
8 they're given in addition to aspirin therapy. And at the
9 initiation of this protocol, we left that decision to the
10 investigators.

11 We then learned some important information from
12 Tom Ryan and the Data and Safety Monitoring Committee from
13 the Beta-cath trial and suggested that there might be a
14 benefit in extended antiplatelet therapy in those patients
15 who received a stent. So on March the 19th, we modified the
16 adjunct antiplatelet regimen, and we recommended at that
17 time that there be a minimum of 90 days of adjunct
18 antiplatelet therapy with the placement of a new stent.

19 The results that you will see today with the 8-
20 month clinical follow-up will comprise 96 percent of
21 patients that were included in the study. This is an
22 updated report from what you have seen in your panel pack.
23 Angiographic follow-up was obtained in over 80 percent of
24 patients, and this is superb for a clinical device trial.

25 I'll just go through some very basic clinical

1 demographics. Let me just say at the outset that they were
2 balanced in the two groups: age, gender, the presence of
3 diabetes, prior myocardial infarction, prior bypass surgery,
4 all balanced without differences between the two groups.

5 The reference vessel diameter was 2.77 millimeters
6 and 2.76 millimeters in the two groups, and there was no
7 significant difference in minimal lumen diameter, the pre-
8 procedural percent diameter stenosis, the lesion length, or
9 the percent of vessels that were treated in the left
10 anterior descending artery.

11 We did use debulking devices frequently to obtain
12 a successful procedure in the study, and rotational
13 atherectomy was used most commonly in approximately 40
14 percent of patients. We also used stents relatively
15 infrequently. Only 20 percent of patients received a new
16 stent. I think it's important to emphasize that we reserved
17 the new stent use for bail-out indications, and those bail-
18 out indications were for a severe residual stenosis or a
19 dissection that occurred after radiation therapy. So the
20 stents that were placed in the study were placed after
21 radiation therapy had been delivered. Again, you got into
22 the study because you were felt to have a successful
23 procedure.

24 Now, I'll just go over just very briefly what
25 antiplatelet therapy the patients actually received if there

1 was not a new stent placed. Forty percent of the patients
2 did not receive antiplatelet therapy or that information was
3 not available to us, which means that 60 percent of patients
4 who did not receive a new stent received subduration of
5 antiplatelet therapy. In patients who had a new stent
6 placed, we did not know or there was no antiplatelet therapy
7 given in 8 percent, but the vast majority of patients
8 received antiplatelet therapy most commonly between 1 to 30
9 days, 11 percent received 30 to 60 days of antiplatelet
10 therapy; another quarter of patients received more than 60
11 days of antiplatelet therapy.

12 I'm going to have to take a moment and explain
13 this slide. You're going to see this slide several times,
14 and I think it's going to be important that I just take a
15 moment and try to set the stage.

16 When restenosis occurs in patients who have a
17 stent, it occurs because there's intimal hyperplasia, tissue
18 growth within the stent. And our conventional method of
19 analyzing that tissue growth is by analyzing the stented
20 segment itself. If tissue grows within the stent, it most
21 commonly occurs within the axial length of the stent, and
22 this is where we would determine our stented segment
23 recurrence rate. So that when we do our conventional
24 analyses, we look at the stented segment itself, and we'll
25 talk about those results, and you'll see percent stenosis

1 and restenosis rates that are specifically confined to the
2 axial length of where the initial stent was placed.

3 You'll also see a total analysis segment, and that
4 total analysis segment includes a lot of things. It
5 includes along the axial length of the vessel where the
6 stent was. In your panel pack, you'll see that we also have
7 the numbers for where the injury was with the balloon, where
8 the radiation was delivered, but the analysis segment that
9 you'll see will include the stented segment, a little bit
10 longer length with the tissue that's injured, a little bit
11 longer length that had the radiation, and then 5 millimeters
12 both proximal and distal to that will be included then in
13 this longer axial length. And we'll talk more about this in
14 just a minute.

15 The restenosis rates were significantly reduced
16 within the stented segment, from 41.2 percent to 14.2
17 percent, a 66-percent reduction. Within the analysis
18 segment, there was also a significant reduction from 45.2
19 percent to 28.8 percent.

20 Now, these are the clinical outcome measures
21 within the study. We'll talk about the primary endpoint,
22 target vessel failure, about major adverse clinical event
23 rates, and then the clinical indices of restenosis of target
24 vessel failure and target lesion revascularization. Those
25 are what we use.

1 There was a 31-percent reduction in target vessel
2 failure, a 31-percent reduction in major adverse clinical
3 event rates, a 34-percent reduction in target vessel
4 revascularization, a 42-percent reduction in target lesion
5 revascularization--again, all the clinical endpoint
6 parameters that we use to demonstrate efficacy in reduction
7 of clinical restenosis.

8 This is a major adverse event-free survival curve.
9 It shows that the two curves are superimposable out to 90
10 days. And then at 90 days, which is the typical time course
11 when restenosis occurs, then there's a separation of these
12 curves showing a benefit of Strontium-90 therapy over
13 conventional--placebo-treated patients, and this is a
14 significant difference out to 360 days.

15 We also need to talk about subacute stent
16 thrombosis, and I want to take a moment and do that. There
17 was one patient in the placebo group who developed subacute
18 stent thrombosis within the first 30 days of the procedure.
19 There was no patient in the treated group that had that
20 event.

21 Between 31 and 240 days, which is the time
22 endpoint, there were no occurrences of subacute stent
23 thrombosis in either group.

24 You will hear in just a moment about one patient
25 who had a clinically adjudicated subacute stent thrombosis

1 that occurred at 244 days, and we'll talk about the details
2 of that patient in just a moment.

3 The angiographic total occlusion rate between the
4 two groups was identical, 3.3 percent versus 3 percent in
5 those treated with placebo. So no differences in the
6 occurrence of late total occlusion at follow-up and no
7 differences between the occurrences of the clinical event of
8 subacute stent thrombosis.

9 Let's go into the details about the patient who
10 developed the event at 244 days. On March 2, 1999, the
11 patient's mid right coronary artery was successfully treated
12 in the radiation group. Following radiation treatment, a
13 new stent was placed because the clinical investigator
14 identified a dissection. Despite the fact that a new stent
15 was placed, there was a 48-percent residual stenosis within
16 the treated area, within the stented segment by the core
17 angiographical laboratory, and that suggests to us that this
18 was a suboptimal initial treatment.

19 On 11/1/99, which was 244 days after the
20 treatment, the patient presented with chest pain and EKG
21 changes which were new in posterior-lateral Q waves. The
22 angiogram showed a total occlusion of the mid right coronary
23 artery. The proximal mid right coronary artery then
24 received a balloon angioplasty and additional stent.

25 Now, there are a lot of issues to something like

1 this, but we were most conservative in how we reported this.
2 And it is not clear to us whether or not this represented
3 progression of disease or a new stent thrombosis, but we're
4 going to classify it in our presentation to you most
5 conservatively as a new stent thrombosis. And to summarize,
6 that means that there was one stent thrombosis in the
7 placebo group that occurred within the first 30 days, and
8 even after the 360 days where we follow up patients now,
9 there's only one additional event that occurred in the
10 Strontium-90-treated group.

11 Next slide.

12 Let's talk about the 8-month safety results.
13 There were four deaths in the study--one death in the
14 placebo group and three deaths in the Strontium-90 group.
15 The overall incidence of about 1 percent death rate is what
16 we'd expect for a randomized clinical trial of this size.
17 The occurrence of myocardial infarction, there were seven in
18 the placebo group, four in the Strontium-90-treated group.
19 There was one patient in your panel pack that was classified
20 as having had an aneurysm. When this was reanalyzed by Dr.
21 Lansky in the angiographic core laboratory, it was found
22 that this aneurysm was present at baseline, and it did not
23 significantly change during the follow-up period.

24 Let's just review very briefly a description of
25 the deaths in this study.

1 The first patient was a 77-year-old patient who
2 was successfully treated with radiation on 12/7/98. He died
3 193 days later of complications following surgical resection
4 of a colonic polyp. The causes of death included his
5 longstanding coronary disease, congestive heart failure, and
6 multi-organ system dysfunction.

7 The next patient was an 83-year-old man who was
8 successfully treated with radiation on 3/4/99. He died 225
9 days after treatment. The cause of death was metastatic
10 prostate and rectal cancer.

11 The third patient in the radiation group was an
12 83-year-old patient who was successfully treated with
13 radiation on 3/5/99. He died 160 days later, two days
14 following a left upper lobectomy for lung cancer. The death
15 was reported as a post-operative myocardial infarction.

16 The final patient was a 69-year-old patient
17 successfully treated in the placebo group on 1/22/99. He
18 died 102 days after treatment with the official cause of
19 death reported as cardiac arrest.

20 We'll talk more about some of the specifics of
21 this in just a moment with Dr. Speiser and Dr. Kuntz's next
22 presentations, but I think what we can take home from the
23 START trial is the following: The 8-month clinical outcome
24 summary shows significant reductions in all outcome
25 parameters, which include target vessel failure, major

1 adverse cardiac event rates, target vessel revascu-
2 larization, target lesion revascularization, and
3 angiographic restenosis. We do not feel that there was an
4 increased risk of thrombosis. There was one subacute
5 thrombosis in the placebo group and conservatively
6 classified one in the Strontium-90-treated patient that
7 occurred at 244 days, and there were no new aneurysm
8 formations found.

9 What I'd like to do is to turn this back over to
10 Dr. Speiser who will discuss device performance and training
11 programs based on what we found in the START trial.

12 DR. SPEISER: Thank you.

13 First, I'd like to briefly go over the device
14 performance. In the START trial, there were 476 patients
15 enrolled. Successful treatment occurred in 467, or 98.1
16 percent. The causes for unsuccessful treatment were the
17 catheter did not optimally cross the lesion in six patients,
18 or 1.3 percent, and the sources could not be sent in three
19 patients, 0.6 percent.

20 The minor device malformations I will refer to as
21 MDMs. There were 89 patients that had successful treatment
22 with MDMs. Now, the total of these different reported MDMs
23 here equals greater than 89 because some of the patients had
24 more than one observation. If, however, you look at the
25 source transit time of greater than 5 seconds and the source

1 marker drift, this accounted for approximately 90 percent of
2 the cases.

3 One of the cases correctly listed here had in your
4 packet been inadvertently categorized as a reported event in
5 the Beta-cath study when, in fact, it was in the START
6 trial. That patient had a reported event but not a mis-
7 administration.

8 Now, with the observations of these MDMs and
9 experience gleaned from the START trial, Novoste has made
10 device modifications to the Beta-cath system, and those
11 modifications are what was submitted to the FDA in this PMA
12 that is being reviewed. In addition, they were able to
13 create an in-depth training program that incorporates
14 experience specifically from the START trial and have
15 modified the user's manual to include detailed instructions
16 on component connections, pressure tests and monitoring, as
17 well as the manual removal procedure, which is a mandatory
18 procedure for any brachytherapy type of treatment.

19 Next, briefly I'd like to discuss the training
20 program, and this will consist of regional training where
21 the individuals and team will train on the device,
22 procedures, both the treatment and safety, and the roles and
23 responsibility of each team member. Specifically, there
24 will be a hands-on session with the devices to familiarize
25 everybody with the device and detailed instructions for the

1 individuals and team with specific experience from the
2 various trials. There's also cross-training of team members
3 on terminology and professional fields so that there's no
4 confusion in the cath lab and radiation safety training for
5 everybody involved.

6 This is then followed by on-site--which means the
7 facility where the device will be used--training and, once
8 again, first will be reinforced the training on the device,
9 the procedures, and roles and responsibilities. Again,
10 detailed instructions will be given a second time. The
11 procedure will be demonstrated, and that in turn will be
12 followed by a mock procedure conducted in the cath lab; and
13 last, but not least, reinforcement of the radiation safety
14 training for all members of the team.

15 And the last phase of the training program will be
16 a proctored program with an estimate of 3 to 5 procedures
17 that will be proctored to assess the team proficiency with
18 the procedure and system and to advise team and individuals
19 on device use and handling.

20 I have one slide I'd like to put in for long-term
21 safety. The BERT trial now has four-year follow-up, and in
22 this particular trial, what I'd like to do is to show that
23 approximately from 12 months to 48 months, the curve is
24 flattened out, which would indicate that there are no long-
25 term problems that are currently being seen with the

1 radiation arm.

2 Thank you.

3 I'd like to next introduce Dr. Kuntz, who is the
4 Director of CDAC, who will be discussing specific clinical
5 topics relevant to the START trial.

6 DR. KUNTZ: My name is Rick Kuntz. I'm an
7 interventional cardiologist at Brigham and Women's Hospital.
8 I'm also the Director of the Academic Contract Research
9 Organization at Harvard Medical School which conducted this
10 trial.

11 I don't have any equity in this company or any
12 other medical device or drug company, and I'm not being paid
13 for my presentation today.

14 Next slide?

15 I'd like to focus on two issues here that I think
16 require further explanation. One is the clinical impact of
17 the minor device malfunctions, and the second is an analysis
18 of the edge effect.

19 Next slide?

20 Starting with the MDM analysis, as Dr. Speiser has
21 shown, 87 percent of those patients classified with an MDM
22 had a radiation-related subcategory of either source drift
23 or a prolonged transit time greater than 5 seconds. The
24 remainder of the MDMs, such as inability to deliver the
25 delivery catheters and others, did not deal with radiation

1 issues.

2 Next slide?

3 We attempted to try to determine whether the
4 classification of an MDM, as written in the protocol, had
5 any substantial clinical consequences. So the first thing
6 that we did was look at one of the primary endpoints, that
7 is, major adverse cardiac event rates at 240 days, and
8 compare them between those patients classified as having
9 either drift or increased transit time compared to those
10 patients who fit within the guidelines.

11 This is an analysis of the placebo patients, and
12 remember, nobody knew if they were placebo or radiation, so
13 drift or increased transit time could occur on both sides.

14 If we look at the placebo group, there was no
15 substantial or significant difference between the incidence
16 of MACE or between those patients classified as having drift
17 or transit time compared to those without.

18 We then evaluated the same endpoint for patients
19 assigned to the active arm. Again, we found no significant
20 difference in the incidence of MACE between those with drift
21 and transit time MDMs versus those without. However, you'll
22 notice there that there may be trend, that is, that there
23 was a tendency for an increased estimate rate associated
24 with patients with drift or transit time MDMs. The
25 difference that was seen here did not reach statistical

1 significance, but we tried to understand what might explain
2 the trend.

3 The first thing we evaluated was the component of
4 the MACE that we felt was the actual safety issue, that is,
5 the occurrence of myocardial infarction or death. What we
6 found was that there was absolutely no evidence of an
7 increased risk of MI or death for patients who had been
8 classified as having drift or increased transit time. There
9 were zero events in those classifications compared to those
10 without drift or those that fit the criteria. So 3 percent
11 of these MIs and deaths overall occurred in patients who met
12 the criteria and protocol and zero occurred in those who
13 were classified as having an MDM. Clearly, there is no
14 significant difference here and no evidence that the drift
15 or transit time issue increased the risk of myocardial
16 infarction or deaths.

17 Next slide?

18 So what explains the difference that we see in
19 this? And the difference is explained by the incidence of
20 target vessel revascularization or the other component of
21 MACE. Now, whether that difference is going to be
22 substantial or clinically important or not leads us to
23 conclude that possibly if there is a difference, it was a
24 difference in the efficacy endpoint of restenosis, and
25 that's it, not in the safety issue of death or myocardial

1 infarction.

2 We can postulate that possibly drift or increased
3 transit time might reduce the overall deliverable radiation
4 effect. It may be--it's something that you can imagine, and
5 maybe that explained the difference in reduction in
6 efficacy. Then, again, there may be no difference at all
7 because the p-value for this is 0.11. But the important
8 point is that when we look at the incidence of major adverse
9 cardiac event rates between these two subsequent
10 applications, we found no evidence that we exposed the
11 patients to an increased risk of a safety endpoint of death
12 and myocardial infarction.

13 Next slide?

14 In order to evaluate whether restenosis
15 differences could have occurred because of the tendency for
16 increased target vessel revascularization was evident in
17 that initial analysis, we looked at another measure of
18 restenosis using quantitative angiography. And in this
19 case, we compared the overall restenosis rates of patients
20 using quantitative angiography between those classified
21 without a source drift or transit increase and those with
22 the MDM. We found absolutely no difference in that measure
23 of restenosis overall.

24 So our conclusion, therefore, is that source drift
25 and source transit greater than 5 seconds were prospectively

1 collected and possibly could have been very conservative
2 measures and potentially ambitious goals, but were still
3 collected as events and identified as primary minor device
4 malfunctions.

5 The clinical impact of this classification of MDMS
6 demonstrated no statistical difference in any safety measure
7 of the Beta-cath system in the treatment of in-stent
8 restenosis; and, moreover, the sponsor has proposed minor
9 device modifications and training measures to reduce the
10 occurrence of these sort of drifts and transit increases.

11 The next issue I'd like to address is the issue of
12 edge effect. As you recall, Dr. Popma's presentation
13 demonstrated that all measures of restenosis showed a
14 significant benefit for patients exposed to radiation
15 therapy compared to those exposed to placebo. But the
16 impact was greatest in the analysis that was confined to the
17 stent area compared to the analysis that was liberated to
18 the wide analysis segment. One of the questions is that if
19 we start to see more restenosis on a wider measure of
20 restenosis compared to one confined within the stent, is
21 there any activity occurring at the areas outside the stent
22 that are measured by the analysis? And this has been
23 brought up by many investigators in the past as a potential
24 problem associated with radiation therapy called edge
25 effect. And there have been a lot of postulated ideas that

1 the radiation therapy can cause narrowing itself in areas
2 that aren't treated, and there may be areas termed
3 geographic miss in which the balloon injures the artery and
4 there's a fall-off of radiation that may cause increased
5 narrowing.

6 So we were very curious to understand whether this
7 analysis was, in fact, accurately depicting edge problems or
8 it was an artifactual result from the limitations of this
9 analysis and may represent nothing at all.

10 Next slide?

11 So in order to approach this, we have to review
12 again how the initial analysis is done. As Dr. Popma showed
13 you earlier, this is a conventional analysis that's been
14 used for all coronary treatments in the last 15 or 20 years;
15 that is, we tend to measure restenosis based on the location
16 of the minimum lumen diameter, and that's how conventional
17 angioplasty has been evaluated. So that if we look at
18 patients and measure restenosis defined by narrowing within
19 the stent segment, a quantitative angiographic technique
20 used by Dr. Lansky at the core laboratory would identify the
21 area of most narrowing and would tell us what that narrowing
22 is and tell us where it's located. If that narrowing is
23 greater than 50 percent of some reference value, we call
24 that restenosis.

25 The other analysis that we can do is to actually

1 measure the minimum lumen diameter across a very wide area.
2 In this case, we called that the analysis segment. Again,
3 we would detect just one single minimum lumen, and we would
4 tell where the location is.

5 What we found is that if you measure the minimum
6 lumen diameter, it tended to be located in the stent for the
7 vast majority of patients assigned to placebo. If we looked
8 at the minimum lumen diameter for those assigned to active
9 arm, a fair substantial minimum lumens were located outside
10 the stent and not in the stent. So we attempted to
11 understand whether there was an increased propensity for
12 radiation therapy to cause more narrowings outside or if, in
13 fact, we were just unmasking disease that was already there.

14 Next slide?

15 So Dr. Lansky re-evaluated the data set with a
16 specific analysis looking at edge effect, and what was done
17 here was that the measurement of restenosis was confined to
18 just the edges of the analysis segment and not the stent
19 itself to specifically address this issue. So at the source
20 end, both proximal and distal, an analysis centered on that
21 end going 5 millimeters on each side was performed in which
22 the minimum lumen diameter was defined, both proximal and
23 distal, and the incidence of restenosis was measured between
24 the two groups. And what she found was that there was no
25 significant difference between restenosis measured both

1 proximally or distally, depending on the assignment of
2 patients to placebo or active arm. And I think this
3 definitively shows that there was no significant increase in
4 edge narrowing seen in patients assigned to radiation
5 therapy.

6 Next slide?

7 So the question is how did we get that disparity
8 between our initial analysis, which showed a 66-percent
9 reduction when we measured restenosis within the stent,
10 compared to 36 percent when we measured restenosis within
11 the analysis segment. We think we can understand how that
12 artifact might have occurred.

13 In a normal vessel that gets treated, most of the
14 narrowing occurs initially within the stent, hence in-stent
15 restenosis. After treatment, we generally clear out the
16 entire lumen with either aggressive balloon angioplasty or
17 debulking device. And then the patient is assigned to
18 either placebo or active therapy with the radiation source
19 train, and then they're followed up six months later to see
20 what happens.

21 Next slide?

22 In patients assigned to placebo, what we found was
23 that there was narrowing that occurred normally on the edges
24 and in the middle, but the vast majority of narrowing and
25 the pattern of narrowing occurred mainly within the stent.

1 And this has been seen in other in-stent trials and even in
2 *de novo* stent lesions; that is, the response of restenosis
3 is a little bit on the edges, and the majority of it occurs
4 within the stent. So when we do analysis of where the
5 minimum lumen diameter is, it tends to occur mainly within
6 the stent for patients assigned to placebo.

7 In the radiation therapy group, we get the same
8 degree of narrowing on the edges, but a profound reduction
9 in restenosis in the middle because that's where the
10 targeted therapy was. This may be familiar to radiation
11 oncologists with the concept of central control, where once
12 you take care of the central control of a tumor, you may
13 start to realize peripheral disease starts to show up. Very
14 similar concept here in that what we see here is that an
15 effective therapy may unmask the already present narrowing
16 that has occurred is identical between the two groups.

17 So as Dr. Lansky showed us, the two narrowings and
18 the size were the same, but if you have an effective therapy
19 in the middle, you may unmask the occasional case of
20 patients who have 50-percent narrowing on the sides. So in
21 that analysis you will see that some of the restenosis
22 occurs in the stent and some of it occurs on the side when
23 you have an effective therapy in the middle.

24 Next slide?

25 So the conclusion of the edge analysis is that

1 significant treatment effects seen between the analysis and
2 stent segments was potentially due in part to injury of the
3 radiation therapy or the masking of progression of normal
4 restenosis seen with any cardiac treatment unmasked by
5 effective radiation therapy.

6 Now, given both of those two possibilities, the
7 data really supports the masking issue of progression of
8 disease rather than an induction of narrowing occurring by
9 radiation therapy.

10 I think now I'd like to turn the podium over to
11 Dr. Popma to make some conclusionary statements.

12 DR. POPMA: Go to the next slide.

13 I think in very rapid sequence you've heard that
14 there is a medical need to treat in-stent restenosis. It is
15 a problem for patients. It's a problem in all of our
16 clinical practices.

17 What we've demonstrated in this trial are the
18 following: First of all, the START trial was the largest
19 trial of its type that used randomization, used triple-
20 masking, used placebo-controlled to demonstrate its
21 conclusions. And the summary of what we've heard so far was
22 that the pre-specified hypotheses were all achieved with
23 statistical significance. Target vessel failure was reduced
24 by 31 percent. Major adverse clinical events were reduced
25 by 31 percent. Target vessel revascularization was reduced

1 by 34 percent, and target lesion revascularization was
2 reduced by 42 percent.

3 On the angiographic analysis, we saw that the pre-
4 specified restenosis hypotheses were achieved with
5 statistical significance. Within the stented segment, there
6 was a 66-percent reduction in angiographic restenosis. And
7 with the analysis segment, there was a 36-percent reduction
8 in angiographic restenosis, both highly significant.

9 We also learned that the treatment with Strontium-
10 90 in the START trial was safe and that there were no
11 differences in the occurrence of death or myocardial
12 infarction between Strontium-90-treated patients and
13 placebo-treated patients.

14 There were no differences in late thromboses,
15 there were no differences in total occlusions, and there
16 were no differences in the occurrence of new aneurysm
17 formation.

18 So what we would conclude from this trial is that
19 the Beta-cath system has been shown to be safe and effective
20 for the treatment of in-stent restenosis.

21 MR. GREEN: That concludes the presentation by
22 Novoste Corporation for the Beta-cath system.

23 ACTING CHAIRPERSON TRACY: Thank you. I'd like to
24 ask the panel members if they have any brief clarifying
25 questions they want to ask. We'll have much more time later

1 for discussion. Anybody? Any brief question?

2 [No response.]

3 ACTING CHAIRPERSON TRACY: Then we'll move on to
4 the FDA presentation.

5 MS. MOYNAHAN: While FDA is setting up, let me
6 mention that we've been joined by Dr. Robert Ayers of the
7 Nuclear Regulatory Commission. Dr. Ayers will be
8 participating as a guest today, and he is identified as a
9 coregulator of the use and licensing of the Novoste Beta-
10 cath system.

11 MS. PETERS: Good morning. My name is Kim Peters,
12 and I'm a biomedical engineer in the Interventional
13 Cardiology Branch of the Office of Device Evaluation. I am
14 also the leader reviewer for the Novoste Beta-cath system
15 PMA submission P000018. Today, Dr. Bram Zuckerman, the
16 medical officer for this submission, and I will present the
17 FDA summary for the Beta-cath system.

18 This presentation will identify the FDA review
19 team members, provide a brief summary of the device
20 description, provide a summary of the non-clinical tests
21 conducted on the Beta-cath system, provide a summary of the
22 clinical investigation of the Beta-cath system, and identify
23 the FDA questions for the panel.

24 Members of the FDA review team include Dr. Sabu
25 Subramanian and Dr. Bram Zuckerman, both from the Office of

1 Device Evaluation; Mr. Tom Heaton from the Office of Science
2 and Technology; Mr. Gary Kamer from the Office of
3 Surveillance and Biometrics; and Ms. Marianne Linde from the
4 Office of Compliance.

5 As described during the sponsor's presentation,
6 the Beta-cath system is comprised of the Beta-cath delivery
7 catheter, the transfer device, the source train, and system
8 accessories. The Beta-cath deliver catheter is a sterile
9 single-use catheter that provides the path through which the
10 source train is delivered to and retrieved from the
11 treatment site. The catheter includes three lumens to allow
12 for the passage of the guide wire, source train, and
13 hydraulic fluid. At the distal end of the catheter, the
14 source train and hydraulic fluid lumens are closed, while
15 the guide wire lumen remains open to the vasculature. The
16 distal end of the catheter also features two radiopaque
17 markers that define the treatment zone of the catheter.

18 The transfer device stores and shields the source
19 train when not in use and controls the hydraulic delivery
20 and return of the source train during the treatment
21 procedure. The transfer device features a series of
22 components intended to protect the health care workers and
23 patient from unnecessary radiation exposure, either by
24 shielding the source train or maintaining proper position of
25 the source train. The transfer device also features a

1 series of components intended to regulate, direct, and
2 manage the hydraulic fluid that controls the delivery and
3 return of the source train.

4 The source train consists of a series of
5 individual, cylindrical, sealed radioactive sources with an
6 inactive gold marker at each end of the train. The
7 radioactive sources are Strontium-90 seeds encapsulated in
8 stainless steel.

9 The system accessories include a procedure
10 accessory pack, an emergency storage container, a response
11 kit, and a medical physicist kit. These components are
12 intended to facilitate the operation of the system during
13 the clinical procedure, permit temporary storage of the
14 system in the event of a disrupted procedure, and facilitate
15 handling of the source train if located outside the system.

16 The clinical investigation for the Beta-cath
17 system, the START trial, was conducted with the Alpha III
18 and Alpha IV models of the transfer device. Approximately
19 83 percent of the clinical data was obtained using the Alpha
20 II model, with the remainder of the data being obtained
21 using the Alpha IV model. Modifications were made to the
22 transfer device in response to reports of device
23 malfunctions during the clinical investigation and to
24 improve the system function.

25 The main difference between the Alpha III and

1 Alpha IV models is the addition of the LED pressure
2 indicators. The pressure indicators provide feedback to the
3 user regarding the pressure necessary to remain the source
4 train at the treatment zone and to send or return the source
5 train to and from the transfer device. The indicators also
6 advise the user when excessive pressure is being
7 administered.

8 The Alpha IV revision two model of the transfer
9 device is the subject of the PMA submission. No clinical
10 data was obtained using this model. As noted in the FDA's
11 summary, the Alpha IV revision two model mainly includes
12 refinement to some of the electronic circuitry and
13 indicators. The Alpha IV revision two model also includes a
14 modification to prevent the transfer device gate from
15 inadvertently being locked prior to the delivery of the
16 source train. FDA believes that these modifications can be
17 evaluated through bench testing.

18 Optional accessories of the Beta-cath system
19 include an introducer sheath and an extension tubing set.
20 Excessive hemostasis valve tightening can restrict the
21 movement of the sources in the Beta-cath system. An
22 optional component, the arrow introducer sheath, may be used
23 to increase the resistance of the catheter to collapse when
24 compressed with the hemostasis valve. The optional
25 extension tubing set provides an additional fluid management

1 system for use during the clinical procedure by allowing two
2 control syringes to be connected to the Beta-cath system.

3 During the START trial, both the 30-millimeter and
4 40-millimeter delivery catheters and source trains were
5 used. The difference between the 30-millimeter and 40-
6 millimeter delivery catheters is the marker spacing at the
7 distal end of the catheters. The spacing identifies the
8 treatment zone of the catheter. The difference between the
9 30-millimeter and 40-millimeter source train is the number
10 of active source seeds. The 30-millimeter source train
11 includes 12 seeds, and the 40-millimeter source train
12 includes 16 seeds.

13 Due to the limited clinical data available for the
14 40-millimeter model, only the 30-millimeter delivery
15 catheter and source train are subject of the PMA.

16 A series of *in vitro* tests were performed to
17 evaluate the mechanical integrity and function of the Beta-
18 cath system and each of the individual components. FDA is
19 working with the sponsor to resolve questions regarding this
20 testing information. The delivery catheter is the only
21 patient-contacting component of the Beta-cath system.
22 Biocompatibility testing completed in accordance with the
23 ISO Standard 10993 demonstrated that the catheter is non-
24 toxic and non-hemolytic.

25 Electrical safety, battery, and electrode magnetic

1 compatibility tests were conducted in accordance with
2 applicable voluntary standards. All test requirements were
3 met.

4 As discussed in the FDA summary, the sponsor has
5 conducted two animal studies using both oversize stent
6 injury and balloon overstretch injury pig models. The
7 results of the animal tests show no difference in restenosis
8 between the Beta-cath system and the control.

9 With regard to the source dosimetry, FDA is
10 working with the sponsor to resolve questions with the
11 dosimetry information and dosimetry labeling
12 recommendations.

13 The sponsor has provided data from three clinical
14 studies: the Beta Energy Restenosis Trial, the Beta
15 Radiation in Europe Trial, and the Stents in Radiation
16 Therapy Trial. The Beta Energy Restenosis Trial was a U.S.
17 feasibility study evaluating the use of beta radiation
18 following PTCA and *de novo* lesions. Eighty-three subjects
19 were enrolled in the study. The Beta Radiation in Europe
20 Trial is a multi-center, non-randomized registry that is
21 studying the use of the Beta-cath system following PTCA or
22 stenting of *de novo* lesions. One hundred fifty patients
23 were enrolled in the study. Summaries of these two trials
24 are provided in the panel pack.

25 The Stents in Radiation Therapy Trial is the

1 pivotal study for the evaluation of safety and effectiveness
2 of the Beta-cath system. Dr. Zuckerman will provide an
3 overview of the trial design and a summary of the results.

4 DR. ZUCKERMAN: Good morning. My name is Bram
5 Zuckerman. I'm a medical officer cardiologist with the Food
6 and Drug Administration.

7 The START trial is the key data set for
8 consideration of this PMA. The sponsor has previously
9 outlined the major elements of the protocol and shown key
10 results. The agency would like to discuss several aspects
11 of this trial prior to presentation of the panel questions.

12 Next slide, please?

13 The START trial was a well-designed trial. A
14 large number of patients were randomized to beta radiation
15 or placebo treatment. Patients, investigators, and core
16 labs were blinded to treatment assignment.

17 Next slide, please?

18 Inclusion criteria indicated that a patient needed
19 a symptomatic case of in-stent restenosis with a reference
20 vessel diameter between 2.7 and 4 millimeters. Visual
21 estimation was performed of reference vessel diameter
22 because this mimics real-world clinical practice. There are
23 few sites in the United States that routinely use online
24 quantitative coronary angiography, QCA, or intravascular
25 ultrasound. But as expected, we saw the discrepancy between

1 our visual reference vessel diameter results and the QCA
2 results.

3 For example, the mean QCA result reported from the
4 core lab for all vessels was 2.76 millimeters. The other
5 point in interpretation of this trial is as noted by Ms.
6 Peters: 95 percent of the data pertains to the 30-
7 millimeter source train.

8 Next slide, please?

9 Vascular brachytherapy represents a new technology
10 for treatment of in-stent restenosis with an unclear
11 risk/benefit profile. As such, a superiority hypothesis was
12 chosen with a primary clinical endpoint: 8-month target
13 vessel failure.

14 Target vessel failure is a conservative endpoint
15 that includes death, non-fatal myocardial infarction, and
16 target vessel revascularization. The angiographic and
17 ultrasound data provided in this panel report should,
18 therefore, be viewed as supporting data that helps to
19 mechanistically explain the effects of vascular
20 brachytherapy.

21 Next slide, please?

22 On this slide, we have acute results presented.
23 You've seen much of this before presented by the sponsor.
24 The only difference with these slides and the next two are
25 that we will have also the 95-percent confidence interval of

1 the differences presented.

2 The key point for acute results was that a
3 respectable post-procedure percent diameter stenosis was
4 obtained, and high rates were reported for both device
5 success and procedure success.

6 Next slide, please?

7 There were nine cases of device failure. These
8 cases have been individually reviewed in your panel pack
9 report.

10 Next slide, please?

11 Please note, however, the definitions used for
12 device failure and procedure success--I'm sorry, for device
13 success and procedure success. Device success was defined
14 as successful placement of the Beta-cath system, and
15 procedure success was defined as a post-procedure percent
16 diameter stenosis less than 50 percent without the
17 occurrence of major adverse cardiac events during the
18 hospitalization. Hence, such problems as initial device
19 failure, minor device malfunction, in all cases where the
20 bail-out box was used emergency, would not necessarily be
21 captured in those two preceding definitions. A balanced
22 assessment of device performance needs to include these
23 variables and results as well as device success and
24 procedure success results.

25 Next slide, please?

1 Eight-month safety results are shown on this
2 slide. At 8 months there was no difference in the rates of
3 death, myocardial infarction, stent thrombosis, site
4 thrombosis, total occlusions, or aneurysms. Two points need
5 to be noted.

6 Firstly, a minority of patients, 20 percent, were
7 restented during this trial. The restented population may
8 be the population at greatest risk for long-term safety
9 problems.

10 The second point is that these are 8-month safety
11 results. As previously noted and as noted by the asterisk
12 on the bottom of the slide, there was already one stent
13 thrombosis reported at greater than 240 days in the beta
14 radiation arm.

15 Next slide, please?

16 Eight-month effectiveness results are reported
17 here. The primary endpoint, target vessel failure, was
18 reduced at 8 months by beta radiation treatment. This was a
19 robust result. Multiple other clinical and angiographic
20 markers of restenosis were reduced by beta radiation
21 treatment, as noted in your panel pack and on this slide.

22 Next slide, please?

23 So, in conclusion, the primary endpoint, target
24 vessel failure, as well as selected clinical and
25 angiographic endpoints, were all reduced by beta radiation

1 treatment. There was no difference at 8 months in the
2 incidence of death, myocardial infarction, stent thrombosis
3 or total occlusion. Device-related malfunctions were
4 observed.

5 MS. PETERS: FDA would like to obtain panel input
6 on the following questions:

7 The original START protocol suggested that the
8 institutional standard of care for antiplatelet therapy
9 after source treatment be utilized for patients who were
10 restented or received PTCA. This regimen was modified based
11 on recommendations from the Data Safety Monitoring Board. A
12 report of the antiplatelet therapy usage during the START
13 trial is provided in the addendum to the START clinical
14 report on page 3. No incidents of stent thrombosis were
15 reported during the START trial.

16 Question 1: Based on this information, please
17 discuss your recommendations for the antiplatelet therapy
18 for patients who receive a new stent and for patients who do
19 not receive a new stent.

20 Table 31 of the START clinical report and the
21 addendum to the START clinical report on pages 13 through 35
22 identify the device failures and malfunctions that occurred
23 during this study.

24 Question 2: Please discuss the clinical
25 importance of the device failure and malfunction events and

1 the evaluation of the safety and effectiveness of the Beta-
2 cath system.

3 As demonstrated by the results included in Table 1
4 of the START clinical report, the incidence of the primary
5 endpoint, target vessel failure, was significantly lower at
6 8 months for the treatment arm compared to the placebo. The
7 incidence of target vessel revascularization, target lesion
8 revascularization, and major cardiac adverse events were
9 also significantly lower over the 8-month follow-up period
10 for the treatment arm compared to the placebo. No incidents
11 of stent thrombosis were detected in the treatment arm, and
12 the frequency of total occlusions was comparable between the
13 treatment and placebo arms.

14 Question 3: Please discuss whether you believe
15 the probable clinical benefit of the radiation treatment
16 outweighs the probable risk of death, myocardial infarction,
17 late total occlusion, and late stent thrombosis posed by the
18 device in the intended patient population.

19 One aspect of the premarket evaluation of a new
20 product is the review of its labeling. The labeling must
21 indicate which patients are appropriate for treatment,
22 identify the product's potential adverse events, and explain
23 how the product should be used to maximize benefits and
24 minimize adverse effects. Please address the following
25 questions regarding the product labeling.

1 Question 4-A: Please comment on the indications
2 for use section as to whether it identifies the appropriate
3 patient population for the treatment with the device.

4 Question 4-B: Please comment on the contra-
5 indications section as to whether it identifies all
6 conditions under which the device should not be used because
7 the risk of use clearly outweighs any possible benefit.

8 Question 4-C: Please comment on the warnings and
9 precautions section as to whether it identifies all
10 potential hazards regarding device use.

11 Question 4-D: Please discuss whether any
12 improvements could be made to the labeling to help minimize
13 the occurrence of device failures and malfunctions as
14 discussed under Question 2.

15 Question 4-E: Please comment on the remainder of
16 the device labeling as to whether it adequately describes
17 how the device should be used to maximize benefits and
18 minimize adverse events.

19 Question 4-F: Does the panel have any other
20 recommendations regarding the labeling of the device?

21 A summary of the physician training program has
22 been provided in Section E of the panel pack and in the
23 addendum to the START clinical report on pages 18 through
24 25.

25 Question 5-A: Please discuss any improvements

1 that could be made to the training program to help minimize
2 the occurrence of device failures and malfunctions as
3 discussed under Question 2.

4 Question 5-B: Please identify any other important
5 elements that should be contained in a physicians' training
6 program for this device.

7 The panel pack includes the available one-year
8 data from the START trial, the available one- to four-year
9 data from the BERT feasibility trial, and the available data
10 from the BRE European trial.

11 Question 6: Based on the clinical data provided
12 in the panel pack, do you believe that additional clinical
13 follow-up data or post-market studies are necessary to
14 evaluate the chronic effects of intravascular radiation
15 administration? If so, how long should patients be
16 followed, and what endpoints and adverse events should be
17 measured?

18 This concludes the FDA's summary presentation.

19 ACTING CHAIRPERSON TRACY: Thank you.

20 At this point we'll move to the open committee
21 discussion, and I'd like to ask Dr. Simmons to begin the
22 discussion with his review, and we'll go around the table
23 after that. And I'd just remind the panel members to
24 restate their names and speak into the microphone when
25 they're asking their questions or making their comments.

1 Dr. Simmons?

2 DR. SIMMONS: Thanks. Well, it's a very nice
3 presentation by the sponsor and the FDA. I have a few
4 questions.

5 The data that the sponsor presented I think
6 presented a better clinical outcome than the data that's in
7 the panel pack or that the FDA presented. I mean, the 8-
8 month stent segment restenosis rate was 14 percent versus 14
9 percent, which is a 27-percent reduction. However, the
10 target lesion revascularization at 240 days was 86 versus
11 76, and you were presenting 31 percent. So it's actually
12 about a 9- to 10-percent reduction in total vessel failure
13 at the 8 months--is that right?--as opposed to the 30
14 percent that I saw on your slides? There's quite a
15 difference between 30 percent versus 9 percent.

16 DR. BAILEY: Is it possible he's talking about a
17 relative reduction?

18 DR. SIMMONS: That's what I'm interested to see.

19 DR. POPMA: I wonder if we could just go back to
20 our slides and the presentation just very briefly, if that
21 will help find out where the discrepancies are.

22 I should also note that it is difficult to take
23 numbers out of the event-free survival curves and then put
24 them back into the rates that are measured. The 240 days is
25 absolutely accurate, but some of the data that's in the pack

1 is extended out further than that and before that. So I
2 want to make sure of the 240-day endpoint for event-free
3 survival curves.

4 Let's go back and discuss those. If you can put
5 up, Richard, just the graph that has the four reductions of
6 TVF, MACE, TVR, and TLR, the one slide with four graphs.

7 ACTING CHAIRPERSON TRACY: Could you also
8 introduce yourself for the--

9 DR. POPMA: Sorry. I'm Jeffrey Popma.

10 Now, we can go by these one by one, if you like.
11 I think that these absolute rates have been relatively
12 consistent and should be consistent to the rates that are
13 reported in the panel pack.

14 DR. SIMMONS: Well, they're not, actually. Why
15 don't you go to page--let's just get on the same page here.
16 Go to your panel pack, page 414, and you've got TVR-free at
17 240 days, 81.4 percent in the treated group versus 72
18 percent in the placebo group, so the difference is only 9
19 percent.

20 DR. POPMA: I got it. I'm going to defer to Rick
21 Kuntz, who did the statistical analysis for this.

22 DR. SIMMONS: And that's different than 32
23 percent.

24 DR. KUNTZ: The numbers that you're referring to
25 is the absolute difference. This is a relative difference.

1 If you look here, the difference in TVF there is 26 minus
2 18, which is 8 percent. And then you're seeing a difference
3 in the MACE-free or TVF-free of 9 percent, which is
4 consistent with the differences between sensoring and non-
5 sensoring survival analysis versus a discrete analysis.
6 Maybe--Dr. Bailey, do you understand what I'm talking about
7 there?

8 DR. BAILEY: In other words, 9 percent is 31
9 percent of 26.

10 DR. KUNTZ: Right. So the difference of 9 percent
11 in the TVF and 8 percent in the event-free survival or vice
12 versa is pretty typical when you're using a sensoring
13 analysis for survival versus one that is a discrete analysis
14 at 240 days. To me that's pretty clear. I'm not quite
15 sure--so 31 percent refers to the relative difference you
16 see here, but the 9 percent you're referring to is the
17 absolute difference. The absolute difference here is 8
18 percent as well.

19 DR. SIMMONS: Okay. So what we're actually
20 talking about is an 8-percent improvement clinically.

21 DR. KUNTZ: Right.

22 DR. SIMMONS: A 9-percent improvement clinically.

23 DR. KUNTZ: Right. Again, if we refer to the MACE
24 or the TVF here, the absolute differences there are in the 8
25 to 9 percent range.

1 DR. SIMMONS: Okay.

2 DR. KUNTZ: Which is similar to the event-free
3 survival we're seeing there.

4 DR. SIMMONS: All right. You know, I understand
5 that most of the device failures did not result in any
6 complications to the patient and any real damage. But, you
7 know, the physicians doing your clinical trials, they're
8 more skilled physicians. They've got more back-up. They've
9 got more interest in what's going on. But we're still
10 talking, I think, a fairly amazing, almost 20-percent
11 incidence of some device failure, either minor or major.
12 And even though it didn't appear to have any clinical impact
13 in the physicians that were performing the study, I'm just
14 wondering what happens when physicians who do two
15 angioplasties a month and don't have any company
16 representative around, whether that's really going to
17 translate into no complications to the patient.

18 In addition, if we look at page 221, where Dr.
19 Zuckerman analyzed the differences between your number three
20 and your number four revision, the incidence of drift and
21 everything else didn't really seem to be affected by putting
22 the lights on the box and changing your connector.

23 So I'm just wondering if one of the clinicians
24 might address this issue, my concerns.

25 DR. SPEISER: Burton Speiser, radiation

1 oncologist. I think the primary problem with both the drift
2 and the transit time is the lack of observation of the
3 observer. It's a very simple process, either by feel before
4 the LED lights were added, or by watching the LED lights.
5 The primary problem is the training of the individuals such
6 that they pay attention to that, and I don't want to put
7 down my colleagues too much, but it isn't a very hard job
8 really to keep the pressure in there.

9 Part of the training process is to use the device
10 and ensure that they know how to keep the pressure constant,
11 which is, in effect, a very easy process.

12 DR. SIMMONS: I don't know. I mean, these are
13 very motivated, highly skilled people. If they can't do it,
14 do you really expect people in other cath labs who aren't as
15 motivated to be able to do it?

16 DR. SPEISER: I think the problem is primarily
17 with the radiation oncologist who is in the cath lab for the
18 first time, feels in a foreign territory, and I think what
19 is necessary is the feeling that they're comfortable there
20 and they know how to use the device. And that will take
21 training. Most radiation oncologists probably have never
22 stepped into a cath lab before, and I think that's why the
23 training and doing the mock procedures is quite important to
24 make sure that each radiation oncologist feels comfortable.

25 If you play with the device, you'll be surprised

1 how simple it really is, so that I do have a little
2 difficulty trying to explain a 20-percent rate of transit
3 problems and drift, when, in fact, if you use the device in
4 a mock session, it's very easy to send the sources out and
5 keep them in station or in place.

6 Now, I know that may not be answering the
7 question, but I do have difficulty in understanding why many
8 of my colleagues had difficulty with it.

9 MR. GREEN: I'd just like to also add to that that
10 the trial did actually include 50 centers, and we believe
11 that was--we tried at least to develop that trial with those
12 centers so that we included both, if you will, the normal
13 usual suspects or the normal usual trial centers that would
14 be seen in these type of trials so we had that component of
15 an understanding of clinical trials, as well as the regional
16 and, if you will, everyday hospitals that would use such a
17 system.

18 We did learn that if you look in your panel pack
19 at the section on the device observations, for instance,
20 with regard to what we call the manual removal procedure or
21 what we term there the bail-out procedure, that as time went
22 on and as enrollment went up, there was a decrease in the
23 rate of the use of that procedure. So, therefore, as people
24 became accustomed to the cath labs, people became--we had
25 improved training from experiences in the trial, we were

1 able to modify not only the device but the user's
2 understanding of the procedure and how they could apply that
3 device and that procedure to obtain the results that we
4 found in the trial.

5 DR. POPMA: If I could just maybe help this a
6 little bit with just a clinical perspective, as the
7 cardiologist, of course, we're not responsible for moving
8 the radiation source, but we actually do have a lot of
9 experience with catheters. And some of the issues with
10 respect to the transit time may well have related to very
11 simple things like having the touie borst (ph) too tightly
12 ratcheted down so that they couldn't move back and forth.

13 We learned these things as we went through, and I
14 can tell you the procedure that I performed today is very
15 different than the procedure that I performed a year ago
16 because I have much more attention to having the touie borst
17 open, it's loose, the catheter being straight, and helping
18 the radiation oncologist deliver the sources more quickly,
19 as well as the fact that we stay on fluoro a lot more to
20 make sure that there's no source drift.

21 All of this is covered in training, and I think
22 it's very important to emphasize that a lot of the things
23 that you're discussing in 20 percent, none of us would want
24 to have that shown prospectively. But they are covered in
25 training, and I do think that we all need--we've all learned

1 from this, and I think that we perform a very much better
2 procedure now than we did before. We just have to, you
3 know, show that, I think, with data.

4 DR. SIMMONS: So even though the radiation
5 oncologist delivers the seeds, the cardiologist should be
6 able to see whether this thing is drifting or not, right?

7 DR. POPMA: Absolutely. And we will do--we do
8 that and have done that, and we're very much aware of that
9 now, by stepping frequently on the fluoro pedal. There's
10 nothing that prevents the--

11 DR. SIMMONS: Well, there is something in the
12 manual about how many times, how often it should be
13 observed, and--

14 DR. POPMA: Exactly.

15 DR. SIMMONS: Do you know what those are?

16 DR. POPMA: Well, right now we--you know, as I say
17 our clinical perspective is--I'll let Drew address what the
18 IFU is, but we do that very, very frequently now in the
19 catheterization lab with a much more heightened awareness
20 about the importance of drift. And by doing that, we
21 actually can catch events. But I'll let Drew--Dr. Green
22 discuss what's in the IFU.

23 MR. GREEN: During the clinical procedure, the
24 protocol described using fluoroscopy, I believe, every 10 to
25 15 seconds to observe to see position and location of the

1 source train and correct--the instructions for use included
2 in your panel pack also reflect that.

3 One thing I would like to also point out is that
4 we did--you know, when we went and designed a clinical
5 trial, we didn't have some of the experience, of course, you
6 go to gain from a clinical trial. So we applied what we
7 believed to be conservative estimates based on bench
8 testing, which, as you know, when you go into a clinical
9 trial is your first available information. We prospectively
10 defined what we believed would be some measures we would
11 want to look at, things like source drift, source transit
12 times, et cetera.

13 The panel pack does now describe its modifications
14 to the instructions for use reflecting source transit, for
15 instance, which also goes into the train. Instead of just
16 having recommendations based on the bench testing, we have
17 recommendations that are based on really three things: one,
18 the bench testing that was conducted over the expected
19 possible pressure ranges by different users; two, experience
20 from the START trial, what was actually reported in terms of
21 time to send or return of source train. If you look in the
22 panel pack there, I think you'll find that the actual--when
23 it was reported, it was reported between 5 and 14 seconds.
24 And then, third, we went to our oncologists and our medical
25 physicists that participated in the trial, and we asked them

1 what would be clinically acceptable based on what the device
2 can do and what is necessary to occur during a clinical
3 procedure. And we put all these things together, and we
4 have a recommendation, which is in the panel pack,
5 instructions for use for 15 seconds for source transit. If
6 you don't see the sources arrive where they're expected to
7 arrive within 15 seconds, you should then perform your
8 manual removal procedure, which I think, again, as we
9 pointed out just a minute ago, as time goes by, as
10 enrollment goes up, as people become more familiar with the
11 system, that is on the decline, indicating that the
12 additional training and the experience they gain in the cath
13 lab as a team--because it is a team approach--is beneficial
14 in changing that.

15 DR. SIMMONS: You know, I've got some questions
16 about the training program, which is actually very much at
17 the end of the thing. Should we--I mean, because I think
18 this addresses part of this problem, but maybe we could put
19 that at the end and I'll just keep going with the clinical
20 stuff now. Okay?

21 ACTING CHAIRPERSON TRACY: Okay.

22 DR. SIMMONS: All right. On page 413 of the
23 submission here, it's interesting--and I think you brought
24 this up--that only 69 out of your 476 patients had--which is
25 15 percent--got 60 to 90 days of antiplatelet therapy, and

1 only 13 of 476 patients, 2.7 percent, got anywhere near 90
2 days of antiplatelet therapy. And yet your recommendations
3 in the labeling are going to be for greater than 90 days of
4 antiplatelet therapy.

5 I mean, I realize that a lot of cardiologists
6 would want their patients on antiplatelet therapy for maybe
7 other reasons. A lot of us keep them on them anyway. But,
8 I mean, is it really necessary that we put in the labeling--
9 I mean, you didn't do it, and you didn't have a problem. So
10 where did the 90 days come from?

11 MR. GREEN: As Dr. Popma presented in his
12 presentation earlier, initially at the beginning of the
13 trial, the initiation of the trial, we had the
14 recommendation for physician discretion. As the trial went
15 on, we had information from the Data Safety Monitoring
16 Board, from the Beta-cath system trial, a *de novo* and
17 restenotic lesion trial that suggested that it may be
18 beneficial to patients to have extended antiplatelet therapy
19 if they receive a new stent. And the recommendation was
20 made by them and submitted to the FDA and the IEE to make
21 that change to the protocol, minimum 90 days antiplatelet
22 therapy for patients receiving a new stent.

23 Therefore, at the end of the trial, we carried
24 that recommendation over into our labeling because it did
25 define, if you will, the experience that we attempted or we

1 implemented to investigate in the trial, and it was a
2 recommendation of the trial. So it's approved protocol and
3 we carried it over.

4 Dr. Benot (ph) is our medical director at Novoste
5 Corporation and an interventional cardiologist at Montreal
6 Heart Institute, and I think he'd like to add something
7 about antiplatelet therapy.

8 DR. BENOT: I think, as was presented either by us
9 or by the FDA, the data that we have concern that type of
10 antiplatelet therapy, that type of adjunct antiplatelet
11 therapy. We have not studied anything else. Why we have
12 the movement from one type of antiplatelet therapy was left
13 to the discretion of the physician at the beginning, because
14 when we start that study in September '98, we don't know, we
15 have not the knowledge of any event related to late stent
16 thrombosis and radiation.

17 At this time we have already the Beta-cath trial
18 in process, which is a different indication, is a treatment
19 of de novo lesion with better radiation. This trial, the
20 Beta-cath trial, has two arms: a PTCA arm, a balloon-only
21 arm, and a stent arm. In the stent arm, we finally find out
22 by the end of October '98--the DSM Committee chaired by Tom
23 Ryan come to us and report some of the complications related
24 to the late stent thrombosis. At this time we started to
25 apply longer adjunct antiplatelet therapy, and at this time

1 we propose two months. By March '99 we were secure, we are
2 sure that the problem was still there, and at this time was
3 implemented the minimum of three months of adjunct
4 antiplatelet therapy.

5 That's the detail we have from the Beta-cath
6 trial. We have never had a problem as reported in the START
7 trial for that as medical officials, and we discussed that,
8 we translate the data we learn from de novo stent and
9 radiation to in-stent restenosis and radiation and apply to
10 ask our investigator to prescribe a minimum of three months
11 when they implant a new stent.

12 DR. SIMMONS: It's just it's interesting that you
13 didn't have any problems in this study and only 2.7 percent
14 of your patients had, you know, anything close to ninety--
15 but I guess there's no harm in putting--

16 DR. BENOT: As the data that we have--and I can
17 have the report from the statistician from the Beta-cath
18 trial, Stuart Pocock, we put on the letter from Dr. Stuart
19 Pocock, and on the Beta-cath system trial, again, different
20 education, de novo lesion, but using a stent as the arm
21 differentiating with the balloon-only. And based on the
22 analysis and listing of the current interim data, the
23 incidence rate of late stent thrombosis, Q wave and non-Q
24 wave MI are all satisfactorily much lower in the patient
25 first randomized and treated in the provisional stent branch

1 following the protocol amendment, which is a minimum of 90
2 days of adjunct antiplatelet therapy (?) .

3 These findings are based on 492 patients
4 randomized in the provisional stent branch before the
5 protocol amendment with a median(?) follow-up of 18 months,
6 and, further, 449 patients randomized in the provisional
7 stent branch after the protocol amendment, which is the
8 minimum 3 months of antiplatelet therapy, and that with a
9 median follow-up of 9 months. That's the detail we have.
10 We have no other data than that to explain why the level of
11 the protocol of adjunct antiplatelet therapy.

12 DR. POPMA: I appreciate your letting my
13 colleagues address the background behind that, but as the
14 principal investigator of the trial, I'm comfortable with
15 the statement of at least 90 days, because we do feel that
16 we only had 50 patients in the study who received new stents
17 and radiation therapy, approximately. So to make a
18 definitive statement that there were no subacute stent
19 thrombosis within the first 242 days, there may be some
20 broad confidence intervals to that statement. In addition,
21 we do have this patient at 244 days that had an episode that
22 could be very consistent with a subacute thrombosis event.

23 So I think I am comfortable for those reasons in
24 saying a minimum of 90 days of antiplatelet therapy.

25 DR. SIMMONS: Now that you are up there--

1 [Laughter.]

2 DR. SIMMONS: So the new stents were discouraged.

3 DR. POPMA: That's correct.

4 DR. SIMMONS: But 20-some percent got new stents.

5 DR. POPMA: I think that's a good point. Let me
6 just again walk through very simplistically exactly how you
7 got into the study and how you got a new stent.

8 First of all, in order to be randomized in the
9 trial, one had to have a successful result. You had to have
10 a 30 percent or lower residual stenosis. So the concept was
11 that in order for you to make the decision for
12 randomization, you didn't want to have a new stent in place.
13 At least we didn't have a new stent in place. So the
14 radiation was then delivered. And time passes during that
15 period of time, and we know some things about stent
16 restenosis.

17 One of the things we've learned from work done at
18 the Washington Hospital Center, Roxanna Mayron and Gary
19 Mintz, is that there is an early recoil that sometimes
20 occurs. The mechanism of treatment for in-stent restenosis,
21 you extrude the tissue outside the stent struts, and then
22 within the first 30 minutes or 40 minutes, there's actually
23 a collapse of that and the tissue comes back within the
24 stent. And there's time that passes as we're delivering the
25 radiation effect.

1 So it's understandable that in some patients the
2 residual stenosis would look a little higher after a time
3 delay than it would be if you just ended the procedure and
4 we called it a successful result.

5 So the new stents that went in went in for two
6 reasons: one was there was a dissection that had to be
7 treated; or, secondly, because there was a residual stenosis
8 that was within the lumen.

9 DR. SIMMONS: Well, wasn't a dissection a contra-
10 indication to giving the radiation therapy in the protocol?

11 DR. POPMA: Yes, that's correct. In order to get
12 into the procedure, one had to have a successful result,
13 which was the absence of dissections. But then, as I say,
14 there is a dynamic change that can occur within the lumen.
15 Sometimes the recoil that occurs from re-extrusion of stent
16 plaque into the vessel wall looks like a dissection. It's a
17 flap that can fall back in. Angiographically it's somewhat
18 difficult to tell those. But what we did is we really
19 lowered what was a prevailing rate of 80 percent new stent
20 use for in-stent restenosis down to 20 percent. And if the
21 truth be known, what we know about the START data now, we'd
22 like to get that even lower. And so some of the things that
23 we would really like to say with this study is that we
24 really want to reserve the use of new stents in the study
25 for bail-out circumstances, some circumstance that happens

1 during the procedure that, after radiation therapy, one has
2 to treat with a stent. That can be a new lesion at a new
3 site or within the site that you're treating initially.

4 DR. SIMMONS: So you would say that these mostly
5 were done--the new stents were put in after the radiation
6 was already given?

7 DR. POPMA: That's correct.

8 DR. SIMMONS: Or where the placebo was given.

9 DR. POPMA: They all were, yes.

10 DR. SIMMONS: And so did--I mean, I guess it could
11 significantly affect the results. Did you look and see were
12 they equally divided on both sides for stents, for the
13 placebo versus the active--

14 DR. POPMA: There was no deleterious effect of
15 radiation on causing more stent in group than another group.

16 DR. SIMMONS: But did one group have more stents
17 than the other group?

18 DR. POPMA: No. They were equally balanced
19 between the two.

20 DR. SIMMONS: How about the breakdown on the
21 diabetic patients? I was especially curious. Did you look
22 at that as far as--

23 DR. POPMA: Diabetic subset? I'm going to let my
24 colleague, Dr. Kuntz, address the subset for diabetes.

25 DR. KUNTZ: Diabetics were evenly distributed, and

1 there was no effect of diabetes on the instance of
2 restenosis in the trial.

3 DR. SIMMONS: Okay. So there was no beneficial
4 effect either.

5 DR. KUNTZ: There was no differential effect.
6 Both groups benefited. It just was no differential effect.
7 That is, the interaction between diabetes and radiation
8 therapy was not positive.

9 DR. SIMMONS: Okay. You know, I guess I'm just
10 going to ask the radiation oncologist--this is exposing my
11 naivete here, I guess, but I just have to ask this. I guess
12 I'm a pessimist by nature and a therapeutic nihilist to a
13 certain extent, and I guess I just don't believe that you
14 can put radiation inside a coronary artery and have a
15 beneficial effect without also having a risk of a negative
16 effect. I mean, there's got to be some downside to this as
17 far as creating aneurysms or scar carcinomas or something.
18 There has got to be a downside.

19 I guess I just want your opinion. You know,
20 what's the downside and how long do we have to wait until we
21 see it?

22 DR. SPEISER: Probably the easiest to dispel is
23 the incidence of cancer from this as an overall problem.
24 About one-tenth of a percent of the total dose is delivered
25 from the Strontium, so that by itself is very insignificant

1 compared to the fluoroscopic dose, which would be a greater
2 concern. The only concern I would have is that it's a very
3 high point dose, so it would be more likely to be a concern
4 as far as increasing fibrosis or aneurysmal formation.

5 At the present time, the data that we do have is
6 the BERT trial that shows neither of those two effects have
7 increased. So the answer is that with the available data,
8 which is very scant, there has not been a late deleterious
9 effect. However, I think most radiation oncologists would
10 agree, because radiation effects are delayed, that we would
11 like to continue looking for it for a longer period of time.

12 DR. SIMMONS: What are talking here? I know some
13 radiation effects, like for lymphomas, can even occur ten
14 years later. Are we looking at something that may all of a
15 sudden show up five or ten years from now with severe
16 scarring in that area?

17 DR. SPEISER: Most late effects, such as scarring
18 or vascular effects, usually occur between 6 and 24 months
19 after completion of radiation. So that I would expect that
20 most of them will show up in that time period. The later
21 effects, as you mentioned, for instance, carcinogenesis, is
22 delayed. Lymphomas are the earliest cancers, about 10 to 20
23 years. Sarcoma is 20 to 30. So that for those we'd have to
24 wait a much longer time. However, I'm not, as a radiation
25 oncologist, concerned about the carcinogenesis, but just the

1 direct immediate effects of the high dose on the vessel
2 wall, and that I anticipate that we should see for the most
3 part between 6 and 24 months from the completion of the
4 procedure.

5 DR. SIMMONS: Okay. Maybe our radiation people
6 will have something more to comment on that. Just one more
7 and then I'll--

8 MR. GREEN: Perhaps, if you'd like, we could have
9 maybe one more opinion from one of our radiation oncologist?

10 DR. SIMMONS: Maybe not. Maybe we'll let our
11 radiation oncologist ask some more questions on that issue
12 since I'm not...

13 As far as your warnings and your contraindications
14 section on your labeling, this study actually did eliminate
15 people with ejection fractions less than 30 percent and it
16 did eliminate people with myocardial infarctions within the
17 last 72 hours. Shouldn't those be--I mean, since those
18 patients weren't studied, I just want to know shouldn't we
19 put some warning or at least some contraindication. I'd
20 like to have your opinion before we discuss it when you're
21 not available to comment.

22 MR. GREEN: When we developed the protocol and we
23 put it together to study this trial, as you do in many
24 trials, there are a lot of things that you put in the trial
25 to try to either limit bias or try to determine what the

1 effect on the patient was. For instance, the one
2 recommendation in the protocol you talked about was
3 myocardial infarction within 72 hours, and that was to be
4 able to delineate the baseline factors for a patient.

5 We are recommending that the patients--the
6 instructions for use reflect what we did in the clinical
7 trial. However, like I said, they were recommendations in
8 the--or exclusion criteria in the protocol that were
9 specifically limited to the ability to evaluate the patients
10 in the follow-up to see if the therapy was effective.

11 DR. SIMMONS: I think that's fine, but for right
12 now what I'd have to say is at least when go to discuss this
13 later on, I would have to say those are things that would
14 have to be added, at least a warning if not a contra-
15 indication, if they aren't there now.

16 DR. POPMA: Respectfully, I'm not a labeling
17 expert, but as a clinician, I don't see that there's any
18 reason to suspect that in a patient who has recurrent
19 refractory in-stent restenosis and ejection fraction plus
20 the 30 percent that this therapy should be contraindicated.
21 And I'd only hope that the trial design construct could be
22 described in the labeling, and then a very careful
23 construction of what was done in the trial and the inclusion
24 sets. But I think at this point in time, I would say from a
25 clinical perspective that I wouldn't know that there'd be

1 data suggesting it should be contraindicated in a patient
2 with an ejection fraction plus the 30 percent.

3 DR. SIMMONS: I've got some other issues on the
4 training session, but maybe we can put those off until
5 everybody else has had a chance.

6 ACTING CHAIRPERSON TRACY: Okay. I'd like to ask
7 our statistician if he has any particular questions, and
8 then we'll go--we'll probably break for lunch before we go
9 around the rest, but if we could have Dr. Bailey ask any
10 questions.

11 DR. BAILEY: Why don't I just list a few
12 questions? Because I tend to get confused when I hear the
13 answers to them.

14 So, in no particular order, well, first of all,
15 I'd like to echo I thought this was a nice study and well
16 reported.

17 With respect to the analysis of these minor--what
18 are they called?--MDMs. That's all I can remember. Device
19 malfunctions. I notice that there was an analysis of the
20 clinical impact, and this probably will just reveal my
21 ignorance. I noted that you pooled the drift and the long
22 transit time. And I was wondering if that's based on a
23 *priori* considerations that those would have the same impact
24 or just--I would have thought, I guess naively, that drift
25 would be a relatively more important issue and perhaps

1 should be analyzed separately. And in that same vein, if
2 you're trying to understand the impact, it would be useful
3 to look at the edge effect with respect to those cases that
4 had drift. I thought that was a nice analysis of the edge
5 effect where you had lots of power for quantitative
6 analysis. That's sort of one set of questions.

7 The second one, actually sort of related: Did you
8 actually look at the cases with these malfunctions to see if
9 there were any patient differences? In hearing the
10 discussion, it sounds like this is more or less a random
11 occurrence, but I didn't know whether you looked at whether
12 there were patient differences in those that had the drift
13 problem.

14 Relative to the recommendation of length of
15 antiplatelet therapy, do these results in the START trial,
16 are they consistent with the earlier--the Beta-cath results?
17 In other words--I'm sort of following up on your comment.
18 If you looked at the results, there's a certain number of
19 patients that did not get what would have been considered
20 the desirable length of antiplatelet therapy and,
21 nevertheless, no events, no thrombosis occurred. And I
22 wondered if it's just too small a sample or if those results
23 could be compared to the earlier results and see if there's
24 anything different about these data. That's question two.

25 The third question has to do with the--I think I

1 saw among the analyses that were in the packet that wasn't
2 presented some modeling done of the effect of lesion length
3 and treatment and an interaction term, which I thought was
4 very interesting and showed that the treatment effect was
5 more pronounced at longer lesion lengths. And I guess I
6 would just ask if this has something to say about the
7 risk/benefit ratio in terms of the labeling aspects, and I
8 thought that analysis should be made more accessible to the
9 user to determine if there's a lesion length that's less
10 than optimal.

11 And then my last question has to do with the
12 heterogeneity between sites, and I saw some analyses in the
13 packet, but I didn't really understand what was being done.
14 And in particular, were analyses done to suggest that sites
15 had different overall restenosis rates, or was there also
16 differences in efficacy rates? So I guess it was more just
17 ignorance that I didn't know what was being presented.

18 ACTING CHAIRPERSON TRACY: I'm not sure how you
19 want to approach that. Maybe one question at a time? You
20 can identify which aspect you're dealing with and also
21 identify yourselves.

22 DR. KUNTZ: Rick Kuntz. I'm a cardiologist and a
23 part-time statistician, I guess, although I'm a little
24 intimidated by Dr. Bailey.

25 Why don't we just start from the top there?

1 DR. BAILEY: Drift and--pooling drift and transit
2 time.

3 DR. KUNTZ: Right. We had a variety of different
4 MDMS. There were four. There was the drift, transit time
5 issues, inability to deliver the catheter, and some other
6 issues that were like one person per category. They
7 represented a variety of different fields that were
8 prospectively collected in the case report form, all
9 classified as potential device malfunctions.

10 In looking at that overall data set, in order to
11 reduce multiplicity and try to deal with, you know,
12 diminishing the alpha to zero, we focused just on the
13 radiation issues. So the decision to pool transit time and
14 drift was an issue of power and reduction of multiplicity.
15 So we haven't looked at the individual events themselves
16 because they were evenly distributed. I think there were
17 something like 80 cases overall between the two groups that
18 were MDMS, and there were about 40 and 40 on each one, drift
19 versus that. So my expectation is that since we generally
20 found no difference in the adverse events with the pooled
21 group that we probably would be completely underpowered to
22 look at the individual groups themselves, and it probably
23 wasn't worth the analysis.

24 DR. BAILEY: What I was thinking, though, is if
25 you looked at the edge effect in a very quantitative way,

1 looking at the delta minimum luminal diameter in that region
2 and separated specifically the drift ones, you might have
3 some power to look at it.

4 DR. KUNTZ: Right. Your second part of the
5 question about did we look at drift effect with the edge
6 analysis specifically, we have not done that, and I agree
7 that would be an interesting analysis because we may have
8 enough power there because they're both continuous measures.
9 And I think that would be an interesting analysis to do.

10 DR. BAILEY: Patient differences in terms of
11 predicting who--was that just a random event?

12 DR. KUNTZ: Yeah, we spent--the question was could
13 we look at--were there any anatomical patient factors that
14 explained patients who had these MDMs. So we spent a lot of
15 time looking at those factors, and we couldn't find, that
16 is, by generally exploratory analysis, that there were any
17 indicators of increased tortuosity, that there were
18 indicators of the distribution of the vessels. For example,
19 was a right coronary artery more likely to drift than a left
20 coronary artery? The amount of calcium that was in the
21 vessel, the age of the patient, I think a variety of
22 different things. I'm thinking off the top of my head. We
23 tried to evaluate whether we could predict who was going to
24 have a drift of source, and we couldn't find them.

25 Practically speaking, the issues of drift to start

1 with and the intransititude to some degree were really issues
2 of the radiation oncologist, an issue about the touie--the
3 hemostatic device more than they were issues of patient
4 factors; that is, the device itself is designed not to have
5 any kinks and is bulky enough that it generally won't be
6 delivered down very, very tortuous vessels to allow the
7 catheter itself to impede delivery. So delivery impedance
8 were issues of maintaining pressure and issues of the
9 hemostatic valve. So we think that those were the things
10 that explained the differences, not issues of patient
11 characteristics, where we can say that this patient is at
12 more risk of a drift than the other patient on initial
13 exploratory analysis. We couldn't identify patient factors.

14 DR. BAILEY: Can you compare the thrombosis rates
15 between patients who in the START trial had new stents
16 placed and did not get 90 days of therapy to the earlier
17 data that were the basis for--

18 DR. KUNTZ: In the Beta-cath trial, we're looking
19 at a 1,500-patient trial compared to a 476-patient trial,
20 the START trial. So the history which was reviewed is
21 important to understand.

22 In that trial, patients were treated initially
23 with balloon angioplasty. Then depending on the result, the
24 physician decided whether they would go down a PTCA branch
25 based on a very, very good result for which the patient was

1 randomized to placebo versus active therapy blinded with no
2 further stent placement, or if the result from balloon
3 angioplasty was suboptimal, they were arbitrarily decided to
4 go down a stent branch and then randomize after that.

5 So we had a fairly large volume of patients
6 initially, as you can imagine, because of the stent interest
7 at that time of patients who had new stents placed, on the
8 order of five or six hundred patients, as opposed to 50
9 patients with new stents in this study. So the opportunity
10 to observe stent thrombosis was greater in the Beta-cath
11 trial than the opportunity to observe stent thrombosis in
12 this trial.

13 So, initially and early on, when the Data Safety
14 Monitoring Committee with its blinded review of the data
15 identified that there were some problems going on when new
16 stents were placed and patients exposed to radiation therapy
17 and expected--and declared that they wanted to extend
18 antiplatelet therapy--this, by the way, was reviewed with
19 the FDA and the protocol was changed. We anticipated that
20 this also might be an issue in the START trial where new
21 stents were placed.

22 However, at the end of the START trial, the 476-
23 patient trial, only 50 patients received a new stent. So we
24 are extrapolating the potential for stent thrombosis, even
25 though we had excellent results in this study, to our

1 experience with over 500 patients early on where there was a
2 higher incidence of stent thrombosis; hence, the interest in
3 potentially having 90 days or more of antiplatelet therapy.

4 DR. BAILEY: The interaction between lesion length
5 and treatment effect.

6 DR. KUNTZ: It's very interesting interaction. We
7 mainly saw it in the restenosis defined by the analysis
8 segment, not by the stent segment. And what we see is that,
9 in general, lesion lengths are associated with a higher risk
10 of restenosis. That's been true with multiple data sets,
11 especially in non-radiation areas. That is, patients who
12 have longer lesions tend to have a higher risk of restenosis
13 than patients with shorter lesions.

14 When we look at the analysis segment, which
15 actually lets us have the opportunity of measuring the
16 minimum lumen over a wide area, we start to see that
17 radiation therapy had an extra effect on patients with
18 longer lesions, and that made sense; that basically the
19 increased risk the patient was exposed to with a longer
20 lesion afforded a more profound treatment effect from
21 radiation therapy than those who had shorter lesions. So
22 the interaction term of longer lesion lengths and radiation
23 therapy made sense to us, understanding the underlying risk
24 the patient had with longer lesions.

25 DR. BAILEY: In fact, based on the coefficients,

1 if you have a lesion length of 8 millimeters, you're at dead
2 even.

3 DR. KUNTZ: Well, right. It's hard to go back and
4 say where the breakpoint is of radiation therapy being
5 ineffective at some level. All we can say is that the
6 continuum shows that longer lesions have more potential for
7 effect than shorter lesions. But these lesions have been
8 linearized in a linear model. We didn't do a lot of non-
9 linear models to see where the breakpoint is. And so I
10 think that it's an interesting extrapolation, maybe the
11 basis of a hypothesis for a new study.

12 DR. BAILEY: I agree it's not very exact, but I
13 think it points to at least an issue if you're a user
14 whether you want to embark on radiation therapy in a shorter
15 lesion.

16 DR. KUNTZ: That's a good point.

17 DR. BAILEY: And, finally, the site heterogeneity.

18 DR. KUNTZ: Right. The site heterogeneity we
19 thought was typical in most of the studies, that is, the
20 overall distribution of treatment effects for a 50-patient
21 trial showed--a majority of patients showed a similar result
22 as the mean effect overall. A couple sites out of the 50
23 had the opposite results, which you typically see in a
24 normally distributed trial.

25 The other heterogeneity issue dealt with--there

1 was one or two sites that tended to use new stents more
2 often than others per se. However, in the overall analysis,
3 the restenosis rates didn't differ, so we didn't see a
4 profound effect on the site.

5 We performed the typical boilerplate pooling
6 analysis for the FDA looking for interactions between
7 treatment site and the overall main effect, and albeit
8 that's another powered analysis usually, we didn't see any
9 deviation from the normal studies that we saw.

10 DR. BAILEY: Thank you.

11 ACTING CHAIRPERSON TRACY: Okay. At this point I
12 think we'll break for lunch, and if we could resume at 1:15.
13 And I would like to remind the panel members not to discuss
14 the contents of this meeting.

15 [Luncheon recess.]

AFTERNOON SESSION

[1:30 p.m.]

ACTING CHAIRPERSON TRACY: We'll resume the open committee discussion, and we'll resume with the panel questions, and I think we'll start down at that end, please. Again, just to remind everybody to identify yourself and to speak into the microphones.

DR. AYERS: Okay. I moved up a little in the order, but I have a couple of questions. One, your presentation indicated that you were giving 18.4 Gray and 23 Gray, depending on the vessel size, but actually since this is non-centered sources, an asymmetric lumen, shadowing by guide wire and stent, what really was the dose range for these studies, minimum, maximum, your estimates to the (?)

MR. GREEN: I think we'll let Dr. Crocker, radiation oncologist and investigator in some of the BERT feasibility studies and who helped in the START trials, answer this question.

DR. CROCKER: My name is Ian Crocker. I'm a radiation oncologist at Emory University. I'm a consultant and a shareholder in Novoste, and, in addition, Novoste has licensed intellectual property from Emory University, and I am co-owner of that intellectual property.

Within the vessel wall, there is a wide range of

1 doses that are delivered, and that's really true of any
2 brachytherapy source. We did prescribe a dose at 2
3 millimeters from the center of the source, and that initial
4 prescription represented a small incremental increase in
5 dose over what was prescribed in the Beta Energy Restenosis
6 Trial based on anticipated shadowing of the source by the
7 stent struts.

8 We had done some measurements which had shown that
9 there was approximately a 10-percent decrement in dose
10 immediately underneath the stent struts, and as a result of
11 that, we recommended increasing the dose that was delivered
12 in the START trial by 2 Gray, which represented an 11- to
13 14-percent increase in dose over what was delivered in the
14 original BERT trial.

15 DR. AYERS: Do you know how much effect the dose
16 varied from the fact you had a 6-millimeter variation in
17 vessel size for the same dose, so I guess that would be 0 to
18 3 millimeters in variation from the vessel wall to the
19 prescription point, and also the fact that it wasn't
20 centered.

21 DR. CROCKER: Right. With the cohort of patients
22 who were treated in this trial, we undertook a retrospective
23 analysis of dosing using intravascular ultrasound images,
24 and information on that has been submitted to the FDA as
25 part of this submission.

1 Basically, the catheter assumes a relatively
2 centered position within the lumen based on these IVUS
3 ultrasound images, and really there are only minor
4 differences in the doses, you know, that are received to the
5 vessel wall with active centering of the catheter within the
6 lumen compared to non-centering of the catheter.

7 DR. AYERS: One other one, I guess just for
8 clarification. It wasn't clear to me. When you added new
9 stents to about--what, 20 percent of the patient population,
10 as I recall.

11 DR. CROCKER: Correct.

12 DR. AYERS: Was that done before or after the
13 radiation therapy or a mixture?

14 DR. CROCKER: Those new stents were added after
15 the radiation therapy was delivered, so that the protocol
16 specifically excluded patients with stent sandwiches or
17 stent within a stent, so that we didn't anticipate that
18 there would be any areas in which there would be stent
19 overlap and more than approximately a 10-percent decrement
20 in dose due to the shadowing effect.

21 You know, I should say that that decrement in dose
22 becomes less important as you get further away from the
23 stent. In other words, there's a relative filling-in of
24 dose at increasing depths beneath the stent.

25 DR. AYERS: Okay. One other thing I noticed, you

1 spent a lot of time on indicating how small the dose to the
2 patient was from--or incremental dose to the patients from
3 the beta therapy, particularly whole-body, which is
4 certainly true. But nowhere in there I saw addressed is how
5 much increase in the dose was due to the fairly substantial
6 additional fluoroscopy in sign, particularly, you know,
7 monitoring the source position every 30 seconds. Do you
8 have any value for the added skin dose for that additional
9 fluoro?

10 DR. CROCKER: I'm not sure that I have any
11 additional information regarding the fluoroscopic dose.
12 Maybe Dr. Popma might want to comment on this.

13 DR. POPMA: These are very short pops of
14 fluoroscopy and not a long length. You really just have a
15 second or less, just to check the position of it, which you
16 can review on your video replay.

17 DR. AYERS: Okay. But every 15 seconds, that
18 would be, what, over a 2-minute treatment time?

19 DR. POPMA: Again, just a second or so each time.

20 DR. AYERS: And I was curious--and the last item I
21 have for right now. I guess I forgot--I'm Robert Ayers,
22 NRC. I didn't identify myself starting this. Novoste
23 introduced later in the study, when problems were uncovered
24 and we investigated some of these and, in fact, generated an
25 information notice on source transport difficulties,

1 particularly through the introducer or touie borst valve or
2 whatever was used for that. It's known that if that's
3 overtightened, it can block the sources going either in or
4 out, and I think well less recognized perhaps by the panel,
5 if you overtighten it too far and go past the elastic limits
6 of the catheter, that blockage stays there even if you
7 loosen the valve. And they introduced this introducer
8 sheath as a corrective measure for that but don't require
9 it, and I wonder how come.

10 Our experience has been that the cardiologists
11 don't use it because it's an extra step, in the one incident
12 we looked at.

13 DR. POPMA: Drew?

14 MR. GREEN: What we found when we--first, of
15 course, you're correct. We did qualify an arrow sheath
16 introducer as an additional accessory that the clinician
17 could use in the procedure. They would place the catheter
18 through the introducer sheath, which is its labeled
19 indication for use for introduction of percutaneous
20 catheters. And, therefore, when they would tighten the
21 touie borst, the hemostasis valve down onto that sheath, it
22 would protect the catheter.

23 What was found in talking with the centers and
24 looking at what was happening, especially the center that
25 you were talking about, was that this had, if you will, a

1 learning--it was part of a learning curve, you know, part of
2 the learning of using the device, and that the
3 interventional cardiologist places the catheter and, you
4 know, it's their job every day to maintain placement of that
5 catheter, angioplasty catheter, guide catheter, what have
6 you. And so that's their job.

7 So now they have a new player, a new team member
8 who's also involved with the interaction of that catheter
9 with the transfer device. So it's part of the training, you
10 know, of the team working together, about moving and then
11 learning when to tighten and how much to tighten on the
12 catheter to allow for passage of the source and, you know,
13 to compensate for another person being in the team.

14 So they wanted the ability to have this as a tool
15 if it was necessary to use in their practice, and if not, or
16 if they felt that they were at--or had demonstrated they
17 were at the learning curve to where they didn't need this
18 tool, they didn't have to use it. And this becomes very
19 important because in the training section of the panel pack,
20 we actually talk about, you know, going through experiences
21 from the trial or all the trials and evaluating what the
22 proper--you know, how a user would use the system. And as
23 part of the hands-on training and the mock training and
24 these pieces of the training, it's important for the
25 clinicians to determine, you know, of the optional

1 accessories such as the arrow sheath and the fluid
2 management system, how they would apply that in their
3 practice; in other words, what works best for them so that
4 they can use this system the most effectively to gain the
5 results that were seen in the START trial.

6 DR. SPEISER: At our institution, the radiation
7 oncologist uses the arrow flex sheath, will flush it and
8 place it over the delivery catheter so that there is no time
9 delay for the cardiologist. And it is my intent in the
10 training program to train the radiation oncologist to do
11 this and to use it all the time, unless the cardiologist
12 specifically says they do not want to use it.

13 DR. AYERS: Well, we're strongly considering
14 making that a mandatory requirement. That's why I wanted to
15 ask the question, at least at our regulatory agency.

16 And one last one was with our upcoming change to
17 our medical regulations--and I'm assuming--I'm not saying
18 that that's, in fact, the way it will work out, but most
19 cases for brachytherapy and particularly high dose rate, we
20 have a mandatory requirement coming now that the licensee,
21 user, medical physicist, you know, the medical institution,
22 however you want to characterize it, is solely responsible
23 for the calibration of the brachytherapy source dose rate.
24 And going over your submission, particularly in Section 2,
25 it is not clear to me that you provide the tools with your

1 system to allow the medical physicist at your customer site
2 to perform proper dose rate calibrations on these sources.

3 MR. GREEN: We understand that some sites do at
4 this time have requirements at a site level, possibly there
5 will be other requirements later for site verification of
6 dose rate, et cetera. And we do have proposals on how to
7 handle that. And I think that Dr. John Lobdel, our Director
8 of Radiation Management, can speak to that and let you--you
9 know, what the proposals we are planning to do at the sites
10 to be able to address that are.

11 DR. LOBDEL: John Lobdel, employee of Novoste. We
12 have a source train that was calibrated at NIST to determine
13 the dose rate at a half millimeter--I'm sorry, at 2
14 millimeters inside the--in water--I'm sorry. Let me go
15 back. We have a source train that was calibrated by NIST to
16 measure the dose rate at 2 millimeters from the center line
17 of the source train in water. This train is our transfer
18 standard. This train is used to calibrate all the trains we
19 send to the clinical sites.

20 Now, during the clinical trials, we had two sites
21 that asked to verify our dose rate. We worked with them on
22 this. The tool we have is a solid water block that
23 positions the center line of the source train at 2
24 millimeters from a film plane. We went to the site,
25 irradiated the trains for the hospitals, then from this film

1 they analyzed the film, determined the dose and dose rate,
2 and also the homogeneity of the train.

3 We actually published a paper on the results. The
4 hospitals were very happy with the results, and it came out
5 in the literature about a year ago.

6 Now, we're also looking into a different method,
7 and that different method is there's another source train in
8 this being calibrated for dose rate and activity. That
9 train will be sent to at least one and probably two
10 accredited dosimetry calibration laboratories, or ADCLs.
11 The ADCL will in turn calibrate their equipment on this
12 train. Then when a hospital wants to know what the dose
13 rate and activity of our trains are, they can simply send
14 their well chamber to the ADCL. It will be calibrated there
15 and returned. And then as often as they wish, they can
16 simply take a source train and put it into the well train to
17 determine activity and dose rate.

18 So we have here one system that has been proven
19 and has been published in the literature. We have another
20 one that we are working on that should be available quite
21 soon.

22 DR. AYERS: The latter was what I was really
23 looking for since our anticipated regulations require our
24 licensees to go to an ADCL or NIST for this type of
25 calibration that you just concluded with. So it sounds

1 good.

2 DR. LOBDEL: Thank you.

3 DR. AYERS: That's all I have.

4 ACTING CHAIRPERSON TRACY: That's it? Okay.

5 Dr. Crittenden, any questions?

6 DR. CRITTENDEN: Yes, I have several. The first
7 question I'm going to direct to Dr. Speiser and Mr. Green.

8 What was the impetus behind the device change--the
9 impetus behind the change in the design for the device going
10 from the Alpha III to Alpha IV? Was it to minimize the
11 problems with source delivery? If so, was there a
12 comparison made between these two devices to see if there
13 was a difference.

14 And then, finally, for this first question, Dr.
15 Speiser stated that the radiation oncologist might feel
16 uncomfortable in the cath lab given that this is a new
17 setting for them and that this may have been a source for
18 some of the source drift or source transit time problems
19 that we saw with the devices.

20 Is it your position, because the analysis showed
21 that there may be no difference when you look at placebo
22 versus the Strontium-treated groups in terms of outcome,
23 whether you have an MDM or not, is it your position that
24 there are no untoward sequelae for source drift or source
25 transit time?

1 MR. GREEN: I'll go first, kind of go in order of
2 the questions. The first question I understood to be, you
3 know, what was the reasoning behind going from the Alpha III
4 transfer device to the Alpha IV device. Basically, we
5 started the clinical trial with the Alpha III transfer
6 device. We were having that transfer device manufactured or
7 built by a subcontractor. We moved to another
8 subcontractor, qualified that subcontractor, and part of the
9 qualification of that subcontractor went through for us was
10 to do an evaluation of the Alpha III transfer device and to
11 propose some improvements in the device that may make it
12 more user-friendly, the user interface a little easier to
13 use. And there were several minor things besides the LED.

14 For instance, the shape of the housing was changed
15 a little better to fit the hand. Some of the graphics were
16 made a little clearer, and the LEDs were added. And the
17 LEDs were added because the subcontractor here determined
18 that they believed that that gave a more accurate feedback
19 than did the mechanism of the Alpha III, which was simply an
20 open window that showed to where you could visualize a
21 spring as part of the pressure relief valve. So they
22 believed this would be a more accurate and calibratable
23 method of providing feedback to the user. And that was the
24 reason we went to that change.

25 So when we implemented that change, one of the