

Premarket Review Considerations for Reprocessed SUDs

**Barbara Zimmerman
Office of Device Evaluation**

Classification

May 1976 premarket process begins

- **Three classes I, II, III**
 - **Class I (least risk)**
 - **Class II**
 - **Class III (greatest risk)**

**All devices grouped by generic type
into one of the three classes**

Important Dates for Reprocessors:

| Premarket Submission Due Dates | | |
|--------------------------------|-----------|---------------------------------|
| | Due By: | Must be Cleared or Approved by: |
| Class III | 2/14/2001 | 2/14/2002 |
| Class II non-exempt | 8/14/2001 | 8/14/2002 |
| Class I non-exempt | 2/14/2002 | 8/14/2002 |

Devices Excluded from Reprocessing Enforcement Strategy

- **Permanently implanted pacemakers (see Compliance Policy Guide, 7124.12)**
- **Hemodialyzers (see specific guidance)**
- **Opened but unused devices**
- **SUDs reprocessed by health care facilities that are NOT hospitals (however, enforcement may be expanded to cover in the future)**

Premarket Requirements

- **Premarket Notification (510(k))**
 - **21 CFR Part 807**
- **Premarket Approval (PMA)**
 - **21 CFR Part 814**

510(k) or PMA?

- **Class III: Submit a premarket approval application (PMA).**
- **Class II: Submit a premarket notification (510(k)), unless device is exempt from premarket requirements.**
- **Class I: no 510(k) or PMA is required except for “reserved” class I devices.**

Examples of Reprocessed SUD's

| Device | Class |
|--------------------------------|-------|
| Surgical Saw Blades | I |
| Electrophysiology Catheters | II |
| Cardiac Ablation Catheters | III |

Additional devices listed at
<http://www.fda.gov/cdrh/reuse/1168a.html>

PMA vs. 510(k)

- **PMA**
 - **valid scientific evidence**
 - **risk/benefit analysis**
- **510(k)**
 - **substantial equivalence**

Reprocessors of SUDs Are Considered to Be Manufacturers

- **All premarket submissions to be handled the same**
 - **No special considerations or allowances, e.g., expedited review or reduced testing requirements**
 - **No extraordinary review criteria or supplemental evaluations necessary**

Specific Considerations for PMA Applications

- **Data to demonstrate safety & effectiveness**
- **Manufacturing processes (including validation)**
- **Pre-approval inspection required**
- **Manufacturing process in compliance with QSR**

510(k) Submission Describes a Device Not a Process

- **Device adequately described & characterized (through testing)**
- **Demonstrate that reprocessed SUD is substantially equivalent to predicate device**

510(k) Submission Describes a Device Not a Process

- **Manufacturing processes assessed through QSR and FDA inspections**
- **Sterility an important device characteristic**

Sterilization and Disinfection Information Should Include

- **Description of decontamination/
cleaning processes**
- **Description of validation method for
cleaning and sterilization/disinfection**
 - **Simulated use protocol should be developed**
 - **Reflect worse case conditions**
 - **Manipulation and simulated use in biological
soil**

Sterilization and Disinfection Information Should Include

- **Endpoints for cleaning and sterilization/disinfection**
- **Sterility assurance level**
- **Conformance with appropriate recognized standards**
- **Description of packaging materials**

Performance Tests

- **Device specific - check for appropriate guidance document**
- **Additional tests may be needed to assess effects of reprocessing**
- **Test devices that have gone through simulated use protocol**
- **Include descriptions of tests performed to determine product release, provide pass/fail criteria**

Biocompatibility Tests

- **If needed for OEM device, will be needed for reprocessed devices**
- **Cleaning and sterilization agents may adversely affect device materials and biocompatibility**

Biocompatibility Tests

- **Testing should be done under worst case conditions (maximum times device is to be reprocessed; include exposure to simulated use).**
- **See specific ODE guidance document on biocompatibility testing.**

Shelf Life Tests

- **If needed for OEM device, will be needed for reprocessed device**
- **Evaluate and document stability under worse case conditions**
 - **should reflect exposure to simulated use protocol and maximum number of times device can be reprocessed**

Labeling

- Labeling (device description, indications for use, adequate directions for use) are required for all devices cleared under 510(k).
- Labeling must include reprocessor's name, address and other contact information.
- Labeling must be included with the device (not on file at the hospital or user facility).

Can One 510(k) Cover Multiple Devices?

- **Yes, if devices would be included in one 510(k) by the OEM**
- **Bundled devices w/same indications for use**
- **Reprocessors could propose specifications for a device that encompasses models from more than one OEM**

To Expedite Review

- **In cover letter identify as a premarket submission for reprocessed SUD.**
- **Consider using special and abbreviated 510(k)s.**
- **Submit master files if information can be used in multiple submissions.**
- **Pre-submission meetings and teleconferences may be useful.**

Office of Device Evaluation



```
graph TD; A[Office of Device Evaluation] --> B[Div. of Dental, Infection Control, & General Hospital Devices]; A --> C[Div. of Reproductive, Abdominal, & Radiological Devices]; A --> D[Div. of General, Restorative & Neurological Devices]; A --> E[Div. of Ophthalmic & ENT Devices]; A --> F[Div. of Cardiovascular & Respiratory Devices]; A --> G[Div. of Clinical Laboratory Devices];
```

The diagram is an organizational chart for the Office of Device Evaluation. At the top is a purple box with the text "Office of Device Evaluation". From the bottom of this box, six white arrows point downwards to six separate blue boxes. The boxes are arranged in two columns. The left column contains three boxes: "Div. of Dental, Infection Control, & General Hospital Devices", "Div. of General, Restorative & Neurological Devices", and "Div. of Ophthalmic & ENT Devices". The right column contains three boxes: "Div. of Reproductive, Abdominal, & Radiological Devices", "Div. of Cardiovascular & Respiratory Devices", and "Div. of Clinical Laboratory Devices".

**Div. of Dental,
Infection Control,
& General
Hospital Devices**

**Div. of Reproductive,
Abdominal, &
Radiological Devices**

**Div. of General,
Restorative &
Neurological
Devices**

**Div. of
Cardiovascular
& Respiratory
Devices**

**Div. of Ophthalmic
& ENT Devices**

**Div. of Clinical
Laboratory Devices**

Contact Persons

- **Coordination of Reuse Policy in ODE
(premarket issues only)**
 - **Philip J. Phillips (301) 594-2022**
 - **Timothy A. Ultatowski (301) 443-8879**
 - **Heather Rosecrans (301) 594-1190**
 - **Barbara Zimmerman (301) 594-2036**

Guidance Documents

- **Specific for Reprocessors**
 - **“Enforcement Priorities for Single-Use Devices Reprocessed by Third Parties and Hospitals”**

www.fda.gov/cdrh/reuse/1168.html

- **“Frequently Asked Questions on Reuse of Single Use Devices”**

www.fda.gov/cdrh/reuse/reuse-faq.shtml

Guidance Documents

- **Specific for Reprocessors**
 - **Letter to Hospitals Re: Reprocessing of Single Use Devices - 4/23/01**
www.fda.gov/cdrh/reuse/042301_reuse.html
 - **Compliance Policy Guide (CPG 7124.16) Section 300.500. Reuse of Medical Disposable Devices - 9/24/87**
www.fda.gov/cdrh/comp/cpgreuse.pdf

Guidance Documents

- **Biocompatibility**
 - **“Use of International Standard ISO-10993, 'Biological Evaluation of Medical Devices Part 1: Evaluation and Testing’”**

www.fda.gov/cdrh/g951.html

Guidance Documents

- **Reusable Devices**
 - **“Labeling Reusable Medical Devices for Reprocessing in Health Care Facilities: FDA Reviewer Guidance”**
www.fda.gov/cdrh/ode/198.pdf

Guidance Documents

- **Reusable Devices**
 - **“Questions & Answers for the FDA Reviewer Guidance: Labeling Reusable Medical Devices for Reprocessing in Health Care Facilities”**

www.fda.gov/cdrh/ode/1198.html

Guidance Documents

- **Use of Standards in 510(k)s**
 - **“Use of Standards in Substantial Equivalence Determinations”**
www.fda.gov/cdrh/ode/guidance/1131.html
 - **“Guidance on the Recognition and Use of Consensus Standards”**
www.fda.gov/cdrh/modact/k982.html

Guidance Documents

- **Use of Standards in 510(k)s**
 - **“Frequently Asked Questions on Recognition of Consensus Standards”**

www.fda.gov/cdrh/modact/faqost.html

Guidance Documents

- **Abbreviated and Special 510(k)s**
 - **“A New 510(k) Paradigm - Alternate Approaches to Demonstrating Substantial Equivalence in Premarket Notifications”**

www.fda.gov/cdrh/ode/parad510.html

- **“Frequently Asked Questions on the New 510(k) Paradigm”**

www.fda.gov/cdrh/ode/92_a.html

Guidance Documents

- **Sterilization**
 - **“Updated 510(k) Sterility Review Guidance K90-1; Final Guidance for Industry and FDA”**

www.fda.gov/cdrh/ode/guidance/361.pdf