

APR - 2 2001

K002856

510(k) Summary  
AT3C70 Transducer  
Sonora Medical Systems

## 510(k) Summary

The following safety and effectiveness summary has been prepared pursuant to requirement for 510(k) summaries specified in 21CFR§807.92(a).

### 807.92(a)(1)

#### Submitter Information

Colleen Hittle, Official Correspondent  
7992 Castleway Drive  
Indianapolis, IN 46250  
Phone: (317) 849-1916  
Facsimile: (317) 577-9070

Contact Person: Colleen Hittle

Date: May 1, 2000

### 807.92(a)(2)

Trade Name: AT3C70 Transducer  
Common Name: Ultrasound Transducer  
Classification Name(s): Transducer, Ultrasonic, Diagnostic  
Classification Number: 90ITX

### 807.92(a)(3)

#### Predicate Device(s)

ATL 3.5CLA/76 K903603

Additional Substantial Equivalence Information is provided in the following substantial Equivalence Comparison Table.

510(k) Summary  
AT3C70 Transducer  
Sonora Medical Systems

807.92(a)(5)

**Intended Use(s)**

The AT3C70 transducer is intended to be used by or under the direction of a physician and in conjunction with standard ultrasound for clinical imaging in fetal and abdominal applications.

510(k) Summary  
AT3C70 Transducer  
Sonora Medical Systems

**Substantial Equivalence Comparison Table**

	Predicate Probe ATL - 3.5CLA/76 (K#903603)	Submitted Probe Sonora - AT3C70
Type:	Curved	Curved
Frequency: Imaging	3.5 MHZ	3.5 MHz
Doppler	3.0 MHZ	3.0 MHz
Operating Mode	B, M, PW, CFM	B, M, PW, CFM
Indications * For Use	Fetal, Abdominal	Fetal, Abdominal

\* Not operable in Fetal Doppler mode



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

APR - 2 2001

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Sonora Medical Systems  
C/O Ms. Colleen Hittle  
The Anson Group  
7992 Castleway Drive  
INDIANAPOLIS IN 46250

Re: K002856  
Trade Name: Model AT3C70 Transducer  
Regulatory Class: Class II/21 CFR 892.1570  
Product Code: 90 ITX  
Dated: February 9, 2001  
Received: February 12, 2001

Dear Ms. Hittle:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

This determination of substantial equivalence applies to the following transducers intended for use with the ATL UM4+ and UM 9 Diagnostic Ultrasound Systems, as described in your premarket notification:

Transducer Model Number

AT3C70

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval) it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895.


A substantially equivalent determination assumes compliance with the Good Manufacturing Practice requirement, as set forth in the Quality System Regulation (QS) for Medical Devices: General (GMP) regulation (21 CFR Part 820) and that, through periodic QS inspections, the FDA will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, the Food and Drug Administration (FDA) may publish further announcements concerning your device in the Federal Register. *Please note:* this response to your premarket notification does not affect any obligation you may have under sections 531 and 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus permits your device to proceed to market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4591. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or at (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

If you have any questions regarding the content of this letter, please contact Rodrigo Perez at (301) 594-1212.

Sincerely yours,



*for* Daniel G. Schultz, M.D.  
Acting Director, Division of Reproductive,  
Abdominal and Radiological Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure(s)

## Diagnostic Ultrasound Indications for Use Form

Fill out one form for each ultrasound system and each transducer.

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application	Mode of Operation									
	A	B	M	PWD (PW)	CWD	Color Doppler (CFM)	Amplitude Doppler (PD)	Color Velocity Imaging	Combined (specify)	Other (specify)
Ophthalmic										
Fetal										
Abdominal		P	P							
Intraoperative (specify)		P	P							
Intraoperative Neurological										
Pediatric										
Small Organ (specify)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac		P	P							
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular		P	P							
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other										

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Additional Comments:

Applicable combined modes:

*Prescription Use*

(Division Sign-Off)

Division of Reproductive, Abdominal, ENT,  
and Radiological Devices510(k) Number K002856

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

## Diagnostic Ultrasound Indications for Use Form

Fill out one form for each ultrasound system and each transducer.

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application	Mode of Operation									
	A	B	M	PWD (PW)	CWD	Color Doppler (CFM)	Amplitude Doppler (PD)	Color Velocity Imaging	Combined (specify)	Other (specify)
Ophthalmic										
Fetal		P	P							
Abdominal		P	P	P		P				
Intraoperative (specify)		P	P	P		P				
Intraoperative Neurological										
Pediatric										
Small Organ (specify)		P	P	P		P				
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular		P	P	P		P				
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other										

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Additional Comments: Small Organ=(specifically, thyroid, testicles and breast)

Applicable combined modes:

***Prescription Use****David A. Syman*  
(Division Sign-Off)

Division of Reproductive, Abdominal and Radiological Devices

510(k) Number **K002851**

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

## Diagnostic Ultrasound Indications for Use Form

Fill out one form for each ultrasound system and each transducer.

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application	Mode of Operation									
	A	B	M	PWD (PW)	CWD	Color Doppler (CFM)	Amplitude Doppler (PD)	Color Velocity Imaging	Combined (specify)	Other (specify)
Ophthalmic										
Fetal		N	N							
Abdominal		N	N	N		N				
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric										
Small Organ (specify)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other										

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Additional Comments:

Applicable combined modes:

*Prescription Use* ✓*David A. Sygum*  
(Division Sign-Off)Division of Reproductive, Abdominal, ENT,  
and Radiological Devices510(k) Number *K002856*

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)