

APR 13 2001

Attachment 4  
510(k) Summary

K 010956

| Category:                              | Comments  |
|--|---|
| Sponsor:                               | Boston Scientific Corporation<br>2710 Orchard Parkway<br>San Jose, CA 95134   |
| Correspondent:                         | Andrea L. Ruth<br>Associate II, Regulatory Affairs<br>2710 Orchard Parkway<br>San Jose, CA 95134                              |
| Contact Information:                   | E-mail: <a href="mailto:rutha@bsci.com">rutha@bsci.com</a><br>Phone: 408.895.3625<br>Pager: 888.509.6375<br>Fax: 408.895.2202 |
| Device Common Name                     | Electrosurgical Probe   |
| Device Proprietary Name                | Cobra® Flex™ Family of Surgical Probe   |
| Device Classification                  | Class II, 21 CFR §878.4400,<br>product code GEI   |
| Predicate Device                       | Electrosurgical Probe   |
| Predicate Device Manufacturer(s)       | BSC/EP Technologies, Inc.   |
| Predicate Device Proprietary Name(s)   | Cobra® Surgical Probe   |
| Predicate Device Classification Number | 21 CFR §878.4400, product code GEI  |
| Predicate Device Classification(s)     | Class II  |

Date Summary  
Was Prepared:

April 11, 2001

Description of  
the Device:

The Boston Scientific Corporation Surgical Probe is a sterile, single use electrosurgical device intended to be used to coagulate soft tissues. The Surgical Probe transmits radiofrequency energy from electrodes which are connected to an Electrosurgical Unit (non-sterile; re-useable) through an Instrument Cable (sterile; re-useable).

**Intended Use:**

The Probe is intended for use only under direct visual control of the physician during open, general surgical procedures to coagulate soft tissues. The Probe may also be used to coagulate blood and soft tissue to produce hemostasis.

**Comparison to  
Predicate  
Device:**

|                       | Predicate Device                      | Modified Device    |
|-----------------------|---------------------------------------|--------------------|
| 510(k) Reference      | K981981                               | Current Submission |
| Intended Use          | Coagulation of Soft Tissue            | Same               |
| Device Description    | Electrosurgical Probe                 | Same               |
| Single Use?           | Yes                                   | Same               |
| EO Sterilized?        | Yes                                   | Same               |
| Manufacturer          | BSC/EP Technologies, Inc.             | Same               |
| Device Classification | 21 CFR §878.4400, product code<br>GEI | Same               |

**Summary of the  
Non-clinical  
Data:**

Where appropriate, testing conformed to the requirements of 21 CFR Part 58 (Good Laboratory Practices (GLP)). Specifically, non-clinical tests conducted for the Device showed the device met its design-input criteria, and was safe & effective for its intended use.



APR 13 2001

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Ms. Andrea L. Ruth  
Associate II, Regulatory Affairs  
Boston Scientific Corporation  
2710 Orchard Parkway  
San Jose, California 95134

Re: K010956

Trade/Device Name: Cobra® Flex Family of Surgical Probes  
Regulation Number: 878.4400  
Regulatory Class: II  
Product Code: GEI  
Dated: March 29, 2001  
Received: March 30, 2001

Dear Ms. Ruth:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

Page 2 - Ms. Andrea L. Ruth

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,

for Celia M. Witten, Ph.D., M.D.

Director  
Division of General, Restorative  
and Neurological Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

**Attachment 2  
Intended Use Statement**

510(k) Number (if known): K010956

Device Name: Cobra® Flex Family of Surgical Probes

**Indication for Use:**

The intended use remains the same as found in K981981 approved September 3, 1998, and reads as follows:

The Probe is intended for use only under direct visual control of the physician during open, general surgical procedures to coagulate soft tissues. The Probe may also be used to coagulate blood and soft tissue to produce hemostasis.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

**Concurrence of CDRH, Office of Device Evaluation (ODE)**

for Mark N. Millerson  
(Division Sign-Off)  
Division of General, Restorative  
and Neurological Devices

Prescription Use ✓  
(Per 21 CFR §801.109)

510(k) Number K010956  
OR Over-the-Counter Use \_\_\_\_\_

(Optional Format 1-2-96)