

1 surgical option down the road in a positive or negative
2 effect?

3 DR. ALVAREZ: Dr. Wright, the answer to your first
4 question is in our initial pilot study, I actually had a
5 patient who had had previous surgery. It was an endoscopic
6 release and he was quite unsuccessful. He was still having
7 quite a bit of pain. That OssaTron treatment on him, in
8 three months, he was basically pain-free.

9 As far as it interfering with future surgery, the
10 thing, as a treating foot and ankle orthopedist, and I can
11 tell you that I see my share of heel pain, as I am sure that
12 Dr. Pfeffer and you see in your practice of foot and ankle,
13 is it actually will give me an opportunity to treat a
14 patient and virtually do minimal or no harm.

15 The question then becomes, well, what is the level
16 of activities of daily living for these people that come in
17 that cannot go back to work. They work on concrete floors.
18 All of a sudden, for my practice, I potentially have
19 something to offer these people.

20 Now, to put this in perspective. When
21 HealthTronics asked me to do the pilot study, I was very
22 skeptical. My first five patients I treated, I did not
23 treat any further than that until I saw that the result was
24 tending to be.

25 What that is going to do for my practice, if the

1 OssaTron is released, it is going to give me a treatment
2 that I can give at about six months for somebody that comes
3 in and tells me that, "I can't stand this pain anymore. You
4 have done night splints. You have done injections. You
5 have done arch support. You have changed my shoes. I just
6 can't continue to move around."

7 All of a sudden, I can have something to offer
8 these people and we hear it. It is very frustrating.

9 Thank you.

10 DR. BOYAN: Dr. Pfeffer?

11 DR. PFEFFER: I have several very brief questions.
12 First, how was the rupture of plantar fascia diagnosed in
13 the two cases that you have?

14 DR. OGDEN: Obviously, initially by history and
15 physical examination. And we documented it with MRI.

16 DR. PFEFFER: Did you look at the effect of prior
17 cortisone injections as a predictor of success?

18 DR. OGDEN: We divided out the treatment groups by
19 each group having had cortisone injections. We didn't
20 specifically look whether or not the prior cortisone
21 injection might have affected the OssaTron treatment,
22 itself. No; we did not do that.

23 DR. PFEFFER: Did you analyze patient weight as a
24 covariate?

25 DR. OGDEN: No.

1 DR. PFEFFER: How did you establish sensitivity
2 and reproducibility of the dolorimeter test?

3 DR. OGDEN: We did not specifically test that with
4 reproducibility. The way it was done, the device was placed
5 against the patient's heel until they experienced the pain
6 that they had in the morning. That number was recorded and,
7 on each successive evaluation, we pushed the palpometer, the
8 dolorimeter, to that number and asked the patient, then, to
9 give us their VAS scale on a 1 to 10.

10 It is a very simple device. I don't think there
11 is an easy way to quantitate it.

12 DR. PFEFFER: I thought that was very fair,
13 personally, but I just wasn't familiar with the device, so I
14 just wanted a little more information.

15 A couple of other quick issues. Would you object,
16 or would the company object, to the term, as the panel has
17 been asked, of plantar fasciitis or, more specific, proximal
18 plantar fasciitis to replace the term heel-pain syndrome in
19 the diagnosis and treatment.

20 DR. OGDEN: I should probably defer to Ms. Marlow
21 on that. I am not sure.

22 MS. MARLOW: When we first submitted the
23 feasibility study IDE, and several supplements afterward,
24 this was one of the most contentious topics we had.
25 Unfortunately for the division we are working with, we also

1 had changing medical officers we were working with. The
2 definition changed every time.

3 We have no objection whatsoever to changing.
4 Actually, as I tried to allude to during my presentation, we
5 finally ended up using your definition, or the definition in
6 the textbook that you wrote the chapter for.

7 DR. PFEFFER: I am honored, but that was ten years
8 ago and I, personally, have changed the way I look at this.
9 Again, it is just a semantic issue. You treated the
10 condition in question. Whether we call it plantar fasciitis
11 or something else--

12 MS. MARLOW: We have no objection.

13 DR. PFEFFER: My last question, very brief; in a
14 warning, in the information that we received, on Page 55, it
15 states, it is 3.12.3 warnings; "ESW treatment with the
16 OssaTron should be performed with experience in the care of
17 patients with foot and ankle disorders." Will that include
18 primary-care doctors, and how are those types of physicians
19 defined?

20 MS. MARLOW: That is a difficult issue for us.
21 Our intent is to target orthopedic surgeons. We believe
22 that this device is best used in the hands of orthopedic
23 surgeons. The way that our training program is set up and
24 the way that our business practices are right now, that is
25 initially what definitely will happen.

1 I can foresee the same scenario that you foresee.
2 I think that if there is a primary-care physician, if I may
3 use the example of one of our clinical-study sites, at our
4 Birmingham site. There are primary-care physicians that
5 work with the orthopedic physicians there who would be
6 appropriate--there would be no problem with training them
7 along with their orthopedic colleagues to use the device.

8 Outside of that example, I can't say that I know
9 what we would do. I know that we plan to hold the training
10 course for everyone. We intend to implement it the same for
11 everyone and, hopefully, that will adjust the situation.

12 DR. PFEFFER: Good. Thank you very much. It is
13 very difficult to study plantar fasciitis because of patient
14 compliance and follow-up issues which we have discussed. I
15 would congratulate you on a fine effort.

16 DR. BOYAN: Dr. Silkaitis?

17 DR. SILKAITIS: Thank you, Dr. Boyan. In an
18 effort to save time, et cetera, et cetera, many of the
19 comments and statements that were made were very
20 appropriate. So, therefore, I don't have anything to add.
21 Other than we all recognize that training is important. It
22 is not necessarily the title of the person who is either
23 doing the training or is receiving the training, but that
24 training is important and somehow that that is documented.

25 That's all.

1 DR. BOYAN: Thank you.

2 We have actually lost our consumer rep, so, while
3 we are waiting for her to come back, why don't we go over to
4 Dr.Aboulafia. You had some questions?

5 DR. ABOULAFIA: I did have just a few questions.
6 First of all, it is not clear how long the patients remain
7 blinded in the--specifically, if you look at all the data,
8 the nonrandomized group dramatically was better than
9 anybody. The best group is the outpatient; here is what we
10 are going to do, and go ahead and do it.

11 While you get great results that way, our job is
12 to look at efficacy. So we have to look at that
13 nonrandomized treatment group differently. So there is
14 clearly a marked placebo effect that we can all agree on,
15 that there is at least a marked placebo effect because the
16 sham group also improves pretty dramatically.

17 Were the patients who filled out their 12-week
18 follow up still blinded to what the treatment was?

19 MS. MARLOW: Yes.

20 DR. ABOULAFIA: How do you blind them? If one
21 group has the liquid, the bag of saline, presumably, between
22 them, are they able to see that there is something between
23 them and the electrode?

24 DR. OGDEN: No. Every patient had a blind put up
25 between them and the treatment device. The only thing that

1 they were able to sense was the machine going off each time.
2 You raised a very interesting question that I brought up and
3 that is the patients who received treatment did not know
4 whether they received treatment.

5 So, in essence, they were a pseudoplacebo patient
6 in that they were as unaware of whether they had received
7 the treatment as the patient that didn't. That may, indeed,
8 have affected their perception of pain. They may have
9 thought that, "Because I am not getting better right away, I
10 didn't receive it," and that may have affected the way they
11 answered some of the subjective questions. There is just no
12 question.

13 I think that is probably a big explanation for the
14 difference between those who received the treatment, not
15 knowing what they got for three months and those who
16 absolutely knew as a training patient they were getting the
17 treatment.

18 DR. ABOULAFIA: Thanks.

19 DR. OGDEN: You're welcome.

20 DR. ABOULAFIA: The other thing is, and I know it
21 is tough to control for these things, but there is probably
22 a huge difference--in your inclusion criteria, you spell out
23 that they had to have undergone certain nonoperative, not
24 necessarily conservative, but nonoperative treatment
25 modalities.

1 There is probably a big difference, though, in
2 giving someone a week's worth of a nonsteroidal
3 antiinflammatory versus a trial of, let's say, four weeks of
4 a nonsteroidal antiinflammatory. There is no
5 quantification, at least that I could get from this, of the
6 duration of nonoperative treatment other than that they were
7 treated for six months.

8 But one group was given, let's say, an orthotic
9 and it didn't work for 5.9 months and the next time the
10 patient was seen, he got a nonsteroidal antiinflammatory
11 medication for three days, that might not be an appropriate
12 inclusion criteria. Does that make sense? Any response on
13 how that could--in other words, I think we, as treating
14 physicians, inherently know what a reasonable nonoperative
15 treatment prescription is. It is not clear to me that each
16 one of those patients fit into that category.

17 MS. MARLOW: I think that is one of the best
18 reasons for doing a randomized trial, because you do have--
19 especially in a situation like ours where there are so many
20 factors, subjective and-objective, that impact a patient's
21 perception of pain.

22 I think that the best answer to that is,
23 hopefully, because we did do a well-run study, that has been
24 controlled for by having a placebo control.

25 DR. OGDEN: Again, that is a very valid point to

1 try to determine and make as cohesive a patient population
2 to do the studies on. We had a minimum of three kinds of
3 treatments, and, on an average, I think it was closer to
4 five different kinds of treatments that patients had.

5 Hopefully, by having that number of treatments
6 prior to doing this, you kind of mellow out the variation in
7 the days that the patient may have taken NSAIDs. The
8 patients who had the orthotic devices were allowed to
9 continue to use those, so we did not stop that, which may
10 have introduced another variable.

11 DR. ABOULAFIA: Okay, great. I don't consider
12 myself a cynical person, but let me give sort of a cynical
13 interpretation of the data and see if you all agree or not,
14 again recognizing this is a difficult problem to treat and
15 you are looking at a group of patients who have already
16 failed, for lack of a better term, on nonoperative
17 treatment.

18 But what you can tell patients who are undergoing
19 this study that--Dr. Alvarez said that people were getting
20 back to work better, they were doing more, after his initial
21 five patients who were the nonrandomized group who had a
22 markedly impressively different response than the patients
23 who were in the randomized group.

24 Looking at the data, I think what you can tell
25 people is that if an investigator examines your foot, we are

1 going to have statistically significantly different results
2 with this treatment than without and that you will have
3 statistically significantly different results in terms of
4 your own assessment of pain on walking early in the morning,
5 but that you will not have any statistically significant
6 difference in activity level or any statistically
7 significant difference in terms of medication use.

8 Is that true or false?

9 DR. DeMUTH: I guess I sort of want to answer that
10 in terms of what we were powered to show and--

11 DR. BOYAN: Before you start, state your name.

12 DR. DeMUTH: George DeMuth. I am a consultant for
13 HealthTronics. I agree, those are the statistically
14 significant endpoints. Actually, pain on walking is
15 marginal and the composite is significant. But I think
16 there is not much more you can say.

17 These other ones, it looks like there is a trend,
18 but we just don't have sample size or significance to say
19 anything about this.

20 DR. BOYAN: One more person has an opportunity to
21 ask any questions she might like to ask.

22 Dr. Butcher, are there any questions you would
23 like to ask either the FDA or the presenters?

24 MS. BUTCHER: Thank you. As the consumer rep, I
25 have paid a little bit more attention to the labeling as

1 opposed to the actual things that you have already been
2 questioned on, so I will go directly to that.

3 My first question is I see that comments have
4 already been made about labeling. There are a couple of
5 pages requesting and asking that other things be done. My
6 question is have they been done.

7 MS. MARLOW: FDA has communicated those to us.

8 MS. BUTCHER: Yes; they have.

9 MS. MARLOW: As part of the requirements to
10 develop final labeling based on your recommendations here
11 today, we have absolutely no problem with finalizing those
12 with FDA.

13 MS. BUTCHER: Okay. Well, I just have a couple of
14 comments about the draft and, not to add anything to what
15 they have already said, but as a consumer, it would appear
16 to me that, perhaps, if you discuss the target patients that
17 you were seeking to address in the first place, in the first
18 instance, that it would give them some relief to know that,
19 "Hey; I am not alone. I need to be in this group of people
20 that have had difficulty with this. It is not like this is
21 the first time we have tried to address this issue." It
22 would give them some degree of comfort in saying, "Let me
23 read on and get more information."

24 The suggestions that were made were valid. I
25 think that, basically, the draft is good. Go for it. I

1 don't have any other questions for you.

2 DR. BOYAN: Are there any other questions from the
3 panel? Seeing none, I would like to invite the FDA forward
4 to read their questions, effectively their charge to us.

5 **Panel Questions**

6 MR. OGDEN: My name is Neil Ogden. I work with
7 the FDA in the General Surgery Devices Branch. The first
8 question we have to the panel; "Although the total number of
9 complications of any type in both the active and control
10 groups were similar, there were some types of complications
11 observed in the active-treatment group that were not
12 observed in the control group. These events included neural
13 injury and irritation, plantar fascial rupture, ecchymoses
14 if the dorsum of the toes. Does this PMA safety profile for
15 the OssaTron treatment compare to the control treatment
16 adequately demonstrate the absence of unreasonable risk of
17 injury?"

18 DR. BOYAN: Go to No. 2.

19 MR. OGDEN: Question No. 2: "Do the data in this
20 PMA demonstrate that there is a reasonable assurance that,
21 in a significant portion of the target population, the use
22 of the OssaTron for its intended use and the conditions of
23 use, when accompanied by adequate directions for use and
24 warnings against unsafe use, will provide clinically
25 significant results?"

1 Question No. 3: "The sponsor proposes an
2 indication of 'The OssaTron is indicated for the use of ESW
3 treatment of chronic heel-pain syndrome in patients who have
4 had symptoms for a minimum of six months and who have failed
5 to respond to conservative treatment.'

6 "Although the term 'heel-pain syndrome' includes
7 the plantar fasciitis diagnosis, it may also be confused
8 with other etiologies like stress fracture of the calcaneus,
9 Achilles tendinitis and tarsal-tunnel syndrome. Please
10 comment on the patient population for which this device
11 should be indicated."

12 Thank you.

13 DR. BOYAN: Thank you.

14 If the person who is handling the slide could go
15 back to Panel Question 1, please. We will leave the
16 question up so everybody can see it while we are discussing.
17 I would like to see if anybody on the panel would like to
18 come forward with an answer to this question and why don't
19 we go--Dr. Pfeffer, would you like to begin the comments?

20 DR. PFEFFER: -I read this information in-detail
21 prior to our meeting and I do not consider these
22 complications consequential to affect approval of this
23 device one way or another.

24 DR. BOYAN: Why don't we come this direction. As
25 each panel member has an opportunity to address this, they

1 don't need to repeat a statement that has already been made.
2 If they agree or disagree, that is all they need to say, or
3 if they have something that they would like to add to the
4 information, add this information at this time.

5 Dr. Wright.

6 DR. WRIGHT: I agree.

7 DR. CHENG: I agree.

8 DR. YASZEMSKI: I will note that I think Dr. Ogden
9 addressed these already and there are no issues remaining
10 related to them.

11 DR. BOYAN: Dr. Finnegan? Oh; she left us. Dr.
12 Larntz?

13 DR. LARNTZ: I have no opinion.

14 DR. LEWIN: I don't have any comments.

15 DR. ROBINSON: Agree.

16 DR. BOYAN: Dr. Goldman; are you reading the
17 question?

18 DR. GOLDMAN: I know my answer. I think that the
19 assessment of nerve injury probably was not adequate to make
20 a determination, in my opinion.

21 DR. ABOULAFIA: I agree with Dr. Yaszemski.

22 DR. BOYAN: Dr. Silkaitis, do you have a comment
23 you would like to make?

24 DR. SILKAITIS: No comment.

25 DR. BOYAN: And none for you, either?

1 MS. BUTCHER: No comment.

2 DR. BOYAN: Dr. Witten, have we addressed this
3 sufficiently to the usefulness of the FDA?

4 DR. WITTEN: Yes; thank you.

5 DR. BOYAN: Okay. Let's go to Panel Question
6 No. 2. Dr. Robinson, would you like to tackle this question
7 first?

8 DR. ROBINSON: The question that reasonable
9 assurance in a significant portion of the target population,
10 its intended use and conditions; my short answer would be
11 yes, and it will provide some clinically significant
12 results.

13 DR. BOYAN: Is there anybody on the panel that
14 would like to make a further comment on this question? Dr.
15 Larntz?

16 DR. LARNTZ: I am not sure if I understand the
17 question totally so what I will say is I think this device
18 will give, has been proven to give, short-term pain relief
19 for a proportion of the individuals probably on the order of
20 less than 50 percent or-around 50 percent of the population.

21 As long as we understand that we are getting
22 short-term pain relief, that is the result of using this,
23 and no indication of duration of relief has been proven and,
24 in fact, no indication that it is a large proportion. It is
25 at least 15 percent, maybe is up to 50 percent, and gives

1 some pain relief. That is what I would say we know about
2 this.

3 DR. BOYAN: Dr. Aboulafia?

4 DR. ABOULAFIA: I would ask the same question
5 again. I hate to come across as the cynic, especially in
6 someone who doesn't routinely treat plantar fasciitis. So I
7 will certainly defer to Dr. Pfeffer and especially Dr. Ogden
8 as well to answer.

9 Convince me that this is a reasonable treatment if
10 the goal of chronic heel pain is long-term lasting results.
11 It seems to me that the effects are statistically
12 significant in one out of four parameters that were used as
13 criteria to evaluate the efficacy of the device at 12 weeks,
14 and that we have no data on anything to suggest that it does
15 provide long-term durable results and that even the 12-week
16 results showed that your heel hurts less if a doctor presses
17 on it, but you still take pain medication and your activity
18 level is unchanged.

19 Am I looking at this cynically? I don't intend
20 to. Dr. Pfeffer, or Dr. Ogden?

21 DR. BOYAN: Actually, they are off line now, but
22 we will give them an opportunity to respond back. They get
23 one more chance. Dr. Pfeffer?

24 DR. PFEFFER: The data speaks for itself. We have
25 all seen it. The only clinical comment I can make is that

1 the study in a group of patients like this longer than 8
2 weeks or 12 weeks is, perhaps, almost impossible. It is
3 just not a group that stays together and is easy to follow.

4 So your data interpretation, I would certainly
5 agree with. To say that our charge to the company would be
6 to go back to study this for six months or one-year follow
7 up, I think would be almost undoable.

8 DR. ABOULAFIA: Let me ask why. We do this with
9 county gunshot-wound patients which I think is a much
10 tougher population, frankly. We do it with intravenous drug
11 abusers with at least 50 percent follow up which I think is
12 a tougher population. So that is item No. 1, the follow up
13 beyond 12 weeks.

14 But I will even ask the question about has
15 industry shown that it is an effective device for 12-week
16 follow up? There are two parts of this. At 12 weeks, can
17 we, as a group, say--I agree that it is a safe device. Can
18 we say that, at 12 weeks, it is effective since there were
19 four parameters tested and one out of four maybe shows a
20 significant difference.-

21 DR. PFEFFER: I have no new comment about the
22 first point. Just the second point is that if someone has a
23 gunshot or someone has a cancer or some serious injury, they
24 are much more likely to stick around their doctor than
25 someone who has had two years of heel pain that has a

1 certain amount of intervention and then may have some
2 decrease in their heel pain.

3 It is a group that is just hard to follow. This
4 is a relatively minor problem, at least in my own personal
5 experience. I think a 6-month study or a year study would
6 be wonderful, but I think it would be very, very hard to do.

7 DR. ABOULAFIA: Except we are selecting for a
8 group of patients who have proven six months of follow up so
9 far, and we are selecting for a group of patients who admit
10 that they have a significant problem because their VAS is
11 5.0 or greater. So I would say, if anything, we are
12 selecting for a group of reliable patients who have proven
13 that they are willing to at least stay with the physician
14 under care and treatment--not maybe one physician, but a
15 physician, for at least a six-month interval.

16 Then there is the third selection, the physicians,
17 themselves, who have not been involved in clinical trials.
18 If a patient doesn't seem like they are going to be able to
19 participate, like they are planning on moving out of town,
20 we don't include them. -

21 DR. BOYAN: I think, to summarize this discussion,
22 we are all in agreement that the 12-week data is adequate.
23 When we get to the final voting and recommendations and
24 comments that we might want to convey to FDA, we certainly
25 are clear, also, that we don't have more long-term data to

1 rely on. So that will come out at that time, I think,
2 pretty adequately.

3 Dr. Cheng, did you have something that you would
4 like to add that is different?

5 DR. CHENG: I was going to add that I am wondering
6 if your concern--statistically, what Kinley said is correct.
7 However, clinically, it may not be that much of a problem.
8 My understanding of this disease is that, once people are
9 better, the likelihood of relapse is pretty low. I defer to
10 my colleagues if I am wrong, but that is my understanding.

11 DR. BOYAN: Thank you for your comment. Is there
12 any other comment that is directly related to Question
13 No. 2? I look around the room. Seeing none, Dr. Pfeffer,
14 if you would like to take a first stab at this one, it would
15 be great.

16 DR. PFEFFER: The separation of these diagnoses
17 clinically is very straightforward. The diagnosis of heel-
18 pain syndrome, if you will, or plantar fasciitis is made
19 easily by maximal focal pain over the medial calcaneal
20 tuberosity that may extend for a centimeter or two distally
21 along the course of the plantar fascia.

22 The most specific, and also well-recognized by the
23 public, diagnosis for this condition is plantar fasciitis
24 or, specifically, proximal plantar fasciitis. That is the
25 term I would recommend be used. It is quite distinct from

1 any of these other diagnoses. That is population that, in
2 fact, was examined by this study.

3 DR. BOYAN: Thank you. Let's just take a quick
4 opportunity to see if anybody else would like to comment on
5 Question No. 3.

6 Seeing none, we are doing very well here. Let me
7 tell you that we are doing better than anticipated, so I
8 have to make it--oh, yes; I have to ask Dr. Witten.

9 Dr. Witten, did we answer Questions 2 and 3
10 adequately for the FDA?

11 DR. WITTEN: Yes; thank you.

12 DR. BOYAN: Now, this is the deal. It is about
13 12:20. I don't want to rush the vote or the discussion
14 after the vote. So I will ask the company one thing. When
15 we come back, after we have lunch, you will be given an
16 opportunity to address anything that you feel you need to
17 clarify.

18 If you feel like you are ready to do that now, I
19 would welcome you to do that. Then, if something comes up
20 in your discussions with each other over lunch that is
21 ground-breaking that you need to bring up, you can.

22 MS. MARLOW: I appreciate that. I think the only
23 thing I would like to do is try to clarify this issue of
24 long-term follow up. I think we got a little bit derailed
25 on that issue.

1 Let me try to explain how we constructed the study
2 protocol. Patients had to stay in the study until 12 weeks.
3 When they signed their consent, they were asked to, "Bear
4 with us, don't go after any other form of treatment, come
5 back for follow up. If you are a failure at 12 weeks, we
6 recognize you are going to want to try to do something else
7 for your heel pain and you may be released from the study at
8 that point in time."

9 At that time, after that 12-week follow-up visit,
10 if the patient says, "Yes; I want to go do something else.
11 I am still having pain. I am unhappy," we said, "Great.
12 Go. But we are going to tell you your options." And then,
13 at that point in time, they were told that they were
14 eligible for a retreatment.

15 Therefore, the patients that did not opt for
16 retreatment could go and have another treatment. There was
17 no point in us following those patients thereafter because
18 we would be studying some other treatment for heel pain. So
19 it is mandatory that the successes were followed. It was
20 voluntary for the failures who did not have other treatments
21 to be followed. We continue to follow those patients.

22 We have an open IDE. We are complying with the
23 requirements of that IDE. The 25 or 26 patients you
24 referred to Dr. Larntz, that were not recorded in the PMA,
25 were not recorded merely because they had not yet reached

1 twelve weeks.

2 We have data for those patients now. We have data
3 for all the patients we continue to follow. I appreciate
4 Dr. Pfeffer's defense of how difficult this patient
5 population is and I will agree. This is a very difficult
6 patient population to follow, but we are continuing to
7 follow them.

8 Our lost-to-follow-up rate is only slightly higher
9 than I reported for the subset of patients that were lost
10 before the 12-week follow up. They are not lost. They are
11 just recorded--

12 DR. LARNTZ: Can I ask, just to make sure I
13 clarify? You are saying that 119 patients, they are the
14 only patients that were eligible for a 12-week follow up,
15 not 130.

16 MS. MARLOW: That's correct.

17 DR. LARNTZ: Okay; I didn't get that.

18 MS. MARLOW: The rest of those patients, we are
19 still following. As a matter of fact, we have enrolled a
20 few more.

21 DR. LARNTZ: Okay; I apologize for that.

22 MS. MARLOW: That's okay. It is absolutely no
23 problem.

24 DR. LARNTZ: With respect to duration, I was only
25 worried about the successes.

1 MS. MARLOW: Right; I understand.

2 DR. LARNTZ: In the successes, or 40 percent of
3 the successes didn't come back at six months.

4 MS. MARLOW: Actually, that is the same story.
5 They were just not eligible. They are out to six months
6 now. We have the data. We had to look at the data briefly
7 to let FDA know whether there were any changes. There is a
8 requirement in the regulations that say threes months before
9 the PMA--or three months after the PMA, excuse me--"Do you
10 have anything new to disclose?"

11 We had the statistician take a quick look. He
12 said, "No; there is nothing new here." So we went about our
13 business. There is nothing new to disclose. What we have
14 taken a look at for the patients we have continued to follow
15 is nothing different than what we have presented here today.

16 When FDA gets the final report on the IDE, that is
17 exactly what will be in there.

18 DR. LARNTZ: But you didn't report how many were
19 eligible at six months. You are saying everyone who was
20 eligible at six months was reported? They all came back?

21 MS. MARLOW: Anyone who was eligible was reported;
22 yes.

23 DR. LARNTZ: When you say "follow up, yes/no," I
24 apologize. I will stop. I am slowing down the process. If
25 it is true, then it certainly is a problem with my reading

1 of your report. And I apologize.

2 MS. MARLOW: It's okay. I'm sorry, but I do want
3 to try to get this cleared up. The other point I would like
4 to make is we tried very hard not to make any claims about
5 duration of results. The reason for that is even if we had
6 follow up on 500 patients, I can't give you any assurance
7 that it is because of the OssaTron treatment.

8 The natural history of the condition is that
9 people are going to get better. It may take them five
10 years. It may take them ten years. But most of them
11 eventually get better.

12 So, even if I had those data to give you, it would
13 probably still be controversial as to whether the continued
14 improvement is the natural history, treatment effect. But
15 we know that, at 12 weeks, it is probably treatment.

16 DR. BOYAN: Thank you very much.

17 We will break for lunch now. Before you get up,
18 we are going to only take a 30-minute lunch. I want to
19 remind everybody that this is a confidential proceeding, and
20 especially remind the panel that we will not be discussing
21 anything that has gone on here this morning, during lunch,
22 so, as we leave the room, this is a--oh; all right.

23 I must clarify that. I am reminded, this is an
24 open session. We want to remind the panel not to discuss
25 the topic while we are at lunch. Now, all the panel will go

at

125

1 down to the restaurant and we will reconvene here in
2 30 minutes.

3 [Whereupon, at 12:30 p.m., the proceedings were
4 recessed to be resumed at 1 o'clock p.m.]

A F T E R N O O N S E S S I O N

[1:12 p.m.]

DR. BOYAN: Right before we broke for lunch, the company was in the process of making their final comments. They have now had 35 minutes to think over their pain. Is there any last thing you would like to say? No? The answer is no.

So we can now move on to the next part of the meeting. We have one brief bit of business. We have to have a second open meeting where the people are invited to address the panel, and it is now until all those people have had an opportunity to speak.

Open Public Hearing

DR. BOYAN: We will now proceed with the open public session of this meeting. I would ask, at this time, that all person addressing the panel come forward and speak clearly into the microphone as the transcriptionist is dependent on this means of providing an accurate record of this meeting.

We are requesting that all persons making statements during the open public session of the meeting disclose which company they represent and whether they have financial interest in any medical-device company before making your presentation to the panel. In addition to stating your name and affiliation, please state the name of

1 your financial interest if any.

2 Is there anybody in the public who is wishing to
3 address the panel? Seeing none, and I have already asked
4 HealthTronics if they had any final comments before the
5 panel proceeds, and they have none. So we will go on to the
6 voting.

7 **Vote**

8 DR. BOYAN: I would like to now ask Mr. Hany
9 Demian to read the voting instructions for the panel.

10 MR. DEMIAN: I will now provide you with the panel
11 recommendation options for premarket approval applications.
12 The Medical Device Amendments to the Federal Food, Drug and
13 Cosmetic Act require that the Food and Drug Administration
14 obtain a recommendation from an outside expert advisory
15 panel on designated medical-device premarket approval
16 applications that are filed with the agency.

17 The PMA must stand on its own merits and the
18 recommendations must be supported by safety and
19 effectiveness data in the application or by applicable
20 publicly available information.

21 Safety is defined in the Act as reasonable
22 assurance based on valid scientific evidence that the
23 probable benefits to health under the conditions of use
24 outweigh any probable risks.

25 Effectiveness is defined as reasonable assurance

1 that, in a significant proportion of the population, the use
2 of the device for its intended uses and conditions of use
3 when labeled will provide clinically significant results.

4 Your recommendation options for the vote are as
5 follows: one, approval. There are no conditions. Two;
6 approvable with conditions. You may recommend that the PMA
7 be found approvable subject to specified conditions such as
8 a resolution of clearly identified deficiencies which have
9 been cited by you, the panel, or FDA staff.

10 All the conditions are discussed by the panel and
11 listed by the panel chair and then voted on one at a time.
12 For example, you may specify what type of follow-up
13 information the panel or FDA should evaluate prior to or
14 after approval. Panel follow up is usually done through
15 homework assignments to one or two primary panel reviewers
16 of the application or to other specified members of this
17 panel. Formal discussion of the application at future panel
18 meetings is usually not held.

19 If you recommend postapproval requirements to be
20 imposed as a condition of approval, then your recommendation
21 should address the following points; the purpose of the
22 requirement, the number of subjects to be evaluated, and the
23 types of reports that should be submitted.

24 The third option is not approvable. Of the five
25 reasons the Act specifies for denial of approval the

1 following three reasons are applicable to your panel
2 deliberations: the data do not provide reasonable assurance
3 that the device is safe under the conditions that are
4 prescribed, recommended or suggested in the proposed
5 labeling; reasonable assurance has not been given that the
6 device is effective under the conditions that are
7 prescribed, recommended or suggested in labeling; and, based
8 on a fair evaluation of all material facts in your
9 discussions, you believe the proposed labeling to be false
10 and misleading.

11 If you recommend that the application is not
12 approvable for any of these stated reasons, then we ask that
13 you identify the measures that you think are necessary for
14 the application to be placed in approvable form.

15 Traditionally, the consumer representative and the
16 industry representative do not vote. Dr. Boyan, as panel
17 chair, would only vote in the case of a tie.

18 Dr. Boyan?

19 DR. BOYAN: Before beginning the voting process, I
20 would like to mention, for both the panel's benefit and for
21 the record, that the votes taken are votes in favor of or
22 against the motion made by the panel. Votes are not for or
23 against the product.

24 Dr. Robinson, are you prepared to make a motion?

25 DR. ROBINSON: Yes; I am. I would make a motion

1 to approve this with no conditions. My rationale for that
2 would be briefly that this is an extremely debilitating
3 disease in some patients. It is extremely frustrating for
4 both patients and physicians in a significant number of
5 instances.

6 Although the group effect was moderate, the
7 individual effect was marked in some individuals. Thus, I
8 think this offers us a new choice for treatment in a
9 frustrating disease.

10 DR. BOYAN: We will have a chance for discussion
11 here in a second. Is there a second for the motion?

12 DR. YASZEMSKI: Second.

13 DR. BOYAN: Yaszemski seconds. Is there any
14 discussion of the motion? Dr. Robinson, do you want to
15 finish your discussion?

16 DR. ROBINSON: I was just going to address there
17 are no safety concerns for me and there are only minor
18 issues concerning labeling that I think can be worked out
19 between the sponsor and FDA.

20 DR. BOYAN: Any other comments? Dr. Finnegan?

21 DR. FINNEGAN: Actually, I think there are some
22 safety concerns. I think that particularly neurologic
23 injury and also the use of the device in certain hands. So
24 I actually cannot support that vote.

25 DR. BOYAN: Any other comments before we vote on

1 the motion?

2 All those in favor of voting for approval without
3 conditions, please raise your hand.

4 [Show of hands.]

5 I count five votes for approval without
6 conditions--six? Raise your hand again.

7 [Show of hands.]

8 Six votes for approval without conditions. All
9 those against approval without conditions, raise your hands.

10 [Show of hands.]

11 I count three votes--four votes. Why am I having
12 trouble counting. You are asking a question? Yes?

13 DR. PFEFFER: Could you outline the types of
14 conditions that might be added?

15 DR. BOYAN: No.

16 DR. PFEFFER: We are just voting on this.

17 DR. BOYAN: That comes next. So you are voting
18 which way? Against approval without conditions?

19 DR. PFEFFER: For there being some conditions.

20 DR. BOYAN: You would like there to be conditions,
21 so you are voting against this motion. So there are four
22 votes against approval without conditions.

23 Are there any abstentions?

24 [No response.]

25 Let me remind you, the motion carries. I just

1 want to remind everybody here, again, that this is a vote
2 for the motion and not for the product. The FDA hears
3 absolutely 100 percent of everything we say, and they will
4 take all of this information back and they will make the
5 final determination, not us.

6 So I think that ends this discussion; right?

7 DR. WITTEN: We need to go around the room.

8 DR. BOYAN: That's right. I forgot the very
9 special part. This is the part that they really listen to
10 so here is where you get to do your thing. We will go
11 around the room, one at a time, and everybody gets to
12 explain why they voted the way that they did.

13 DR. WITTEN: And state your vote, too.

14 DR. BOYAN: And state your vote out loud for the
15 record. Dr. Aboulafia, would you start, please.

16 DR. ABOULAFIA: I think everyone knows what I am
17 going to say because I already said it, so I will just
18 outline it very briefly. I am not concerned about safety
19 issues. I think all those things have been appropriately
20 addressed. I would add that I thought it was a well-
21 designed study and the integrity of the data is not in
22 question at all.

23 I thought labeling concerns were well addressed.
24 I think the question that was raised about who can use it
25 and who can't use it is an impossible question to answer. I

1 am licensed to administer general anesthesia. I would never
2 do it. You can't label a device for one doctor to use and
3 not another if they are a licensed, practicing physician.
4 So that was not a concern for me.

5 The only concern I had was whether the data
6 supports that it is an effective product. I have said that
7 before and there is no point repeating it.

8 DR. BOYAN: Dr. Goldman?

9 DR. GOLDMAN: Although I raised some issues
10 regarding the process of the data collection, I think that
11 it was a well-done study. I also think that it also was
12 effective at the time point with the primary endpoint.
13 Although there are trends, it is not clear that it has long-
14 term effects, although it probably does.

15 My concern, and I did approve this without
16 conditions, is that in the labeling, which would be, I
17 guess, minor issues to clarify the labeling, is that it
18 should include precautions for people that are smokers, that
19 may have microvascular disease, people with diabetes who
20 might have both microvascular disease and peripheral
21 neuropathy, and anyone else who might have a peripheral
22 neuropathy such as people with a long history of alcohol
23 use.

24 So my concerns only involve labeling.

25 DR. WITTEN: Could you just state whether you

1 voted for the motion or against it?

2 DR. BOYAN: He did. He was for.

3 DR. WITTEN: Thank you.

4 DR. ROBINSON: I voted for approval without
5 conditions and then blurted out my rationale before I should
6 have. My minor concerns are exactly what Dr. Goldman is
7 mentioning plus the fact I think the FDA and sponsor need to
8 talk about just assurance of an adequate training program.

9 DR. LEWIN: I voted for approval without
10 conditions. As I mentioned before, I primarily looked into
11 the technical specs and I was very impressed with the solid
12 and very complete documentation which the company provided.
13 They definitely know what they are doing

14 The company strongly supports training of the MDs
15 or whoever will be performing the treatment. The device
16 offers pain relief when all other treatments fail to do
17 this. Overall, I haven't seen any serious contraindication.
18 So I am convinced that they will do a good job.

19 DR. LARNTZ: I voted against approval without
20 conditions. I would have voted for approval with -
21 conditions. The conditions would have involved finishing
22 the duration analysis and making sure that no danger or
23 problem came in long-term use of the product.

24 DR. FINNEGAN: I voted against the approval
25 without conditions. I would have voted for approval with

1 conditions. My concerns are several. There is a long
2 history of instruments getting into the hands of people who
3 have not been properly trained and this actually usually
4 comes back to haunt everyone involved in the instrument.

5 I think mandatory training is essential. Also, I
6 think that there is enough data from gunshot wounds, in
7 particular, to show that shock-wave injury to nerves is a
8 problem. Around the foot, either loss of sensation or motor
9 function causing deformities is a significant problem and I
10 think there should probably be some postmarket surveillance
11 on this.

12 DR. YASZEMSKI: I voted for approval and I have
13 nothing to add to what has already been said.

14 DR. CHENG: This is Cheng. I voted for approval
15 without conditions. My only concern is dealing with
16 inappropriate usage for, perhaps, acute disease or I think
17 every practitioner who takes care of foot problems, M.D. or
18 otherwise, will use this as reimbursement issues will drive
19 them to use it more frequently. But I think that is
20 difficult to enforce through any type of condition-or,
21 perhaps, outside the FDA purview.

22 DR. WRIGHT: I voted for approval.

23 DR. BOYAN: Do you want to add any other comments?

24 DR. PFEFFER: Glenn Pfeffer. I voted against the
25 blanket approval. I certainly would have supported this

1 with conditions. The conditions that I would have are,
2 perhaps, out of the purview, however, of this panel and the
3 FDA. The conditions I would like to see is that everyone
4 who uses this device has an appropriate training course and
5 that the device not be used for the treatment of plantar
6 fasciitis or heel-pain syndrome in a patient who has
7 symptoms for less than six months.

8 Otherwise, I completely support this product.

9 DR. BOYAN: Although you didn't vote, is there any
10 comment that you would like to make, Dr. Silkaitis?

11 DR. SILKAITIS: No. I have no comment.

12 DR. BOYAN: Any comment?

13 MS. BUTCHER: No.

14 DR. BOYAN: Dr. Witten, have you received enough
15 information?

16 DR. WITTEN: Yes; and I would like to thank the
17 panel and the sponsor and the FDA presenters here.

18 DR. BOYAN: Okay; so this part of the panel
19 meeting is now--let me just make sure I am covering the
20 territory here. Ah; I have to state it over.

21 The panel is recommending that the premarket
22 approval application for HealthTronics OssaTron be approved
23 without conditions.

24 I now bring this meeting to an end.

25 DR. WITTEN: This portion. We have another

1 product.

2 DR. BOYAN: Well, yes; I know. But HealthTronics
3 is freed from captivity. It is the panel that is not free.

4 **Open Public Hearing**

5 DR. BOYAN: We will go ahead and open up the open
6 public hearing. This is an open public hearing session. I
7 would like to ask at this time that all person addressing
8 the panel come forward and speak clearly into the microphone
9 as the transcriptionist is dependent on this means of
10 providing an accurate record of this meeting.

11 We are requesting that all persons making
12 statements during the open public session of the meeting
13 disclose which company they represent and whether they have
14 financial interest in any medical-device company before
15 making your presentation to the panel. In addition to
16 stating your name and affiliation, please state the name of
17 your financial interest if any.

18 Is there anybody in the public who is wishing to
19 address the panel?

20 **Session 2: Howmedica Osteonics PMA-P000013**

21 DR. BOYAN: We will now proceed with the second
22 PMA for a ceramic-on-ceramic total hip arthroplasty. This
23 will be application P000013, Howmedica Osteonics Corporation
24 ABC/Trident Systems. We will now consider the premarket
25 approval application for the Howmedica Osteonics Corporation

1 ABC/Trident Systems.

2 I would like to remind public observers that,
3 while this portion of the meeting is open to public
4 observation, public attendees may not participate except at
5 the specific request of the panel.

6 We are now ready to begin with the sponsor's
7 presentation, followed by the FDA presentation. We have
8 several new members of the panel that we need to introduce.
9 I think what we should do is ask the Executive Secretary, is
10 there anything else you need to do before we start this
11 panel?

12 MR. DEMIAN: Actually, do you have copies of your
13 presentation? Do you have copies for the panel members?

14 MS. STAUB: We can make some but we don't
15 currently have them, no.

16 DR. BOYAN: We can handle that. What I would like
17 to do is introduce Dr. Steve Li. Steve, why don't you say
18 who you are, where you are from and what you do?

19 DR. LI: Steve Li, Senior Scientist, Hospital for
20 Special Surgery, Department of Biomechanics and Biomolecular
21 Design, New York City.

22 DR. BOYAN: I am thinking, Executive Secretary,
23 since we have a whole new company here we should probably go
24 once more around the room and have everybody introduce
25 themselves for the record.

1 MR. DEMIAN: I agree.

2 DR. BOYAN: So, let's start. Dr. Aboulafia?

3 DR. ABOULAFIA: I am still Albert Aboulafia --

4 [Laughter]

5 -- I am working at Sinai Hospital and the
6 University of Maryland, both in Baltimore.

7 DR. LARNTZ: Kinley Larntz. I am a statistician.
8 I am Professor Emeritus of Statistics at the University of
9 Minnesota and my research interests are in clinical design
10 and analysis of data.

11 DR. FINNEGAN: Maureen Finnegan, Associate
12 Professor, U.T. Southwestern. My areas of interest are
13 trauma and sports.

14 DR. YASZEMSKI: Michael Yaszemski, Departments of
15 Orthopedics and Bioengineering, Mayo Clinic; clinical and
16 adult reconstruction and spine surgery and research and
17 tissue engineering.

18 DR. BOYAN: I am Barbara Boyan. I am Professor
19 and Director of Orthopedic Research at the University of
20 Texas Health Science Center at San Antonio, and my specialty
21 is bone and cartilage biology.

22 DR. CHENG: My name is Edward Cheng. I am with
23 the Department of Orthopedic Surgery at the University of
24 Minnesota, and my interest is in reconstructive surgery,
25 muscle cell oncology and osteonecrosis.

1 DR. LYONS: I am John Lyons. I am from Erie,
2 Pennsylvania. I am a private practice orthopedist. Adult
3 reconstruction is my area of interest and I am also a
4 biomedical engineer.

5 DR. SILKAITIS: My name is Raymond Silkaitis. I
6 am the industry rep., a non-voting member of the panel. I
7 have a Ph.D. in pharmacology and I am a registered
8 pharmacist.

9 MS. BUTCHER: My name is Vicky Butcher. I am --
10 what am I? I am the consumer rep., also a non-voting
11 member. My background is in teaching the law, and I have
12 served as the consumer consortium member for the FDA.

13 DR. BOYAN: Thank you. If I have already read
14 this, that is how it is! I would like to remind public
15 observers that while this portion of the meeting is open to
16 public observation, public attendees may not participate
17 except at the specific request of the panel.

18 We are now ready to begin with the sponsor's
19 presentation, followed by the FDA presentation. I would
20 like to ask that each speaker state his or her name and
21 affiliation to the firm before beginning the presentation.

22 The sponsor's presentation will include an
23 introduction by Beth Staub and Michael Manley; product
24 description by Thomas McCarthy; laboratory testing by
25 Michael Bushelow; implantation technique by James D'Antonio;

1 clinical data by Michael Manley and, finally, summary and
2 conclusions again by Michael Manley.

3 **Sponsor Presentation**

4 **Introduction**

5 MS. STAUB: Good afternoon.

6 [Slide]

7 My name is Beth Staub, and I am the vice president
8 of quality assurance, regulatory affairs and clinical
9 research for the Howmedica Osteonics Corp.

10 [Slide]

11 On behalf of Howmedica Osteonics, we are all
12 pleased to be here this afternoon to present safety and
13 efficacy data that we have collected demonstrating the
14 safety and efficacy of two alumina on alumina ceramic
15 bearing surfaces, ABC and Trident.

16 Since their introduction in the '70s, the wear
17 resistance and biocompatibility of ceramic couplings has
18 been widely reported, and significant improvements have been
19 made to materials and manufacturing processes. The ABC and
20 Trident Systems incorporate these improvements, as-well as
21 design elements that are important in ceramic/ceramic total
22 hip arthroplasty.

23 [Slide]

24 In 1996, the Osteonics Corporation initiated a
25 prospective, controlled, randomized multi-center trial to

1 collect clinical data on the ABC alumina bearing system.

2 This study included two styles of acetabular shells, an HA,
3 or hydroxyapatite-coated femoral stem, and a polyethylene
4 control group.

5 [Slide]

6 In 1999, Howmedica Osteonics obtained FDA approval
7 to begin a supplement to the original ABC study, evaluating
8 the Trident bearing design. The trident bearing differs
9 from the ABC in its locking mechanism. This enhanced
10 locking mechanism affords the surgeons more revision
11 options, and helps protect the ceramic component from
12 chipping during insertion. The surgeons with the greatest
13 number of implants in the original ABC study were selected
14 for the Trident arm, and the same polyethylene control group
15 was used.

16 [Slide]

17 This table provides an overview of the components
18 used in the study. All are commercially available by the
19 510(k) process with the exception of the two we are asking
20 the panel to recommend for approval today, the ABC-and
21 Trident alumina inserts.

22 ABC System I used a microstructured or porous
23 coated shell, the ABC alumina insert with the aluminum head
24 and the HA hip stem. The ABC System II shell featured a
25 Secur-Fit coating, titanium Arc-Deposited with an HA-coated

1 surface. All other components in System II are identical to
2 those in System I.

3 The Trident arm of the study was similar to the
4 ABC System II, using an Arc-Deposited HA-coated shell but
5 substituting the Trident locking mechanism on the acetabular
6 components.

7 The control group for both the ABE and Trident
8 arms received a microstructured shell, along with the
9 standard polyethylene liner and cobalt chrome femoral head.

10 [Slide]

11 Today we will be presenting a summary of our
12 clinical and non-clinical data. Dr. Michael Manley,
13 Howmedica Osteonics chief scientific advisor, will moderate
14 our program, and I will now turn the agenda over to Mike.

15 Introduction

16 DR. MANLEY: Good afternoon. I am Michael Manley,
17 chief scientific advisor for Howmedica Osteonics.

18 This presentation is in five parts. First, we
19 will discuss the specific design of the devices used;
20 secondly, the lab testing that was done on those devices;
21 thirdly, the implantation technique for the devices; fourth,
22 a summary of the clinical data, and, finally, within the
23 discussion we will address the questions raised by FDA.

24 [Slide]

25 To describe the design of the devices used, let me

1 call on Tom McCarthy, project engineer of the acetabular
2 team at Howmedica Osteonics. Tom?

3 **Product Description**

4 MR. MCCARTHY: Thank you, Mike. I would like to
5 take the next few minutes to describe the designs of the
6 implants associated with this study.

7 [Slide]

8 We call the ceramic-on-ceramic design the ABC
9 System. ABC is an acronym that stands for alumina bearing
10 couple. The study consisted of three separate systems
11 referred to as Systems I, II and III.

12 Systems I and II had ceramic bearing surfaces,
13 while System III, the study control, had a polyethylene
14 bearing surface. All implant combinations in this study
15 used the same femoral stem, an Osteonics Omnifit HA-coated
16 stem. And, all acetabular shells had the option of using
17 the same screws.

18 [Slide]

19 I will start out by describing the study control.
20 It consisted of an Osteonics PSL microstructured titanium
21 shell with porous coated beads, a polyethylene liner, and a
22 cobalt chrome femoral head. Only neutral poly liners gamma
23 irradiated in an inert atmosphere were used.

24 [Slide]

25 System I consisted of an ABC PSL microstructured

1 titanium shell with porous coated beads, an alumina
2 acetabular insert, and an alumina femoral head.

3 I want to note that the ceramic inserts are the
4 only components under IDE investigations in this study as
5 all other components have been 510(k) cleared.

6 [Slide]

7 In System II, an ABC Secur-Fit Arc-Deposited
8 coated titanium shell with HA was used. The alumina femoral
9 head and insert bearing combination is the same as in System
10 I. [Slide]

11 A unique feature with the ABC shells is the
12 ceramic protection rim. The lip of the shell extends beyond
13 the surface of the ceramic insert in order to protect the
14 ceramic insert from neck impingement and possible damage.

15 [Slide]

16 A cementable poly liner was also offered as a
17 revision option for Systems I and II which used the ceramic
18 insert. It is cemented directly into the shell.

19 [Slide]

20 The study arm-used ceramic components from the
21 Trident System. The design goal of Trident was to enhance
22 the ABC System with design advantages beyond the ceramic
23 insert.

24 The Trident System design allows for ease of
25 ceramic liner insertion, intraoperative flexibility; and

1 more revision options. The shell allows for independent
2 locking of ceramic and poly inserts, both of which can be
3 chosen intraoperatively or at time of revision.

4 [Slide]

5 The Trident study arm, again, used the same
6 femoral stem, ceramic heads, and screws as with ABC. In
7 addition, a 36 mm femoral head was used. The shell was an
8 Arc-Deposited coated titanium shell with HA. The shell
9 outside geometry and coating were identical to that of the
10 ABC System II. The ceramic insert has a permanently
11 assembled titanium sleeve on the outside with the same
12 ceramic insert protection rim feature.

13 [Slide]

14 On the left you see a cross-section of an ABC cup,
15 and on the right, a Trident cup. The metal sleeve has a
16 taper to taper fit within the shell, as is shown right
17 there. We were able to add a metal sleeve essentially
18 without decreasing the thickness of the ceramic insert. The
19 most important features remain the same. Both ABC and
20 Trident have identical bearing surface dimensions and
21 tolerances for both the ceramic head and the inserts;
22 identical range of motion and ceramic protection rim; and
23 identical shell contact with the bone. Thank you.

24 [Slide]

25 DR. MANLEY: Michael Manley again. Thank you,

1 Tom. I would now like to introduce Michael Bushelow, who is
2 assistant director of device evaluation at Howmedica
3 Osteonics, who will discuss the laboratory testing behind
4 these devices. Mike?

5 **Laboratory Testing**

6 MR. BUSHELOW: Thank you, Mike.

7 [Slide]

8 I will be spending a few minutes describing some
9 of the mechanical testing and analyses that were preformed
10 to ensure safety of both the ABC and Trident acetabular cup
11 systems.

12 Specifically, I will describe the strength testing
13 that was performed on both system designs based upon ISO,
14 ASTM and FDA procedures and guidelines, as well as the
15 finite element analysis that was performed to look at bone
16 stresses at the fixation interfaces.

17 [Slide]

18 Presently, there are no standard test methods for
19 evaluation of ceramic inserts used in acetabular cups.

20 Therefore, standard test methods for evaluation of ceramic
21 femoral heads were modified and used for these evaluations.

22 The test methods used include ultimate compression
23 strength testing, axial fatigue strength testing, and post-
24 fatigue ultimate compression strength testing.

25 The figure on the slide shows the general test

1 set-up used for all three test methods. The methods are
2 based upon ISO standard 7206-5 and the present ASTM draft
3 standard for evaluation of ceramic femoral heads.

4 Based upon the FDA guidance document for ceramic
5 femoral heads, the performance requirements for the ceramic
6 inserts were established. Inserts must have an average
7 ultimate compression strength value of 46 kN with no single
8 insert having a strength lower than 20 kN. Inserts must
9 survive 10 million cycles of fatigue at loads between 1.4
10 and 14 kN. Finally, inserts that have been axially fatigued
11 must have UCS values greater than 20 kN.

12 [Slide]

13 This slide shows the results of the ultimate
14 compression strength testing. Please note that the ABC
15 System is available in only 28 mm and 32 mm sizes, while the
16 Trident has an additional 36 mm size insert.

17 The recommended FDA guidance performance standard
18 is indicated by the yellow line on the graph. It can be
19 seen that all components exceed this standard with average
20 strength values ranging from approximately 56 kN to 67 kN.
21 It should be noted that in all cases individual inserts had
22 strengths greater than 43 kN.

23 [Slide]

24 This slide shows the results of the axial fatigue
25 testing and the post-fatigue ultimate compression strength

1 testing. All had survived 10 million cycles without failure
2 at the 1.4 kN to 14 kN load range.

3 The recommended FDA guidance performance standard
4 for the post-fatigue ultimate compression strength testing
5 is indicated by the orange line on the graph. Results show
6 that the inserts post-fatigue ultimate compression strength
7 significantly exceed the performance standard, yielding
8 average strength values ranging between 52 kN and 62 kN.

9 [Slide]

10 This slide summarizes the data from the previous
11 two slides, showing all of the pre- and post-fatigue UCS
12 data for the 32 mm size inserts. This graph shows that the
13 ultimate compression strength for both the ABC and Trident
14 inserts is minimally, if at all, affected by the fatigue
15 loading. Note that minimal differences are shown when
16 comparing pre- and post-fatigue ultimate compression
17 strength values.

18 Secondly, it should be noted that all testing,
19 both pre- and post-fatigue showed average ultimate
20 compression strength values that exceeded the 46 kN limit,
21 with no single head showing a value less than 20 kN. In
22 other words, all tested inserts from both the ABC and
23 Trident Systems exceeded the FDA guidance document pre-
24 fatigue performance criteria.

25 [Slide]

1 Finally, this graph is shown to put some
2 perspective on the performance standards used to evaluate
3 these ceramic inserts. Showed in red and yellow are
4 performance standards presently used when evaluating other
5 components within the total hip system. ISO 7206-8 is the
6 accepted performance standard, 2300 N, for the body region
7 of a total hip arthroplasty femoral component. While no
8 standards are available for the neck region, CeramTec has
9 proposed a value of 4300 N, while at Howmedica Osteonics we
10 use a minimum value of 5300 N.

11 The next two bars represent the FDA guidance
12 performance standards for ceramic femoral heads, once again
13 a fatigue strength of 14 kN and an average ultimate
14 compression strength of 46 kN.

15 Finally, the last three green bars show results of
16 the static testing of components, including the alumina
17 femoral head used in both the ABC and Trident acetabular cup
18 systems.

19 [Slide]

20 I am going to switch gears here and discuss
21 analyses that were performed to look at the fixation
22 interfaces and compares stresses in the bone due to
23 implantation of the ABC and Trident acetabular cups.

24 All analyses were performed using the ANSYS finite
25 element software. All models were 2-dimensional

1 axisymmetric models and, therefore, represent full 3-
2 dimensional geometry. The models included surface-to-
3 surface contact elements at the bone/shell interface. The
4 models looked at stresses in the bone due to initial press-
5 fit of the acetabular cups.

6 Multiple analyses were performed to look at the
7 effect of overall stiffness of the components and stresses
8 at the fixation interface. Shell stiffness was analyzed as
9 a function of shell design, ABC and Trident, and within each
10 shell design by modifying the insert and shell materials.

11 The ABC and Trident inserts were evaluated with
12 both alumina and polyethylene inserts. Additionally,
13 analyses modeling the shells as cobalt chrome alloy were
14 performed to look at a case where the shell is as stiff as
15 the presently available components.

16 [Slide]

17 Shown here are the two models used for the
18 evaluations. The upper left model shows the ABC
19 configuration and the lower right model shows the Trident
20 configuration. Material properties of the various-
21 components were modified to look at the effect of shell and
22 insert material on both stress.

23 [Slide]

24 A typical maximum principal stress contour plot
25 for the ABC and Trident cup systems are shown in this slide.

1 Please note that the stress patterns and magnitudes are
2 similar for both designs. These results remained consistent
3 for all design and material variations analyzed.

4 [Slide]

5 The bar chart shown here shows the maximum
6 principal stress calculated for each analysis performed.
7 Results indicate that bone stresses at the shell fixation
8 interface are similar at both the ABC and Trident cup
9 systems. Additionally, the insert material and shell
10 material had little effect on the stress magnitude and
11 distribution at the fixation interfaces.

12 [Slide]

13 In summary, variables investigated included shell
14 design, shell material and insert material. Results
15 indicated a maximum 5 percent difference in bone stress due
16 to these design variables, and it can be concluded that for
17 these cup designs overall there was minimal effect on bone
18 stress magnitude and distribution at the fixation
19 interfaces.

20 [Slide]

21 Presented today are results from several of the
22 tests and analyses that were performed to ensure the safety
23 of both the ABC and Trident acetabular cup systems.
24 Mechanical testing, recommended by ISO, ASTM and the FDA for
25 ceramic heads, and performance standards recommended by the

1 FDA were modified and used to evaluate the alumina inserts
2 of these two systems.

3 Results indicate that the inserts easily exceed
4 the recommended parameters standards. Finally, analyses to
5 look at the bone/shell fixation interface showed similar
6 bone stress magnitudes and distribution when comparing
7 between the ABC and Trident Systems. Thank you.

8 DR. MANLEY: Michael Manley again. Thank you,
9 Mike.

10 [Slide]

11 Dr. James D'Antonio, an orthopedic surgeon from
12 Sewickley, Pennsylvania, will now describe the surgical
13 technique used with these ceramic/ceramic bearings. Dr.
14 D'Antonio has an academic appointment to the University of
15 Pittsburgh, and is a principal investigator in the ABC
16 study. Jim?

17 **Implantation Technique**

18 DR. D'ANTONIO: Thank you, Mike. My name is James
19 D'Antonio, and I have disclosed financial arrangements with
20 Howmedica Osteonics, which do include some compensation for
21 consulting services. I also own stock in the parent
22 corporation, Stryker.

23 At this time, I would like to outline for you the
24 surgical technique that is used for the implantation of both
25 the ABC and Trident hip systems.

1 [Video presentation]

2 The preparation of the acetabulum is by standard
3 technique using hemispherical reamers. This first case is a
4 38-year old active man, disabled from osteoarthritis. The
5 hemispherical reamers prepare the acetabulum in an under-
6 sized fashion so that when this implant is placed it will
7 fit snugly.

8 This implant is from the System II, the Secur-Fit
9 implant with an HA coating. It will now be impacted into
10 the prepared cavity within the acetabulum, and following
11 impaction, again, will achieve a snug, tight fit. The
12 positioning and placement is important.

13 Shown now is the insertion of the ABC alumina
14 shell. It is first softly placed within the metal liner. A
15 finger is passed around the rim to ensure that it is fully
16 seated before impaction. Once that is ascertained, then
17 impaction to secure the peripheral taper lock occurs.

18 We then irrigate and flush and examine to make
19 sure it is seated and if, indeed, there were any problems or
20 defects or chips that would easily be identified by staining
21 of the body fluids. One sees there the implantation of the
22 shell and the alumina ceramic liner.

23 The femoral side is prepared with tapered reamers
24 and broaches. It is machined in a fashion to receive this
25 tapered titanium stem that is coated with hydroxyapatite in

1 the proximal third. This stem has a 13-plus year track
2 record now, with a 99-plus percent fixation success.

3 The alumina ceramic femoral head is then impacted,
4 again, to secure the taper lock. This shows the system
5 fully implanted. This is the ABC System.

6 [Slide]

7 This is the postoperative x-ray in this individual
8 patient.

9 The second case is that of a 31-year old woman, a
10 large woman, 6 ft., 190 lbs., college basketball player who
11 now is a high school basketball coach. She had a slipped
12 capital epiphysis when she was a teenager and now has
13 disabling osteoarthritis.

14 [Video presentation]

15 The acetabulum was prepared in the same fashion.
16 The socket here has the same outside geometry and surface as
17 the previous one that I showed you, the Secur-Fit titanium
18 Arc-Deposited with an HA coating. It also will then be
19 placed in a slightly under-sized hemispherical cavity to
20 create a very tight, snug fit once it is placed and impacted
21 into its proper orientation.

22 Once that is accomplished, the inside of the shell
23 will be irrigated, cleansed, and cleared of any debris, as
24 with the ABC System.

25 Here, the Trident insert, the alumina insert which

1 shows the metal backing as you have seen, will now be placed
2 first by hand. One of the advantages of the Trident is that
3 this is technically easier to do than with the ABC System.
4 The locking tabs will be lined up. It will gently be seated
5 by hand, and then the impaction device placed to fully seat
6 and impact to engage and secure the peripheral taper lock.

7 Here, again, after irrigation by inspection one
8 assures himself of full seating and lack of any injury to
9 the implants.

10 This shows the fully seated Trident System in
11 place after reduction of the hip.

12 [Slide]

13 This is the postoperative x-ray for the Trident
14 System.

15 [Video presentation]

16 Finally, I would just present the case
17 presentation. This is a 55-year old, 6'3", 230 lb. dentist.
18 He received two total hips with alumina on alumina ceramic
19 bearings at the age of 52. He works daily, walks the golf
20 course two or three times per week. He is typical of the
21 active and physically demanding patients in my practice who
22 have received the alumina on alumina ceramic bearings.
23 Thank you, Mike.

24 **Clinical Data Summary**

25 DR. MANLEY: Michael Manley again.

1 [Slide]

2 I will now present a summary of the clinical data.
3 These data are based on the three-month update document sent
4 to FDA. Panel members may refer to Amendment 4, Volumes 1
5 and 2 of their documents for reference. As the panel has
6 the full data set, I will go through these data fairly
7 rapidly.

8 [Slide]

9 To remind you, these are the three systems under
10 test. System I has an alumina/alumina bearing. System II
11 has an alumina/alumina bearing and the control system has a
12 polyethylene liner with a cobalt chromium alloy head.

13 [Slide]

14 The randomization procedure is outlined fully in
15 your documentation. The important point to note is the
16 surgical site called the sponsor with patient identification
17 etc., and day of surgery, and the randomly selected implant
18 was then sent to the site.

19 [Slide]

20 the inclusion-criteria are also in your document.
21 The important point here is that all of the patients in this
22 study were non-inflammatory joint disease cases.

23 [Slide]

24 Investigational sites are listed. You will note
25 there are only two of them with less than ten implantations.

1 A total of 515 implants were placed in each of the study
2 arms.

3 [Slide]

4 Here is the number of cases with two-year follow-
5 up. These are the data that we are presenting today. In
6 System I there were 131 patients with two years of follow-
7 up; in System II, 129; and in the control group, 119.

8 [Slide]

9 I have put some of these data as pie charts simply
10 because it is easier to visually compare the data sets.
11 This is diagnosis by system. You will note that in System I
12 the majority of the patients have osteoarthritis. This is
13 also true of System II and also of the control group.

14 [Slide]

15 Here is the demographic breakdown by age of the
16 three systems. You will notice that the randomization
17 procedure did, in fact, assign patients to every age decade
18 from 21 to 75. The mean age of this study population is
19 quite young for total hips, around 55 or a little less than
20 55.

21 [Slide]

22 Here are the demographics by gender for the three
23 systems. You will note that in all three systems a majority
24 of the patients were male. This is consistent with the
25 literature for young groups of patients where the majority

1 of total hip patients happen to be male.

2 [Slide]

3 Let's quickly go through the results, firstly the
4 clinical findings.

5 [Slide]

6 Here is the follow-up available for the study.
7 You will note that in all three systems more than 90 percent
8 of the patients are available for follow-up two years after
9 surgery.

10 [Slide]

11 When we look at this group of patients and look at
12 their Harris Hip Scores -- the preoperative score is shown
13 here and the two-year data on the far right -- we see that
14 there is a statistically significant increase in Harris Hip
15 Score over that two-year period, and at the two-year period
16 all three systems are equivalent.

17 [Slide]

18 This data set is the patients' own rating of their
19 satisfaction of the hip system from the Hip Society scoring
20 system. If you look at System II, you see that the patients
21 are completely satisfied with that system. But when we scan
22 through all of the data, we find equivalence between System
23 I, System II and the control for increased function,
24 decreased pain, less pain medication and satisfied patient.

25 [Slide]

1 A different type of patient satisfaction measure,
2 this is the HSQ-12 scoring system from the Medical Outcomes
3 Institute. Here are the preop scores for the three groups,
4 and here are the postop scores at 24 months. You will note
5 that there is a statistically significant increase in score.
6 The three systems, however, are equivalent to one another at
7 two years follow-up.

8 [Slide]

9 This is a very busy slide. This is operative site
10 adverse events by category and by system. Let me take you
11 through the top line first. The top line is alumina insert
12 chips. We note that in System I two inserts were chipped at
13 the time of insertion; two inserts in System II and, of
14 course, the controls do not suffer from this problem. I
15 will come back to the issue of insert chips again later.

16 If we look at the other adverse events, we find
17 that Systems I, II and control are equivalent to one
18 another, and these data are fairly typical of a total hip
19 series in this sort of patient population.

20 [Slide]

21 The time course distribution of operative site
22 adverse events shows that the majority of them in all three
23 systems occurred within the first three months post-surgery.

24 [Slide]

25 The definitive endpoint, of course, of hip

1 replacements is revisions or reoperations. These data are
2 shown on this slide. You note that for revisions and
3 reoperations there seem to be less revisions and
4 reoperations in Systems I and II compared to the control.
5 For reoperations the data sets are equivalent to one
6 another.

7 [Slide]

8 Here is a breakdown of these revisions. In System
9 I, there was 1/140 revised; in System II, 2/140; and in the
10 controls, 5/133, for a total number of 8 revisions.

11 [Slide]

12 Here is the reason for those revisions in System
13 I. The one implant that was revised was a postoperative
14 femoral fracture at nine months post-surgery. In System II,
15 there were two revisions, one for deep joint infection at
16 ten months post-surgery, and one for recurrent dislocation
17 soon after surgery, at five days. The five revisions in the
18 control group are, one for postoperative femoral fracture;
19 one for leveling discrepancy; one for deep joint infection;
20 and two for recurrent dislocation.

21 [Slide]

22 Now I would like to turn to the radiographic
23 findings. The reviewer of the radiographic data was Dr.
24 Peter Bonutti, from Effingham, Illinois. Dr. Bonutti has an
25 academic appointment at the University of Arkansas. I

1 should note here that he reviewed all these radiographs in a
2 blinded fashion.

3 [Slide]

4 Here is a sample radiograph from the ABC System I.
5 This is in a 30-year old female. Here is the seven-week
6 film and here is her film at two years follow-up. You
7 notice that there are no lucent lines; no implant migration.
8 This is consistent with a stable hip.

9 [Slide]

10 Here is a 32-year old male in System II, here at
11 seven weeks follow-up and here at two years follow-up. We
12 note in this particular instance that the surgeon chose to
13 use bone screws with the acetabular components. Again,
14 there are no lucent lines; no implant migration -- another
15 stable result.

16 [Slide]

17 Finally, a control implant at seven weeks follow-
18 up, and the two-year follow-up on the right. This is a 55-
19 year old male. Again, the findings are no lucent lines, no
20 implant migration and a-stable result.

21 [Slide]

22 This is the number of radiographs available for
23 follow up. This is the percentage of radiographs from those
24 patients that returned for their two-year clinical
25 evaluation. So, more than 90 percent of those patients had

1 films that were evaluable for the radiographic review.

2 [Slide]

3 One of the failure criteria outlined in the study
4 is femoral radiolucencies greater than or equal to 2 mm at
5 two years follow-up. We see that in System I there is one
6 patient with a radiolucency in zone 1 of the Gruen zones; in
7 System II, one implant had a distal radiolucent line; and in
8 the control group one implant also had a distal radiolucent
9 line.

10 [Slide]

11 Cortical erosion is related probably to release of
12 debris from the articulation. We see that in System I one
13 patient was diagnosed as having cortical erosion in zone 1.
14 There were no cases in System II. In the controls, one case
15 had cortical erosion in zones 1 and 7.

16 [Slide]

17 Here is that particular case from the control
18 group. Here is the seven-week film and here is the
19 patient's two-year film. We see here an area of cortical
20 erosion close to the resection level, and here also in the
21 greater trochanter. This is the only case of this type in
22 the control group.

23 [Slide]

24 Here is the case from System I. The cortical
25 erosion was read to be here, within the greater trochanter,

1 at six-months follow-up although the reviewer said that at
2 three years follow-up that now seems to be a radiolucency
3 and not an area of cortical erosion at all. It is, however,
4 reported as cortical erosion in your documents.

5 [Slide]

6 Turning to the acetabulum, again radiolucencies
7 greater than or equal to 2 mm at two years follow-up, in
8 System I there was one in the dome of the implant. There
9 were none in system II and in the control group one patient
10 had radiolucency superior to the cup.

11 [Slide]

12 If we look at stability using the criteria
13 outlined in your document, there were no unstable implants
14 at two years follow-up.

15 [Slide]

16 So, in summary for the ABC study at two years
17 follow-up, here are the success/failure criteria for
18 revision. There was one case in System I, two in System II
19 and five in the control group.

20 Patients with Harris Hip Scores of less than 70
21 were also taken to be failures. There were two in System I,
22 two in System II and three in the control group.

23 The other failure causes are equivalent to each
24 other, and there were no findings, except there was one
25 patient in System II who had a femoral component subsidence.

1 This was after a traumatic event following surgery.

2 [Slide]

3 The radiographic criteria we have presented
4 suggest that these hip components utilizing the ABC and the
5 control bearings are consistent with stable, pain-free hips.

6 [Slide]

7 The clinical criteria -- the Harris Hip Scores
8 indicate that at two years follow-up System I and System II
9 is at least equivalent to the control. For the overall
10 success/failure rates these also suggest that System I and
11 System II are at least equivalent to the control implants.

12 [Slide]

13 Here are the Kaplan-Meier survivorship curves for
14 these three different systems. From these data, it seems
15 that Systems I and II are at least equivalent to the
16 control.

17 [Slide]

18 So, in conclusion from the ABC study, we believe
19 the study demonstrates an equivalent performance of hips
20 with the ABC bearings to that of hips with control-cobalt
21 chromium polyethylene bearings at two years post-
22 implantation.

23 [Slide]

24 I would like to now turn to the Trident arm of the
25 study.

1 [Slide]

2 Just to remind you, this is the Trident System.
3 The ceramic insert is backed with a shrunk-fit metal shell
4 which is permanently fitted to the articulation, and then
5 the metal shell is implanted into the acetabular component.

6 [Slide]

7 Here is a simple radiograph of a Trident case at
8 six months post-surgery.

9 [Slide]

10 The investigators in this portion of the study
11 were those investigators who had done the most ABC cases.
12 This is the total number of ABC cases each investigator
13 performed. Here are the controls and here are the numbers
14 for Trident.

15 [Slide]

16 This is the total number of cases implanted
17 throughout the time of the study. For the ABC System there
18 were 173; System II, 177; Trident, 159; and 165 controls.

19 [Slide]

20 Here is the diagnosis breakdown for Trident. We
21 see that the majority of the cases are osteoarthritis, as
22 were the other three systems.

23 [Slide]

24 When we look at demographics by age of the Trident
25 superimposed on the Systems I, II and control demographic

1 breakdown, we see that they are similar to one another,
2 although no Trident implants have been implanted in the
3 earliest age decade.

4 [Slide]

5 Breaking down gender demographics by system, we
6 see once again that there is a majority of males in the
7 Trident group, as there were in the other systems.

8 [Slide]

9 If we look at revision/reoperation by study system
10 within 75 days of surgery, we see that for Trident there
11 were no revisions or removals within this time period, as
12 there were no System I removals within this early time
13 period. The different systems are equivalent to one
14 another. The same is true of reoperations. At 75 days
15 follow-up the systems are equivalent to one another.

16 [Slide]

17 Operative site adverse events by study system,
18 with Trident superimposed here on the other systems, we see
19 that for Trident the intraoperative adverse events are less
20 than they were for Systems I and II. The reason for this is
21 that there are no cracks of the acetabular liner with
22 Trident. The Trident articulation is completely protected
23 by the metal backing.

24 [Slide]

25 So, here is the data for the operative site

1 adverse events by category for Systems I, II, control and
2 Trident. If we go back to intraoperative insert chips, four
3 of those occurred with System I in the entire group; four
4 with System II. This does not happen with polyethylene.
5 But none have occurred with Trident. For the other
6 operative site adverse events, the different systems are
7 equivalent to one another.

8 [Slide]

9 So, in summary for Trident, if we compare Trident
10 to the ABC System, we have shown that the demographics are
11 similar to one another, patient demographics for Trident and
12 ABC.

13 Adverse events, there were less for trident than
14 there were for ABC. Revisions within 75 days of surgery --
15 there have been none for trident; there was 0.3 percent for
16 ABC, and reoperations -- there have been none for Trident
17 versus 1.2 percent for ABC.

18 [Slide]

19 When we continue the comparison of the systems,
20 the articulating bearing surfaces between Trident and ABC
21 are identical to one another. We have shown that the
22 stresses on the bone for Trident and ABC are equivalent, and
23 both Trident and ABC meet FDA's standard for alumina femoral
24 heads of 46 kN.

25 [Slide]

1 We looked at the risk/benefit analysis of Trident
2 and ABC. One of the potential risks is breakage of the
3 alumina insert. There have been none for ABC and none for
4 Trident.

5 Disassembly of modular components -- none for ABC,
6 none for Trident.

7 Revision options -- for ABC the revision option is
8 the cementable polyethylene insert into the metal shell; for
9 Trident the option is to use either a polyethylene insert
10 which fits the shell or an alumina insert which fits the
11 shell. This gives the surgeon greater scope if revision
12 surgery is needed. As far as intraoperative chipping is
13 concerned, there were 3.4 percent for ABC and none for
14 Trident.

15 [Slide]

16 So, in summary, there are minimal risks, we
17 believe, with ABC. We believe also there are fewer risks
18 with Trident. The advantage of Trident is that the titanium
19 sleeve protects the ceramic insert. The dual locking
20 mechanism adds greater versatility to the system -
21 intraoperatively, and the system has multiple revision
22 options.

23 [Slide]

24 I would like to now turn to the three questions
25 raised by FDA with regard to these systems.

1 [Slide]

2 The first question relates to the issue of insert
3 chipping.

4 [Slide]

5 There are four issues here that we have to address
6 with insert chipping. The first is clinical consequences if
7 a chip occurs. The second is implantation technique -- how
8 do these chips occur and how are they prevented? The third
9 is labeling, and the fourth is training and education. For
10 the first two, clinical consequences and implantation
11 technique, I would like to ask Dr. D'Antonio to address
12 those two issues.

13 DR. D'ANTONIO: Thank you, Mike.

14 [Slide]

15 This slide shows an example of a typical case.
16 There is a peripheral chip here. The fragment is laid
17 inside, there, and just taped so that you can see the
18 fragment inside.

19 Mike talked about 4 chips in that group of
20 patients who had a minimum 2-year follow-up. In fact, there
21 have been a total of 16 chips in the ABC study, 9 of which
22 occurred in the study group and 7 of which occurred in the
23 continued access group.

24 In this group of 16, 3 of the liners were left in
25 place in the patient and continue to remain in the patient.

1 All of the others were replaced. Now, on x-ray review of
2 all of these cases there is no evidence of any retained
3 ceramic fragments. All of the implants are stable without
4 signs of migration, reactive lines or lysis.

5 On a clinical evaluation, all but one of the
6 patients is doing well. The one patient with a low Harris
7 Hip Score has diffuse pain, including lower extremity pain
8 and including pain in the operative side. The pain is of
9 unknown etiology and x-ray review by four orthopedic
10 surgeons has shown that the implants are well positioned and
11 appear to be secure, without any adverse findings.

12 These cases have not resulted in an increased
13 operative time, an increase in either the total
14 intraoperative or postoperative complication rates of any of
15 the study groups, and it is my opinion that if chipping does
16 occur the additional risk to the patient is minimal as long
17 as the ceramic liner is not left in a canted position within
18 the metal shell. Chipping occurs at a very low rate and
19 appropriate physician warning, as well as training, should
20 address this issue in the future when this device is used.

21 I would like now just to go over the insertion
22 technique to give you an idea of what we are talking about
23 and how these occur.

24 [Video presentation]

25 This illustrates placing the ABC liner within the

1 shell, passing your finger around to make sure that it is
2 seated by about 2 mm circumferentially all the way around.
3 Then, the impactor can be placed and the taper lock secured.

4 I will now position this in a canted position.
5 This canted position is extreme and, if you could imagine,
6 it could be a lot less than this and create a small chip.
7 That is canted and not fully seated. If you now apply
8 pressure and try and force that into place, then a little
9 peripheral chip will occur in this area. If, indeed, these
10 do get canted, then simply tapping on the metal rim loosens
11 them and you can then softly seat them with your index
12 finger and then secure them with the peripheral taper lock.

13 DR. MANLEY: Michael Manley. I would now ask Beth
14 Staub to address the issue of labeling and training and
15 education.

16 MS. STAUB: While our initial instructions for use
17 in the surgical protocol that accompanied the study devices
18 did touch on the issue of proper alignment of the liner in
19 the shell, there was no specific reference to chipping,
20 which we were not aware-would be an issue at the time. We
21 believe we can develop a program of labeling and surgeon
22 education that will discuss the possibility of chips and
23 provide surgical techniques, such as those discussed by Dr.
24 D'Antonio, on how to avoid them.

25 DR. MANLEY: Thank you, Beth.

1 [Slide]

2 Let's turn now to question number two. This
3 question relates to the Trident data and whether the seven-
4 week data available from Trident could be representative of
5 the two-year data from the ABC System. Again, Beth Staub
6 should answer this question.

7 [Slide]

8 MS. STAUB: The FDA has published draft guidance
9 identifying the least burdensome approach to premarket
10 approval. We believe that the combination of mechanical and
11 clinical data that we have provided on Trident adequately
12 addresses any potential risks.

13 [Slide]

14 ABC and Trident are similar designs, use similar
15 components and have identical bearings.

16 [Slide]

17 ABC and Trident have shown comparable mechanical
18 performance to ultimate compression strength, axial fatigues
19 and post-fatigue testing, off-axis fatigue testing, fretting
20 testing evaluating the metal/metal interface, axial
21 distraction and the bone/shell interface analysis.

22 [Slide]

23 Additionally, we have provided clinical results
24 demonstrating equivalency of demographics, adverse events
25 and 75-day revision and reoperation data between ABC and

1 Trident.

2 [Slide]

3 So, to echo Mike's comments from before, we
4 believe there are minimal risks with ABC, even fewer risks
5 with Trident. The titanium sleeve protects the ceramic
6 insert. The dual locking mechanism adds versatility and
7 provides intraoperative flexibility to the surgeons, and we
8 have now multiple revision options with this system.

9 DR. MANLEY: Michael Manley.

10 [Slide]

11 Finally, I would like to address panel question
12 number three, which is the issue of postmarket surveillance
13 with the systems and, again, I would like to ask Beth Staub
14 to address the issue.

15 [Slide]

16 MS. STAUB: Howmedica Osteonics proposes to
17 continue to follow the subjects who have been involved in
18 this study annually until the patients have obtained two-
19 year follow-up. That would include following the 515 ABC
20 and 114 continued access patients until the last patients
21 from the original cohort reaches two years, giving us four-
22 year follow-up on the early patients. We also intend to
23 follow the 213 Trident cases until the last patient reaches
24 two years, providing three-year follow-up on the early
25 cases.

1 [Slide]

2 DR. MANLEY: So, in final summary, this radiograph
3 shows a single patient who has a Trident bearing on the left
4 and an ABC bearing on the right. We believe that data has
5 shown that ABC bearings are safe and effective as compared
6 to the control bearings, and we believe that the mechanical
7 testing and the early intraoperative data with Trident also
8 shows that these bearings will be safe and effective.

9 Thank you. That concludes our presentation. I
10 would like to turn this back to Dr. Boyan.

11 DR. BOYAN: Thank you. I am going to ask the FDA
12 to make their presentation now, and ask Peter Allen, the
13 lead reviewer to come up and give his analysis of the
14 preclinical and clinical application.

15 **FDA Presentation**

16 **Preclinical and Clinical Information**

17 MR. ALLEN: Good afternoon.

18 [Slide]

19 My name is Peter Allen, and I am a biomedical
20 engineer in the Orthopedic Devices Branch of the Office of
21 Device Evaluation at FDA.

22 DR. BOYAN: Mr. Allen, before you start. Sponsor,
23 actually it is time for you to go back and be in the
24 audience. Thanks. Okay.

25 MR. ALLEN: I am also the lead reviewer for this

1 PMA.

2 I would like to thank Howmedica Osteonics for
3 their presentation this afternoon, and the panel for your
4 attendance here today.

5 We are here to discuss the premarket approval
6 application for the Osteonics ABC and Trident ceramic-on-
7 ceramic hip systems. I will provide a brief review of the
8 preclinical and clinical information in the PMA and Dr.
9 Harry Bushar, of the Division of Biostatistics, will provide
10 a review of the statistical data.

11 [Slide]

12 These hip systems are intended for use in patients
13 requiring primary total hip replacement who are diagnosed
14 with non-inflammatory degenerative joint disease, which is
15 defined by the indications listed here.

16 [Slide]

17 The ABC System is available in two versions,
18 referred to here as System I and System II. Both versions
19 feature a ceramic-on-ceramic bearing couple. The bearing
20 couple consists of an alumina ceramic femoral head and an
21 alumina ceramic acetabular insert. It is this ceramic
22 bearing couple that makes these systems investigational.

23 Both systems use commercially available Omnifit
24 hydroxyapatite-coated hip stems, and all components of both
25 systems are intended to be implanted without cement.

1 The primary design difference between these two
2 systems involves the exterior coating on the acetabular
3 shells. System I features a titanium shell with an Arc-
4 Deposited titanium coating beneath a plasma-sprayed
5 hydroxyapatite coating. System II features a titanium shell
6 with an Arc-Deposited titanium coating beneath a plasma-
7 sprayed hydroxyapatite coating.

8 [Slide]

9 Like the ABC Systems, the Trident System also
10 features a ceramic-on-ceramic bearing couple. It uses the
11 same ceramic femoral head and Omnifit hip stem as the ABC
12 Systems. All components of the Trident System are also
13 intended to be implanted without cement.

14 The Trident represents the latest design iteration
15 of the ABC Systems. The primary design difference between
16 the Trident and ABC Systems involves the acetabular
17 insert/shell interface. The Trident incorporates
18 modifications to the insert locking mechanism that helps to
19 eliminate intraoperative chipping of the ceramic insert, and
20 improve the use and revisability of the device.

21 The Trident alumina ceramic insert is pre-
22 assembled to a titanium alloy sleeve at the factory. This
23 insert and sleeve assembly mates with the Trident acetabular
24 shell via a taper lock fit. This metal-to-metal
25 interference fit eliminates the potential for chipping of

1 the ceramic insert that can occur with the ceramic-to-metal
2 interference fit of the ABC Systems.

3 The Trident acetabular shell is manufactured from
4 titanium alloy and has an Arc-Deposited titanium coating
5 beneath a plasma-sprayed hydroxyapatite coating, similar to
6 the coating on the ABC System II shell.

7 [Slide]

8 The sponsor performed these preclinical mechanical
9 tests in support of the ABC and Trident Systems. A detailed
10 description of these tests and results was provided in the
11 PMA. In addition, wear test data on the ceramic bearings
12 was provided in a master file from the ceramic bearing
13 supplier, CeramTec of Germany. FDA believes that the
14 preclinical testing is adequate and has no further issues
15 with it.

16 [Slide]

17 The next four slides depict the criteria under
18 which the clinical data was collected and analyzed for the
19 purposes of supporting this PMA.

20 [Slide]

21 The following primary safety and effectiveness
22 data were to be collected at the designated follow-up
23 evaluations until all patients reached the two-year study
24 endpoint. Efficacy was to be based on Harris Hip Score,
25 which includes pain and function components, and on

1 radiographic assessment. Safety was to be based on
2 component revision events and overall adverse events.

3 [Slide]

4 Patients were considered a failure if they met any
5 one of the following criteria at the two-year study
6 endpoint. That is, a total Harris Hip Score of less than
7 70; a radiographic failure; or a revision of any of the
8 device components.

9 [Slide]

10 A radiographic failure was defined as meeting any
11 one of the criteria defined here.

12 [Slide]

13 Study success was defined as not detecting, as
14 statistically significant, an increase of greater than or
15 equal to 7.5 percentage points in the 2-year patient failure
16 rates for Systems I or II over the 2-year failure rate for
17 the control, and complication rates that are statistically
18 no worse than the control.

19 [Slide]

20 I will now focus my review on the clinical data
21 collected and analyzed for the ABC Systems. This will
22 include the updated information provided in Amendment 4 to
23 the PMA. I will then discuss the Trident data a little bit
24 later.

25 [Slide]

1 The ABC System study was a randomized,
2 prospective, controlled, multi-center trial with 515 cases
3 enrolled at 16 sites. The PMA contains data on 413 cases
4 who have either reached the 2-year study endpoint or were
5 revised prior to their 2-year follow-up.

6 This PMA also contains data on 102 cases who have
7 not yet reached the 2-year study endpoint. As a result,
8 only their 1-year safety data has been examined by the FDA
9 for the purposes of evaluating study success.

10 [Slide]

11 The 413 cases serve as a primary analysis group
12 that supports the safety and efficacy analysis for this PMA.
13 Of the 413 cases in this group, 140 cases were implanted
14 with System I, 140 cases with System II, and 133 cases with
15 the control system.

16 The control system is a standard metal-on-
17 polyethylene hip that consists of the Howmedica Osteonics
18 components listed here. All four components are
19 commercially available for use in the U.S. The acetabular
20 shell of the control system has the same titanium porous
21 coating as the ABC System I shell.

22 [Slide]

23 With respect to the results obtained for the
24 primary efficacy measures, here we have the Harris Hip Score
25 results taken from the primary analysis group. A score

1 greater than 90 is considered excellent, 80-90 good, 70-80
2 fair, and below 70 poor.

3 As you can see from the first line, the mean
4 preoperative scores were virtually identical for all three
5 groups. At two years there is no significant difference in
6 mean scores as all three groups are in the excellent range.

7 You will also note that there is no significant
8 difference in the percentage of cases with scores less than
9 70. Remember that a score less than 70 is one of the
10 patient failure criteria. Two cases from System I, two
11 cases from System II, and two cases from the control system
12 had scores below 70 at the two-year follow-up.

13 [Slide]

14 With respect to the previously defined
15 radiographic failure criteria, one case from System II was
16 defined as a radiographic failure due to progressive
17 subsidence of the femoral component. No radiographic
18 failures were detected for System I or the control system.

19 [Slide]

20 With respect to primary safety measures, here we
21 have a summary of the revision and adverse event rates.
22 This table includes all 515 patients enrolled in the
23 clinical study including those 102 cases with less than 2-
24 year data.

25 With regards to revision rates, the control system

1 demonstrated a slightly higher revision rate than the ABC
2 Systems. Specifically, System I had one revision, System II
3 had two revisions and the control system had five revisions.

4 Operative site and systemic complications appear
5 comparable for all three groups. However, within the
6 operative site interoperative complications we do find that
7 among the ABC devices there were a few occurrences of
8 intraoperative chipping of the ceramic inserts which
9 contributed to their slightly higher intraoperative event
10 rate.

11 [Slide]

12 Of the 172 cases implanted with System I, there
13 were 5 reports of chipping of the ceramic insert during
14 insertion of the device. That is, 2.9 percent of the cases
15 experienced this event. There were 4 chipping events
16 reported for the System II components, for an occurrence
17 rate of 2.3 percent.

18 This chipping complication is unique to ceramic
19 inserts due to the brittle nature of ceramic materials. The
20 chipping is a potential-concern because chipping of the
21 ceramic insert, if undetected, could lead to catastrophic
22 fracture of the insert postoperatively.

23 It should be noted that the chipped inserts
24 reported on here were replaced intraoperatively with no
25 further complications, and these patients were all doing

1 fine at their last evaluation.

2 [Slide]

3 Here we have the overall failure rates for each
4 system based on the number of patients who have met at least
5 one of the three patient failure criteria.

6 If we look at the overall patient failure rates
7 for the three systems we see that the results are not
8 significantly different based on the defined study success
9 criteria. ABC System I had a 2.1 percent failure rate,
10 System II a 3.6 percent failure rate, and the control system
11 a 6 percent failure rate at the 2-year endpoint. Based on
12 the failure rates and adverse event rates, both ABC Systems
13 meet the defined study success criteria.

14 [Slide]

15 In addition to the primary and safety efficacy
16 measures, the sponsor provided secondary data based on two
17 patient satisfaction assessment tools.

18 The Hip Society Patient Satisfaction Assessment is
19 directly related to the total hip process. The questions
20 are taken from the Hip Society Clinical Evaluation. The
21 percentages recorded here are the percentage of patients
22 responding "yes" to these questions. As you can see, the
23 results are comparable for all three systems.

24 The Health Status Questionnaire, (HSQ) -12 is a
25 measurement of patient's general health which includes

1 physical and mental health components. It is based on a
2 100-point scale, with a higher score representing an
3 improvement in health. The preoperative and 24-month mean
4 scores are provided here and, again, you see that the
5 results are comparable for the three systems both
6 preoperatively and postoperatively.

7 [Slide]

8 After enrollment of the original 515 cases was
9 completed, FDA approved a continued access study for an
10 additional 336 ABC System cases. These cases were followed
11 to provide additional safety information.

12 To date, 116 cases have been implanted. Data from
13 114 is included in the PMA. All 114 cases are out past 7
14 weeks, and 86 cases are out to their 1-year postoperative
15 time point. The vast majority of these cases, 113, were
16 implanted with System II. Only 3 cases were implanted with
17 System I, as it appears the study surgeons have a strong
18 preference for the hydroxyapatite-coated acetabular shells
19 of System II.

20 Since most of the cases received System II, the
21 adverse event safety information was pooled together. The
22 adverse event rates were unremarkable in that they were
23 comparable to the rates previously discussed for both
24 Systems I and II and the control system. There were no
25 reported revisions. However, 6 chipping incidents were

1 reported for the 116 cases.

2 Combining the continued access cases with the
3 original group of cases receiving Systems I and II, the
4 overall chipping rate for the ABC System inserts is 3.4
5 percent.

6 [Slide]

7 Now I would like to comment on the clinical study
8 for the Trident system. The Trident System was added as an
9 additional study arm to the ABC System IDE last year.
10 Remember that the Trident is an updated design to the ABC
11 System, with the main modification involving the locking
12 mechanism between insert and shell.

13 [Slide]

14 The Trident arm is a non-randomized, prospective,
15 controlled, multi-center trial for 213 cases, conducted at 6
16 of the original 16 ABC System study sites.

17 Patients implanted with the Trident were compared
18 to the control system data collected in the original ABC
19 System IDE. Trident study patients were evaluated using the
20 same clinical protocol as the ABC System patients.

21 To date, 159 cases have been enrolled in this
22 study arm, 157 of which are included in the PMA data
23 analysis, and 135 cases have reached the 7-week
24 postoperative evaluation time point and 27 cases have
25 reached their 6-month evaluation time point. No cases have

1 yet reached the 1-year evaluation point.

2 [Slide]

3 As a result, clinical data on the Trident is
4 preliminary at this tie as only a handful of cases have
5 completed their 6-month evaluation.

6 In addition, no radiographic data has been
7 provided for these cases due to the short postoperative
8 follow-up times. If we look at the available 6-month
9 results for both Trident and the control, we see that the
10 average Harris Hip Scores are comparable for both systems.

11 The adverse event rate for the Trident is lower
12 than the adverse event rate for the control at 6 months. Of
13 particular note, there were no revisions reported for the
14 Trident, and no occurrences of chipping of the Trident
15 ceramic insert. In addition, the mean HSQ-12 score at 6
16 months is slightly higher for the Trident.

17 It is the sponsor's contention that, in addition
18 to the preclinical mechanical testing and short-term safety
19 data from the Trident System, the clinical data for the ABC
20 Systems may be used to support the safety and efficacy of
21 the Trident System. This is based on the use of the
22 identical ceramic bearing surfaces, exterior shell
23 geometries, and femoral components in both the ABC and
24 Trident Systems.

25 [Slide]

1 I will now turn the floor over to Dr. Harry Bushar
2 to discuss the statistical analysis. After Dr. Bushar's
3 presentation I will provide a slide which summarizes the
4 issues we would like you to think about during your
5 discussion. I will then provide the specific panel question
6 slides once we get into that part of your discussion.
7 Harry?

8 **Statistical Analysis**

9 DR. BUSHAR: Thank you, Peter.

10 [Slide]

11 My name is Harry Bushar. I am the statistical
12 reviewer for the Howmedica Osteonics ABC/Trident Systems
13 PMA.

14 [Slide]

15 I am going to discuss the ABC System clinical
16 trial which, of course, you have already heard about from
17 Peter and the sponsor. I am going to focus on a few
18 statistical points. The original study was prospective,
19 controlled by Osteonics ABC System III, which is a standard
20 hip system, and they do-have 2-year follow-up on that. It
21 was randomized between this control and 2 concurrent study
22 arms, each getting about the same number of patients. These
23 were Osteonics ABC Systems I and II which are, of course,
24 experimental, and they each have 2-year follow-up. This was
25 a multi-center study with 16 investigational sites.

1 [Slide]

2 The Trident System clinical trial was also
3 prospective but, since it started late, it was historically
4 controlled. Instead of doing a randomization, they decided
5 to simply borrow the control from the previously completed
6 Osteonics ABC System III. So the control was used for three
7 different purposes, to compare to System I, System II and
8 also Trident. This late study arm only has operative
9 follow-up to speak of. In other words, they do not have any
10 2-year follow-up. This multi-center study used 6 of the 16
11 previous investigational sites to keep the results
12 comparable.

13 [Slide]

14 In terms of follow-up, what we now have for System
15 I is 172 hips, referred to as cases in most of the other
16 presentation, with 140 hips out to 2 years. System II is
17 similar, 177 hips with operative follow-up and 140 out to 2
18 years. With Trident there are 157 hips with operative
19 follow-up and none out to 2 years. System III or control is
20 165 hips and 133 of these are out to 2 years. -

21 [Slide]

22 I am going to look at safety first and I am going
23 to focus in on the operative site intraoperative adverse
24 event rates. The reason for doing this is that this is the
25 category in which chipping occurred. Of course, you can see

1 the effect of that, even though I haven't broken it out as
2 such. The event rates for System I and II are comparable,
3 11.6 in once case and 9.6 in the other. If you combine
4 these, and there is no reason not to, there are similar
5 results, and the demographics are very similar for System I
6 and II. You get 10.6 percent. Then, if you look at
7 Trident, you see that this is quite low, 2.5 percent and the
8 control rate is 6.7 percent, somewhere in between the two.

9 [Slide]

10 What I have done here, I have looked at the
11 difference of binomial proportions. So, I am comparing the
12 results that were shown on the previous slide. What I am
13 constructing are 90 percent confidence intervals. The
14 reason for using 90 percent is I am going to focus in on one
15 end of the other to try to make a statement as to what one
16 can say when one compares Systems I and II combined to
17 control.

18 What you could say from that interval is that
19 Systems I and II combined are no worse than 9.1 percentage
20 points higher than the control, and you could make that
21 statement at the 95 percent confidence interval because you
22 are just borrowing the upper end of the interval.

23 One can do a similar comparison of Trident to
24 control, and there you can see that the worst that can occur
25 is a 1.1 percentage point increase of Trident over control

1 at the 95 percent confidence interval.

2 The best way to look at the bottom interval is to
3 look at the bottom part of the bottom interval, namely the
4 3.2 percent, and that shows that System I and II combined is
5 at least 3.2 percentage points worse than the Trident, and
6 that can be said with 95 percent confidence. So, it does
7 appear that the Trident in the operative site intraoperative
8 adverse event rates is doing something that has certainly
9 reduced the chipping to zero.

10 [Slide]

11 As far as effectiveness goes, I am going to focus
12 in on the failure, and failure was defined as revision,
13 total Harris Hip Score less than 70, or radiographic
14 failure. There was only one radiographic failure, and that
15 was due to a progressive femoral component, subsidence which
16 was greater than or equal to 5 mm.

17 [Slide]

18 The effectiveness failure rates are, for System I
19 2.1 percent, for System II 3.6 percent. Again, these are
20 close and if you combine them you get 2.9 percent. The
21 control rate was 6 percent.

22 [Slide]

23 I have done a very similar thing here. I have
24 looked at the difference of binomial proportions,
25 constructing 90 percent confidence intervals, and here one

1 would want to focus on the upper end of the interval. You
2 can see that at worse the failure rate for System I is 1.9
3 percentage points greater than the control, with 95 percent
4 confidence. System II is no worse than 3.6 percentage
5 points greater than control, with 95 percent confidence.
6 Then, if you combine the two and compare it to control you
7 can see that the upper limit now drops to 1.3 percent so you
8 are no worse than 1.3 percentage points greater for System I
9 and II combined compared to the control, with 95 percent
10 confidence.

11 That is it. I have finished. Thank you very
12 much.

13 MR. ALLEN: Finally, we have our last slide.

14 [Slide]

15 This is a list of our discussion topics. I have
16 individual slides for each of these that go into a little
17 bit more detail, and we can go through those one by one when
18 you are ready to address them. You also have a draft hard
19 copy of these questions that was provided to you earlier, in
20 Tab C of your white binder.

21 Anyway, here is a condensed version of the issues
22 for which we are seeking your input: Number one is with
23 regard to the chipping events reported with the ABC ceramic
24 inserts; number two the short-term clinical data on the
25 Trident System, and number three, the possible need for

1 longer-term clinical information for either or both device
2 systems.

3 Thank you very much, and this concludes FDA's
4 presentation.

5 DR. BOYAN: Thank you very much. I think now what
6 we will do is have the lead reviewers from the panel give
7 their assessments and then we will begin our official
8 discussion. So, Dr. Li, would you please present your
9 review of the preclinical data?

10 **Panel Reviews**

11 DR. LI: Sure.

12 [Slide]

13 These are just some of my comments that I have on
14 the application. First, to kind of explain where I will end
15 up, I thought I would tell you about how I got there.

16 One is to kind of look at the ceramic-on-ceramic
17 historical overview prior to this particular device which
18 Osteonics supplied in their application. The previous
19 problems with ceramic-on-ceramic devices were probably in
20 three categories: frank-fracture of the ceramic itself; a
21 loosening usually of the acetabulum but often the stem as
22 well; and impingement which led to other problems directly
23 with the ceramic.

24 [Slide]

25 As far as fracture goes, after reading the

1 application and all the data, it appears that the problem is
2 solved, and it is solved in a way that is with a scientific
3 basis, basically in control of the grain boundaries by
4 CeramTec of the ceramic, and followed by kind of an every
5 product testing protocol where every particular product is
6 actually load tested prior to sale.

7 This appears to have lowered the incidence of
8 fracture of the ceramic to something probably less than the
9 fracture of the femoral stems. So, I think the issue of
10 fracture of the ceramic appears to be behind us as far as
11 all the laboratory data goes.

12 [Slide]

13 However, interestingly enough, although the
14 driving force for ceramic-on-ceramic is the reduction of
15 wear, especially the avoidance or prevention of osteolysis,
16 and although all the lab wear tests for ceramic-on-ceramic
17 have always been essentially zero, there is osteolysis
18 reported in the literature. The first case that I could
19 find is back in 1991, which is a single case report in Acta
20 Scandinavica. More recently, in '94, Shih reported 8/134
21 ceramic-on-ceramic devices that had osteolysis. Then in
22 JBJS, January of '98, Yoon, in Korea, reported, amazingly,
23 66/103 ceramic-on-ceramic had osteolysis.

24 Now, it should be pointed out that these devices
25 were a completely different design. The ceramics were

1 previous generations of ceramics. But I think it is
2 important to note that just being a ceramic-on-ceramic
3 device does not guarantee solutions of low wear and
4 osteolysis.

5 [Slide]

6 The issue of loosening probably is the one that
7 has the question mark in my head. The loosening in ceramic-
8 on-ceramic devices, by and large, is aseptic loosening in
9 the absence of osteolysis. Dr. Laurence Sedel, in France,
10 probably has the longest history of this. In one particular
11 case, out of 401 ceramic-on-ceramic devices, 44 of them
12 became loose at the 15-year time period. In none of these
13 cases was there osteolysis.

14 Again, I wish to point out that these devices were
15 of a different design and different materials and after a
16 15-year follow-up, and really not a direct reflection,
17 again, of the current device that we are talking about, but,
18 again, wear has never really been the problem with ceramic-
19 on-ceramic, with a few exceptions, it has always been
20 loosening.

21 [Slide]

22 Kind of a fallout out of this is why do these
23 patients loosen with these devices? Kind of a general
24 consensus is it may have something to do with the design of
25 the socket, many of which were threaded ceramic sockets.

1 Without a lot of analysis, we have always kind of pointed
2 the finger at that, saying that must be the problem although
3 there really hasn't been much follow-up data to actually
4 demonstrate that that is it.

5 An interesting feature out of most of the ceramic-
6 on-ceramic data is that if you break it out by age, the
7 younger patients always do better in ceramic-on-ceramic
8 devices. Of the 401 cases that we talked about of Dr.
9 Sedel's, the green arrow indicates the survivorship at 15
10 years of those patients who were less than 50, and the white
11 arrow indicates the survivorship of those patients that were
12 older than 50. So, it is kind of the opposite of metal-on-
13 polyethylene. In Metal-on-polyethylene, the younger you
14 are, the worse you are. In ceramic-on-ceramic, it appears
15 the younger you are, the better off you are.

16 [Slide]

17 So the question is although the fracture problem
18 may be solved and osteolysis appears to be less frequent but
19 certainly not necessarily at zero, the question is, is
20 loosening solved? I think perhaps that is the feature that
21 probably is going to make or break this device, and it is
22 hopeful that the use of the cementless metal back liners
23 will have addressed the loosening issue, although that is
24 yet to be demonstrated in a long-term series.

25 [Slide]

1 This has been presented a couple of times but I
2 want to talk a little bit about can you compare the
3 different groups. Just as a quick reminder, perhaps for the
4 fourth time that you have seen these slides, but I want to
5 throw these up just as a kind of a reminder of how they are
6 different either by material or by design. The ABC I and II
7 essentially differ by whether or not there is an HA coating
8 on the outside of the metal shell. The Trident is a little
9 different. It has an HA coating but the key issue on the
10 Trident, in my mind, is that the ceramic liner has
11 essentially a shrink-fit titanium alloy sleeve that goes
12 around the outside of it so you end up with a metal-on-metal
13 junction in the shell rather than metal-on-polyethylene.

14 [Slide]

15 The other issue is that the Trident comes in one
16 additional size, the 36 mm. The ABC I and II do not. I
17 think the issue on the expanded sleeve is that, one, I think
18 it was put in there to be a little more forgiving, if you
19 will interface and if the installation was supposed to be
20 easier, although I didn't really see anything that
21 demonstrated that. Perhaps Dr. D'Antonio can actually
22 comment on whether that is true and how they actually
23 documented that it is easier.

24 I think the materials issue on the sleeve is that
25 the way the sleeve is put on there, it is inductively heated

1 to kind of make the metal soft and expanded, and then they
2 cool it and on cooling, it essentially clamps itself around
3 the alumina ceramic, as I understand. Is that right? So,
4 the titanium alloy sleeve is essentially in tension around
5 the ceramic cup.

6 The only issue there that I could think of that
7 might be a long-term issue is corrosion. Titanium is known
8 to corrode, you know, at implant time of five or seven
9 years, for crevice corrosion if you have like a mixed metal
10 head on it. But, certainly, in this case where you have a
11 piece of titanium that is under tension, the tendency to
12 corrode is essentially higher depending on how much tension
13 is applied. But, I didn't see any number in there for how
14 much residual tension there is in the sleeve. It is not
15 unrealistic if you have really high tension to increase the
16 corrosion rate by a factor of five or ten quite easily.
17 And, corrosion is not going to be seen in a short time
18 period. When it occurs on stems, it usually comes at least
19 at the five to seven year mark.

20 [Slide]

21 So, the Trident System is supposed to facilitate
22 alignment, although I am not clear how the facilitation was
23 documented, and it is supposed to be more forgiving in
24 putting it in. Dr. D'Antonio called it canting. When I
25 have handled these devices, it has always been kind of an

1 issue with me. It is possible, as Dr. D'Antonio showed, to
2 put the liner in kind of off alignment. So, it is kind of
3 stuck in there but it is actually not seated correctly and I
4 have never been able to actually tap it out like he did in
5 the movie. Sometimes the thing is wickedly in there and you
6 have to work pretty hard to get it out. So, I am not sure
7 how a general surgeon, or one that does these not so often,
8 can actually guarantee that the liner is actually, in fact,
9 aligned each and every time. And, if the Trident System
10 actually helps to do that, it would have been helpful to
11 have documented the benefit of that.

12 [Slide]

13 For the mechanical testing, especially in the
14 absence of a guidance document, you have to compliment the
15 applicants for doing extensive and actually very well done
16 tests, and I really have no issues with those tests. I
17 think the tests are appropriate and you passed them all
18 well.

19 [Slide]

20 As a reviewer-you can't say that and stop. I
21 think the only thing that I am a little surprised about,
22 quite frankly, is that there is no independent wear testing
23 done. Now, I did not have access to the actual CeramTec
24 document that the applicants referred to but I guess my
25 question, in the absence of seeing that data, is I am not

1 sure that comparison is appropriate. In other words, were
2 the tolerances and the CeramTec thing the exact same
3 tolerances that were used here, and what were the loading
4 conditions and variations?

5 I am a little concerned also about testing of the
6 devices -- if I had to comment on the testing, it is that
7 they are done essentially under a single condition and not
8 really addressing the wide variety of conditions that might
9 be encountered surgically. For instance, on the metal and
10 polyethylene case, things like abduction angle and
11 anteversion, at least in extreme cases have been known to
12 affect the wear and the question is did those same rules of
13 thumb that govern metal-polyethylene surgical procedures,
14 did those exact same rules hold for ceramic-on-ceramic?
15 Maybe they do; maybe they don't. I just don't know.

16 [Slide]

17 So, in my opinion the load wear must be directly
18 evaluated. In general, I think wear is a relatively poorly
19 understood phenomenon. It is unclear to me how over
20 probably the 25 or 30-year history of ceramic-on-ceramic in
21 some people's hands the earlier ceramic-on-ceramic devices
22 had wear; in other cases they don't. So, clearly, our
23 understanding of wear in general is not all that well
24 understood, and although I don't expect the results to be
25 anything but sterling out of the laboratory given the