

1 DR. McNALLY: But then we looked at lower
2 levels and we thought, at lower levels, this is
3 still an okay clinical level for it.

4 DR. PULIDO: I could see why you chose,
5 what was the reasoning behind the 5 percent. But
6 now have you have a 3 percent incidence with the
7 Acuvue lens and you have a 6 percent incidence with
8 this SEE3. So there is a doubling there. What
9 would happen to that null hypothesis if you used
10 greater than or equal to 4 percent?

11 DR. McNALLY: I have to say I am not sure.

12 DR. PULIDO: I would like to know that,
13 though. Did the statistician do that?

14 DR. SUGAR: Dr. Cutter, why don't you come
15 to the podium. You might be more comfortable.

16 DR. CUTTER: I doubt it.

17 DR. SUGAR: You notice that we are
18 interested in your comfort.

19 DR. CUTTER: I think, in a way, one could
20 do that calculation and the numbers would come out.
21 You take a percent difference and divide it out. I
22 haven't done it. I suspect I know where it is
23 going to come out. You specify the hypothesis for
24 decision-making in advance. You end up making a
25 decision on the basis of the a priori evidence that

1 you set up for this -- after you have observed the
2 results to kind of go in. I think you are in an
3 estimation procedure.

4 The absolute difference is less than 3
5 percent. So I think you have to look at what was
6 observed and maybe put some confidence intervals on
7 the difference and look at the difference for the
8 magnitude and the size of the difference rather
9 than really going back to a hypothesis-testing
10 mode.

11 I am not trying to be evasive, but I think
12 it has to do with conceptualization. We plan
13 trials with the best information that is available
14 years in advance of when we actually do the
15 analysis. We set up, and we sort of live or die by
16 that proposal.

17 When you have the data in hand and you can
18 see whether or not your assumptions were correct or
19 whatever, I think it is appropriate to look at the
20 size and the magnitude of the absolute difference
21 or, if you want, proportional difference, whatever
22 you are looking at, and discuss it in those
23 contexts in terms of what size confidence interval
24 you have.

25 DR. SUGAR: Dr. Bandeen-Roche?

1 DR. BANDEEN-ROCHE: Dr. Bandeen Roche.
2 You did mention a confidence interval on the
3 difference in rates. I did a truly back-of-the-
4 envelope calculation. That calculation showed that
5 a 95 percent confidence interval did exclude 0 so
6 that the rates were significantly different. Does
7 this agree with your calculation

8 DR. CUTTER: Yes.

9 DR. BANDEEN-ROCHE: Because mine really
10 was back-of-the-envelope.

11 DR. CUTTER: Yes.

12 DR. PULIDO: So, again, there is a
13 difference between the rate with the SEE3 and the
14 Acuvue, because, from my quickie thing, too, it
15 looked to me like the confidence intervals had a
16 difference in overlap

17 DR. CUTTER: Again, not to split hairs,
18 the only other thing I would do is that would be --
19 you would adjust the p-value for the multiple-
20 hypothesis tests you are doing because the primary
21 hypothesis was a noninferiority test.

22 DR. McNALLY: If I could add a few
23 comments. When we looked at the unadjusted rates,
24 and, by that, is the number of patients with an
25 endpoint infiltrate divided by the number of

1 patients dispensed, and when we do the sort of
2 classical comparison of those, we find no
3 statistical difference between the unadjusted
4 rates.

5 Then, when you perform the life-table
6 analysis to take account of all the patients who
7 discontinued and the other things that happened.
8 Then you come up with these other rates and you get
9 into the discussion you were just getting into.

10 But, remember, we excluded from that the
11 two peripheral ulcers, so those were two of the
12 more serious of these endpoints, from the control
13 group. This was a conservative thing because one
14 statistically throws us out -- you know, it makes
15 the rate go to 5.7 percent if we included that last
16 endpoint for Acuvue.

17 We thought, you know, this really isn't
18 representative for Acuvue to go from a 3.3 to a 5.7
19 because of a statistical foible, I will say, for
20 the life-table analysis. Secondly, we included the
21 second ulcer -- we didn't include the second ulcer
22 at six months. I think if you include those in the
23 analysis, I think that the conclusions then change
24 and you find that there is overlap with 0.

25 DR. CUTTER: That is correct. And the

1 foible that Dr. McNally was talking about is since
2 the number of patients who are observed beyond one
3 year, the actual date of their exam starts
4 diminishing, the life-table rate is based upon the
5 event rate in the interval where the event occurs,
6 the number of patients that are still around.

7 This really makes the rate probably not
8 representative of the control group so we chose to
9 use a conservative -- leaving out those other
10 events. So one could split hairs about whether or
11 not it gets significant or not.

12 If you include the other events that are
13 known and occurred, but they occurred outside the
14 interval and you are using a life-table estimate
15 for it -- but we have done the analysis looking at
16 if that had occurred at day 365 and what did that
17 do.

18 DR. MATOBA: This is pertaining to --

19 DR. PULIDO: It is pertaining -- so,
20 again, you are speculating. There was some
21 speculation that is going on. You had mentioned
22 before, Dr. Cutter, that you could speculate what
23 would happen if the null hypothesis had been
24 changed to greater than or equal to 4 percent.
25 What is your speculation on that?

1 DR. CUTTER: I haven't done it with all
2 the events included. Again, I think, obviously,
3 someone would want me to include all the events. I
4 don't know who that might be. But I haven't done
5 the calculation where you take all the events and
6 then look at it relative to a 4 percent difference.

7 DR. PULIDO: You had speculated before?
8 Give me a speculation, p less than 0.05 or not

9 DR. CUTTER: The absolute difference is
10 slightly over 2 percent. You have got your
11 envelope. The standard error doubles. I think it
12 wouldn't be significant, actually, with 4 percent
13 but, again, if you are including at least one of
14 the two events that are left out.

15 If you include both events, I am almost
16 certain that they are not.

17 DR. PULIDO: If you are going to exclude
18 certain events, you also excluded a severe red eye
19 as a problem with the SEE3 lens. Nothing like that
20 ever happened with the Acuvue lens. So what was
21 this severe red eye?

22 DR. McNALLY: I think I can address your
23 question, contact-lens acute red eye, perhaps, as
24 we explained contact-lens acute red eye. In the
25 contact-lens industry, we tend to put things into

1 little definitions and boxes which maybe not
2 everybody agrees with. Contact-lens acute red eye
3 is, as defined in our protocol, an acute event
4 involving infiltrates overnight when this occurs
5 and they have pain in the morning and redness and
6 so forth.

7 The critical part of that definition is
8 that there are infiltrates. However, what happened
9 in the study is that if anybody's eye became
10 injected in an acute way, the investigators often
11 marked acute contact-lens red eye. However, there
12 were no infiltrates and so it really didn't fit
13 that definition.

14 So we removed that definition but we
15 included them, then. If there were infiltrates, we
16 included them in the endpoint if it met the
17 endpoint criteria. We included them then in the
18 adverse events under infiltrative keratitis if,
19 indeed, there were infiltrates. If all there was
20 was injection overnight, then they were included in
21 the appropriate place which would be, if it was
22 grade 4, it would be biomicroscopy greater than
23 grade 4 in the table.

24 So they were not excluded. They were in
25 the table under a more descriptive definition of

1 the findings and symptoms that occurred.

2 DR. MATOBA: Alice Matoba. The patients
3 who were discontinued for biomicroscopic findings,
4 they were not included in your final analysis; is
5 that correct?

6 DR. McNALLY: No; all patients were
7 included in the final analysis. The life-table
8 particularly takes into account patients who are
9 discontinued.

10 DR. MATOBA: So the 3.1 percent versus 5
11 percent, that first incidence of endpoint
12 infiltrates, that includes those five patients in
13 the SEE3 and the one patient in the Acuvue group
14 who were discontinued and had endpoint infiltrate
15 as a biomicroscopic finding?

16 DR. McNALLY: It does. That doesn't
17 include the one corneal ulcer that was seen by the
18 ophthalmologist because we didn't have infiltrate
19 data provided, although we had diagnosis of corneal
20 ulcer as well as a scar later. But that wasn't
21 included in that unadjusted 3.1 percent rate.

22 DR. SUGAR: I have a couple of questions.
23 One, what was your instruction to the investigator
24 concerning infiltrates; that is, were they told to
25 remove the lens, treat them with a specific

1 medication or do whatever was standard therapy and
2 under what circumstances, if any, were they
3 instructed to culture?

4 DR. McNALLY: Dr. McNally. I haven't
5 identified myself. It was standard of care was
6 what we were using and so if culturing was felt to
7 be needed by the investigator, then filtering was
8 done.

9 DR. SUGAR: Do you know what the specific
10 treatment was in terms of medications for the
11 infiltrative keratitises?

12 DR. McNALLY: We have included a table
13 showing the different treatments. You can refer to
14 the table for the different endpoint infiltrates.
15 This was included in the report. Normally,
16 actually, in the SEE3 there were 9 percent that
17 were just treated by removal and with the control
18 lens, it was a little less than that. I think it
19 was maybe 5, while we are looking for the number.

20 The rest were either antibiotic steroid
21 combinations or antibiotics, mostly siloxane and
22 this type of antibiotic was used. But we have
23 listed this in the table.

24 DR. SUGAR: And you also mention a patient
25 who had an adverse response to Tobradex implying

1 that steroid antibiotic combinations were used in
2 some of these "nonmicrobial keratitises." We don't
3 have a good definition of microbial keratitis
4 industrywide, but it says that some of these
5 patients may, indeed, have had infiltrative
6 keratitis if they were not instructed to culture
7 these patients and they were treated with the
8 steroid antibiotic combination.

9 DR. McNALLY: Dr. McNally, again. We
10 tried in the report to list everything we could
11 know in terms of how fast they resolved, was there
12 any outcome that was negative for the wearer. We
13 did list, in table 12 on page 58 of 85, the various
14 pharmaceutical agents or other treatment for each
15 of events, the endpoint-event infiltrates that
16 occurred.

17 DR. SUGAR: Thank you.

18 Dr. Zadnik?

19 DR. ZADNIK: Karla Zadnik. You reported
20 that you did not find an association between number
21 of consecutive nights wear before an event;
22 correct. With that small number of events, what
23 statistical power did you have to find that, if it
24 existed? Do you know?

25 DR. McNALLY: Actually, I don't know if

1 Gary can answer that one, but we didn't set up the
2 study to test that. So I am not really sure what
3 the power was.

4 DR. ZADNIK: My concern is that you are
5 saying there is no association and really we are
6 strongly behind that in recommending that language
7 be removed from the product labeling. What I want
8 to know is what power you got to report that there
9 is no association before something as bold as
10 deleting that from the labeling would happen. So I
11 think that is a fairly important number to find
12 out.

13 DR. SUGAR: Dr. Weissman?

14 DR. WEISSMAN: You reported one case of
15 Thygeson's and one case of herpes keratitis. Were
16 those patients in which -- I think it was mentioned
17 that the Thygeson's was a second episode for that
18 particular patient. Was there any reason why that
19 patient got into the protocol? Wouldn't they have
20 failed protocol by having had no history of
21 previous eye disease?

22 DR. McNALLY: We didn't eliminate previous
23 history of eye disease as one of the -- I don't
24 know the word here, but they could get in the study
25 if there was a previous history as long as the

1 investigator felt that the eyes were quiet and that
2 they were suitable candidates for contact-lens
3 wear.

4 So it was left to the investigator's
5 decision as to whether they felt this was an
6 appropriate candidate. However, the exclusion
7 criteria -- I thought of the word -- if there was
8 any active corneal inflammation or other things
9 like that at the time of enrollment, they were not
10 able to be enrolled in the study.

11 But a previous history did not exclude
12 them from participating.

13 DR. SUGAR: Dr. McMahon?

14 DR. McMAHON: Tim McMahon. Clarify
15 something for me, Dr. McNally. The postapproval
16 study, the primary endpoint item that will be
17 looked at again will be infiltrative keratitis or
18 it will be infectious keratitis?

19 DR. McNALLY: The goal is to determine the
20 rate with the proper sample size of infectious
21 keratitis. Because of this problem in terms of is
22 it or isn't it, anything that starts with an
23 infiltrate, the data will be collected and
24 presented and have an independent review board to
25 then determine, by a definition set up in advance

1 which we haven't done yet because that would be a
2 role of the independent board to say this is the
3 definition, these are the criteria, that will call
4 this a microbial keratitis.

5 Then we will be able to collect the
6 information. We took infiltrate as an endpoint for
7 collecting data because that is pretty clear when
8 there is an infiltrate. It is just not clear in
9 terms of how you would diagnosis or what you would
10 call that entity.

11 So that is the entry criterion to collect
12 the data and then that data can be evaluated by the
13 independent board to determine is it a microbial
14 keratitis or is it not.

15 DR. McMAHON: Then, as follow-up question,
16 as Dr. Holden showed, there seems to be some
17 cumulative risk for microbial keratitis in
18 conventional, nonsilicon hydrogels. Why did you
19 select a follow-up period of only a year?

20 DR. McNALLY: My first answer to that is
21 that this was the recommendation in the discussions
22 with the FDA.

23 DR. ROSENTHAL: Could I just make a
24 comment?

25 DR. SUGAR: Please.

1 DR. ROSENTHAL: Dr. Rosenthal. We do
2 establish, for the new panel members and for the
3 old panel members, these criteria that they have
4 set up for the clinical trial, often many years
5 before the clinical trial comes to you. So, some
6 of these issues like the hypotheses and so forth
7 are based on the best information at the time.

8 One of the biggest problems that the
9 agency gets into is when the company decides to
10 change their hypotheses during the course of the
11 clinical trial. I would add that I think that
12 panel should try to accept the hypothesis since it
13 was accepted by us at the time the clinical trial
14 was designed based on the best information
15 available.

16 DR. SUGAR: Dr. Weiss, Dr. Grimmett. Then
17 I have a question.

18 DR. WEISS: Dr. Jayne Weiss. I wanted
19 some clarification of the rates of the CLPC. I
20 understand from the data that the patients who had
21 the SEE3 lens versus the Acuvue, the SEE3 category
22 had a much higher rate of having preexistent CLPC.
23 But, if we remove those, then you have a SEE3
24 incidence of 3.2 percent for CLPC versus 0.9
25 percent for the Acuvue.

1 I wanted to find out if you separated out
2 that group and looked at the onset of CLPC, in the
3 total group, you indicated the onset of CLPC was 70
4 percent within the first three months in SEE3 group
5 versus after three months in the Acuvue group for
6 75 percent.

7 In other words, SEE3 had a much earlier
8 onset. But if you take out those who had
9 preexistent CLPC, was there still an earlier onset
10 in the SEE3 group which, to me, might imply that
11 the polymer, itself, would give you a better chance
12 of getting CLPC or were the onsets, then, similar.

13 I am referring to page 41 of 58 in table
14 18. I think it is in part 2.

15 DR. McNALLY: This is Dr. McNally. I
16 don't recollect the answer directly to your
17 question but I think, when we looked at our data
18 overall, we did say that a number of these people
19 had previous CLPC. We don't have enough
20 explanation, perhaps, at this time to say whether,
21 if they did not have that, would they have had a
22 lesser rate.

23 So we wanted, in our labeling, to say that
24 there is a potential increased risk, particularly
25 if you have had this in the past. We had several

1 hypotheses in the report, but these are just
2 hypotheses.

3 One hypothesis was potentially mechanical
4 and that was the timeframe and the drawings that
5 some of the investigators gave as well as some
6 experience we have had from our international
7 trials where, instead of being a generalized
8 capillary conjunctivitis, it was localized in a
9 particular place like you might have seen with the
10 stitches that used to be used, if they still are; I
11 don't know.

12 But it looked like it was of mechanical
13 origin. Then when you saw when I showed the one
14 fitting where the edge was lifted a little bit off
15 the edge, that the inferior -- and we are thinking,
16 with no proof to bring to you today, but we are
17 thinking that if the lens can lift at the upper
18 part, you might get some irritation up there as
19 well.

20 But these were hypotheses. We found no
21 correlation with deposits, filming or dirty lenses
22 which is often the other thing blamed for CLPC.

23 The other difference in the factors is 30
24 versus 6, and whether that makes a difference, we
25 are unable to answer at this time. So we hoped to

1 address this in the labeling because we feel this
2 is not a sight-threatening event. It is an
3 irritating event. But we hope to address that in
4 the label.

5 DR. WEISS: I think it would be, because
6 this is a new polymer which is basically why it can
7 be used under the basis it is, I think it would be
8 interesting and probably easy for the company to
9 look at for the 3.2 percent who had no preexistent
10 condition in the SEE3 versus the 0.9 percent in the
11 Acuvue group to see if the onset was at similar
12 times or much earlier in the SEE3 because, if it is
13 much earlier in the SEE3 group and it is
14 statistically significant, then that would imply
15 that, perhaps, this is going to give you a higher
16 chance of having this condition in these patients.
17 That might need to be indicated in the labeling as
18 well.

19 DR. SUGAR: Dr. Grimmatt?

20 DR. GRIMMETT: Michael Grimmatt. I first
21 just wanted to congratulate the sponsor for a very
22 thorough presentation and detailed booklet brought
23 to panel after the study was completed and not in
24 progress. I thought it was a very nice job.

25 I have one observation. In prior reviews

1 that I have participated in, generally speaking,
2 the prescribing range matches the testing range. I
3 just wanted to point out that the testing range for
4 the lens here is +6.00 to -6.00 and the prescribing
5 range is certainly much greater.

6 I wanted clarification on the exact
7 prescribing range sought because, in different
8 places in the notebook, I was seeing different
9 numbers. In the summary of safety and
10 effectiveness the range is listed from -20.00 to
11 +20.00. In the package insert, it is listed from -
12 20.00 to +10.00 and, in the handout of the slide
13 copies we received today, it is listed as -20.00 to
14 +10.00. What is the exact range the sponsor is
15 seeking?

16 DR. ROBIRDS: This is Scott Robirds. The
17 approval range that we are seeking is +20.00 to -
18 20.00. That that will be available for dispensing
19 is the +6.00 to -10.00, initially. The approval
20 range would be +20.00 to -20.00.

21 DR. SUGAR: Can I ask for clarification
22 from the agency? If we approve the lens in a given
23 range, but the guidance is that it can be
24 manufactured and distributed in a broader range; is
25 that correct?

1 DR. ROSENTHAL: I must say I will have to
2 defer to one of my staff.

3 DR. SUGAR: Dr. Lepri?

4 DR. LEPRI: Dr. Lepri. The agency has
5 established a policy over the years, based on their
6 experience and the maturity of the contact-lens
7 technology, that, during the investigation, they do
8 not need to investigate all the available powers.

9 However, at some point during the approval
10 process, they will have to submit to us the effects
11 of varying thicknesses of the contact lenses in the
12 whole range of contact lenses available to evaluate
13 the safety issue of the oxygen permeability. At
14 that point, the agency makes the determination of
15 the final approval range for the lenses.

16 DR. SUGAR: Thank you.

17 Dr. Grimmett.

18 DR. GRIMMETT: I have two more questions,
19 just of clarification mostly. At one point, it was
20 indicated that about 2 percent, 1.9 percent, of the
21 SEE3 eyes lost two lines or greater wearing the
22 contact lenses. I assume they were correctable
23 with overrefraction or spectacle correction; is
24 that correct?

25 DR. McNALLY: Dr. McNally here. That is

1 correct. There was no loss of two lines of acuity
2 in any patient through the study.

3 DR. GRIMMETT: Just one more housekeeping
4 point. Under adverse device effects, it was listed
5 at one point that a patient had optic neuritis. I
6 am assuming there is no implication that the lens
7 had anything to do with that. It was listed as a
8 matter of all complications seen in all these
9 patients; is that correct?

10 DR. McNALLY: Dr. McNally. That was in
11 the control group as well and it was listed for
12 completeness.

13 DR. GRIMMETT: Thank you.

14 DR. SUGAR: Quick question. You included
15 in your study that patients that were pregnant or
16 lactating could be entered. The draft guidance
17 suggests that those patients be excluded. How many
18 such patients were entered and were there any
19 adverse events in those patients?

20 DR. McNALLY: As far as reported to us at
21 entry, there were none. Over the course of a year,
22 there were ten subjects in the control group who
23 became pregnant and two in the SEE3 group. So
24 there was a differential there.

25 DR. SUGAR: So it does have a birth-

1 control effect. I guess it depends on where you
2 put the lens.

3 DR. McNALLY: That's right. We will cover
4 that in the labeling. The two in the SEE3 group
5 completed the trial without problem. Of the ten,
6 eight completed the trial, one with an adverse
7 event which was just a grade 3 staining which
8 resolved just with removing the lens. Two of the
9 eight in the control group were discontinued, one
10 because she was confined to bed rest at some point
11 in the pregnancy and the other just because she
12 felt like she didn't want to wear the lenses
13 anymore.

14 But there were no adverse events of any
15 significance related to that very small group of
16 patients.

17 DR. SUGAR: Dr. Jurkus?

18 DR. JURKUS: I had some questions for
19 clarification regarding acuity and poor vision. In
20 your discontinuation rate, you had indicated and
21 you had showed the slide of the defect, but I was
22 still wondering, do you have any information about
23 the number of people who discontinued because of
24 poor acuity who did not have the lens defect.

25 Second, sort of going along with that, was

1 there any correlation between poor lens fit and
2 poor acuity? Were they combined in the statistics
3 or were they separated out specifically?

4 DR. McNALLY: Dr. McNally, again. I don't
5 remember the exact number of fits, how that
6 distributed. It is in one of the tables, but in
7 terms of the rating by the investigators, the vast
8 majority of the fits were rated as optimal. Then,
9 on SEE3, they tended a bit towards the acceptably
10 loose side. We found no correlation with vision in
11 that group.

12 In those who discontinued for lens fit, I
13 actually didn't look at that data to see if that
14 dropped the vision. We didn't examine every lens
15 that was for a patient who discontinued for lens
16 acuity, but I tried to look at the data in terms of
17 when did it occur. It all did occur in the very
18 beginning, so if you had the first of the lenses
19 you were wearing for 30 days or your second, this
20 was something where you would say, "Well, this is
21 unusual," I wouldn't want to continue.

22 But if you have cycled through a few
23 lenses and you get one that you can't see with, you
24 say, "Well, let me get another one." So it was all
25 in the first three months that we had these acuity

1 discontinuations.

2 Again, there was no loss of best corrected
3 acuity. Most of the ratings, as I showed, I think
4 98 percent were the same at baseline and 83 percent
5 were 20/20. So there may have been an occasional
6 patient where they didn't get good acuity, but I
7 think it would be no different than any soft
8 contact lens.

9 DR. SUGAR: Dr. Matoba and then Dr.
10 Bandeen-Roche.

11 DR. MATOBA: You had a slide in your
12 presentation, Dr. McNally, that was not in the
13 original report and that is the average wearing
14 time for the people in the SEE3 group approached 27
15 days by twelve months. I wanted to know if that
16 graph was generated from the same tele-diary raw
17 data that was used for the second table which shows
18 that only 67 percent of time were the patients in
19 the SEE3 group wearing their lenses for 20 to 31
20 days.

21 DR. McNALLY: This graph was new graph for
22 you, but, because of the questions, I thought I
23 should show it. But it was in the trend analysis
24 profile, table 13. So that is directly taken from
25 there. This was taken from the report and the

1 case-report forms when, at each visit, they asked
2 what had been your continuous nights in a row that
3 you had worn the lens. That data came from there.

4 DR. MATOBA: So that is the tele-diary
5 system? It is the same raw data that generated the
6 second table? Is that what you are saying?

7 DR. McNALLY: No; the second table is from
8 tele-diaries. The first is from the case-report
9 forms.

10 DR. MATOBA: Okay, because eyeballing the
11 two, they seem disparate to me because they
12 achieved average wearing time of almost 25 days
13 within one month and they stay at that range, 25 to
14 27. It seems very different from what the tele-
15 diary data reveals.

16 DR. McNALLY: We looked and we included
17 it. It wasn't included in the panel packet but it
18 was included in the PMA application. We looked at
19 the correlation between case-report form and the
20 tele-diary and they matched very closely.

21 Then, to maybe address your question here,
22 the tele-diary graph shown there includes all
23 visits. That is the reports for all visits
24 including the first month and the first week and
25 the whole thing where, as you see in the graph

1 before, the wearing time averages during that first
2 month were lower.

3 So they do correlate because we looked
4 very closely to see is there a difference in the
5 reporting with the tele-diary versus the case-
6 report form and presented in the PMA packet that
7 they correlate.

8 DR. SUGAR: Can I, just to understand. If
9 they took the lens out and cleaned it and put it
10 back in, didn't leave it out overnight, that would
11 still shorten their wearing time or --

12 DR. McNALLY: No; that did not. If they
13 left it out overnight --

14 DR. SUGAR: If they left it out overnight,
15 it did. Okay.

16 Dr. Bandeen-Roche?

17 DR. BANDEEN-ROCHE: First, I would like to
18 add my congratulations to sponsor for their study
19 and their presentation. I especially appreciated
20 the matched design and the wide variety of
21 investigators and the really good-faith attempt to
22 provide adequate power. So thank you very much.

23 I have three questions, one of which is
24 pretty general and the other two are statistical.
25 The general one first may follow up on Dr. Weiss's

1 comment about this being a new polymer. So,
2 certainly, you cited the decreased dryness symptoms
3 in SEE3 but, in conjunction, there was an increase
4 in burning and tearing and I think lens awareness
5 symptoms.

6 It just made we wonder whether there could
7 be a subgroup of patients who don't well tolerate
8 the material. I wonder if you could comment on
9 that.

10 DR. McNALLY: Dr. McNally. There may be a
11 subgroup that doesn't tolerate it. This is why we
12 have stress the first month because we did see --
13 we were surprised that we had a number of events
14 happening in the first month. We tried to look
15 through the data to come to, what can we find about
16 that. We found a few things. The few things we
17 found, the lens fit. I think these were the
18 discomfort and the awareness, and these things I
19 think are easily explained by the lens fit.

20 The burning and stinging, I can't explain
21 directly. So there may be some, I think it is like
22 most contact lenses. There are lenses that
23 patients don't like. In this case, we can't say,
24 here is a patient that may not like this lens.

25 So we really tried to emphasize in the

1 fitting guides, and if there are suggestions on how
2 to better emphasize, we welcome these -- the first
3 month follow-up time to determine, first of all,
4 are you comfortable with it and are you suitable
5 for a 30-night indication.

6 DR. BANDEEN-ROCHE: Thank you. And the
7 two statistical questions. The first is that this
8 was a matched design within investigators. So, to
9 my reading, the analyses did not account for the
10 matching or for correlation within investigators in
11 any explicit way; is that right?

12 DR. McNALLY: Yes.

13 DR. BANDEEN-ROCHE: So FDA, when it
14 reviews the ultimate materials, I think, should
15 look at analyses that do account for that because
16 it has to do with the believability of the
17 confidence intervals and estimates of incident
18 differences within provider rather than across
19 providers.

20 The second question has to do with the
21 adverse-event table, all adverse events. I believe
22 you referred to adjusted versus unadjusted. My
23 reading is that those were not life-table
24 estimates. They were just -- I think it is
25 important to provide life-table estimates for those

1 as well because there was such a difference in the
2 time at risk due to the differential dropout of
3 SEE3 subjects early on.

4 DR. SUGAR: Dr. Pulido. Then I think we
5 are going to be ready for our break and the FDA's
6 presentation after lunch. Go ahead, Jose.

7 DR. PULIDO: Following up on Dr. Bandeen-
8 Roche, page 1251, which is product labeling,
9 Package Insert, rather. It has the adverse device
10 effects were reported at the following annual rate.
11 When I add those numbers up, it is 4.63 percent.
12 That is less than the annualized rate estimation
13 for the primary safety endpoint which was
14 6.1 percent. So did I add up improperly?

15 DR. ROBIRDS: This is Scott Robirds. What
16 we selected in the labeling were just the corneal
17 inflammatory event. The subset of the table that
18 is one page 12 of 21 in your summary of safety and
19 effectiveness, the very first section, where it is
20 a comprehensive list of adverse device effects.

21 We elected to focus on just those corneal
22 inflammatory events which totals -- that total that
23 you mentioned. But there were other events,
24 obviously.

25 DR. PULIDO: But the primary safety

1 endpoint was infiltrative keratitis; correct?

2 DR. McNALLY: It was infiltrates grade 3
3 or greater or with overlying staining.

4 DR. PULIDO: Right. And you have an
5 annualized rate estimation of 6 percent. So how do
6 you justify saying, later on, that the annual rate
7 is 4.63 percent.

8 DR. McNALLY: This is Dr. McNally. First,
9 in the proposed labeling, I will make the comment
10 that any recommendations the panel will make in
11 terms of what you include in here, we are very fine
12 with that. These rates here, they don't include --
13 if you look at them, they don't include the ones
14 that were under serious adverse device effect and
15 perhaps they should have.

16 But there were a couple of cases in there
17 with anterior-chamber reaction. So we pulled these
18 as a first proposal directly off the table. We
19 didn't include everything in there and we would be
20 very happy to include whichever the panel thinks is
21 important.

22 DR. SUGAR: Sally has some comments.

23 MS. THORNTON: I just wanted to let the
24 panel know that the lunches that you have ordered
25 are here. We have reserved room 20G for panel

at

130

1 folks to eat in, and the sponsors, we have reserved
2 room 20H.

3 We would like to advise everyone to leave
4 the room. We have to clear the room completely
5 during the lunch break for security purposes.

6 DR. SUGAR: Everyone, please try to be
7 back here by 1 o'clock.

8 [Whereupon, at 12:05 p.m., the proceedings
9 were recessed to reconvene, at 1:00 p.m., this same
10 day.]

A F T E R N O O N P R O C E E D I N G S

[1:05 p.m.]

DR. SUGAR: We will now proceed with the
FDA presentation on PMA P010019.

FDA Presentation

DR. SAVIOLA: Thank you, Dr. Sugar. At
this time, I will introduce Myra Smith who is a
microbiologist in our branch and the project leader
for this review group. Any additional comments I
will reserve until after we present the questions,
if you have any questions for us.

MS. SMITH: I am Myra Smith. The primary
panel reviewers for this PMA were Dr. Matoba and
Dr. Jurkus. The FDA team responsible for review of
this PMA included Dr. Bernard Lepri, clinical
review, Dr. Gene Hilmantel, statistical review, Dr.
Daniel Brown, toxicology review, Dr. Jimmy Chen,
chemistry review and myself for the microbiology
review.

Dr. Lepri will now present the clinical
issues.

DR. LEPRI: Good afternoon, members of the
panel, sponsors and other guests.

[Slide.]

I am about to present to you just some key

1 elements upon which I believe your review and
2 recommendations should focus regarding this device
3 for the SEE3 Focus contact lens.

4 [Slide.]

5 The history of extended-wear contact
6 lenses is one of low patient satisfaction,
7 unfavorable rates of complications and higher risks
8 of complications.

9 [Slide.]

10 The primary complication of concern, both
11 historically and here today is that of corneal
12 ulcers. The relationship of hypoxia and the
13 development of complications, namely infiltrates
14 and ulcer development, is well known as reflected
15 in this slide.

16 [Slide.]

17 The sponsor believes that the development
18 of the SEE3 lens, lotrafilcon A, addresses these
19 issues. Lotrafilcon A is a very high Dk lens. The
20 Dk of SEE3 is 140. The characteristics of their
21 device and the sponsor's presentation emphasizes
22 the role of oxygen permeability in its performance.

23 [Slide.]

24 The range of power of lenses studied in
25 this investigation were from -6.00 to +6.00

1 diopters with a mean of -3.05. In the next slide,
2 one can see that there is a notable difference in
3 the range of lens powers tested in the clinical
4 trial as compared to the ranges available by the
5 sponsor.

6 [Slide.]

7 Unlike other refractive devices whose
8 ranges of affectedness are limited to those
9 studied, FDA has established policy over the years
10 to deal with this technical discrepancy. This
11 policy addresses the issue of the safety with
12 respect to lens thickness and higher powers as
13 related to oxygen permeability.

14 [Slide.]

15 This has been established by FDA based
16 upon the maturity of contact-lens technology and
17 FDA's experience in dealing with this issue. FDA
18 determines the appropriate range of power approval
19 for extended-wear lenses based upon the effects of
20 lens thickness on lens permeability. The sponsor
21 will have to demonstrate these data to FDA before
22 final approval.

23 [Slide.]

24 In order to achieve their goal of
25 marketing the SEE3 lens for 30-day extended wear,

1 the sponsor, in communication with FDA, designed a
2 prospective, randomized open-label clinical trial
3 for a determination of noninferiority to the
4 control device.

5 A note I would like to add is that, in
6 conversations and communication with the company,
7 in the preparation of this IDE several years ago,
8 the wide range of rates reported in the literature
9 were what contributed to their selection of the 8.6
10 percent infiltrate rate to use as a benchmark for
11 targeting a sample size that would yield sufficient
12 number of patients to provide some reasonable
13 assurance of safety and effectiveness when combined
14 with the postapproval study. So it was the attempt
15 to not have an overly burdensome investigation and
16 yet not have one that produced so few patients that
17 we had absolutely no confidence in the data.

18 [Slide.]

19 Based upon the reported Acuvue infiltrate
20 of 8.6 percent as reported in the literature, this
21 surrogate endpoint was selected utilizing the
22 criteria presented in this slide.

23 [Slide.]

24 This surrogate endpoint was chosen because
25 of its effects upon sample size and due to the fact

1 that most infiltrates are not infectious.
2 Infiltrate development usually precedes ulcer
3 development. This endpoint would provide an
4 estimate of safety upon which a postapproval study
5 would be conducted to attempt to determine the true
6 rate of microbial keratitis for this device.

7 [Slide.]

8 The design was based upon an enrollment
9 number and endpoints determined by a per-patient
10 perspective. 697 test patients were enrolled and
11 this translates to an enrollment of 1,394 eyes for
12 SEE3 which provided reasonable sampling to achieve
13 an estimate of the rate of infiltrates.

14 [Slide.]

15 It is unreasonable to speculate that
16 everyone who is fit with extended-wear lenses could
17 or should wear them for 30 days. Special
18 consideration was given to this fact in the design
19 of this study. It was intended that this study
20 would determine the proportion of patients that
21 could safely wear this type of contact lens for 30
22 days. The endpoints in the study were tailored
23 according to this consideration.

24 [Slide.]

25 I am now going to present to you some key

1 clinical-trial observations that FDA believes
2 should be taken into consideration in your
3 recommendations regarding the SEE3 30-day extended-
4 wear lens.

5 [Slide.]

6 One interesting observation in this study
7 was that an infiltrate event in one eye carries a
8 six-times greater risk of a second event in the
9 same or fellow eye as compared to having a first
10 event.

11 [Slide.]

12 Another finding is that SEE3 infiltrate
13 endpoints occurred earlier in the study than did
14 Acuvue endpoints. Standard contact-lens labeling
15 generally states that the incidence of ulcers
16 increases with the length of wear time. FDA
17 requests that the panel's discussion of labeling
18 will address whether this general warning about
19 ulcers regarding wear time should be kept in the
20 labeling or should labeling reflect the findings of
21 this specific study for SEE3.

22 [Slide.]

23 42.4 percent of the 33 SEE3 subjects who
24 developed infiltrate endpoints experienced them at
25 one month whereas only 23.8 percent of the Acuvue

1 subjects did at this time. From the second month
2 on, the number of endpoint events was similar in
3 number and timing of occurrence with 19 occurring
4 for SEE3 and 16 for Acuvue.

5 [Slide.]

6 For subjects that experienced more than
7 one endpoint event, which were 10 in number for
8 SEE3 and 4 for Acuvue, 70 percent of SEE3 subjects
9 experienced the endpoint in the first month as
10 compared to 25 percent for Acuvue. This can be
11 inferred to mean that SEE3 events occur early on in
12 wear when patients are most closely monitored.

13 [Slide.]

14 The study results also revealed that there
15 were no differences in gender or age for
16 infiltrates. For this study, infiltrates were not
17 restricted to daily wear or new lens wearers. 9
18 out of 13, or 69.23 percent of SEE3 and 100 percent
19 of the 9 Acuvue subjects who experience infiltrates
20 had worn extended-wear lenses on a 7-day basis
21 prior to participation in this study.

22 [Slide.]

23 Adverse events as related to wear time are
24 a major issue in the evaluation of extended-wear
25 contact lenses. The average wear time in this

1 study for all completed patients was 27 days at 12
2 months. This was achieved by 67.2 percent of the
3 dispensed cohort who had completed the study. Even
4 28.9 percent of discontinued patients wore the
5 lenses for an average of 27 days.

6 [Slide.]

7 Of the discontinued patients, only 2.4
8 percent were discontinued for positive
9 biomicroscopy findings. The majority were
10 discontinued for lens-fit discomfort and acuity
11 followed by lost-to-follow-up. All of these issues
12 have been addressed in the sponsor's presentation
13 this morning.

14 [Slide.]

15 Some of the most important aspects of a
16 clinical trial are those that occur between
17 scheduled visits. In order to attempt to obtain
18 some of this information, the sponsor included a
19 patient-managed daily diary in this investigation.
20 Review of this information by the sponsor revealed
21 that SEE3 patients had fewer complaints of dryness
22 than the Acuvue patients. 19.8 percent of SEE3, as
23 compared to 24.2 percent of Acuvue patients,
24 reported dryness.

25 [Slide.]

1 Statistical analyses of these subjective
2 reports were found to be significant by the
3 sponsor. The sponsor proposes that the labeling
4 claim that SEE3 lenses reduce dryness symptoms
5 associated to wearing hydrogel lenses. The panel
6 should address the issue of this finding and its
7 clinical significance in the labeling discussion
8 here today.

9 [Slide.]

10 Question No. 1: do the data presented in
11 PMA P010019 provide reasonable assurance of safety
12 and effectiveness for the proposed indication for
13 use?

14 [Slide.]

15 This is the indication statement
16 concerning the first two issues which we believe
17 are the focus of this PMA discussion today
18 regarding the general indication for refractive
19 conditions and length of wear.

20 [Slide.]

21 Question No. 2: does the panel recommend
22 any modification of the proposed wording of the
23 indication statement?

24 [Slide.]

25 Question No. 3: please discuss the merits

1 of including the maximum 30-day time period in the
2 indication statement. Does the panel recommend
3 that it be included in other sections of the
4 product labeling rather than the indication
5 section?

6 [Slide.]

7 Question No. 4: does the panel have any
8 specific recommendations for the proposed product
9 labeling in terms of warnings, precautions,
10 clinical data outcomes or practitioner-directed or
11 patient-directed labeling?

12 [Slide.]

13 Question No. 5: does the panel recommend
14 that the sponsor conduct a prospective postapproval
15 study within the U.S. population to gather
16 information on the incidence of microbial
17 keratitis?

18 [Slide.]

19 Following that question is No. 6 in topic
20 and in number: in consideration of the potential
21 differences in the standard of care and device-
22 usage patterns outside of the United States, does
23 the panel have any recommendation concerning the
24 use of foreign data in the postapproval study?

25 Thank you for your time.

1 DR. SUGAR: Thank you. Does that end the
2 FDA presentation or do you have more?

3 DR. LEPRI: That is pretty much it.

4 DR. SUGAR: Are there questions for FDA?
5 Jose?

6 DR. PULIDO: Jose Pulido. Dr. Lepri,
7 again, as I try to resolve in my mind what is my
8 biggest concern, knowing what I had discussed this
9 morning, do you feel that there is a 1.5-fold or
10 greater risk of infiltrative keratitis for this
11 lens versus a 7-day-wear lens?

12 DR. LEPRI: The data show that the rates
13 are definitely higher than they are for 7-day
14 lenses. But then, again, that was expected. In
15 fact, at the panel meeting when we discussed these
16 issues last November, Dr. Hilmantel's presentation
17 was asking the panel to conjecture on what X amount
18 of fold increase would the panel find acceptable
19 for marketing a new 30-day lens.

20 Those numbers that were recommended by the
21 panel were much higher, 2, 3 and 4 times, when he
22 presented those data. This is actually much lower
23 than I would have expected to see. But it is
24 definitely higher than 7-day. It stands to reason.
25 The longer you wear it, the longer the cornea is

1 stressed.

2 I would like to make one more comment that
3 I forgot to make.

4 DR. SUGAR: Go ahead.

5 DR. LEPRI: That was to thank and commend
6 the sponsors for providing me with a very concise,
7 succinct and fluent document to review and for
8 their extreme cooperation and helpfulness in
9 working through this entire process in the past
10 five years that I have been with FDA.

11 DR. SUGAR: Thank you. Are there other
12 questions for FDA? If not, the sponsor, if they so
13 choose, can make comments. We have ten minutes for
14 that. Do you wish to retake the floor? Seeing no
15 desire, we will then proceed -- I think we will
16 reserve the right to question both the agency and
17 the sponsor if the need arises in our
18 deliberations.

19 We will now move on to the deliberations.

20 **Committee Deliberations**

21 DR. SUGAR: We are going to begin with the
22 primary reviews, the first of which is Dr. Jurkus.

23 DR. JURKUS: This is Jan Jurkus, the
24 primary reviewer. I would like to start out my
25 review by saying thank you very much to Dr. Lepri

1 and Dr. Hilmantel for their excellent, excellent
2 reviews that were given to me and also to Ciba for
3 a very readable and sort of straightforward report.

4 Much of my review has been already talked
5 about so I will try to make it brief in terms of
6 the highlights. Things that I find to be of major
7 interest in this report include the lens material,
8 itself, this being a lens that has a
9 transmissibility of 175 times 10^{-9} in the -3.00
10 power. This certainly, as a high oxygen-
11 transmissible lens, does, indeed, exceed the
12 criteria that was set forth by Holden and Mertz of
13 '87 as well as the more current criteria proposed
14 by Lehood of 125.

15 So the actual oxygen transmission is
16 something that I think, as a practitioner, we
17 certainly look forward to. In reviewing the study,
18 itself, some things that I found of interest,
19 starting with the number of lenses that were not
20 dispensed in terms of the trial lenses. There were
21 39 subjects who did not get lenses dispensed to
22 them as part of the study.

23 Well over -- or actually about 50 percent
24 were due to the inadequate fit. So, changing or
25 adding an additional base curve would, at least in

1 theory, take care of 50 percent of the people that
2 were unable to have this lens prescribed for them.

3 In the number that had been discontinued,
4 again, poor vision, discomfort and lens fit made up
5 the majority of the reasons for the lens to be
6 discontinued. There, too, the statement that the
7 sponsor makes that a flat fit with the SEE3 may
8 also result in small amounts of edge lift, that may
9 be judged better by subjective reports of lens
10 awareness or discomfort than biomicroscopy findings
11 is one that is very interesting and I think needs
12 to be highlighted very carefully in the
13 practitioner manual.

14 In the past, as practitioners, we were
15 always looking to fit the loosest lens that was
16 stable on the eye where here the loosest lens may,
17 indeed, not be the most appropriate for a
18 particular patient or a particular group of
19 patients. I think that should be certainly
20 addressed in the labeling portion of the
21 practitioner guide.

22 When it comes to the safety endpoints, I
23 agree with Dr. Hilmantel's assessment that it can
24 be concluded that the SEE3 lens is not inferior
25 within a tolerance of 0.05 to the Acuvue lens with

1 regard to the primary safety endpoint of corneal
2 infiltrates with staining or grade 3 infiltrates.

3 The timing, again as illustrated in all of
4 the reviews thus far, was something that I did find
5 to be very interesting and, again, should be
6 indicated very much in labeling that the infiltrate
7 existence was much sooner with the SEE3 than with
8 the control lens. I think that should certainly be
9 highlighted and stressed to practitioners.

10 When it comes to the percentages of
11 serious significant adverse events, nonsignificant
12 adverse events, the study did show that they were
13 really remarkably similar between the SEE3 and the
14 Acuvue lenses.

15 The thing that did stand out, as commented
16 earlier, was the development of CLPC, the contact-
17 lens-induced papillary conjunctivitis. This is
18 something that I think practitioners, again, need
19 to be made very much aware because GPC as a whole,
20 or CLPC, had been sort of dwindling in clinical
21 practice and now this may be a resurgence to be
22 checking for, although the incidence rate of 4.6
23 percent certainly did fall within the percentages
24 for extended wear that are included in the
25 literature. Those that I could find were between 2

1 and 16 percent.

2 So it wasn't outrageous but, certainly,
3 within the portions that are currently available.

4 When it comes to subjective symptoms, the
5 report of dryness was 19.8 for the SEE3 group and
6 24.2 for the control group. Although statistically
7 that has been shown to be significant, when you
8 think of it in clinical terms, it sort of breaks
9 down to, with the SEE3 group, one out of five
10 people is going to tell you that they experienced
11 dryness. With the Acuvue, one out of four people
12 is going to tell you that they experienced dryness.

13 From a clinical standpoint, when you are
14 in a busy practice, that one-person difference
15 doesn't seem to make a huge influence on a
16 practitioner's selection of choice. So I still
17 have a hard time with the statement that they had
18 proposed for labeling regarding dryness.

19 When it comes to the visual outcome, I was
20 very pleased that the 98.1 percent of the test
21 group maintained acuity within two lines of
22 dispensing as well as the efficacy outcome that
23 95.5 were able to wear to lenses for 22 to 31 days.

24 So I guess, to sort of summarize my review
25 of this, when I looked at the whole thing, putting

1 it in very simplistic terms, what I was hoping this
2 document would answer would be two things; one is
3 does the lens work and, secondly, does the lens do
4 any harm.

5 To answer those two specific questions, I
6 can say that the answer to, does the lens work,
7 does it do what they say it is going to do, I would
8 have to say the answer to that has shown to be yes,
9 that people certainly can see with this lens on and
10 that it does provide extended-wear capabilities.

11 The vision measurement to be 20/20 was
12 achieved by 83 percent of the subjects while
13 maintaining Snelling contact-lens acuity within two
14 lines of dispensing was achieved by 98.1 percent.
15 This I thought was a very remarkable and very
16 laudatory achievement.

17 Continuing with that answer, can people
18 wear this contact lens on an extended-wear basis
19 for up to 30 days, again, the numbers were a little
20 bit confusing between the 67 percent and the 95
21 percent, but I would certainly say that it is safe
22 up to about that 30-day for most people.

23 Looking at the second question, does this
24 do any harm, for that part, we are not really sure.
25 At this point, the study had showed that there was

1 no significant harm done and that the 5 percent
2 endpoint infiltrate rate does not seem to be
3 totally different than what is currently available.

4 I do believe very strongly that the
5 postmarket surveillance study will give us a much
6 better answer to that particular question.

7 So, at this point, in my opinion -- I am
8 not supposed to give my opinion yet, until we have
9 completely discussed this?

10 DR. SUGAR: You may.

11 DR. JURKUS: I can give my opinion? Okay.
12 In my opinion, I think labeling certainly can
13 address some of the issues that we have discussed
14 but this has been shown to be a safe and effective
15 lens.

16 DR. SUGAR: Thank you.

17 Dr. Matoba?

18 DR. MATOBA: Thank you. I just had a few
19 relatively minor points. Since many of the points
20 actually pertain directly to the questions, I
21 thought maybe I would just go down and discuss each
22 of the questions that the FDA has posed.

23 The first question was do the data
24 presented in PMA P010019 provide reasonable
25 assurance of safety and effectiveness for the

1 proposed indications for use. Initially, I was
2 troubled by the fact that only 67 percent of
3 subjects in the SEE3 study had worn the lenses for
4 the 22 to 31-day period whereas 92 percent of
5 subjects in the Acuvue study had worn the lenses
6 for the 5 to 7-day period. So for a basis of
7 comparing the incidence of the endpoint
8 infiltrates, it seemed to me that they really have
9 not compared 7-day wear versus 30-day wear and yet
10 wanted approval for 30-day wear.

11 On the other hand, based on the data
12 presented today looked at it another way, the
13 average wearing time goes up to 25 days and
14 approaches 28 to 27 days by the end of the 12-month
15 period. I am no longer as bothered by that
16 discrepancy. In terms of other questions I had
17 regarding the nature of the infiltrates that were
18 seen in patients who were discontinued from the
19 study, the sponsor has addressed my questions from
20 the initial review.

21 So my answer to No. 1 would be yes.

22 For 2, would I recommend any modifications
23 in the proposed wording of the indications
24 statement, I still have a problem with the dryness
25 symptom as an indication for use of this

1 potentially 30-day extended-wear lens.

2 If a practitioner were to look at that
3 indication only and then not go down to look at
4 contraindications, it may be construed as
5 recommending that a patient with potentially
6 aqueous-tear deficiency and related ocular-surface
7 disease may be an appropriate patient for
8 dispensing of the 30-day lens. I would be very
9 concerned about that possibility.

10 So rewording of this indication or some
11 other modification, as the sponsor has already
12 indicated they may be willing to consider, would be
13 appropriate, I think.

14 In terms of the third question, I had no
15 problem with the maximum 30-day indication in the
16 statement.

17 In terms of the fourth question, proposed
18 labeling changes, I think that the fact that this
19 lens did have a statistically significant increased
20 incidence of GPC in their study patients should be
21 included and sponsor has already indicated that
22 they would include that in the labeling.

23 The second thing I would like to suggest
24 is that labeling include the fact that once a
25 patient has had one infiltrate, they are at greatly

1 increased risk for a second infiltrate. So the
2 practitioner should use extra caution in monitoring
3 those patients who have had at least one
4 infiltrate.

5 In terms of the fifth question, does that
6 panel recommend that the sponsor conduct a
7 prospective postapproval study, I would. The
8 sponsor has already indicated that they have plans
9 to proceed with that study.

10 The sixth question was are there any
11 special considerations for the study outside the
12 U.S. My answer would be no.

13 DR. SUGAR: Go ahead Ralph.

14 DR. ROSENTHAL: Dr. Rosenthal. Could I
15 just make two comments about the questions. I
16 don't know whether you want me to make them now or
17 whether you want me to make them before you start
18 to discuss them specifically.

19 DR. SUGAR: I think it is fine. Go ahead.

20 DR. ROSENTHAL: Thank you. I think I am
21 getting this right. Jim, correct me if I am not.

22 DR. SAVIOLA: I am listening carefully.

23 DR. ROSENTHAL: In the past, we have not
24 prepared patient-directed labeling with contact
25 lenses. Question 4 specifically asks should we

1 require the company to prepare patient-directed
2 labeling as opposed to just practitioner-directed
3 labeling.

4 DR. SUGAR: Could you define that for us?
5 What does patient-directed labeling mean?

6 DR. ROSENTHAL: As with other devices, the
7 patient-information booklets have to be made up by
8 the company and provided to the practitioner who
9 dispenses the lenses. Sorry; that would be if they
10 were required.

11 So Dr. Matoba just brushed by that
12 question rather quickly and I want to be sure you
13 discuss that issue because I think that brings up
14 the second issue is that we really -- if a company
15 is proposing a 30-day approval study and, of
16 course, we are agreeing that it should be included
17 as a condition of approval, there are still
18 questions out there about the safety of the lens
19 over a 30-day period.

20 The past has shown that, long before I
21 came to the FDA, lenses were approved for a certain
22 period of time and then, because of problems out in
23 practice, they had to reign in the time of
24 approval. I am not sure you are all aware of that.
25 I wasn't in this country when it all happened, but

1 apparently it happened in the '80's.

2 So there have been issues in the office
3 about whether or not the incidence should include
4 the 30 days when we may, in fact, have to reign it
5 in after a postapproval study. That is why we
6 specifically asked about that issue, if the
7 postapproval study shows a very incidence of
8 microbial keratitis.

9 DR. SUGAR: Could I ask what our options
10 are? If we feel that this is demonstrated safe and
11 effective for 30-day use but there is a concern
12 about the postapproval study, what middle ground do
13 we have to approve the lens but reserve the option
14 which you always have, of course, to change the
15 indications in the future.

16 DR. ROSENTHAL: Correct. Dr. Saviola will
17 be happy to answer that for you.

18 DR. SAVIOLA: Let me take a step back
19 first before I answer that direct question. What
20 Ralph is alluding to is a couple of points. On
21 Question 4 regarding panel-specific recommendations
22 for labeling and patient-directed labeling, one of
23 the things, as you brought up, the patient booklets
24 are passed out to the different doctor offices by
25 the different account managers, detail people,

1 whoever.

2 They are not necessarily always passed out
3 to the patient. So the thought is, is there some
4 other vehicle that people should perhaps get risk
5 information about this device and, along those
6 lines, the traditional ones, as Ralph described,
7 with the package insert that is directed to the
8 practitioner, the practitioner fitting guide and
9 the patient information booklet, what we normally
10 have seen, but the concept, perhaps, of a patient
11 package insert might be something that you want to
12 think about or talk about in the context of your
13 discussion similar to what you have seen often in
14 different pharmaceutical advertisement, the back
15 page of an ad will have the patient package insert,
16 essentially, which has some information in it that
17 talks about the fundamental information that is
18 found in the regular package insert which is
19 warnings, precautions, contraindications, et
20 cetera, but not using technical terms or technical
21 terminology to that degree. That was the first
22 thing.

23 The second thing Ralph was talking about
24 was the idea that while the proposed indication has
25 lots of different elements in it, wear time being

1 one of them, the idea that -- and this is, again,
2 in the context of our questions that we posed to
3 you for discussion, the idea that the wearing
4 period, or the recommended wearing period, may or
5 may not have to be part of that specific
6 indication.

7 As it is proposed now, it is written up to
8 30 days as recommended by your eye-care
9 practitioner. Well, that second part, as
10 recommended by your practitioner sort of gets
11 forgotten and it becomes a 30-day lens.

12 So in the context of a failed postapproval
13 study where the rate is significantly higher, if we
14 were to have to make an adjustment later on in the
15 maximum wearing period, the indication, if it just
16 said for correction of refractive error would
17 remain the same, and the modification would occur
18 in a different part of the labeling, such as
19 prescribing information, wearing time, what have
20 you.

21 So, in essence, our questions are getting
22 to the discussion of your clinical viewpoints, pros
23 and cons of having the maximum 30-day period in the
24 indications statement as opposed to some other part
25 of the labeling such as prescribing information.

1 DR. SUGAR: Dr. Pulido?

2 DR. ROBIRDS: Jose Pulido. I don't
3 understand. You are saying that, fine, comes out
4 the postapproval study and it shows there may be an
5 increased risk and you want to back off a little
6 bit, so you have to change it down in your scheme
7 of things, down where it says length of time, but
8 it wouldn't be in the indications.

9 The way it is set up now, you would just
10 change it in the indications, so what is the
11 difference?

12 DR. SAVIOLA: The other element that is --

13 DR. ROSENTHAL: Let me just clarify
14 something. It is not our scheme of things. We are
15 asking the panel's recommendation. So we have not
16 proposed either. The company has proposed an
17 indication statement including a 30-day lens. We
18 have just raised the specter of another
19 possibility.

20 DR. PULIDO: The question, still, is
21 what's this difference.

22 DR. SAVIOLA: Far be it for us to lead you
23 in your determinations. We were asked to bring the
24 idea to you for your comments and, during the
25 course of the discussion, perhaps alternatives.

1 While we do know that there is reasonable assurance
2 of safety and effectiveness up to 30 days in the
3 population that has been studied, we know from
4 experience that these preapproval studies don't
5 translate into the general population.

6 So the true incidence of microbial
7 keratitis in the general population really hasn't
8 been studied or established. Having gone down this
9 road before, we have some experience here so, in
10 the internal discussions in the office, there is a
11 mix of opinions one of which is that there might be
12 some merit in not including the length of wear time
13 specifically in the indication statement because,
14 A, that might push people to wear it as a 30-day
15 lens and forget the second part, "as directed by
16 your eye-care practitioner."

17 Two, because we don't really have the full
18 picture -- we have a preliminary picture at this
19 point in time and we know from experience that the
20 preliminary picture didn't translate to the general
21 population.

22 DR. SUGAR: Dr. Matoba?

23 DR. MATOBA: Alice Matoba. It seems to
24 me, whether it is in the indication or not, if it
25 is anywhere in the labeling or the advertising, it

1 is going to be considered 30-day lens no matter
2 where you put it. My question is, actually, I was
3 surprised that you allowed data from the patients
4 who hadn't worn the lenses for 30 days to be
5 included in your study.

6 If you had problems before, why would you
7 not have told the sponsor up front that they would
8 have to design a study that would strictly compare
9 7-day versus 30-day wear?

10 DR. SAVIOLA: The initial brilliant idea
11 we had in how to deal with these devices the second
12 time around was to allow subjects in the study to
13 wear the lens for whatever period of time that they
14 would tolerate so we would have a distribution of
15 7-day, 14-day, 21 to 30-day, a strata of outcomes
16 to get some sense for how often people could really
17 tolerate this lens because, again, our sense is
18 going to be that not everybody is going to wear
19 this lens for a month.

20 That idea didn't really pan out because,
21 in the course of the study, they ramped everybody
22 up to 30 days. So, it didn't really disturb us
23 much that there were people who completed the study
24 in less than 30 days because that was sort of one
25 of our original expectations of the outcomes, that

1 there be a natural demarkation of people who could
2 tolerate different periods of wear.

3 We are seeing here that that really played
4 out. So, to say you must have a 30-day wear period
5 of be discontinued from study wasn't really
6 consistent with the way we were trying to get some
7 information about how this would translate into the
8 total population.

9 DR. SUGAR: Dr. Zadnik? I'm sorry; Alice.
10 Did you have a follow up to that?

11 DR. MATOBA: No.

12 DR. SUGAR: Dr. Zadnik?

13 DR. ZADNIK: I am not sure I understand
14 who the postapproval study that is proposed is
15 going to resolve any of this. If you enroll 2,000
16 people for a year and this lens at 30 days, or 20
17 days, or however long the people end up wearing it
18 or the practitioners recommend it is just as good
19 as or bad as our experience has been so far, you
20 are going to get 4; right?

21 Does that mean we are going to be sitting
22 here and saying, well, we got 6, so it is a lot
23 worse. Is it going to have to be that we got 40 so
24 it is a lot worse? Or we got 2, so it is a lot
25 better? Or we got none out of those 2,000 people

1 in a year so it is a ton better?

2 You know, I realize there has to be a
3 limited scope to that to make it even feasible, but
4 I am not sure doing something of a limited scope
5 for the sake of doing it if it is not going to
6 really answer the study question -- that is, is the
7 annualized microbial-keratitis rate greater than 20
8 per 10,000 in Focus Night and Day wearers, I don't
9 get the point of the postmarket study other than us
10 --

11 DR. SAVIOLA: You are correct that, at
12 that scope, it won't answer the question. We ask
13 the questions to you for discussion purposes of
14 your opinions about a postapproval study. We did
15 not want to get into a whole discussion of the
16 specifics of that simply because it would get
17 really convoluted very quickly.

18 We had a discussion of this at our
19 November meeting last fall and got some sense for
20 it. We have already had discussions with the firm.
21 We had discussions with the industry in general.
22 There are some considerations in terms of how much
23 it is going to cost companies to do these studies.

24 We are not going to be satisfied with
25 2,000 people.

1 DR. ZADNIK: I just want to --

2 DR. ROSENTHAL: Excuse me, Dr. Zadnik.
3 This is Dr. Rosenthal. This is proposed by the
4 company.

5 DR. ZADNIK: I understand.

6 DR. ROSENTHAL: This has not been agreed
7 to by the agency.

8 DR. ZADNIK: I just want to sort of enter
9 a cautionary note that if, as we have these
10 discussions, we fall back in, "It's okay; there is
11 this postmarket study and some of these questions
12 will be answered." I am just not sure in my head I
13 could design a feasible study that would answer
14 some of these questions until this lens is out in
15 the hands of practitioners.

16 DR. SAVIOLA: Right. Our initial goal was
17 a 10,000 to 15,000-patient study which is
18 significantly expensive to conduct. It will be
19 somewhere between 10,000 and 2,000 as an initial
20 study. Then, depending on the outcomes, we go from
21 there. If it is something that shows consistent
22 with the preapproval data, then we are all set. If
23 it is something that shows it is questionable,
24 maybe there is need for additional studies after
25 the first one. Who knows?

1 But, as it is presented to you in this
2 context, this is the company's initial offer, so to
3 speak.

4 DR. SUGAR: The answer to question 6 has
5 an impact on that, also, in terms of what other
6 data they can recruit for dealing with the
7 question.

8 Did you have something else you wanted to
9 say, Jim? What we will do is have a little bit
10 more of this general discussion and then we will go
11 specifically question through question.

12 Sally wants to know if you want to sit
13 down, Jim.

14 DR. SAVIOLA: Unless you have any other
15 questions, I will sit down.

16 DR. SUGAR: She is very interested in
17 everybody's comfort today. That is nice to see.

18 Dr. Weissman?

19 DR. WEISSMAN: This is Weissman. I had a
20 question specifically about the indication for
21 aphakic use. As far as I know, none of the
22 subjects in the initial study were aphakic and many
23 of us who have seen aphakic patients have a bias
24 that aphakic patients may have a higher rate of
25 infection.

1 So the question I have is why is that on
2 the table?

3 DR. SUGAR: I think we can ask the sponsor
4 to make a brief comment on that.

5 DR. McNALLY: This is John McNally with
6 Ciba Vision. That we did not, indeed, study in
7 aphakic patients. The labeling, that is the
8 standard labeling for most contact-lens approvals.
9 So we did not study it. We think it might be an
10 interesting thing to study and proceed with, but we
11 put that in because that is the standard labeling
12 for contact lenses. That is very much open for
13 discussion.

14 DR. WEISSMAN: It might be impossible to
15 do because there are not many adult aphaks running
16 around. I just wondered why it was there. It
17 might be something that the agency might want to
18 consider.

19 DR. SUGAR: It may be something that we
20 may want to, at the end, with -- a change in the
21 labeling.

22 Dr. Bandeen-Roche? She is just agreeing.
23 Okay. Other general comments? Dr. Zadnik?

24 DR. ZADNIK: Dr. Karla Zadnik. Dr.
25 Matoba, you mentioned, I think, the papillary

1 conjunctivitis sort of label warning kind of thing.
2 One of the problems I see in this dataset for
3 making that kind of statement, that either this
4 lens is riskier in terms of that developing or less
5 so, is that the randomization didn't work for a
6 previous history of giant papillary conjunctivitis.

7 There are a lot more patients in the SEE3
8 group than in the Acuvue group who had a previous
9 history of contact-lens papillary conjunctivitis, I
10 think. Isn't that what the data say? So I think
11 that to then say this lens is at increased risk, I
12 think is impossibly confounded, perhaps, by that
13 historical risk factor.

14 DR. MATOBA: Alice Matoba. I think you
15 might say something like, "in the study, a greater
16 incidence of GPC was found."

17 DR. ZADNIK: Could it say something, that,
18 specifically in people who had a previous history
19 of or could it mention -- in other words, if you
20 are a previous GPC sufferer, this might not be the
21 lens for you.

22 DR. MATOBA: That would be fine. But I
23 don't think you can just throw it out because you
24 found a way to explain it, because that is the
25 study that was done and that is what it showed.

1 DR. ZADNIK: I think you want to advise
2 patients who should and who should not try this
3 lens

4 DR. MATOBA: Yes; that would be fine.

5 DR. ZADNIK: Maybe that is who shouldn't
6 try it.

7 DR. MATOBA: Right.

8 DR. SUGAR: Go ahead, Dr. Weiss?

9 DR. WEISS: Jayne Weiss. That was the
10 comment that I was addressing my question to the
11 sponsor before is that when they separated out
12 those patients who had not had previous GPC, and I
13 will call it GPC because it is just so much easier
14 for me, they had an approximately 3 percent rate in
15 the SEE3 category of GPC in those who did not have
16 a previous history of this, but a 0.9 percent in
17 the Acuvue.

18 I don't know if that is statistically
19 significant and the sponsor evidently didn't have
20 that data. I don't know if the onset was earlier.
21 So it is a question that I think the sponsor should
22 go back and answer. I don't know if and how we
23 should address that particular thing in the
24 labeling.

25 DR. ROSENTHAL: Rosenthal. May I just say

1 that it is obvious that it is important and if you
2 tell us you would like it to be addressed in the
3 labeling, based upon the comments that you have put
4 forward, we will insure that the company does the
5 appropriate analysis to insure that the labeling
6 reflects the various issues.

7 DR. SUGAR: At this point, I would like to
8 organize our discussion around the six questions
9 that the agency presented us with, and begin with
10 the first question. Dr. Matoba, do you want to --
11 you have already, but go ahead and just make a --

12 DR. ROSENTHAL: Excuse me; Rosenthal.
13 Could you take the first question last because,
14 essentially, it is what you are going to be asked
15 to vote upon. So I would rather you -- well, you
16 can do as you wish.

17 DR. SUGAR: I am not sure that that is the
18 case. I think that we can deal with the issue and
19 then deal with the details. That is what the
20 subsequent questions are.

21 Go ahead, Alice. Just restate your stance
22 on the question.

23 DR. MATOBA: All right. Alice Matoba. I
24 am going to restate my stance on the question.

25 DR. SUGAR: Thank you.

1 DR. MATOBA: Word for word?

2 DR. SUGAR: Any way you want.

3 DR. MATOBA: Initially, I was bothered by
4 the fact that only 67 percent of the subjects in
5 the SEE3 group wore their lenses for 22 to 31 days
6 whereas 92 percent of the subjects in the Acuvue
7 group wore the lenses for 5 to 7 days, the upper
8 range of the wearing time.

9 So it seemed to me that there was a bias,
10 possibly due to this discrepancy. But,
11 subsequently, sponsor did show other data
12 indicating that the average wearing time went up to
13 25 days within 1 month and approached 27, 28 days
14 over the next 11 months.

15 So I believe that there was a fair
16 comparison of approximately 30-day wearing time
17 versus 5 to 7 days for the Acuvue group. The
18 incidence of the endpoint infiltrates, the
19 surrogate for a microbial keratitis, was lower than
20 expected and is an acceptably low range for both
21 groups.

22 Other concerns I had regarding the nature
23 of the peripheral ulcerations and infiltrates in
24 the subjects who were discontinued for
25 biomicroscopic findings were addressed by the

1 sponsor. So, at this time, I feel that there has
2 been provided reasonable assurance of safety and
3 effectiveness for the proposed indications for use.

4 DR. SUGAR: Is there anyone who feels
5 otherwise and would like to discuss it? Please.
6 Dr. Bandeen-Roche.

7 DR. BANDEEN-ROCHE: I am Dr. Bandeen
8 Roche. I am not saying that I feel otherwise, but
9 I did want to make just a couple of statements
10 about my view of the data. So I have to rely on my
11 panel associates' judgment to some degree, their
12 clinical judgment.

13 The first issue is how good of a surrogate
14 are corneal infiltrates for the outcome that we
15 ultimately care about, microbial keratitis. I
16 would be very interested to see what the corneal
17 infiltrate rate at the grade that has been defined
18 in this study was in the old time extended-wear
19 studies because it would be a real cautionary tale
20 if the rates were similar and yet things ultimately
21 didn't turn out well.

22 Secondly, I would like to reiterate Dr.
23 Pulido's concerns about noninferiority, the
24 tolerance chosen. It is not unreasonable but, to
25 some extent, it is arbitrary. I think that Dr.

1 Matoba did highlight the more important thing which
2 is what is the rate, is it acceptable rather than
3 does it fall within a 5 percent tolerance.

4 The adverse-event rates are likely
5 understated in the SEE3 group because a life-table
6 analysis was not used to develop those rates and
7 the record was such an appreciably higher early
8 drop out in the SEE3 group than in the Acuvue
9 group, so this is something that should be taken
10 into account in evaluating whether those rates are
11 acceptable or not.

12 I am talking about the other adverse-even
13 rates at this point, and then what other we decide,
14 I just feel that it is important that patients
15 understand what we mean by safety and
16 effectiveness, including the sorts of outcomes that
17 this study has not established and was not intended
18 to establish.

19 DR. SUGAR: Other comments? Thank you for
20 those wise comments. I am not supposed to make
21 judgments, but -- okay. We are not going to vote
22 on the answers to these questions, but we are
23 getting a sense of the panel for the agency's sake.

24 The next, and I think important, issue is
25 do we recommend any modification of the proposed

1 wording of the indication statement. Janice, do
2 you want to comment on changes --

3 DR. JURKUS: In the proposed wording, one
4 of the -- or, actually, there were a couple of
5 things that were not included that I would like to
6 have included in the indication statement. In the
7 alternative practices and procedures section, I
8 think it would be important that we include the use
9 of daily-wear contact lenses and also a different
10 alternative to this would be refractive surgery
11 LASIK.

12 DR. SUGAR: Is that for the indication for
13 the labeling? I think that is more a labeling
14 issue.

15 DR. JURKUS: That is more a labeling
16 issue. Okay. Then, in the indications statements,
17 as stated right up here --

18 DR. MATOBA: That is not the whole
19 indications statement, is it?

20 DR. JURKUS: Yes; that is what was in the
21 book.

22 DR. WEISS: It is on page 1348 of 1314 of
23 the sponsor's manual, if you are looking for it.
24 There are two statements that are missing from that
25 on the screen.

1 DR. ZADNIK: Karla Zadnik. I think the
2 one that Dr. Jurkus mentioned in her comments was
3 on this expanded version of it. It is the dryness
4 issue.

5 DR. JURKUS: Right. Dr. Jan Jurkus,
6 again. The indications for use where they do have
7 dryness symptoms, that the Night and Day contact
8 lens may reduce dryness symptoms that are present
9 with regular hydrogen soft contact lenses. I
10 object to that. I think it should be eliminated.
11 They did not truly study what I would consider to
12 be regular hydrogel contact lenses. They looked at
13 one specific type and there are many other types
14 that had not had any indications for study.

15 So, at this point, I would exclude that
16 statement.

17 DR. SUGAR: Are there other agreements,
18 disagreements? I agree. Jayne?

19 DR. WEISS: Jayne Weiss. I would agree
20 with that. I also think Dr. Weissman's comment was
21 an important one is that the lens was not studied
22 in aphakic patients. So I am not sure that should
23 be included as an indication, although I think
24 perhaps, later on in the labeling, we can address
25 the fact that it may be useful in aphakic patients

1 although it was not studied. I feel a little
2 uncomfortable saying it is indicated for aphakic
3 persons when there wasn't one patient in the study
4 who was aphakic.

5 DR. SUGAR: Dr. Yaross?

6 DR. YAROSS: In the context of that, if
7 that is, in fact, what is referred to as the class
8 labeling indication, typically, if something is an
9 across-the-board indication for a class, industry
10 looks to see if there is a specific reason to
11 exclude a specific product from the class.

12 So I think the question there is is there
13 some special reason to believe that this product is
14 specifically inappropriate to aphaks. You might
15 want to consider that as part of this class
16 indication issue because I would expect that many
17 of the other products that carry this indication
18 also have not been specifically studied in aphaks.

19 DR. SUGAR: Dr. McMahon?

20 DR. McMAHON: Actually, responding to that
21 issue, Dr. Weissman is correct in that there is
22 past data with the older form of hydrogels and
23 extended wear did show a higher complication rate
24 in the aphakic population, particularly in the
25 vantage groups. There has been no evidence to

1 controvert that with this particular material, so I
2 would support considering removing aphakic. The
3 next question would be about pseudophakia. I am
4 less worried about that, in addition to the dryness
5 issue.

6 DR. PULIDO: Just a question. How is
7 diabetes taken care of in the contraindications.
8 They have any systemic disease which may be
9 exacerbated by or interferes with contact-lens
10 wear. Then, they have before that, corneal
11 episthesia.

12 Do we need to worry about the effects of
13 diabetes on corneal surface and the ability to use
14 these lenses? Was that even evaluated? I am a
15 retina person, so I am just asking the panel.

16 DR. WEISS: Jayne Weiss. They have, soon
17 after that, a contraindication, with any systemic
18 disease which may be exacerbated or interferes with
19 contact-lens wear. I think that would be fairly
20 global to go through various conditions that could
21 cause, let's say, decreased immunity or increased
22 sensitivity of infection.

23 We can go through -- there are multiple
24 diseases, aside from diabetes. But I think they
25 have good will in terms of trying to indicate that

1 there may be other diseases that a practitioner
2 might want to consider not using the lens.

3 DR. SUGAR: Let's stay with the
4 indications. Dr. Weissman?

5 DR. WEISSMAN: This is Weissman. I agree
6 with Dr. Weiss that there are an awful lot other
7 immune diseases. I think covering with a global
8 statement is appropriate. In the aphakia think, I
9 want to make it plain that I don't consider,
10 necessarily, this lens to be a problem for aphakic
11 patients, that aphakic patients often have a lot of
12 comorbidities that is what maybe has caused the
13 problem in the past. But the data was while not
14 absolutely convincing, given the old statistics and
15 the few numbers, quite compelling at the time that
16 aphaks did run into an awful lot more trouble
17 attempting extended wear than phakiks did.

18 DR. SUGAR: Go ahead, Dr. Jurkus.

19 DR. JURKUS: Jan Jurkus. One of the
20 reasons, from my understanding that most of the
21 aphaks did have more difficulty, could have also
22 been with the oxygen transmission through the older
23 types of lenses where, indeed, this lens, having a
24 much higher oxygen transmission, may actually
25 benefit that aphakic population as opposed to

1 saying that we shouldn't use it for them.

2 DR. WEISSMAN: I don't disagree, Jan, but
3 I think that needs to be shown and then the
4 indication added. That is what I would like to
5 see, if you can find enough aphaks to study. That
6 is other thing.

7 DR. SUGAR: Dr. Weiss?

8 DR. WEISS: Jayne Weiss. I would pose
9 this question to Dr. Saviola in terms of the lens
10 studies coming through here, what percentage have
11 indicated that the lens is used for, or can be used
12 for, aphakia when no aphakic patients have been
13 included in the study. If, as you are commenting,
14 most of the studies have not included aphakic
15 patients but have included aphakia as an
16 indication, then we shouldn't have any higher
17 requirements for this sponsor than anyone else.

18 DR. SAVIOLA: As you saw in the sponsor's
19 presentation, they are only making the lens in low-
20 plus powers. In the protocols that we have seen,
21 there have been a limited power range of people who
22 were enrolled. The historical perspective, from
23 our standpoint -- Dr. Lepri gave you some
24 information about how we look at power ranges and
25 permeabilities and things like that based on lens

1 thickness.

2 We have applied that, certainly in daily
3 wear, across the plus, minus-20 power range. If
4 you can do a +12.00 to +20.00, it should be for a
5 aphakic population. Even though it is really hard
6 to find new aphaks, there are certainly aphaks out
7 there who need different contact lenses. Other
8 than rigid lenses, those get harder and harder to
9 find in the soft-lens arena.

10 For this particular device, the thing to
11 consider is that they want to stay up to 30 days.
12 If you feel strongly as a panel that, in an aphakic
13 population, there are some different
14 considerations, you might be fine with this up to 7
15 days for aphakic wear but if, up to 30 days, you
16 have reservations, well then we should hear about
17 that.

18 Generally speaking, though, we apply the
19 permeability analysis and decide how high they can
20 go based on safe levels of oxygen.

21 DR. WEISS: So just as a follow up, does
22 the FDA have any concerns that a +10.00 lens or a
23 +15.00 lens would have any higher -- just because
24 of the lens makeup, have any higher risk than a
25 +5.00 lens.

1 DR. SAVIOLA: I would have to see their
2 thickness analysis and look at the difference in
3 permeability before I could answer that question.

4 DR. WEISS: It sounds like, from the
5 clinician's standpoint and the FDA standpoint, it
6 is a big question as to whether this is as safe in
7 aphakia.

8 DR. SAVIOLA: Again, in light of a 140 Dk
9 material compared to materials out there that are
10 in 18 and 28 Dks that are currently approved for
11 extended wear, they are going to have a pretty
12 thick lens in order to raise some transmissibility
13 concerns with us.

14 DR. SUGAR: Dr. Weissman and then I would
15 like to --

16 DR. WEISSMAN: I don't mean to monopolize
17 it, but, as a clinician, I would like to see this
18 lens available particularly for aphakic infants.
19 But I just think that possibly a different wording
20 at some point in the labeling might be appropriate.

21 DR. SUGAR: So is it correct that the
22 sense of the panel is that the word "aphakia"
23 should be removed from the first indication. Is
24 there agreement by head nodding or something? We
25 are not allowed to vote on this. Jayne?

1 DR. WEISS: Jayne Weiss. I am reluctant
2 to remove it but I would like to indicate that the
3 study that was performed was not performed on
4 aphakic patients. I wouldn't want to make it such
5 that a clinician could not use this for an aphakic
6 patient because of our stringent criteria in the
7 indications statement. I would like to give the
8 clinician some leeway at the same time as
9 indicating that we don't have any data.

10 DR. SUGAR: Ralph is bristling there.

11 DR. ROSENTHAL: Rosenthal. A clinician in
12 the practice of medicine can use an approved device
13 as they see fit if there is nothing in the -- you
14 know, unless there is something in the labeling
15 that warns them they better not use it. Then, of
16 course, even then they can still use it as the
17 practice of medicine.

18 So whether it is in the indications
19 statement or not does preclude whether or not a
20 physician or an eye-care practitioner can use a
21 lens in a certain population. So it is just not in
22 the indications statement. It is a regulatory
23 issue.

24 DR. WEISS: Jayne Weiss, again.
25 Clinicians from legal aspects may be a little bit

1 more reluctant when it is not included in the
2 indications.

3 DR. SUGAR: The sense that I am getting is
4 that we want to remove the word "aphakic" from this
5 but add to the labeling a statement that
6 information on the performance of this lens in
7 aphakia is not yet available, or something to that
8 effect; is that correct, or am I misstating it?

9 Go ahead, Jose.

10 DR. PULIDO: Jose Pulido. After hearing
11 this discussion, I feel, in my mind, at least, that
12 for me I would feel more comfortable leaving the
13 aphakic there and then, later on in the warnings
14 section, put, "the study did not involve patients
15 that were aphakic so the results in these patients
16 should be looked at very carefully," something to
17 that effect.

18 Does that satisfy Dr. Weissman and Dr.
19 Weiss?

20 DR. WEISS: Jayne Weiss. I would agree
21 with Jose's recommendation.

22 DR. SUGAR: And Dr. Grimmett? There are
23 enough nods that I think we can proceed. Are there
24 other modifications of the proposed wording for the
25 indications statement? One was, then, to eliminate

1 the fourth bullet which is Focus Night and Day
2 lenses may reduce dryness symptoms that are present
3 with regular hydrogen soft contact lenses.

4 That was what Dr. Jurkus proposed. The
5 issue is -- I think the issue is, one, does this
6 imply that the lens is indicated more for dry-eye
7 patients. The other is, is this just a statement
8 that this lens performed better than another lens.

9 In the earlier discussion, the issue came
10 up that we are talking about the Acuvue lens, not
11 all hydrogel soft contact lenses. So one
12 modification would be to make it specific. The
13 other would be to eliminate this and have it in the
14 discussion in the labeling. The other would be to
15 just eliminate this. I think those are the
16 options.

17 Would someone like to champion one of
18 those?

19 DR. EDRINGTON: This is Edrington. I
20 would recommend eliminating the statement.

21 DR. SUGAR: Are there those who feel
22 otherwise? Dr. Grimmer?

23 DR. GRIMMETT: I would agree with that
24 because it may be statistically significant, but I
25 don't think it is clinically relevant. So I agree.

1 DR. SUGAR: So the sense of the panel is
2 that we would eliminate that fourth bullet.

3 The other indication, the lenses may be
4 prescribed for daily wear or extended wear for up
5 to 30 nights of continuous wear as recommended by
6 the eye-care professional. I guess that that gets
7 discussed in our third question. Anything else on
8 the second question? The third question is really
9 still dealing with the indications; that is, does
10 the panel recommend that the 30-day statement be
11 included in only other sections of the product
12 labeling rather than the indications statement,
13 with the agency discussing the option of removing
14 the 30 days from the indication and putting it
15 elsewhere in the labeling, assuming that we require
16 a package insert to be presented to the patient
17 receiving the lens.

18 Dr. Grimmett?

19 DR. GRIMMETT: Dr. Grimmett. I just want
20 to point out that, at least by my review, I did see
21 the statement in the package insert, tab 8, part 7,
22 in the professional fitting guide, tab 9, part 7
23 and in the patient booklet, tab 10, part 7. So, as
24 it stands now, at least as per my review, I saw the
25 statement at least in four locations. So it seems

1 like it is all over the place.

2 DR. SUGAR: I think that the issue that
3 the agency has is with it being up front in the
4 indications statement and perhaps being a marketing
5 issue that this is marketed as indicated for 30-day
6 wear, for wear up to 30 days. This is for us to
7 discuss.

8 Karla?

9 DR. ZADNIK: Karla Zadnik. I guess I
10 would ask Dr. Saviola would the alternative be to
11 make the second bullet say, "the lenses may be
12 prescribed for daily or extended wear as
13 recommended by the eye-care practitioner," because
14 then they could be 60 days, or 120 days, or years?

15 DR. SAVIOLA: That is the other side of
16 the coin; yes. It is an extended-wear lens and the
17 doctor decides. In the other parts of the
18 labeling, as Dr. Grimmett said, the 30-day wear
19 period still remains in those sections of the
20 labeling. It is just not in the indications
21 section.

22 DR. ZADNIK: Or would you recommend an
23 alternative that said, "for up to 7 nights of
24 continuous wear?" I mean, I am trying to get sort
25 of if you reject this, what is the alternative

1 option.

2 DR. ROSENTHAL: Excuse me. Really,
3 though, that is for you. If you reject this, that
4 is for the panel to recommend. Please.

5 DR. SAVIOLA: I am not going to recommend
6 to you what you should do.

7 DR. ZADNIK: Okay.

8 DR. SAVIOLA: But, again, the reason this
9 came up and the time period in the indications, if
10 you say to us, "we don't think it should be an
11 indication but it should be as it stands now in the
12 other parts of the labeling," well, yes; it still
13 could be promoted and sold and whatever as a 30-day
14 lens. But still, technically it is not indicated
15 for that.

16 We had a rate number before for the
17 hydrogels based on epidemiological data that got
18 published. And we said, "That is too high a rate
19 number. We can't live with that rate number, 30."
20 Everything went back to 7 days.

21 Okay; we don't have the rate number now so
22 that is part of the problem with saying, "Yes; we
23 can go with 30," because we have a missing piece of
24 the puzzle which you won't have until later on.

25 Again, it is up to you to discuss the pros

1 and cons of the idea that someone would run it up
2 to 2 or 3 months certainly is an issue. If you
3 feel like that is a strong enough issue to say,
4 "You guys work out the regulatory details. We
5 think it should be 30 days," well, then, say that
6 to us.

7 DR. SUGAR: This is presented -- we are
8 reviewing this a 30-day lens. To eliminate it from
9 the indications, I think, is game-playing and
10 really eliminates the basic issue. So I,
11 personally -- am I not supposed to say "I,
12 personally," anything?

13 I am supposed to vote in ties, but I
14 personally think that we ought to leave it in.

15 Go ahead, Mike.

16 DR. GRIMMETT: Mike Grimmett. My belief
17 would be to leave it in the indications statement
18 as well as in the other sections that I already
19 mentioned, the reason being that, for all practical
20 aspects, the manufacturer-sponsor would still
21 advertise it up to 30 days. I don't see the
22 difference in practical terms to the clinician if
23 you somehow hide it out of the indications
24 statement.

25 The sponsor did do a study up to 30 days.

1 So I would leave it in the indications statement.

2 DR. SUGAR: Is there anyone that feels
3 otherwise? Do you want to comment, Janice?

4 DR. JURKUS: Just a possible
5 consideration. We could change that, instead of up
6 to 30 days to include use of the lens from 1 to 30
7 days. That way, people would not get the idea that
8 you have to use it for 30 days. It can be used 1,
9 2, 3, 4, 5, you know, any number of days at a time.

10 DR. SUGAR: I am not sure that "up to"
11 says anything different than that.

12 Jayne?

13 DR. WEISS: Jayne Weiss. It is almost as
14 if there is an elephant in the room and we just
15 want to ignore the fact that the elephant is
16 sitting next to us.

17 DR. SUGAR: That is not Mike that you are
18 talking about.

19 DR. GRIMMETT: I have lost some weight
20 recently.

21 DR. ROSENTHAL: So have I.

22 DR. WEISS: Diplomacy has never been my
23 strong suit, as you can tell.

24 The sponsor did an excellent study to show
25 that this can be used in many patients successfully

1 up to 30 days. Let's give them credit for that.
2 Let's put it in the indications. If postmarket
3 shows something else, then we will change it.

4 DR. SUGAR: Other comments? We will move
5 on to Question 4; does the panel have any specific
6 recommendations for the proposed product labeling
7 in terms of warnings, precautions, clinical data,
8 outcomes or practitioner-directed or patient-
9 directed labeling?

10 Dr. Saviola specifically pointed out the
11 option of patient-directed labeling which -- I need
12 to understand this. People, when they get contact
13 lenses now, do not have a package insert with the
14 lens? Is that correct?

15 DR. ROSENTHAL: This is Rosenthal. Have
16 you seen the advertisements that are published in
17 papers?

18 DR. SUGAR: Sure; on the back of the page,
19 they have listed the --

20 DR. ROSENTHAL: Ad infinitum, all the
21 issues. I think that is a patient-directed
22 advertisement. So, up to now, there has not been
23 that type of requirement.

24 DR. SUGAR: For contact lenses.

25 DR. ROSENTHAL: For contact lenses.

1 DR. SUGAR: But there is a requirement
2 that a package insert be given or not?

3 DR. ROSENTHAL: Oh, yes; of course a
4 package insert is included in all --

5 DR. YAROSS: A package insert and
6 advertising issues are really quite distinct. One
7 falls under the restricted device regulation and
8 that is distinct from labeling that is disseminated
9 to the practitioner to then distribute to the
10 patient.

11 I guess the question is does any other
12 contact lens at this time have patient brochures
13 that are provided to the practitioner to provide to
14 the patient.

15 DR. ROSENTHAL: There is, apparently, a
16 patient brochure required -- not required --

17 DR. SUGAR: But is not in the package.
18 When you open the box of your Acuvue lens, it is
19 not there.

20 DR. ROSENTHAL: It is not in the package.
21 It is required, but I understand for many years, it
22 is sort of made up but no one ever uses it.

23 DR. YAROSS: Sponsors do have no control
24 over what the practitioners do in that respect.

25 DR. ROSENTHAL: That's correct.

1 DR. SUGAR: But if it was required to be
2 in the box, this would be dispensed, I presume, in
3 clusters of six lenses or whatever. Then we could
4 request that.

5 DR. SAVIOLA: Let me just, again, recap.
6 For current lenses, as the sponsor described
7 before, there is labeling guidance out there. It
8 talks about a package insert that is directed to
9 the practitioner with information to review with
10 the patient. There is a practitioner fitting guide
11 and there is also a patient information booklet.

12 Those are all elements of labeling that we
13 review as part of the approval for either daily or
14 extended-wear lenses.

15 The words "patient-directed labeling" that
16 we put into the question bring up to idea, do you
17 think there should be something else besides those
18 three elements currently, such as a patient-
19 directed package insert.

20 Whether or not that patient-directed
21 package insert gets printed on the back of an
22 advertisement is a restriction issue and that is
23 something that we don't need to discuss within the
24 context of panel because we make the decision
25 whether or not we are going to restrict it under

1 502(q) and (r). So don't get the two confused. We
2 are basically saying, do you think -- and, of
3 course, we have contact-lens consultants and panel
4 members so they are quite familiar with the
5 labeling that is out there, hopefully, that you
6 have seen it.

7 Do you think, in the context of this
8 thing, a new class of lens device, that there
9 should be something else besides what is currently
10 out there.

11 DR. SUGAR: We are going to have to break
12 up this discussion into how we are going to change
13 the labeling and then what we are going to do with
14 the labeling.

15 Go ahead, Dr. McMahon.

16 DR. McMAHON: Tim McMahon. Currently, we
17 have extended-wear lenses that are approved out
18 there already for 7 days. We don't require
19 corporations to provide this patient-specific
20 instructions, if you will.

21 This proposal is being held, basically, to
22 a comparison to those already approved lens
23 designs. I don't feel that we need to add an
24 additional burden to them after we have already
25 pretty much come to a pretty close consensus that

1 they have adhered to that.

2 So I don't see the rationale, despite the
3 fact that this is a new class of lenses at this
4 point, to add that patient-specific label.

5 DR. SUGAR: Dr. Bandeen-Roche?

6 DR. BANDEEN-ROCHE: I may be
7 misunderstanding, but I have to respectfully
8 disagree. I do feel that there should be patient-
9 directed labeling. Again, maybe I am just
10 misunderstanding semantics, but I think the patient
11 absolutely must receive certain information given
12 the history of continuous-wear lenses, that there
13 has been an unfortunate history with them.

14 We have a promising new product before us,
15 but I do think that patients absolutely need to see
16 the data in some understandable form that we have
17 seen here today and they also need to understand
18 that the ultimate endpoint has yet to be evaluated.
19 At least, that is my opinion.

20 DR. SUGAR: Dr. Pulido?

21 DR. PULIDO: I think that intraocular
22 lenses are extended wear also. Patients don't
23 receive -- do they receive? They don't receive a
24 patient --

25 DR. SUGAR: The physician, with the

1 package with the intraocular lens, receives the
2 intraocular lens and the package insert that
3 details the indications, the contraindications.

4 DR. PULIDO: But nothing for the patient.
5 It is the doctors --

6 DR. SUGAR: The doctor chooses to either
7 give it to the patient or not give it to the
8 patient.

9 DR. PULIDO: So it is the doctor's
10 responsibility to let the patient know of the
11 adverse events, et cetera. I think that there is
12 nothing more extended wear than an intraocular lens
13 and there is nothing special done in that
14 situation.

15 DR. SUGAR: Dr. Edrington?

16 DR. EDRINGTON: Edrington. Intraocular
17 lenses, I assume there is informed consent where
18 the majority of extended-wear patients are not
19 having informed consent.

20 DR. SUGAR: Dr. Zadnik?

21 DR. ZADNIK: Karla Zadnik. I really think
22 this patient education has to happen when the
23 patient initially starts with these lenses. That,
24 almost by definition, has to come from the doctor
25 as it should with all other extended-wear products

1 that are already approved.

2 Remember, these patients are going to
3 start getting six packs or two packs or something
4 of these and they are going to have that same
5 monster patient-directed labeling inside every one
6 for them to pitch over their shoulder because they
7 are on their third year of wear and they are on
8 their sixth box of lenses.

9 So I think for the education to be
10 meaningful and for the communication to do what it
11 is supposed to do, it really has to come from the
12 practitioner and be directed in that way rather
13 than directed at the patient each and every time.

14 DR. SUGAR: Although we cannot control
15 that.

16 DR. ZADNIK: Of course not.

17 DR. SUGAR: Dr. Bandeen-Roche, and then
18 Dr. Jurkus.

19 DR. BANDEEN-ROCHE: That would just be my
20 question. That sounds absolutely right but how do
21 we insure that it happens. I guess the answer is
22 that we can't.

23 DR. JURKUS: Jan Jurkus. It seems that
24 the thing that we are looking for is to try to
25 protect the patient, to let the patient know that

at

193

1 not everything is absolutely perfect with this type
2 of lens or could not be absolutely perfect.

3 So something to consider might be, on the
4 package, on the box that it comes with, giving the
5 patient directions on what to do if they have any
6 signs of irritation, redness, change in vision or
7 lens discomfort, "Take out the lens and see your
8 practitioner."

9 If we put that on every box that the
10 patient gets, hopefully, it would prevent the more
11 serious complication from happening and it is
12 something that the patient might actually read
13 instead of little tiny pieces of paper that they
14 would throw out.

15 DR. SUGAR: Dr. Pulido.

16 DR. PULIDO: Jose Pulido. On page 1247,
17 it says, the first thing after it says package
18 insert, Focus Night and Day extended-wear soft
19 contact lenses, prescription only. So, whenever
20 you give a prescription to a patient, you always
21 tell them -- I mean, you are legally bound to tell
22 them the risks and benefits.

23 Whenever I prescribe Timoptic, I don't
24 give them the package insert of the risks and
25 benefits of the Timoptic. It is my duty to have

1 already told them that. So we are making this much
2 different than we are doing anything else that
3 comes by prescription only.

4 DR. SUGAR: Dr. Weiss?

5 DR. WEISS: Jayne Weiss. Just to play a
6 little bit of devil's advocate, I agree with Dr.
7 Pulido that we shouldn't make the rules and
8 regulations for this any more stringent than
9 anything else we do and that would be unfair. But,
10 in my own practice, I see more and more patients
11 ordering lenses by phone or mail and maybe never
12 even interacting with an eye-care professional.

13 In that sense, is this lens going to be a
14 higher risk than the other lens and maybe we should
15 have a different set of criteria. I don't have an
16 answer for that. I am just throwing that out.

17 DR. ZADNIK: But they are not ordering
18 this different modality of lenses the first time
19 they get them from Lens Express and Linda Carter;
20 right? They can't get these the first time, I am
21 assuming.

22 DR. PULIDO: It is prescription only.

23 DR. ZADNIK: It is by prescription so they
24 have got to have a prescription for this type of
25 lens initially. That is where the education from

at

1 the doc comes in and then only when they get their
2 replacement lenses through and 800 number would
3 they not be receiving this additional labeling.

4 But I think Jan's idea of something that
5 really gets the message across; "if your eyes hurt
6 or you can't see, take your lens out and call your
7 eye doctor." That is really what you want the
8 message to be. Why not have the sponsor think
9 about delivering it in a way that the patient might
10 actually see it.

11 DR. SUGAR: Is that within our purview to
12 suggest that to the agency? Ralph?

13 DR. ROSENTHAL: The panel may suggest
14 whatever they like, Dr. Sugar.

15 DR. SUGAR: I think we are suggesting
16 that. Is that the sense that there is not strong
17 support for patient-directed labeling but there is
18 support for a warning label on the package that
19 states, "this is prescription only." It already
20 says that on the package, I think, for all lenses,
21 for all soft lenses, and that in the event of pain,
22 redness, discharge, you should seek attention from
23 an eye-care practitioner.

24 Dr. McMahon?

25 DR. McMAHON: I think it is a good idea

at

196

1 but I have some qualms about specifically directing
2 this particular sponsor with this particular
3 product only doing this. This should be something
4 that goes across the spectrum which I don't know is
5 within the purview of this panel.

6 I have this sort of grumbling feeling that
7 this is being somewhat unfair.

8 [Many panel members in agreement.]

9 DR. JURKUS: Jan Jurkus. I happen to
10 disagree. We are looking at a totally new modality
11 that could possibly have more patient
12 noncompliance. We are giving our stamp of approval
13 to something that is different than what is already
14 out there. So I would think that, from whatever we
15 have from this point forward, if it falls into the
16 same grouping, should have pretty much the same
17 requirements.

18 But this is different than anything else
19 that we currently have on the market.

20 DR. SUGAR: What I guess I would like to
21 do is straw poll the panel in terms of how many
22 feel that there should be special specific package
23 labeling as we just described; that is, the warning
24 labeling on the package. All those who would like
25 to suggest that that be done, signify by raising

1 your hand.

2 [Show of hands.]

3 DR. SUGAR: There are four. Again, I am
4 not counting and we are not voting. Those who
5 would not like to see that done?

6 [Show of hands.]

7 DR. SUGAR: Six. So you have a sense that
8 that issue was raised but there is not overwhelming
9 support for it.

10 Go ahead, Dr. Edrington.

11 DR. EDRINGTON: I would just like to add,
12 what Karla said about the way patients access
13 lenses these days. We are starting, in a sense, a
14 new modality. In the '80's, when we ran into
15 extended-wear problems, the practitioner was
16 primarily the one delivering the lenses to
17 patients.

18 So, I think when you look at those two
19 things together, I think the extra labeling is not
20 a bad idea.

21 DR. SUGAR: Dr. McMahon, we are going to
22 move from this into the specific wording of the
23 labeling.

24 DR. McMAHON: Tim McMahon. That straw
25 vote should be recognized as that is pertaining

1 just to this. I think it would be interesting to
2 know how the panel feels about adding that
3 particular type of warning to all lenses.

4 DR. SUGAR: Because airplanes leave at
5 7:00 tonight and earlier for some people, I would
6 like to leave it with the issue at hand and not get
7 more global. But if someone wants to usurp my
8 ability to do that, go ahead.

9 I would like to now ask for suggestions
10 for specific changes to the labeling. The labeling
11 is in the back. I have three sections, the
12 labeling for the physician, for the patient and I
13 don't know what the third one is.

14 Things that were brought up include --

15 DR. GRIMMETT: Dr. Sugar, I wrote them all
16 down as each doctor made recommendations. So, even
17 though they are not voted on yet, I did keep a
18 record of that.

19 DR. SUGAR: Why don't you just go through
20 those and then we can discuss them one by one.

21 DR. GRIMMETT: Mike Grimmatt. Jan Jurkus
22 initially suggested that a looser fit may not be
23 best. I hope I summarized that correctly. Did
24 anyone else have concerns in the labeling that we
25 should address the issue of making this particular

1 lens tighter in select patients?

2 DR. SUGAR: Why don't we deal with these
3 one at a time, then. So, Janice is suggesting that
4 the specific fitting recommendations be made in the
5 practitioner labeling that discuss the issue of the
6 fact that this lens may be appropriately fit
7 differently than standard lenses.

8 DR. GRIMMETT: Due to the discomfort, I
9 assume, that was experienced with this lens.

10 DR. JURKUS: Right.

11 DR. SUGAR: Is there a sense that this
12 would be an appropriate addition to the labeling?

13 DR. EDRINGTON: Dr. Edrington. I would
14 probably stay away from that, telling the
15 practitioner, in a sense, to go tighter on the fit.
16 I think they will, either by word of mouth or at
17 meetings, whatever, determine that or by patients'
18 symptoms of discomfort.

19 I think if you erred so that they were
20 fitting them too tight and went to the steeper base
21 curve as their default system that, perhaps, other
22 complications that we are not currently aware of
23 could occur.

24 DR. SUGAR: Other comments on that? Dr.
25 McMahon?

1 DR. McMAHON: Tim McMahon. This trial was
2 done with one base curve and all the patients are
3 basically forced into one particular lens shape.
4 The sponsor has recognized that they need a
5 different base curve and the data that we are
6 presented did not give us the opportunity to
7 determine whether an alternative base curve would
8 have had some influence on those things.

9 So I agree with Dr. Edrington, we should
10 probably stay away from it.

11 DR. SUGAR: There are statements in the
12 practitioner guide that talk about how you measure
13 the fit and the push-up test and that kind of
14 thing. So you are suggesting that something
15 different be added.

16 DR. JURKUS: Jan Jurkus, again. What I am
17 suggesting is that a statement regarding the
18 patient response of, "The lens doesn't feel right,"
19 or some discomfort be, in some way, highlighted
20 because of the fact that the lens may look perfect
21 on the eye where the sponsor said that it may all
22 look okay but if the patient says that it is
23 initially not comfortable, that you might want to
24 go to the steeper design wear. Patient
25 symptomatology plays an important part in the