

DEPARTMENT OF HEALTH AND HUMAN SERVICES  
FOOD AND DRUG ADMINISTRATION  
CENTER FOR DEVICES AND RADIOLOGICAL HEALTH

CIRCULATORY SYSTEM DEVICES PANEL

Wednesday, October 23, 2001

8:39 a.m.

Walker/Whetstone Room  
Gaithersburg Holiday Inn  
2 Montgomery Avenue  
Gaithersburg, Maryland

## C O N T E N T S

AGENDA ITEM:	PAGE
Call to Order	
Dr. Cynthia M. Tracy, Chairperson	3
Presentation: Dr. Dale R. Tavris, Office of Surveillance and Biometrics, FDA "Adverse Events and Deaths Associated with Hemostasis Following Cardiac Catheterization: Comparison of Manual Compression versus Collagen Plug and Suture Hemostasis Devices"	3
Conflict of Interest Statement	23
Open Public Session	28
Sponsor Presentation	
Ms. Jean Chang, EMBOL-X, Mountain View, California	29
Dr. Nicholas T. Kouchoukos, Study Co- Principal Investigator	32
Dr. Richard E. Kuntz, Study Trial Design	41
Dr. Keith B. Allen, Site Principal Investigator	54
Questions and Answers	82
FDA Presentation	
Catherine Wentz - Engineering Review	96
Dr. Julie Swain - Clinical Review	100
Dr. Gerry Gray - Statistical Review	110
Questions and Answers	123
Open Committee Discussion	134
Open Committee Discussion - Continued	157
Questions for Panel	212
Adjournment	268

## P R O C E E D I N G S

## Call to Order

DR. TRACY: Good morning. I'd like to call to order this meeting of the Circulatory System Devices Panel.

We are going to start with a presentation from the Office of Surveillance and Biometrics on "Adverse Events and Deaths Associated with Hemostasis Following Cardiac Catheterization: Comparison of Manual Compression versus Collagen Plug and Suture Hemostasis Devices."

If the presenter is present, would he please take the podium?

Presentation by Office of Surveillance  
and Biometrics

DR. TAVRIS: Thank you.

This study was a collaborative effort between the FDA and the American College of Cardiology.

[Slide.]

Over the 5-year period between 1996 and 2000, 1,880 reports of serious injuries and 36 reports of deaths associated with the use of hemostasis devices used to prevent bleeding from the femoral artery following cardiac

1 catheterization were reported to the FDA through  
2 its medical device reporting system.

3 Most of these serious injuries and deaths  
4 involved hemorrhagic complications. It was also of  
5 interest that a large majority of injury and death  
6 reports involved women, even though cardiac  
7 catheterization is more common in men than in  
8 women.

9 This study considered the earliest two  
10 types of hemostasis devices--the suture device,  
11 Perclose, and the two collagen plug devices,  
12 VasoSeal and AngioSeal.

13 Because of the continued receipt of  
14 adverse event reports involving injuries and  
15 deaths, the FDA was concerned about the safety of  
16 these devices. Of course, these reports themselves  
17 do not causally implicate these devices in the  
18 injuries and deaths, since these can also occur  
19 following manual compression, in which case the  
20 events would not be reported to the FDA.

21 But the medical literature on this subject  
22 also gave us cause for concern. Of 13 studies that  
23 we could find that utilized manual compression  
24 control groups to assess risk of serious injury  
25 associated with hemostasis device use, 9 showed no

1 difference in the rate of adverse events between  
2 device users and controls, and 4 demonstrated a  
3 higher rate in device users. None demonstrated a  
4 higher rate in controls than in device users.

5 Two important weaknesses of the studies  
6 found in the medical literature are small sample  
7 sizes and the use of a single or small number of  
8 institutions. For example, the 13 studies referred  
9 to above utilized a total of 19,582 procedures,  
10 including a little over 15,000 controls and a  
11 little over 4,000 device users. Most involved a  
12 single institution.

13 By contrast, this study, which utilized  
14 the American College of Cardiology's National  
15 Cardiovascular Data Registry, involved 214  
16 participating institutions and 166,680 procedures,  
17 including over 113,000 controls and over 53,000  
18 device users--more than 7 times as many controls  
19 and 13 times as many device users than all the  
20 other studies combined.

21 This study included information from all  
22 cardio catheterization lab admissions representing  
23 the 214 institutions included in the ACC's data  
24 registry from the year 2001. Excluded from the  
25 analysis were outpatients and any patient for whom

1 critical data was not available.

2 Outcomes assessed in these analyses  
3 included hemorrhage, arterial occlusion and loss of  
4 arterial pulses, artery dissection, the development  
5 of an AV fistula or pseudo-aneurysm, and death  
6 associated with any of these events. Hemorrhage,  
7 by far the most common of these events, was defined  
8 as "blood loss requiring transfusion or prolonging  
9 the hospital stay or causing a drop in hemoglobin  
10 of greater than 3."

11 Stepwise backward multiple logistic  
12 regression analysis was performed in order to  
13 control for the effects of potential confounding  
14 variables. The main independent variable of  
15 interest was hemostasis device use.

16 Potential confounding variables that were  
17 assessed included demographic variables, comorbid  
18 conditions, type of procedure--that is  
19 interventional versus diagnostic cardiac  
20 catheterization--presence of left main coronary  
21 artery stenosis, and indications for the procedure.

22 Of the more than 160,000 subjects in this  
23 analysis, by far the most frequent complication was  
24 hemorrhage. There were 1,756 episodes of this, a  
25 rate of 1.1 percent. This represented 73 percent

1 of the subjects who were characterized by any  
2 adverse event in this study.

3 The mortality rate for the 2,418 subjects  
4 in this study reported within the adverse event was  
5 5.5 percent.

6 First, I will present only the  
7 multivariate results that used the reporting of any  
8 adverse event as the outcome. Female gender and  
9 the use of interventional as compared to diagnostic  
10 cardiac catheters were found to be the biggest risk  
11 factors, both with odds ratios of 2.3.

12 Several comorbid conditions and  
13 indications for the procedure were also found to be  
14 associated with adverse events. These included an  
15 emergency indication for the procedure, plaque as  
16 an indication for the procedure, acute myocardial  
17 infarction, history of renal failure, New York  
18 Heart Association class, and peripheral vascular  
19 disease, all with odds ratios between 1.2 and 1.8.

20 The use of a hemostasis device was found  
21 to be protective as compared with the use of no  
22 device, especially the use of the collagen plug  
23 devices, which demonstrated an odds ratio of 0.79,  
24 which was highly statistically significant.

25 [Slide.]

1           This table depicts odds ratios for those  
2 associations that were statistically significant at  
3 the P less than .05 level. To a large extent, the  
4 risks associated with any complication pertained to  
5 many of the specific complications as well.

6           Female gender was associated with five of  
7 the seven specific complications, all with odds  
8 ratios of greater than 2 and P values of .0002 or  
9 less.

10           The use of interventional cardiac  
11 catheterization was statistically associated with  
12 four of the seven specific complications.

13           Of the nine comorbid conditions or  
14 indications for the procedure that were assessed in  
15 this analysis, six of them were associated with  
16 three to five of the seven specific outcomes in the  
17 multivariate analysis. Those that weren't  
18 associated with these outcomes were probably  
19 precluded from this by their close association with  
20 the other comorbid conditions.

21           [Slide.]

22           As for the hemostasis devices, the most  
23 pronounced protective effect pertained to  
24 pseudo-aneurysms. Both types of hemostasis devices  
25 were characterized by odds ratios of approximately



1 one-half with respect to this outcome.

2 The hemostasis devices as a group  
3 demonstrated a protective effect with regard to  
4 hemorrhage, with an odds ratio of 0.89, although  
5 this protective effect was not statistically  
6 significant.

7 When the collagen plug devices alone were  
8 compared with manual compression, the odds ratio  
9 was 0.85 with a P value of .035.

10 Neither hemostasis device demonstrated a  
11 statistically significant protective effect with  
12 regard to vascular complication-related deaths.  
13 The collagen plug devices demonstrated an odds  
14 ratio of 0.56 with respect to this outcome, but  
15 because there was only a total of 144 of these  
16 deaths, the study was not powerful enough to attain  
17 statistical significance even with that low odds  
18 ratio.

19 The risks associated with female gender,  
20 interventional cardiac cath, and several of the  
21 comorbid conditions were not surprising, as these  
22 had been demonstrated previously. But we were  
23 surprised to note the protective effect of  
24 hemostasis devices given that concern over their  
25 safety was the main reason for conducting the study

1 and that the medical literature had generally shown  
2 them to be associated either with a greater than or  
3 equal risk of adverse events compared to manual  
4 compression controls.

5 Possible explanations for this apparent  
6 discrepancy are the following. First, confounding  
7 variables that were not controlled for in this  
8 analysis. However, this explanation seems unlikely  
9 to us given the large range of comorbid conditions  
10 that were controlled for in this study, although we  
11 did not control for coagulation status.

12 Second, it could be that the medical  
13 providers who participated in this study were more  
14 skilled in the use of these devices than most other  
15 medical providers. This is a complicated issue.  
16 On the one hand, that explanation is made less  
17 convincing by the large number of participating  
18 institutions in this study; but on the other hand,  
19 the protective effect seen in this study was not  
20 very great--only a 21 percent decrease in total  
21 complication rate for the collagen plug devices and  
22 a 15 percent decrease for Perclose.

23 None of the studies that we found in the  
24 medical literature was large enough to detect that  
25 small of a protective effect, even if one assumes

1 that the skill of the providers was comparable to  
2 that found in this study. But that would still not  
3 explain why a minority of studies in the medical  
4 literature have demonstrated a harmful effect  
5 associated with hemostasis device use.

6 Of course, since these studies generally  
7 involved a single institution, it is possible that  
8 the physicians involved in one or more of these  
9 studies were less experienced or skilled in the use  
10 of the hemostasis devices than average.

11 Finally, a third possible explanation for  
12 the protective effect found in this study compared  
13 with other studies is that over time, the users of  
14 these devices have become more skilled and thus  
15 more likely to produce better results than those  
16 seen in other studies.

17 That concludes my presentation, and I  
18 would be happy to answer questions.

19 DR. TRACY: Dr. Krucoff?

20 DR. KRUCOFF: Dale, with great respect for  
21 the noble intention here, I am really concerned  
22 that maybe we're just pouring worms into a can of  
23 worms.

24 In a nonrandomized platform there is an  
25 enormous bias involved in how you and when you

1 choose as an operator to use these things. Most of  
2 us will actually do a femoral injection before you  
3 position one of these, so anybody, where you are  
4 involved in a plaque or at a bifurcation point or  
5 where you have had multiple sticks in a lesion, you  
6 don't even deploy these things.

7 So immediately in this kind of registry,  
8 there is an intrinsic bias just by case selection.  
9 And I would certainly list that amongst your  
10 possible explanations for the findings that you are  
11 looking at.

12 The other question I have is could you  
13 detail for us what you are aware of as far as any  
14 kind of quality control at the sites? Was the  
15 operator or somebody related to the operator who  
16 placed the device also the one who followed up on  
17 the patient and reported on complications, which of  
18 course is another source of bias--if you put one of  
19 these things in, you may tend to look at a little  
20 oozing or bleeding as just something that is going  
21 on in the tract because the patient is  
22 anticoagulated, whereas if you with manual  
23 compression have subsequent bleeding, you may  
24 report it differently.

25 How as the data quality-controlled at the

1 site level?

2 DR. TAVRIS: There were a number of  
3 quality control procedure that mostly included  
4 educational interventions to train the sites on how  
5 to record the data, and also overview of the data  
6 for completeness.

7 As far as potential bias in the recording  
8 of data, I'm not sure that any of the quality  
9 controls could have favorably influenced that.

10 The earlier comment about potential  
11 selection bias sets--part of what I meant when I  
12 talked about potential confounding variables that  
13 could have affected this--do you feel that the  
14 selection of patients would work in a way that  
15 would make the patients who received the devices a  
16 group that was less likely to experience  
17 complications?

18 DR. KRUCOFF: Definitely, absolutely,  
19 positively.

20 Lastly, you called this a "study." Did  
21 patients whose data were recorded in this registry  
22 provide informed consent?

23 DR. TAVRIS: I'm pretty sure they did.  
24 This was the American College of Cardiology's  
25 registry, and I'm pretty sure they did.

1 DR. KRUCOFF: We do not consent patients  
2 for registry data for the ACC.

3 DR. TAVRIS: Okay.

4 DR. TRACY: Dr. White?

5 DR. WHITE: Thank you.

6 I'd like to echo Dr. Krucoff's statement.

7 As a user of these devices, I think you cannot  
8 underestimate the selection bias that goes into  
9 this compared to a randomized device trial. The  
10 published literature--our hands are tied about how  
11 patients are treated. You absolutely enroll the  
12 patients, they get one or the other, and you work  
13 very hard to make the patient fit the trial.

14 In our regular practice, however, we don't  
15 poke skunks with sticks. If that groin doesn't  
16 look good, it doesn't get a device, and we don't  
17 look for that trouble, and I think that is a major  
18 reason to explain why it appears that devices are  
19 safer.

20 The other thing is I think you may not  
21 completely understand the ACC Data Registry, which  
22 is a very voluntary, self-selected population. It  
23 certainly isn't a widely-adopted process. It  
24 doesn't even represent a majority of the  
25 catheterization laboratories in the U.S.--not

1 because it's not a good thing, it's just that it is  
2 expensive. For example, in my institution, I would  
3 love to participate, and my institution cannot  
4 afford to participate. So there is a lot of data  
5 missed, I think, if that is the database.

6 Then, finally, I think that that database  
7 is used for quality control, not for scientific  
8 generally reporting, so I think patients generally  
9 are not consented that their data would be  
10 collected and used. And I'm not sure what the  
11 implications are for that in terms of informed  
12 consent and the use of the data. That would be  
13 something that would be of concern, I think, to my  
14 local IRB that release of that kind of information.

15 DR. TRACY: Dr. Laskey?

16 DR. LASKEY: I would be very wary now and  
17 for the foreseeable future in using NCDR data. It  
18 is non-quality-assured. It is non-verified. There  
19 is no routine auditing as far as I can tell, as far  
20 as I know. And as Chris mentioned, it is an  
21 entirely voluntary registry which, maybe on one  
22 side of the coin bespeaks honesty on the part of  
23 the reporting sites, but it is extremely spotty,  
24 and--

25 DR. TAVRIS: When you say "voluntary,"

1 doesn't that mean that it is voluntary as far as  
2 the institution, but the institutions that are  
3 participating in it do routinely collect data from  
4 all procedures?

5 DR. LASKEY: Yes, but the rigor with which  
6 that data is collected cannot be vouched for, and  
7 it will vary from site to site.

8 DR. WHITE: And the discipline should not  
9 be confused with the discipline for a randomized  
10 trial, which has audits, and you're pretty sure  
11 about that.

12 DR. LASKEY: It has become very  
13 fashionable in the last few years to do things like  
14 propensity scores and all that to try to adjust  
15 away for some of these biases in these  
16 observational trials. I would bet if you did that,  
17 you would still be left with the same answer, but  
18 it is probably worthwhile going through the  
19 exercise since you have a number of people who  
20 didn't get the device.

21 It has been our experience that it goes to  
22 more than the groin. If the groin or the artery  
23 doesn't look right, or if the patient doesn't look  
24 right at the end of the procedure, they do not get  
25 a closure device, either.



1 DR. TAVRIS: What data would you use to  
2 calculate the propensity scores?

3 DR. LASKEY: You have folks who got the  
4 device and folks who didn't get the device. That  
5 is your endpoint for the propensity score, and then  
6 you would put into the soup all the other variables  
7 that you just looked at as step one. And then,  
8 step two is to put the propensity scores into your  
9 final analysis. That is generally the way it is  
10 done. I am not supporting that as a way of  
11 verifying this data, but I think it is an  
12 interesting exercise.

13 DR. TRACY: You mentioned that there was  
14 an exclusion of outpatients. Does that mean  
15 outpatient caths, and procedures and interventions  
16 were excluded from analysis--because I think that  
17 in many institutions, most procedure are done as  
18 outpatient.

19 DR. TAVRIS: In this database, most of  
20 them were not outpatients, but yes, the data that I  
21 showed excluded outpatients although before we  
22 excluded them, we did analyze the data with  
23 outpatients in it, and we got very similar results.

24 The reason we excluded the outpatient  
25 afterward was that we felt that there might be some

1 potential for additional bias because there might  
2 have been disproportionate followup in one group or  
3 the other. We couldn't be sure how good the  
4 followup was.

5 DR. LASKEY: Dale, if available, just one  
6 other key variable which may not be in the database  
7 but that has clearly been associated with  
8 hemorrhagic complications is the extent of  
9 anticoagulation; what kinds of ACTs are in this  
10 patient population in a nonrandomized format.

11 Again, that may actually affect whether or  
12 not an operator would deploy a device, so you may  
13 have your higher ACTs on IIb/IIIa's in the wrong  
14 group.

15 DR. TAVRIS: We wanted to do that, but  
16 that data wasn't available. We do intend to do  
17 that in our next study. But from review of the  
18 medical literature, there is only one article, I  
19 believe, that did control for coagulation status,  
20 and in that article, those with hemostasis device  
21 use were more anticoagulated than those without  
22 use, and I would think that that would tend to make  
23 them look worse given that hemorrhage was the main  
24 complication here, and that group would have been  
25 more likely to experience higher hemorrhage rates.

1           So that would have contributed bias, I  
2 would think, in a direction that would have made  
3 our results, the protective effects seen for the  
4 hemostasis devices even less likely.

5           DR. TRACY: Dr. Zuckerman?

6           DR. ZUCKERMAN: Just by way of background,  
7 the Agency is always interested in ways of  
8 potentially addressing the pre- and post-market  
9 balance in terms of what we require pre-market  
10 before PMA approval of a hemostasis device versus  
11 post-market. And certainly this has been one  
12 attempt at looking at post-market datasets given  
13 the controversial nature of some of these devices.

14          Panel members have expressed problems with  
15 this particular registry, so the first question is  
16 are there any other datasets that might be useful  
17 to explore.

18          The second question refers to some of the  
19 implications of the panel discussion for our  
20 pre-market approval data requirements for these PMA  
21 hemostasis devices. Because there are large  
22 opportunities for selection bias, et cetera, our  
23 general standard has been to require a randomized  
24 trial versus manual compression.

25          Many sponsors have indicated that there is

1 a large historical database of manual compression  
2 results and have suggested other trial designs. I  
3 would be interested in any comments on how to  
4 evaluate these devices and with minimization of  
5 bias.

6 DR. KRUCOFF: I think one key question  
7 would be whether your issue is safety or efficacy.  
8 When you started the presentation, it seemed to me  
9 that the concern, because of the reporting  
10 mechanism that kicked all this in, was safety. And  
11 ultimately what your conclusions are leaning toward  
12 is are you demonstrating some kind of efficacy  
13 impact.

14 I would at least start by being clear on  
15 what the question that is being addressed is, and  
16 if complication rates associated with these devices  
17 are what you need to learn about, then I think you  
18 need to make sure that the data is collected in a  
19 way that you can understand whether the  
20 complication rates are higher than your target.

21 Now, in terms of manual compression, which  
22 I think has got to be the target, whether you could  
23 do a comprehensive job of characterizing an  
24 historic control where you understood that it was  
25 matched across important parameters to a study

1 population where that device is deployed, I don't  
2 think an historic control would be out of the  
3 question, but I think with this kind of registry,  
4 the trouble is you are absorbing a selection bias  
5 that probably has as much or more to do with any  
6 observed results. At least an historic control,  
7 you could try to structure to a population so you  
8 wouldn't have that kind of implicit bias.

9 DR. TRACY: It seems that maybe part of  
10 the problem is that there is selection bias for the  
11 type of device that is used in any given  
12 individual. It is going to be different, and the  
13 operator experience tremendously influences whether  
14 they do or do not use a particular device.

15 So I think that if there were a new device  
16 that was coming along, the only control that you  
17 could use would either be manual compression in  
18 perhaps the ACC database versus one of the approved  
19 device studies.

20 But I think it is just very difficult  
21 given the amount of bias that is inherent in this  
22 type of device to come up with a clean comparison.

23 DR. WHITE: I think the true value of this  
24 trial is not whether or not your operators are more  
25 skilled than the PMA published papers are, because

1 you can bet the people who are doing those trials  
2 are good at it. They wouldn't be doing those  
3 trials if they weren't fairly skilled, particularly  
4 at the randomized level. Those of us who do these  
5 trials get really good at these devices.

6 I am reassured that in real use, people  
7 aren't getting hurt with these devices, and they  
8 are able to select patients and perhaps use these  
9 devices to some optimum. I mean, the fact that you  
10 found some benefit here, or some lower risk,  
11 reassures me that people know how to use these  
12 devices to their optimal ability. In the  
13 randomized trials, you have really skilled  
14 operators doing the best they can, and it is sort  
15 of an even ground.

16 So I am reassured by your data, and I just  
17 think it means that out there in the real world,  
18 they are being used pretty well.

19 I don't think you can use historical  
20 controls, because the patient populations are so  
21 highly variable. Whether we are talking about  
22 diagnostic catheterization, interventional  
23 catheterization, the level of anti-coagulation are  
24 huge impacts, and you can make things look better  
25 or make things look worse depending upon how you

1 select. So the randomization here becomes a key.  
2 Now, if you want to accept randomization  
3 against another closure device, that would be an  
4 interesting model if there is enough data to make  
5 you convinced that you know--if you want to  
6 randomize PercuSurge against the next level of  
7 device, that might be something I would be willing  
8 to do, but I would want to make sure that it was  
9 randomized so that the risks were evenly  
10 distributed in both populations.

11 DR. TAVRIS: I certainly agree that  
12 randomization would be by far a preferable way of  
13 looking at this. The problem is that what  
14 randomized data we have is very, very small and  
15 would not be able to detect small differences.

16 DR. WHITE: I agree.

17 DR. TRACY: Are there any other comments?

18 [No response.]

19 DR. TRACY: If not, thank you very much  
20 for that presentation.

21 We'll move on to the discussion of  
22 premarket notification of the Embol-X aortic  
23 filter.

24 MS. WOOD: The following announcement  
25 addresses conflict of interest issues associated

1 with this meeting and is made part of the record to  
2 preclude even the appearance of an impropriety.

3 To determine if any conflict existed, the  
4 Agency reviewed the submitted agenda for this  
5 meeting and all financial interests reported by the  
6 Committee participants. The conflict of interest  
7 statutes prohibit Special Government Employees from  
8 participating in matters that could affect their or  
9 their employers' financial interest.

10 The Agency has determined, however, that  
11 the participation of certain members and  
12 consultants, the need for whose services outweighs  
13 the potential conflict of interest involved, is in  
14 the best interest of the Government.

15 Therefore, a waiver has been granted for  
16 Dr. Thomas Ferguson for his interest in a firm that  
17 could be affected by the Panel's recommendation.  
18 The waiver involves a grant to his institution for  
19 the sponsor's product study in which he had no  
20 involvement and for which funding was less than  
21 \$100,000 per year. Copies of this waiver may be  
22 obtained from the Agency's Freedom of Information  
23 Office, Room 12A-15, in the Parklawn Building.

24 In the event that the discussions involve  
25 any other products or firms not already on the



1 agenda for which an FDA participant has a financial  
2 interest, the participant should excuse him or  
3 herself from such involvement, and the exclusion  
4 will be noted for the record.

5 With respect to all other participants, we  
6 ask in the interest of fairness that all persons  
7 making statements or presentations disclose any  
8 current or previous financial involvement with any  
9 firm whose products they may wish to comment upon.

10 DR. TRACY: Thank you.

11 I'd like to ask the panel members to  
12 introduce themselves, starting with Mr. Morton.

13 MR. MORTON: My name is Michael Morton. I  
14 am the industry representative, and I am employed  
15 by Soren Kolb [phonetic] Cardiovascular.

16 DR. WHITE: Good morning. My name is  
17 Chris White. I am an interventional cardiologist  
18 from the Ochsner Clinic in New Orleans.

19 DR. LASKEY: I am Warren Laskey, an  
20 interventional cardiologist from the National Naval  
21 Medical Center in Bethesda.

22 DR. KRUCOFF: I am Mitch Krucoff, an  
23 interventional cardiologist from Duke University.

24 DR. AZIZ: I am Samil Aziz, adult cardiac  
25 surgeon in Denver and associate clinical professor

1 at the University of Colorado.  
2 DR. DeMETS: I am David DeMets. I am a  
3 biostatistician at the University of Wisconsin in  
4 Madison.  
5 DR. TRACY: I am Cindy Tracy. I am an  
6 electrophysiologist at Georgetown University  
7 Hospital.  
8 MS. WOOD: Geretta Wood, Executive  
9 Secretary.  
10 DR. EDMUNDS: I am Hank Edmunds,  
11 University of Pennsylvania, a surgeon.  
12 DR. MARLER: I am John Marler, Associate  
13 Director for Clinical Trials at the National  
14 Institute of Neurological Disorders and Stroke, and  
15 I am a neurologist.  
16 DR. FERGUSON: Tom Ferguson, a cardiac  
17 surgeon at Washington University in St. Louis.  
18 DR. PINA: Ileana Pina, heart failure  
19 transplant cardiologist, Case Western Reserve in  
20 Cleveland.  
21 MR. DACEY: Robert Dacey, Consumer  
22 Representative, from Boulder County, Colorado.  
23 DR. ZUCKERMAN: Bram Zuckerman, Director,  
24 Division of Cardiovascular Devices, Food and Drug  
25 Administration.

1 DR. TRACY: Thank you.

2 MS. WOOD: Pursuant to the authority  
3 granted under the Medical Devices Advisory  
4 Committee Charter dated October 27, 1990, and as  
5 amended August 18, 1999, I appoint the following  
6 individuals as voting members of the Circulatory  
7 System Devices Panel for this meeting on October  
8 23, 2002: Christopher White, M.D.; L. Henry  
9 Edmunds, Jr., M.D.; Mitchell W. Krucoff, M.D.; John  
10 Marler, M.D.; Thomas B. Ferguson, M.D.; David L.  
11 DeMets, Ph.D.

12 For the record, these people are Special  
13 Government Employees and are consultants to this  
14 panel and other panels under the Medical Devices  
15 Advisory Committee. They have undergone the  
16 customary conflict of interest review and have  
17 reviewed the material to be considered at this  
18 meeting.

19 This is signed by David W. Feigel, Jr.,  
20 M.D., M.P.H., Director, Center for Devices and  
21 Radiological Health, and dated October 10, 2002.

22 Pursuant to the authority granted under  
23 the Medical Devices Advisory Committee Charter of  
24 the Center for Devices and Radiological Health  
25 dated October 27, 1990 and as amended August 18,

1 1999, I appoint the following individual as a  
2 voting member of the Circulatory System Devices  
3 Panel for the meeting on October 23, 2002: Ileana  
4 L. Pina, M.D.

5 For the record, Dr. Pina is a consultant  
6 to the Cardiovascular and Renal Drugs Advisory  
7 Committee of the Center for Drug Evaluation and  
8 Research. She is a Special Government Employee who  
9 has undergone the customary conflict of interest  
10 review and has reviewed the material to be  
11 considered at this meeting.

12 This is signed by William K. Hubbard,  
13 Senior Associate Commissioner for Quality Planning  
14 and Legislation, and it is dated October 18, 2002.

15 DR. TRACY: Thank you.

16 At this point, we'll move to the open  
17 public hearing. There were no scheduled speakers,  
18 but is there anyone in the audience who wishes to  
19 address the panel on today's topic or any other  
20 topic?

21 [No response.]

22 DR. TRACY: If not, we will close the open  
23 public hearing and move on to the presentation.

24 MS. WOOD: I would just like to remind the  
25 speakers to introduce yourself and state your

1 conflict of interest.

2 Sponsor Presentation

3 EMBOL-X, Inc.

4 K022071, EMBOL-X Aortic Filter

5 MS. CHANG: Thank you.

6 My name is Jean Chang. I am the Chief

7 Operating Officer for EMBOL-X, and I would like to

8 thank the FDA, our panel reviewers, and all panel

9 members for the opportunity to present our clinical

10 results today.

11 [Slide.]

12 This is the presentation that we have

13 planned. After I do a company overview, Dr.

14 Nicholas Kouchoukos, Missouri Baptist Medical

15 Center, the co-principal investigator, will give an

16 overview of atheroembolism in cardiac surgery.

17 Dr. Richard Kuntz, from Brigham and

18 Women's, will present our EMBOL-X clinical trial

19 design.

20 And finally, Dr. Keith Allen, who is the

21 site PI at Saint Vincent, will present the clinical

22 study results.

23 [Slide.]

24 EMBOL-X is a small, privately-funded

25 company that was founded in 1996 by two physicians.

1 It is in Northern California, and has less than 50  
2 employees.

3 The product focus from the start has been  
4 intra-aortic filtration utilized during cardiac  
5 surgical procedures.

6 The first clinicals were done at the end  
7 of 1997; CE marked the product in the end of 1998.  
8 And the product has been commercially available in  
9 Europe since 1999, with over 2,000 documented  
10 cases.

11 [Slide.]

12 What we show here are the two devices that  
13 make up the EMBOL-X intra-aortic filtration system.  
14 The top device there is the EMBOL-X aortic cannula,  
15 which has the premarket [inaudible] this past  
16 September is a modified standard cannula.

17 The bottom device is the subject of our  
18 presentation, which is the EMBOL-X intra-aortic  
19 filter. The distal filter basket there is composed  
20 of two primary components. It is a polyester mesh  
21 with 120 microns, and the polyester mesh is heparin  
22 coated with a duraflow heparin coating, which is  
23 the same heparin coating that is used [inaudible]  
24 filters.

25 What this demonstrates is the principle of

1 use of the filter in cardiac surgery. It has been  
2 inserted through the sideport in the cannula, and  
3 again, as you will note here, the filter captures  
4 particulates that arise from the heart up to and  
5 proximal to the arterial cannula. It does not  
6 capture particulate that is distal to the cannula,  
7 including arterial flow.

8 What you see on the right there is a  
9 representative sample of the particulates that are  
10 captured. The grid marks there are 3 mm, and Dr.  
11 Allen will talk more about particulate capture as  
12 we go forward.

13 [Slide.]

14 Embolic protection, embolite capture, is  
15 not a new technology, and the EMBOL-X intra-aortic  
16 filter follows along the same lines as existing  
17 devices that are either currently approved or under  
18 investigation. And as we discussed earlier, the  
19 extracorporeal filter is standard in cardiac CPB  
20 surgery. The PercuSurge distal protection device  
21 is used for SVG, and Dr. Kuntz will talk a little  
22 about that.

23 And finally, in other arterial beds, there  
24 are other distal filtration devices that filter  
25 particulate emboli in the [inaudible] vein area

1 graft, which is for carotid intervention.

2 And finally, in a different area, there  
3 are vena cava filters.

4 [Slide.]

5 All these devices capture particulate  
6 emboli, and it is the basis for the indications for  
7 our device. "The EMBOL-X aortic filter is  
8 indicated for use with the EMBOL-X aortic cannula  
9 in cardiac surgery procedures to contain and remove  
10 particulate emboli." This is the basis for our  
11 clinical study design and for the clinical study  
12 results which you will hear later on today.

13 I would now like to present Dr. Nicholas  
14 Kouchoukos from Missouri Baptist.

15 DR. KOUCHOUKOS: Thank you very much,  
16 Madam Chairman, members of the panel.

17 I am Nicholas Kouchoukos from the Missouri  
18 Baptist Medical Center in St. Louis, Missouri. I  
19 served as the co-principal investigator in this  
20 trial and was the principal investigator at the  
21 Missouri Baptist Medical Center.

22 I have no financial interest in the  
23 company or any equity investment in the company. I  
24 have been reimbursed for my travel expenses, and a  
25 grant on my behalf for services rendered has been



1 made to the Educational Research Foundation of the  
2 Missouri Baptist Medical Center.

3 [Slide.]

4 Since the beginning of open heart surgery,  
5 employing cardiopulmonary bypass has been  
6 recognized, and embolization of atheromatous debris  
7 from the atherosclerotic aorta is a cause of stroke  
8 and other embolic-related complications.

9 Until the last decade, there were  
10 scattered case reports implicating atherosclerosis.  
11 In 1992, in a landmark study published from the  
12 Cleveland Clinic by Christopher Blauth and  
13 colleagues, they autopsied 221 patients who had  
14 died following cardiac surgical procedures, and  
15 they observed a high prevalence of atheroembolism  
16 in these patients and were able to correlate the  
17 presence of atheroembolism with increasing age and  
18 aortic atherosclerosis, as well as the presence of  
19 peripheral vascular disease.

20 [Slide.]

21 Among the patients with atheroembolism who  
22 had atherosclerosis of the ascending aorta, the  
23 presence of atheroembolism to various organs was 37  
24 percent; among the patients who had no significant  
25 atherosclerosis, the prevalence was 2 percent.

1 This was a highly significant difference.

2 [Slide.]

3 There was a high correlation with  
4 increasing age and the presence of atherosclerotic  
5 disease in the ascending aorta.

6 In a clinical trial we conducted at  
7 Washington University in St. Louis using epiaortic  
8 scanning to determine the severity of  
9 atherosclerosis in the ascending aorta, we also  
10 observed a substantial correlation between  
11 increasing age and the prevalence of severe  
12 atherosclerosis in the ascending aorta.

13 Among the patients over the age of 80, for  
14 example, 33 percent of the patients had moderate or  
15 severe atherosclerosis.

16 [Slide.]

17 The prevalence of other risk factors for  
18 increased mortality and morbidity in patients  
19 undergoing cardiac surgical procedures such as  
20 coronary bypass grafting has increased over time.

21 This is a study from the Society of  
22 Thoracic Surgeons Database looking at a subset of  
23 Medicare patients, that is, those over the age of  
24 65, and looking at the prevalence of important risk  
25 factors for mortality and morbidity over a decade

1    between 1990 and 1999. This involved over 620,000  
2    patients.  
3        In this analysis, there was a substantial  
4    increase in many of the important risk factors  
5    associated with mortality and morbidity.  
6        [Slide.]  
7        This is an example of atherosclerosis in  
8    the ascending aorta with very friable material  
9    located circumferentially in this aorta, and this  
10   is the material that is at risk for dislodgement  
11   during cardiac surgical procedures where  
12   manipulation of the aorta with interventions such  
13   as cannulation or clamping is prone to dislodge  
14   this material.  
15       [Slide.]  
16       In the study by Blauth and colleagues  
17   looking at the organs that were affected with  
18   atheroembolism, the most common site was in the  
19   brain, and this was followed by the spleen and the  
20   kidney. This is not surprising because  
21   approximately 40 percent of the cardiac output is  
22   delivered to these two organs.  
23       [Slide.]  
24       These are some examples of small  
25   atheroemboli in the cerebral circulation, and

1 below, in the presence of a large atheroembolism in  
2 a medium-sized artery. These are cortical infarcts  
3 in a patient following cardiac surgical procedure,  
4 and again, a large embolus of atheromatous material  
5 present in one of the renal arteries.

6 [Slide.]

7 We also looked in the early 1990s at the  
8 association of atherosclerosis as a predictor and  
9 the development of postoperative renal dysfunction.  
10 The index of renal dysfunction was the elevation of  
11 the keratinine to a level above 2.0, or an increase  
12 of 50 percent from baseline. And we correlated  
13 these changes in renal function with the presence  
14 of ascending atherosclerosis determined by  
15 epi-aortic scanning.

16 There was a correlation with the severity  
17 of disease and the prevalence of renal dysfunction.

18 [Slide.]

19 The analyses were performed on Day 1 and  
20 Day 6, and using multivariate analysis on the first  
21 postoperative day, ascending atherosclerosis was  
22 the only independent predictor of renal  
23 dysfunction. On Day 6, it was one of three  
24 predictors of renal dysfunction, along with low  
25 postoperative cardiac output and preoperative left

1 ventricular dysfunction.

2 [Slide.]

3 There are a number of interventions which  
4 have been designed and implemented in an attempt to  
5 reduce the frequency and severity of embolization  
6 from the ascending aorta.

7 Extracorporeal filtration will remove  
8 embolic material, but it is not likely to remove  
9 any material from the ascending aorta.

10 The interventions that are commonly used  
11 are those that involve minimal manipulation of the  
12 ascending atherosclerotic aorta. The use of a  
13 single cross-clamp rather than placement of  
14 multiple clamps reduces the frequency of  
15 manipulation of the aorta and, presumably, the  
16 dislodgement of atheromatous debris.

17 The use of proximal anastomotic devices to  
18 avoid the placement of clamps on the aorta may have  
19 a protective effect.

20 The use of off-pump surgery avoids  
21 placement of clamps on the ascending aorta, and  
22 other techniques such as hypothermic fibrillation  
23 and circulatory arrest have been utilized, again,  
24 to avoid clamping and other manipulation of the  
25 ascending aorta.

1           Although these techniques may be  
2 effective, they have certain limitations. All of  
3 the atheroembolism cannot be eliminated with these  
4 techniques. For example, the proximal and  
5 anastomotic devices do involve manipulation of the  
6 ascending aorta with the potential for  
7 dislodgement. Off-pump surgery also involves  
8 manipulation of the heart and the aorta, despite  
9 the fact that no clamps are placed on the aorta.

10           Furthermore, there are other sources of  
11 emboli. The left atrial appendage can release  
12 thrombus; there can be neural thrombus in the left  
13 ventricular cavity that can be released, and also  
14 debris from diseased mitral and aortic valves, and  
15 also surgical debris.

16           The intra-aortic filter has the capacity  
17 to capture this debris as well.

18           [Slide.]

19           In a study by Dr. Denise Barbut and her  
20 colleagues, looking at embolization of particulate  
21 matter, they utilized transesophageal  
22 ecocardiography and transcranial doppler and  
23 identified the release of emboli during cardiac  
24 surgical procedures.

25           This is just the distribution of the

1 particle size of these emboli that were released in  
2 a study in 10 patients. And above is shown the  
3 diameters of various vessels in the intracranial  
4 circulatory system--the leptomeningeal vessel, the  
5 small cortical arteries, the posterior cerebral  
6 artery, the branches of the middle cerebral, and  
7 here, the larger middle cerebral artery and the  
8 internal carotid artery.

9 The diameters of these particles  
10 corresponds to the diameters of these arteries.

11 [Slide.]

12 In another study, Dr. Barbut and her  
13 colleagues looked at the temporal sequence of  
14 release of emboli from the aorta during the conduct  
15 of a cardiac surgical procedure. They found that  
16 the majority of these particles were released at  
17 the time of the release of the aortic cross-clamp  
18 from the ascending aorta.

19 In fact, over 70 percent of the emboli  
20 were released in 20-second interval following the  
21 release of the clamp. This is the rationale for  
22 inserting the intra-aortic filter just before  
23 release of the aortic clamp during the cardiac  
24 procedure.

25 [Slide.]

1           As I have indicated, there are important  
2 complications that can result from atheroembolism.

3           Stroke is the one that has caused the  
4 greatest concern because of the important mortality  
5 and morbidity that results from stroke. And it is  
6 now clearly recognized that atheroembolism is the  
7 principle cause of stroke following cardiac  
8 surgical procedures.

9           There is also evidence for renal and other  
10 organ system dysfunction. Pathologic and clinical  
11 studies that we have presented suggest that  
12 embolization may be an important cause of renal  
13 dysfunction postoperatively.

14          Embolization is also a possible  
15 contributing factor to postoperative neurocognitive  
16 dysfunction.

17          There have been strategies employed to  
18 reduce serious embolic-related complications.  
19 Obviously, prevention would be the best option, and  
20 this would involve minimal or no manipulation of  
21 the aorta, but this is not 100 percent effective in  
22 eliminating embolization.

23          Resection of the diseased aorta is a way  
24 to eliminate the emboli, but is only applicable to  
25 a very small percentage of patients.



1           Reduction of the embolic load is an  
2 attractive way to reduce this embolization.

3           Diversion is one option, to divert the  
4 material away from the central nervous system, but  
5 this would merely disseminate this material to  
6 other organs.

7           And capture, using an intra-aortic filter,  
8 is an attractive method for capturing this embolic  
9 material.

10          [Slide.]

11          In the subsequent presentations, Dr. Kuntz  
12 and Dr. Allen will discuss the design and execution  
13 of a large, randomized clinical trial evaluating  
14 the safety and efficacy of the EMBOL-X intra-aortic  
15 filter. This study is well-designed, in my  
16 opinion, and clearly demonstrates that use of the  
17 EMBOL-X filter is a rational, safe, and beneficial  
18 intervention for removal of atheromatous and other  
19 embolic material from the ascending aorta of  
20 cardiac surgical patients.

21          Thank you.

22          DR. KUNTZ: Good morning. My name is Rick  
23 Kuntz. I am a cardiologist at Brigham and Women's  
24 Hospital in Boston.

25          I got involved with this group about 3

1 years ago because of my interest in designing and  
2 working with trials on embolic protection in the  
3 heart, in the brain, and in the kidney. This  
4 afforded me an opportunity to work with a company  
5 who was looking at another way to impact on the  
6 reduction of embolic problems associated with  
7 cardiac surgery.

8 My interest in this study is mainly  
9 academic. I have no financial interest in the  
10 company. I have no equity. I am being reimbursed  
11 for my travel, and a small grant was made to the  
12 Department of Medicine on behalf of this  
13 consultation.

14 [Slide.]

15 The purpose of this study--and I am going  
16 to talk about the rationale as to how we came up  
17 with the design for the study--was from the outset  
18 to demonstrate the ability of this device to safely  
19 and effectively remove visible particulate emboli  
20 during cardiac procedures.

21 So at the outset, there was an assumption  
22 that these particles were bad--that they floated  
23 around in the bloodstream, that they would be  
24 released with the cross-clamp, and that they  
25 probably don't do good things if they go around the

1 arch of the aorta.  
2 So from the beginning, it was important to  
3 understand that we were trying to remove these  
4 things, and how to measure them clinically was the  
5 biggest challenge in trying to design the trial.  
6 [Slide.]  
7 So the goal was to design a clinical trial  
8 that evaluates the utility of the device aimed to  
9 prevent the dissemination of released emboli  
10 following cross-clamping.  
11 Now, here is the dilemma. We have  
12 particles that we can pull out, but the question is  
13 going to be what will these particles mean--is it  
14 really important to take them out or not. So the  
15 bet way to correlate that is to try to find hard  
16 clinical endpoints that could be collated overall.  
17 And we were struck with trying to design a trial to  
18 demonstrate that, because we were dealing with a  
19 problem of embolic showers that might not manifest  
20 themselves as frank organ infarctions.  
21 So for example, if you have small  
22 particulate emboli that cause microvascular  
23 injury--organs such as the brain, the kidney, the  
24 spleen, and others--it might not be demonstrated as  
25 a frank, say, NIH-level major stroke or as a kidney

1 infarction.

2           And there was some stuff that we had  
3 learned from heart trials which I will talk about  
4 in a second, but ultimately, this was the biggest  
5 issue we had to deal with.

6           One of the potential roles in looking at  
7 the impact of shower emboli was to measure  
8 cognitive dysfunction, and I will talk about the  
9 availability of instruments at the time of the  
10 trial design and whether there was consensus in the  
11 surgical community as to whether that could be  
12 applied or not.

13           So these various study designs were  
14 explored and discussed with the FDA, and  
15 ultimately, in multiple discussions at which I was  
16 present, the focus was to demonstrate safety of  
17 this device with removal of particulate emboli as  
18 some demonstration of efficacy.

19           Let me give you a parallel about the  
20 importance of shower emboli and how you can measure  
21 it in an organ that actually does give you a  
22 clinical outcome with shower emboli.

23           [Slide.]

24           There is a device on the market to protect  
25 emboli from intervention on vein grafts through the

1 heart. This is a picture of the typical amount of  
2 emboli that is removed in the vein graft  
3 intervention, and they are manifested mainly by  
4 outcomes measured by cardiac enzyme elevation.

5 So the heart is a nice organ, because it  
6 actually can show the impact of shower emboli,  
7 mainly manifested by elevations of cardiac enzymes.

8 [Slide.]

9 If we look at the primary endpoint of this  
10 study, which was 8 percent in those patients  
11 randomized to protection, where we actually removed  
12 particles, compared to nothing at all where  
13 particles were not removed, there was a 50 percent  
14 reduction in the major endpoint of the trial.

15 [Slide.]

16 But if we look at an index like frank  
17 organ infarction, which would be QMI, something we  
18 could pick up clinically, such as a change in the  
19 EKG or Q-wave, it only represented about 10 percent  
20 of the outcome. The majority of the outcome of  
21 this endpoint was measured by enzyme elevation  
22 which didn't manifest itself as a frank organ  
23 infarction.

24 So it is important to understand that at  
25 least in the heart, shower emboli do have an impact

1 that led to approval of the device, but wasn't  
2 mainly manifested in anything else other than a  
3 cardiac enzyme elevation.

4 [Slide.]

5 So if we look at this issue, the shower  
6 emboli from vein grant intervention, it does not  
7 usually manifest as a frank MI, and in the SAFER  
8 trial, the availability of cardiac enzyme rise is  
9 essentially used to show utility of the device.

10 This reduction in myocardial infarction  
11 led to the approval of the device, and I think  
12 there was a general consensus across the  
13 interventional community that this was a good thing  
14 to use and now is considered to be a standard of  
15 care for vein grant intervention.

16 Now, the same endpoints are also used in  
17 the heart to approve the whole classification of  
18 IIb/IIIa inhibitors. That is, another valuable,  
19 considered standard therapy, across our area was  
20 based on the reduction, mainly in emboli, that  
21 manifested themselves as cardiac enzyme elevations  
22 but not frank organ infarction.

23 This was not applicable to the EMBOL-X  
24 system because it was north of the heart; this  
25 device wasn't designed to protect emboli in the

1 coronary arteries following bypass surgery. And  
2 there were few parallel sensitive measures that  
3 were available for noncardiac end organs. That is,  
4 we don't have an enzyme elevation for the brain or  
5 for the kidney like we do for the heart to measure  
6 the impact of these emboli overall.

7 [Slide.]

8 So if we look at the distribution of  
9 organs that Dr. Kouchoukos showed that were targets  
10 for emboli from previous studies, we have a lot of  
11 important organs that we don't want in embolism,  
12 obviously, but we don't have good, readily  
13 available measures to demonstrate their injury  
14 pattern from shower emboli. And this was a  
15 conundrum that we were stuck with in trying to come  
16 up with an effective endpoint ultimately to  
17 demonstrate utility of this device by a clinical  
18 signal.

19 [Slide.]

20 So the issue raised by the FDA in our  
21 meeting was that there are few sensitive and  
22 specific measures available to look at the  
23 noncardiac end organs and their impact from shower  
24 emboli.

25 Neurological assessment was obviously a

1 very important one to evaluate, and there was much  
2 time spent look at all the available ways of  
3 looking at the neurological outcomes, because after  
4 all, the brain does receive approximately 20  
5 percent of the circulation of the cardiac output  
6 and was obviously a target that we wanted to  
7 reduce.

8 Well, if we wanted to look at frank  
9 reduction in stroke, as Dr. Kouchoukos showed--that  
10 stroke is likely involved with emboli per se--the  
11 instance of stroke following cardiac surgery was  
12 large enough that this would have to be a very,  
13 very large sample size in order for us to  
14 demonstrate a reduction. Now, a 20 percent  
15 reduction is pretty small, but still, if you are  
16 looking for 30 or 40 percent, we are talking about  
17 5,000 to 10,000 patients minimum to demonstrate a  
18 reduction in the 2 to 3 percent stroke rate seen  
19 postoperatively.

20 So what about measuring the cognitive  
21 function per se. Well, there are a lot of issues  
22 raised regarding using cognitive function as an  
23 endpoint in this study, and it is very  
24 controversial. First of all, the cell deficits may  
25 be due to diffuse small vessel embolism, to be



1 sure, but there are other multifactorial causes of  
2 cognitive dysfunction after surgery that may  
3 involve general anesthesia.

4 And even though it will be important to  
5 look at that at the time of the study and, I would  
6 even argue today to some degree, there is still no  
7 great consensus about instruments available for  
8 psychometric or neurological outcomes that measure  
9 cognitive dysfunction that has been accepted in the  
10 cardiovascular community, and at the time of this  
11 study, we couldn't get consensus along the lines of  
12 understanding whether to apply a battery of tests,  
13 most of which still have not been validated.

14 [Slide.]

15 So the practical approach was that the  
16 huge sample size to show a reduction of frank  
17 infarction such as stroke was just not feasible,  
18 because this is a large, randomized cardiac  
19 surgical trial, and it was unlikely that we could  
20 do a 5,000 to 10,000 patient study. Cognitive  
21 dysfunction could not be readily measured with  
22 mature, validated instruments was the conclusion  
23 that we reached in discussions with the FDA, and  
24 the proof of safety plus demonstration of captured  
25 emboli seemed to be the most feasible and logical

1 approach to go forward.

2 So when you look at this study to say why  
3 wasn't there a clear clinical measure of efficacy  
4 of this endpoint per se, it is because we wrestled  
5 with endpoints that had consensus to demonstrate  
6 true efficacy from a clinical perspective.

7 And the final conclusion was that we would  
8 demonstrate safety by the safety endpoints to  
9 demonstrate this didn't cause any increases in  
10 those elements in the safety endpoints. And if we  
11 demonstrated that with removal of actual  
12 particulate emboli, at least there would be some  
13 measure of utility. Now, whether the utility would  
14 be enough for product approval, I think will be the  
15 discussion of this panel.

16 Therefore, the approved IDE study design  
17 was safety equivalency for the composite primary  
18 endpoint and effectiveness through demonstration of  
19 particulate capture.

20 [Slide.]

21 So given that, there was a prospective  
22 study design; multi-center, 21 sites; a sample size  
23 of 1,289 patients was calculated using a Blackwater  
24 [phonetic] formula for equivalency of an expected  
25 outcome of 15 percent plus a 5 percent delta; the

1     EMBOL-X aortic catheter was randomized through  
2     Standard J tip cannula; and the primary endpoints  
3     again were effectiveness with the demonstration of  
4     particulate emboli capture and safety with  
5     equivalence of the safety profile using the current  
6     standard procedures.

7             [Slide.]

8             So this safety endpoint, which might be  
9     viewed as also a measure of efficacy, was  
10    necessarily not refined enough to demonstrate  
11    efficacy based on this sample size. It was mainly  
12    used to demonstrate that there would be no increase  
13    in problems associated with the instance of death,  
14    myocardial infarction, renal insufficiency, GI  
15    complications, limb-threatening embolisms, or  
16    neurologic deficit, either mild or severe, using  
17    the NIH Stroke Scale and other stroke measures.

18            The safety endpoint was designed to  
19    demonstrate freedom from device complications.

20            Now, it is important to point out that  
21    this is mainly a safety endpoint, and for example,  
22    the inclusion of myocardial infarction is important  
23    to have in a study looking for safety, but we  
24    wouldn't aim to actually improve myocardial  
25    infarctions, because the device is north of the

1 heart, as it were.  
2 [Slide.]  
3 The effectiveness endpoint and hypothesis,  
4 therefore, was successful capture of the emboli,  
5 and this was defined as retrieved particles  
6 observed at 10X power at the operating table.  
7 And the hypothesis was that we would  
8 capture greater than 75 percent of the cases that  
9 would have emboli that was evident.  
10 [Slide.]  
11 The sample size was driven by the safety  
12 endpoint--small sample size needed to demonstrate  
13 primary effectiveness endpoint of particulate  
14 capture, and the 1,286 patients were used to  
15 demonstrate safety, and there was a calculation for  
16 one interim analysis using a boundary condition  
17 under Bryant-Fleming [phonetic] for the Blackwater  
18 test.  
19 [Slide.]  
20 Safety was monitored by blinded,  
21 independent Clinical Events Adjudication Committees  
22 and the independent [inaudible] Monitoring  
23 Committee. There was an independent medical  
24 monitor. The Core Laboratories were blinded.  
25 Randomization was performed just prior to

1 cannulation in the operating room, and a  
2 neurological examiner and the patients themselves  
3 were blinded to treatment assignment.

4 [Slide.]

5 We had independent EKG Core Laboratories  
6 and histological laboratories to evaluate the  
7 emboli.

8 [Slide.]

9 There was a separate ecocardiographic  
10 imaging core laboratory as well for the epi-aortic  
11 as well as TEE endpoints.

12 [Slide.]

13 And to put this into perspective, as large  
14 randomized trials in surgery are difficult to do,  
15 this ranks among the top enrolling randomized  
16 studies in the history of randomized trials in  
17 cardiac surgery. So this was quite an effort to do  
18 this well-designed trial in order to demonstrate  
19 the endpoints that Dr. Allen will review.

20 [Slide.]

21 So if we summarize this, safety was to be  
22 demonstrated under an equivalence endpoint in the  
23 agreed-upon IDE using a broad net composite safety  
24 endpoint chose, which included myocardial  
25 infarction, for example. The safety endpoint was

1 not optimized to demonstrate clinical efficacy or  
2 superiority.

3 Now, there is no question that this could  
4 be used if we had a huge sample size to demonstrate  
5 some reductions in embolic injury, but we didn't  
6 want to fool ourselves by thinking this initially  
7 would be the primary viewpoint of this endpoint  
8 overall to demonstrate utility.

9 Therefore, the utility was focused on  
10 demonstrating safety first, followed by efficacy to  
11 show frequency of actual particulate removal from  
12 the operating room.

13 Now I'll turn it over to Dr. Allen.

14 DR. ALLEN: Thank you, Madam Chairman and  
15 members of the panel.

16 My name is Keith Allen, and I was a site  
17 principal investigator. I practice as a  
18 cardiovascular and thoracic surgeon out of Saint  
19 Vincent Hospital in Indianapolis.

20 [Slide.]

21 From a financial disclosure standpoint, I  
22 have no financial interest in the country and  
23 certainly no equity investment in the company. I  
24 was reimbursed for my travel and time expenses to  
25 come today.

1 [Slide.]

2 On behalf of the 88 investigators at 22  
3 centers across the U.S. and one Canadian site, I  
4 thank the panel for the opportunity to present our  
5 clinical results. I think, as you can see from our  
6 centers that were utilized in this study, they  
7 represent a broad spectrum of cardiac surgery in  
8 North America, involving both private, academic,  
9 and community centers across the board.

10 [Slide.]

11 The inclusion and exclusion criteria are  
12 summarized. Obviously, as with any large study  
13 like this, particularly when you are looking at  
14 safety as an endpoint, there are a number of  
15 exclusion criteria to confine your sample size to  
16 patients who are going to demonstrate safety for  
17 you.

18 The inclusion criteria were confined to  
19 patients who were 60 years and older who either had  
20 primary CABG or primary valve procedure.

21 Some of our exclusion criteria that we  
22 feel are important were dialysis dependent, a  
23 patient who had a previous stroke who had a  
24 residual deficit, or previous surgery or damage to  
25 the aorta.

1 Obviously, the filter has various sizes  
2 and is part of the randomization process. To be  
3 included in the study, you had to have an internal  
4 diameter of the ascending aorta that would  
5 appropriately fit a filter that you could put in  
6 the patient.

7 [Slide.]

8 About 15 percent of patients screened for  
9 this study ultimately met inclusion and exclusion  
10 criteria, resulting in 1,394 patients available for  
11 the study.

12 As is common with any study of this size  
13 and nature in which a new device is being placed in  
14 a clinician's hands, we had as a component of our  
15 study a roll-in phase. Each investigator was  
16 required to do at least one nonrandomized patient  
17 to gain familiarity with the device and understand  
18 how it could be used and inserted appropriately.

19 While this does not impact the study  
20 results, I will concentrate the rest of our data  
21 analysis on patients who were actually randomized  
22 between filtered and nonfiltered arms.

23 We ended up with 1,289 patients who were  
24 evenly distributed between filter with standard  
25 cannula or simply receiving the standard cannula



1 alone, without intra-aortic filtration.

2 It is important to understand our  
3 randomization stratification stream, and it really  
4 involved three components. Patients were  
5 stratified based on whether they were a valve or a  
6 primary CABG and, importantly, we randomized  
7 patients based on injection fraction.

8 [Slide.]

9 There were a number of key baseline and  
10 medical variables that were obviously evaluated in  
11 this study. There were four variables that  
12 differed between groups. One variable that favored  
13 a control arm was a patient-given history of aortic  
14 disease.

15 There were three variables--atrial  
16 fibrillation, valvular dysfunction and severe  
17 carotid disease--which all favored the filter arm.  
18 It is important, though, when a multivariable  
19 analysis was done on these discrepant variables at  
20 the end of the study, there was no interaction or  
21 impact on our results.

22 [Slide.]

23 It is amazing, as different as cardiac  
24 surgeons are across the board how uniformly this  
25 operation was done across the centers that were

1 involved in this study. There was a good  
2 distribution between CABG and valve patients.  
3 There was a good distribution between whether a  
4 partial clamp was used or whether a single-clamp  
5 technique was used.

6 We tried to look at things like whether  
7 the aortic cross-clamp was repositioned as an  
8 impact on embolic release, and that was similar  
9 between arms. Things like cross-clamp time and  
10 cardiopulmonary bypass time were also familiar and  
11 similar between groups. And, importantly, we  
12 wanted to look at nuances like were the number of  
13 proximal anastomoses done between groups similar,  
14 because obviously, you are manipulating the  
15 ascending aorta, and we wanted to ensure that one  
16 group wasn't having more proximal anastomoses done  
17 than the other. And in fact, they were identical  
18 between the two arms.

19 When we looked at filter dwell time, which  
20 obviously is not applicable to the control, the  
21 filter dwell time in our patients was approximately  
22 21 minutes.

23 [Slide.]

24 As outlined very nicely by Dr. Kuntz our  
25 primary composite endpoint was a safety endpoint.

1 It is important to reemphasize the fact that this  
2 was an equivalence safety endpoint, and I think we  
3 did achieve that safety endpoint. 17.1 percent of  
4 the treatment group compared to 18.9 percent of the  
5 control group had a composite event that was a  
6 priori defined. And once again, it is a safety  
7 endpoint, and it was not intended to capture  
8 clinical effectiveness of the device. It is  
9 important that the panel understand this and that  
10 they don't confuse this safety endpoint as a  
11 surrogate for clinical efficacy, because as Dr.  
12 Kuntz pointed out, if you were designing a  
13 composite endpoint to demonstrate clinical  
14 efficacy, you certainly would not have included  
15 myocardial infarction which the device can have no  
16 impact on and that occurred and represented  
17 approximately half of the events in our composite  
18 endpoint.

19 [Slide.]

20 Any time somebody presents or uses a  
21 composite endpoint, as an investigator, I always  
22 want to see the details of all the components that  
23 were involved in creating that composite endpoint  
24 to ensure that there are not trends favoring one or  
25 the other that even out when you do just the

1 composite endpoint.

2       When you look across-the-board at the  
3 components of our composite safety endpoint, we  
4 don't make claims of superiority in any area, but  
5 what you clearly see is equivalence and safety of  
6 the device in a very large prospective randomized  
7 trial, not only with the composite endpoint, but  
8 with the individual components of that composite  
9 endpoint.

10       [Slide.]

11       Clearly, in a trial like this where you  
12 are presenting people with a major surgical  
13 operation, you look at other serious adverse  
14 events, and once again, it is striking how evenly  
15 distributed these are across centers in this very  
16 large study. When you look at serious adverse  
17 events across the board, there was no statistical  
18 difference between either arm. And when you come  
19 down to actually tallying up whether or not  
20 patients had a serious adverse event in this very  
21 large study, the bottom line down at the bottom is  
22 that they were absolutely identical between both  
23 the control and the filter group.

24       [Slide.]

25       As an investigator who was asked to

1 participate in this trial, one thing that we were  
2 interested in is that obviously, you are placing  
3 something inside the ascending aorta. So one of  
4 the key adverse events that I was interested in was  
5 does placing this device inside the ascending aorta  
6 potentially cause harm to the patient. And we  
7 captured that using, out of 18 of the 22 centers,  
8 sophisticated imaging--primarily epi-aortic  
9 scanning, but also transesophageal  
10 ecocardiography--to try to capture whether the  
11 device was leaving some type of footprint within  
12 the ascending aorta.

13 And we looked at it as both did it cause  
14 ascending aortic dissections, and did it have  
15 ecocardiographic or imaging abnormalities that we  
16 might term aortic wall or intimal changes.

17 [Slide.]

18 I think that for the surgeons on the  
19 panel, this picture really doesn't need much  
20 explanation. As an investigator, this is what I  
21 was most concerned about when I was going to put  
22 this device in a patient--was I going to cause a  
23 clinically significant and relevant aortic  
24 dissection?

25 Clearly, on the left, the blue, engorged

1 aorta that makes your heart skip as a cardiac  
2 surgeon when it occurs is a dreadful complication  
3 that has serious clinical implications.

4 As a surgeon, this is very apparent, and  
5 while the TEE is dramatic in showing it, epi-aortic  
6 scanning or transesophageal ecocardiography aren't  
7 necessary for me to make this diagnosis.

8 And it is interesting when we look at this  
9 serious clinical event, there were two ascending  
10 aortic dissections seen in the control arm, and  
11 there were no ascending aortic dissections in the  
12 EMBOL-X filter arm.

13 [Slide.]

14 What we did see was a footprint that may  
15 be left by the device. As I told you and as Dr.  
16 Kuntz pointed out in his study design, 18 of our 22  
17 centers utilized either epi-aortic scanning or  
18 transesophageal ecocardiography peri-procedurally  
19 to look at the ascending aorta. And Dr. Weismann  
20 at the core lab for ecoocardiography did that very  
21 detailed blinded review. And, as will be pointed  
22 out later by the FDA, there was an incidence of  
23 endothelial disruptions or what I call intimal  
24 abnormalities seen more frequently in the filter  
25 group compared to the control arm.

1           What were these endothelial disruptions,  
2 and what clinical context can we put them in?

3           [Slide.]

4           I think a series of images will hopefully  
5 clarify that.

6           On your left is an epi-aortic scanning of  
7 the first endothelial disruption identified very  
8 early on in the study by a surgeon. There were  
9 three endothelial disruptions that were identified  
10 by surgeons early on in the study and that were  
11 elected to be repaired. One of those endothelial  
12 disruptions was an inadvertent stab from an  
13 11-blade knife to the posterior wall of the  
14 ascending aorta. The filter certainly didn't cause  
15 that.

16           But there were two endothelial  
17 disruptions, both occurring in the first four  
18 months of the study, both at the same center, in  
19 which surgeons elected to repair them. There was  
20 no historical basis for these. They weren't in the  
21 setting of an acute ascending dissection. But the  
22 surgeon had no background about what these  
23 endothelial disruptions are, and what you see--and  
24 it is hard to see unless you turn the lights down  
25 and so forth--is this small disruption or intimal

1 flap that is right there.  
2 The patient had a 0.1-centimeter fibrinous  
3 strand removed after the ascending aorta was  
4 opened, sent for pathology, closed the ascending  
5 aorta, and the patient suffered no sequelae.  
6 DR. EDMUNDS: Could you point out the arch  
7 vessel?  
8 DR. ALLEN: This is actually mid-ascending  
9 aorta, so it is beyond; the arch vessels wouldn't  
10 be seen in this particular vein. It is not  
11 scanning farther on down there.  
12 DR. EDMUNDS: That is the pulmonary artery  
13 going across.  
14 DR. ALLEN: The pulmonary artery is right  
15 here.  
16 [Slide.]  
17 Here is another example. Once again, the  
18 surgeon identified this endothelial disruption, and  
19 you see it right here. It is a little easier to  
20 see it than on the last one. This occurred a  
21 little later on in the study after we had  
22 experience from the core lab telling us that we  
23 were seeing these ultrasound abnormalities, and in  
24 this case once again, there was no clinical  
25 dissection. The patient was doing fine. And this



1 surgeon because he had been provided with some of  
2 this historical information didn't repair it. And  
3 in fact in 10 out of the 13 endothelial disruptions  
4 that were identified by surgeons, those surgeons  
5 decided not to repair it, and those patients didn't  
6 suffer sequelae from it.

7 DR. TRACY: Before we leave this slide, I  
8 think Dr. Edmunds wants some clarification on a  
9 couple of things.

10 DR. EDMUNDS: Can you point out the  
11 location or probable location of the deployed  
12 filter in relation to these so-called injuries?

13 DR. ALLEN: I can tell you that we did an  
14 analysis on the location of the endothelial  
15 disruptions. In both of these cases, the  
16 endothelial disruptions were in the mid-ascending  
17 aorta. The filter was downstream from these  
18 devices, so they weren't in the area of--

19 DR. EDMUNDS: Where are you  
20 cannulating--the sinuses of the falsalva  
21 [phonetic]? Where was the cannula, then?

22 DR. ALLEN: The cannula was approximately  
23 right at the innominate [phonetic] artery.

24 DR. EDMUNDS: So the filter is deployed  
25 upstream to the cannula tip?

1 DR. ALLEN: No, sir.  
2 DR. EDMUNDS: I am totally confused.  
3 DR. ALLEN: The filter is deployed just  
4 proximal to the cannulation. It is part of the  
5 filtering process of the cannulation itself. It is  
6 part of the cannula that goes into the aorta.  
7 DR. EDMUNDS: Yes, but when it is  
8 deployed, like a parachute, where is that in  
9 relation to the nozzle of the bypass cannula?  
10 DR. ALLEN: It is posterior to it.  
11 DR. EDMUNDS: It is proximal to the aortic  
12 cannula spigot?  
13 DR. ALLEN: Yes, and it is distal to the  
14 cross-clamp.  
15 DR. TRACY: Maybe at the end, we will  
16 review your Slide 6.  
17 DR. ALLEN: I can show you another slide  
18 of that.  
19 DR. TRACY: Okay, but let's go ahead with  
20 your presentation.  
21 [Slide.]  
22 DR. ALLEN: This is an example of an  
23 interoperative photograph that I borrowed from Dr.  
24 Banberry [phonetic] at the Cleveland Clinic in a  
25 patient who was undergoing a routine aortic valve

1 replacement in which we hypothesized what these  
2 endothelial disruptions might look like based on  
3 the one pathologic specimen that was sent and that  
4 was resected by a surgeon.

5 [Slide.]

6 These are two examples of epi-aortic  
7 scans, once again done at Dr. Kouchoukos' center,  
8 one involving a filter patient, one involving a  
9 control patient. These are based on--we asked the  
10 core lab to provide us with representative slides,  
11 and once again, without the lights turned down, it  
12 is difficult in this patient to see this, but the  
13 slight endothelial disruption here, and in a  
14 similar area, right here, that are tagged as being  
15 what we are calling these intimal injuries or  
16 endothelial disruptions.

17 Once again, these were not identified by  
18 the surgeons at the time of the operation even in a  
19 center that has a vast experience with this  
20 technology and were not repaired by the surgeons  
21 and had no clinical sequelae because of that.

22 [Slide.]

23 Well, simply telling you that they weren't  
24 repaired and that they might not impact things  
25 isn't enough for me, and I certainly wouldn't think

1 it is enough for you. And we asked that we do an  
2 analysis on was there a correlation between  
3 endothelial disruptions and adverse events.

4 If you looked at, for example, patients in  
5 the filter arm and compared those patients who had  
6 EDS with those who didn't have EDS, there certainly  
7 was not a correlation with acute adverse composite  
8 events.

9 If you similarly looked at the control  
10 patients, and one of the control patients had an  
11 endothelial disruption, or nine control patients  
12 had endothelial disruptions, and compared those  
13 with EDS to those without EDS, there certainly  
14 wasn't a correlation with EDS to an adverse event.

15 And, more importantly, then, if you just  
16 forgot whether they were randomized or not and  
17 looked at all patients who had EDS and compared  
18 adverse composite events to those that didn't have  
19 EDS, there clearly is not a correlation to EDS with  
20 adverse acute events.

21 [Slide.]

22 Are there long-term consequences of EDS?  
23 And I think, as part of the presentation, it is  
24 important to understand that long-term followup in  
25 this study was not part of the protocol. But

1 surgeons and investigators were interested in that,  
2 and we needed to know that what we were doing to  
3 our patients wasn't going to hurt them.

4 So we developed a methodology to try to  
5 follow these patients and assess the long-term  
6 impact of EDS on composite and individual event  
7 rates. As I said, this is not part of the original  
8 protocol. But in order to obtain appropriate  
9 followup in the image patients, we targeted centers  
10 that had the imaging and centers that had EDS; we  
11 looked at centers that were high enrollers in order  
12 to have less variability between arms, and we also  
13 needed to be able to get timely IRB approval for  
14 this longer-term followup.

15 [Slide.]

16 We ultimately looked at four  
17 high-enrolling centers in which 90 percent followup  
18 was obtained. We wanted, though, to look  
19 specifically at EDS, and obviously, there were some  
20 EDSs occurring outside of those four high-enrolling  
21 centers. So we wanted to get followup on all EDS,  
22 even patients who were roll-in. So we ended up  
23 trying to find followup on 58 patients.

24 Seven patients couldn't have followup.  
25 Six of those were simply because we could not get

1 appropriate IRB approval. We know the patients  
2 were alive, but we just couldn't get followup.  
3 There was one patient lost to followup.

4 [Slide.]

5 We ended up having a 360-day or almost  
6 one-year mean followup. And when you looked at  
7 composite event rates between filter and control  
8 patients out to 360 days mean, there were events  
9 occurring, but they were occurring absolutely  
10 identical at 6.1 percent between both arms.

11 [Slide.]

12 Once again, I asked the question--well,  
13 that's great, but I want to know what about the  
14 patients who have EDS. And once again, when you  
15 look at the filter patients who had EDS compared to  
16 those who did not during long-term followup, there  
17 was no correlation to an adverse outcome. And when  
18 you look at all patients, once again, there wasn't  
19 a correlation during long-term followup.

20 [Slide.]

21 I come back to the issue of aortic  
22 dissections and the development of aneurysms,  
23 because while I tell you that long-term followup  
24 did not correlate with an acute composite event,  
25 what about the development of a late dissection or

1 the development of an aneurysm.

2 I told you that in the study, two  
3 dissections occurred acutely in the control arm and  
4 none in the filter. During followup, no patients  
5 were operated on for the development of acute  
6 dissections in either arm. There were three  
7 additional aneurysms that were seen in control  
8 patients--two thoracic aortic aneurysms, one that  
9 was repaired, and one abdominal aortic aneurysm  
10 that was also repaired. Obviously, none of these  
11 were in areas where EDSs were identified, and in  
12 fact none of these three patients even had EDS.

13 [Slide.]

14 So from a summary standpoint--and I think  
15 it is an important safety issue, and that is why we  
16 have spent time on this--this was primarily a  
17 finding on aortic imaging. Seventy-eight percent  
18 of surgeons, despite using sophisticated epi-aortic  
19 scanning, were not able to identify these  
20 endothelial disruptions. They should not be  
21 classified as clinically significant aortic  
22 injuries, and while they were seen more frequently  
23 in the filter arm, they were seen in both arms.

24 These are not aortic dissections, and I  
25 think our acute data demonstrates no correlation

1 acutely to composite events, and I think our  
2 long-term data and due diligence in collecting that  
3 also does not demonstrate a safety issue.

4 [Slide.]

5 We had a second primary endpoint which, as  
6 Dr. Kuntz pointed out, was an effectiveness  
7 endpoint. The hypothesis was that we could capture  
8 particulate emboli in greater than 75 percent of  
9 the EMBOL-X aortic filters. And in that case,  
10 successful emboli capture was defined as retrieved  
11 particles observed at 10 times power before  
12 histologic processing.

13 [Slide.]

14 And I think as these photographs  
15 demonstrate, our primary effectiveness endpoint was  
16 indeed met. Ninety-six-point-eight percent of  
17 filters prior to histologic processing visualized,  
18 documented, and photographed captured particles.

19 [Slide.]

20 As surgeons, we were interested in what  
21 these particles might be composed of, and as part  
22 of the trial, we employed a pathologic core lab  
23 that would analyze this data.

24 We anticipated, based on what the core lab  
25 was telling us, that because this small amount of



1 tissue was going to go through extensive histologic  
2 processing and handling that there were going to be  
3 some specimens lost, displaced or dissolved and not  
4 available for analysis. And indeed that is what we  
5 found.

6 We found that approximately 21 percent of  
7 previously photographed and visualized specimens  
8 were not ultimately available for histologic  
9 analysis. However, more than 85 percent of  
10 specimens that were available for analysis  
11 demonstrated that the material atheromatous in  
12 nature.

13 [Slide.]

14 There were various other things captured  
15 by the filters, and I show as this example RBC  
16 thrombus or clot, this polyploid structure in both  
17 specimens that, based on its organized nature on  
18 pathology, likely came from an intercavitary  
19 source, as Dr. Kouchoukos mentioned in his talk.

20 [Slide.]

21 It is important that we assure the panel  
22 that this device isn't causing what it captures and  
23 that it is not thrombogenic. And we did extensive  
24 bench-testing to demonstrate that with filters  
25 having the 95 percent heparin bonding with a 2-hour

1 dwell time, they were not thrombogenic. But what  
2 about in a human?

3 [Slide.]

4 And to assess for thrombogenicity, we used  
5 scanning electron microscopy, and it is important  
6 to note when the SEM was done.

7 Unlike the histologic data, which was done  
8 obviously after processing, we looked at filters  
9 prior to histologic processing. The original  
10 intent of the study was to evaluate 10 percent of  
11 the filters, but after 5.6 percent of filters had  
12 been examined and presented to the FDA, the FDA  
13 agreed in a letter on October 12, 2001 that the  
14 scanning electron microscopy demonstrated no  
15 significant platelet thrombus formation, and we  
16 could discontinue doing additional scanning EMs.

17 [Slide.]

18 So from a summary standpoint, there  
19 certainly were captured particles that were  
20 documented and visualized as part of our a priori  
21 defined effectiveness endpoint that were not  
22 available for histologic analysis. But I contend  
23 that they are not available because of the  
24 extensive histologic processing that went on and  
25 the small amount of material that they represented,

1 and that scanning electron microscopy does not  
2 demonstrate that filters are thrombogenic.

3 [Slide.]

4 What about the number of particles  
5 captured? I think Dr. Kuntz touched on that by  
6 what does it mean whether we capture one particle,  
7 five particles, or 20 particles.

8 The study data demonstrate that there was  
9 a mean number of particles captured of 5.6.

10 [Slide.]

11 I think it is more interesting to look at  
12 the quantity of the particles that were captured,  
13 and here, you see a slide looking at the  
14 distribution of sizes of particles captured, and we  
15 superimposed the previous slide that Dr. Kouchoukos  
16 showed you of representative arteries such as  
17 middle cerebral artery branch or posterior  
18 circulation artery. And we also put in renal  
19 interlobular arteries versus the size of renal  
20 arcuate arteries, or intralobular renal arteries.

21 And you see that the vast majority of the  
22 size of the particles that we captured are filtered  
23 or would be filtered by small arcuate or cortical  
24 cerebral arteries.

25 [Slide.]

1           So from a clinical study overview, we feel  
2   that we have met both of our primary endpoints  
3   successfully--that effectiveness was demonstrated  
4   in that 96.8 percent of filters did capture emboli  
5   as documented under 10 times magnification; and we  
6   certainly feel that our safety endpoint of  
7   equivalence was met and that no clinical adverse  
8   events were associated with the findings of  
9   epi-aortic scanning.

10           I think I reemphasized the point that our  
11   study was not powered nor was it designed to  
12   demonstrate superiority in this low-risk patient  
13   which was specifically selected to demonstrate  
14   safety.

15           [Slide.]

16           Following the completion of the study, the  
17   FDA, because this was a safety equivalency study,  
18   asked the company as well as the investigators  
19   whether study data could be extrapolated to  
20   clinical efficacy, and this was one of the  
21   questions that they had actually posed to the  
22   panel.

23           And I think it is important that we in an  
24   attempt to answer this question did some additional  
25   analysis. But this additional analysis is in no

1 way claims for labeling, and we don't make claims  
2 of superiority. It is simply an attempt to answer  
3 those questions raised by the FDA.

4 From a surgeon's standpoint, we know that  
5 clinical outcomes are influenced by preoperative  
6 risk variables, so we felt that if you looked at  
7 the high-risk patients in our population, could we  
8 extrapolate some clinical benefit or efficacy.

9 [Slide.]

10 We utilized the Cleveland Clinic score, as  
11 published by Dr. Higgins in JAMA in 1992, to assess  
12 for preoperative risk. We specifically chose the  
13 Cleveland Clinic score because it looks at both  
14 morbidity and mortality unlike, for example, the  
15 STS or the New York State Index, which only look at  
16 mortality.

17 The Cleveland Clinic score which is  
18 utilized at our center is a validated preoperative  
19 risk score that has been validated in over 9,000  
20 patients, and a score of 5 or higher has been  
21 validated for increased morbidity and mortality.

22 Eighteen-point-seven percent of the 1,289  
23 patients randomized in this study met the criteria  
24 for moderate to high risk.

25 [Slide.]

1           When we looked at whether high- or  
2 low-risk patients had composite events, you will  
3 see that when we control filter to control in  
4 low-risk patients, a Cleveland Clinic score of zero  
5 to 4, there was absolutely no difference between  
6 the two groups.

7           But when we looked at patients who had  
8 Cleveland Clinic scores greater than 5 as defined  
9 in that Higgins paper, we saw that the trend  
10 certainly favored statistically patients who  
11 received the filter.

12           [Slide.]

13           When we looked at components of that  
14 composite event and compared those to patients who  
15 had Cleveland Clinic scores greater than 5, across  
16 a broad range of all unselected components of the  
17 composite endpoints, there weren't statistical  
18 differences between groups except for renal  
19 insufficiency.

20           I think it is important that we break it  
21 out as dialysis patients, patients without  
22 dialysis, or all patients who had renal  
23 insufficiency as defined in our study. And when  
24 you look at all patients with renal insufficiency,  
25 it was significantly less in the filter patients.

1 [Slide.]

2 Now, renal insufficiency is not as sexy  
3 and glamorous as preventing frank stroke, but I  
4 think renal insufficiency has a significant  
5 clinical impact. We know that based on the STS  
6 data base, renal insufficiency is a predictor of  
7 increased morbidity and mortality postoperatively,  
8 but when you look at patients who had renal events,  
9 their length of stay was a little under 15 days.  
10 For those patients who did not have renal  
11 insufficiency in this study, their length of stay  
12 was 7.2 days.

13 And this looks at all patients with renal  
14 insufficiency. The data isn't simply driven by  
15 those who had dialysis. If you take out the  
16 dialysis patients, which you would expect to have  
17 an even longer length of stay, the length of stay  
18 only drops by one day.

19 [Slide.]

20 I think this slide has put it in  
21 perspective for a lot of the investigators in the  
22 study, because it helps us evaluate the  
23 risk-benefit of this device. It is an odds ratio  
24 comparison of baseline variables that, as surgeons,  
25 we know are predictive of increased postoperative

1 morbidity and mortality.

2       You can see the middle line demonstrating  
3 no benefit one way or the other; to the left, the  
4 filter arm tends to do better; to the right, the  
5 control patients tend to do better.

6       Once again, we don't make claims of  
7 superiority, but it is interesting when you look at  
8 this odds ratio table that the vast majority of  
9 events are mitigated by placement of a filter.

10       [Slide.]

11       This additional analysis to address the  
12 possible clinical efficacy I think does demonstrate  
13 that at least in moderate- to high-risk patients,  
14 there may be a benefit--and I underline "may."  
15 Captured particles were predominantly of the size  
16 associated with cerebral and renal cortical  
17 arteries. The reduction in renal insufficiency  
18 events can be demonstrated with a sensitive marker  
19 in those high-risk patients. But the study was not  
20 designed to assess for neurocognitive dysfunction,  
21 and we make no claims for that.

22       [Slide.]

23       In summary, our study objective was to  
24 demonstrate that particular captured could be  
25 safely accomplished in lower-risk populations



1 against the backdrop of known detrimental effects  
2 of particulate emboli.

3 The risks of this device I think have been  
4 safely assured by the equivalency of our composite  
5 safety endpoint. Across this large prospective  
6 study, serious adverse events were identical  
7 between groups.

8 Certainly, using epi-aortic imaging, we  
9 did demonstrate an increased incidence of these  
10 endothelial or intimal disruptions, but I think  
11 that certainly there was no acute correlation to  
12 adverse events, and our due diligence to try to  
13 provide you with some long-term clinical followup  
14 hopefully provides that there is not a correlation  
15 with long-term events.

16 The benefits of this are that particulate  
17 capture was clearly demonstrated in 97 percent of  
18 filters, and we feel that the additional analysis  
19 asked for by the FDA does demonstrate that clinical  
20 efficacy can be reasonably extrapolated from  
21 particulate capture.

22 MS. CHANG: I would like to thank Dr.  
23 Allen, Kouchoukos, and Dr. Kuntz for the  
24 presentation today.

25 We believe that the study design and

1 clinical [inaudible] supports the following  
2 indications for the study, and, Dr. Edmunds, if you  
3 would like me to answer your earlier question, I  
4 can do that for you, as to location of the  
5 epi-aortic imaging.

6 DR. EDMUNDS: I think you have done that.

7 MS. CHANG: Okay.

8 DR. EDMUNDS: I have to say I missed  
9 it--it was in the writeup, but your diagram was a  
10 little misleading, or at least I didn't interpret  
11 it right.

12 MS. CHANG: Sorry.

13 Questions and Answers

14 DR. TRACY: Okay. I would like to ask the  
15 panel if they have any brief clarifying questions  
16 for the sponsor at this time. This is not the open  
17 committee discussion, but just clarifying  
18 questions.

19 Yes?

20 DR. MARLER: I was interested in the  
21 details of the discussion that led to the  
22 conclusion that there were not cognitive outcomes  
23 that could be used, and who participated.

24 DR. KUNTZ: In the discussion with whom?  
25 With the FDA?

1 DR. MARLER: You mentioned that you  
2 reached the conclusion that there was no outcome  
3 that could be agreed upon, and I was just wanting  
4 to hear some more details about that and who  
5 couldn't agree.

6 DR. KUNTZ: All right. I may refer this  
7 to [inaudible] who is an expert in this area, but  
8 we reviewed--I think the Stump [phonetic] battery  
9 of criteria at that time, which I think was  
10 probably the best candidate overall, but they were  
11 a collection of approximately 6 to 12 instruments  
12 and batteries, and at that time, we didn't think  
13 there was a consensus or a study that had  
14 demonstrated or validated that those outcomes could  
15 be correlated with changes in cognitive dysfunction  
16 in 1999, and we were not aware of any validations  
17 of that emerging battery of tests, which I think is  
18 being refined and is probably a good set, but at  
19 that time, it was difficult to say if there was  
20 consensus.

21 Maybe I could call up two other people who  
22 may want to make some comments.

23 Dr. Gold?

24 DR. GOLD: Thank you.

25 My name is Jeff Gold, and I am a cardiac

1 surgeon from New York, and I do indeed own some  
2 stock options in EMBOL-X and have been involved in  
3 a study of neurologic and cognitive function  
4 associated with cardiac surgery for perhaps longer  
5 than I care to remember sometimes.

6 The definition of cognitive changes  
7 associated with cardiac surgery has been a complex  
8 and moving target for an extremely long time.  
9 There have been, as I am sure you are all aware, at  
10 least one consensus panel and several others  
11 looking at development of a battery of tests.

12 However, at the time that this study was  
13 conceived, not only was there not a defined battery  
14 of tests that cardiac surgeons across the board  
15 could agree upon, let alone psychometricians and  
16 neurologists, but the etiology of cognitive  
17 abnormalities was also highly controversial at the  
18 time.

19 You might recall a very interesting study  
20 published by a Dr. Rousseau, who looked at the  
21 incidence of cognitive function abnormalities in  
22 patients undergoing total knee replacement surgery  
23 under local anesthesia. The incidence of cognitive  
24 abnormalities in that 524-patient cohort was  
25 exactly equal to an equivalent study done in

1 patients undergoing coronary bypass surgery.

2 So perhaps the art has improved. Perhaps  
3 if we were to redo this at a time in the future, we  
4 could agree upon a panel of tests. But if you were  
5 to ask about the significance of cognitive  
6 abnormalities among practicing cardiac surgeons  
7 today, and our ability to reliably demonstrate  
8 them, I would say they are poor.

9 DR. TRACY: Dr. Aziz?

10 DR. AZIZ: From what I understand, the way  
11 that you place this catheter, this could not  
12 protect against the sandblast effect of emboli  
13 being dispersed; is that right?

14 DR. ALLEN: Yes, sir.

15 DR. AZIZ: Okay. And secondly, so this  
16 cannula has to be inserted proximal to the enormid  
17 [phonetic] artery; is that right?

18 DR. ALLEN: The cannula is inserted  
19 identically. There is no difference in what you  
20 do. If you were to, for example, cannulate the  
21 arch, as we sometimes have to, yes, you wouldn't  
22 use this cannula in somebody, for example, that you  
23 were going to cannulate the mid-arch.

24 DR. TRACY: Dr. Laskey?

25 DR. LASKEY: Dr. Allen, these

1 ECO-abnormalities detected in the core lab, I am  
2 assume that there are pre-/post-placement ECOs. is  
3 this a serial analysis where the ecocardiogram is  
4 obtained only at one time, or was there some  
5 protocol that people adhered to where you had  
6 baseline ECOs and then ECOs during place and ECOs  
7 after placement? How was that done?

8 DR. ALLEN: That's a very good question.  
9 It was quite laborious to do this. Actually, our  
10 center did not do imaging. I will let Dr.  
11 Kouchoukos' center answer that, because his study  
12 was actually doing epi-aortic scanning.

13 DR. KOUCHOUKOS: The epi-aortic scanning  
14 was performed after the pericardium was opened,  
15 before instituting cardiopulmonary bypass, and at  
16 the completion of the procedure, after removal of  
17 all the cannulas, administration of protamine,  
18 another scan was performed in the longitudinal and  
19 transverse planes for the whole ascending aorta.

20 DR. EDMUNDS: Nick, could you explain how  
21 this was deployed? You say you didn't leave the  
22 filter up for 60 minutes, and you deployed it when  
23 you first started to manipulate and got the cannula  
24 in so you could deploy it.

25 DR. KOUCHOUKOS: Yes. The cannula was

1 deployed immediately before removal of the aortic  
2 cross-clamp--in other words, after completion of  
3 the proximal and distal anastomosis and the bypass  
4 operation or closure of the aorta or the left  
5 mitral valve replacement--

6 DR. EDMUNDS: It was removed.

7 DR. KOUCHOUKOS: --it was removed. And  
8 then it was left in place generally until the  
9 protamine was administered. And the safety  
10 analyses indicate it could be left for an hour, but  
11 it was left on average for 20 minutes.

12 DR. EDMUNDS: So you put the aortic  
13 cross-clamp on without the filter deployed?

14 DR. KOUCHOUKOS: That's correct. The  
15 filter was deployed into the aortic cannula  
16 immediately before removal of the cross-clamp.

17 DR. EDMUNDS: So you are not claiming that  
18 you got all the emboli; you just got what you got.

19 DR. KOUCHOUKOS: Well, based on the  
20 analyses that I presented, the majority of these  
21 emboli are released at the time of removal of the  
22 cross-clamp, and that was the logic for deploying  
23 the filter immediately before release of the clamp.

24 DR. EDMUNDS: Well, you did the studies,  
25 but on Dick Clark's study, I thought you got a

1 shower of emboli when you put the clamp on, also;  
2 is that not correct?

3 DR. KOUCHOUKOS: Well, I think I showed  
4 you a slide showing that there are emboli released  
5 at various times during the course of a cardiac  
6 surgical procedure, but that the majority of the  
7 emboli are released at the time of release of the  
8 cross-clamp.

9 DR. EDMUNDS: Oh, I agree with that, yes.

10 DR. TRACY: Okay. Are there any other  
11 clarifying questions? We'll have the open  
12 committee discussion in a few minutes after the FDA  
13 presentation, but are there any other clarifying  
14 questions?

15 DR. KRUCOFF: Just a plumbing question.  
16 Relative to, say, any other commercial cannula that  
17 you would use routinely, is this cannula different  
18 from a flexibility or a dimensional perspective?

19 MS. CHANG: Let me show a picture of the  
20 product again.

21 [Slide.]

22 MS. CHANG: Actually, this main body is  
23 what a standard cannula looks like, and we have  
24 modified it so that we have added a sideport there.  
25 There is still only one hole, so it is actually



1 virtually identical to existing commercial  
2 cannulas. So the big difference is the sideport.  
3 DR. KRUCOFF: Okay. So my question is  
4 does the presence of the sideport affect in any way  
5 the portion of the cannula that actually goes  
6 through the aorta relative to a commercial, in  
7 either dimension, or just how it feels?  
8 DR. ALLEN: The short answer is "No," and  
9 I think Dr. Kouchoukos would concur with that.  
10 DR. TRACY: Okay. One more question.  
11 DR. DeMETS: I would like to ask what  
12 about the randomization process. You didn't  
13 describe it in your presentation, but your writeup,  
14 as I understand it, there were some patients who  
15 did not get treated as randomized.  
16 Could you walk me through that process so  
17 I can understand exactly what happened and what you  
18 did about it?  
19 DR. ALLEN: The specific details--the  
20 patient was obviously net inclusion and exclusion  
21 criteria. The one exclusion criterion that  
22 couldn't be determined until you actually got in  
23 the operating room was whether his ascending aorta  
24 was of an appropriate size. So the patient had a  
25 stronotomy [phonetic], and then you measured the

1 ascending aorta approximately 1.5 centimeters below  
2 the nominid [phonetic] artery, which is about where  
3 the filter would be deployed, and if it fell into  
4 an appropriate range, then a randomization card was  
5 opened, and the patient then was either randomized  
6 to have a standard J-tip cannula or a modified  
7 J-tip cannula inserted.

8         There were nine patients who, when the  
9 card was opened early on in the study, the way the  
10 card read at the top was "The EMBOL-X Study," and  
11 actually, the very first patient that I randomized,  
12 our coordinators read the patient was in the  
13 EMBOL-X study, and they assumed he was a filter  
14 patient when indeed, you had to read the line below  
15 it, which said whether he was a control or a filter  
16 patient.

17         Those happened very early on in the study,  
18 and once they were educated about it, they ceased  
19 to happen. And in fact, in three of those nine, we  
20 actually caught the mistake before we actually put  
21 one device or the other in.

22         Does that clarify that for you?

23         DR. DeMETS: That clarifies the first  
24 part. The second part is given that that happened,  
25 which I now understand how it happened, in which

1 group were the patients left?  
2 DR. ALLEN: Intent to treat.  
3 DR. DeMETS: So they were left in the  
4 group that they should have been randomized to?  
5 DR. ALLEN: Correct. And they represent a  
6 very small proportion of the number of patients  
7 that we put in the study, but we did an intent to  
8 treat analysis.  
9 DR. DeMETS: Thank you.  
10 DR. TRACY: Dr. Ferguson?  
11 DR. FERGUSON: I missed one point about  
12 the cannula. I thought from my reading that you  
13 used the cannula with the sideport even in the  
14 control group. That is not the case?  
15 DR. ALLEN: No. The standard cannula is  
16 the cannula that I trained with in Chicago, which  
17 is just a standard J-tip cannula. So essentially,  
18 if you take the--  
19 DR. FERGUSON: Is that a cannula  
20 manufactured by this company?  
21 DR. ALLEN: No. That's a standard J-tip  
22 cannula.  
23 DR. FERGUSON: Did everybody use--  
24 DR. ALLEN: Everybody.  
25 DR. FERGUSON: --were they instructed to

1 use that same cannula?  
2 DR. ALLEN: Yes. Everybody used the same  
3 cannula. So all control patients, regardless of  
4 what your standard cannula was at your site, you  
5 had to use the same standardized cannula.  
6 DR. FERGUSON: And I would ask again, if I  
7 may, when the standard cannula and then this  
8 cannula are affixed to the aorta, the aortic size  
9 and so forth at that point of entry into the aorta  
10 are both the same size?  
11 DR. ALLEN: It is identical. The only  
12 part that is inside the ascending aorta is right  
13 there.  
14 DR. FERGUSON: I understand that, but I  
15 just want to be sure that the impact of the sidebar  
16 does not enlarge that--  
17 DR. ALLEN: No, sir, it doesn't. That's a  
18 great question, but no, sir, it doesn't.  
19 DR. TRACY: Mr. Morton?  
20 MR. MORTON: Madam Chair, does the sponsor  
21 have an example available, and would you mind if  
22 the panel could see it?  
23 DR. TRACY: No. Are we allowed to do  
24 that?  
25 DR. ZUCKERMAN: Yes. We can take a look

1 at an example.  
2 DR. TRACY: Okay. Do you have that  
3 available?  
4 DR. ALLEN: [Handing.] I'll give it to  
5 the cardiac surgeons first.  
6 [Laughter.]  
7 DR. EDMUNDS: Can I ask a question? Is  
8 there any connection between the flow path from the  
9 pump and the deployment path of this filter?  
10 DR. ALLEN: If I'm understanding you--  
11 DR. EDMUNDS: In other words, this sort of  
12 filter deployment apparatus is just riding shotgun  
13 on the cannula. There is really no hole between  
14 the two.  
15 DR. ALLEN: Yes.  
16 DR. EDMUNDS: You vent the air out by  
17 blood coming around the wire that is around the  
18 filter.  
19 DR. ALLEN: Actually, that's a great  
20 question, and it involves a safety issue with how  
21 the air is vented.  
22 The filter--I don't know if we actually  
23 have a filter to show you--there is a plug, much  
24 like is on the standard cannula, which allows air  
25 to be vented. So when you put the filter in, you

1 see the white plug turn red, indicating that blood  
2 has come up and evacuated in the air. And  
3 obviously, if you don't see that, you need to  
4 change filters or do something differently.

5 DR. MARLER: I had a question on page 68,  
6 Tables 7-18 and 7-19--and I'm sure there is an  
7 explanation, but I just didn't understand why there  
8 were adverse events, NIH Stroke Scale greater than  
9 4, 13 in the control group, with 644 patients--

10 DR. TRACY: I think this may actually be  
11 more appropriate for the open committee discussion.  
12 Unless there are some very brief clarification  
13 questions, I'd like to stop at this point for a  
14 break.

15 DR. EDMUNDS: Where does the damn thing  
16 come out? Does it come out this hole or some other  
17 hole?

18 [Dr. Allen handing sample to Dr. Edmunds.]

19 DR. EDMUNDS: Why don't you show  
20 everybody, because I can't be the only one  
21 confused.

22 DR. ALLEN: The cannula is inserted as you  
23 would any other cannula; the cannulation is no  
24 different. Once the cannula is inserted, patients  
25 are put on cardiopulmonary bypass, everything is

1 done--you manipulate the heart, you do your  
2 proximals, you do your distals, and so forth. Just  
3 before you release the cross-clamp, you take out  
4 this plug, and you insert the filter, and the  
5 filter goes in like this.

6 The air venting that you alluded to, for  
7 those of you--and I will pass this around--there is  
8 a white hemostatic pump that allows fluid to vent  
9 out and push the air out just like our cannulas do  
10 today. Once you confirm that it is vented  
11 appropriately, the device is deployed, just like  
12 this.

13 DR. EDMUNDS: Why don't you pass that  
14 around?

15 DR. ALLEN: And then the cross-clamp is  
16 released.

17 DR. TRACY: While that thing is making its  
18 way around--Dr. Aziz?

19 DR. AZIZ: So during the time that you  
20 take the stop off to put the actual filter in,  
21 could air get in, or is there a one-way valve that  
22 is--

23 DR. ALLEN: No. It's one-way. It is a  
24 one-way valve.

25 DR. TRACY: Okay. I think at this point,

1 while that is working its way around, we'll take a  
2 15-minute break and resume at approximately 20 of  
3 11.

4 DR. ALLEN: Thank you.  
5 [Break.]

6 DR. TRACY: I'd like to reconvene the  
7 meeting at this point and ask the FDA to begin  
8 their presentation.

9 FDA Presentation

10 MS. WENTZ: Good morning. My name is  
11 Catherine Wentz, and I'll be opening up the FDA  
12 presentation for the EMBOL-X aortic filter.

13 This will be done in four parts. I will  
14 do an introduction. Dr. Julie Swain will follow up  
15 with her clinical summary. Dr. Gerry Gray will  
16 then do his statistical summary, and I will then  
17 close with the questions to the panel.

18 [Slide.]

19 I'll start with a brief description.  
20 EMBOL-X gave you a better one than I can, but this  
21 is just a short reiteration.

22 The EMBOL-X aortic filter is used in  
23 conjunction with the EMBOL-X aortic cannula which  
24 was cleared this past September and is "intended to  
25 contain and remove particulate emboli from the



1 ascending aorta during and following cross-clamp  
2 removal and as the heart resumes ejection."

3 The heparin-coated filter has a pore size  
4 of 120 microns and is mounted on a nitinol frame.  
5 The filter is inserted into the ascending aorta via  
6 a sideport on the EMBOL-X cannula. The flexible  
7 wire filter frame expands upon insertion into the  
8 vessel and is available in five sizes. The filter  
9 is then retracted back through the same sideport at  
10 the end of the procedure.

11 [Slide.]

12 In the next three slides, I would just  
13 like to reiterate briefly some regulatory  
14 information that I think you all received in your  
15 training this morning.

16 I would also like to reiterate that this  
17 is just for your information and should not enter  
18 into the discussion of the EMBOL-X study. It will  
19 be FDA's responsibility to take the recommendations  
20 made today at the panel meeting to make a final  
21 decision within the 510(k) realm.

22 So just to reiterate some definitions, the  
23 510(k) requires a manufacturer to demonstrate  
24 substantial equivalence or SE to a legally marketed  
25 predicate device.

1           To further define substantial equivalence,  
2 substantial equivalence basically means that the  
3 two devices have the same intended use, similar  
4 technology, and if the technology is not similar,  
5 there are means by which to demonstrate that the  
6 new technology does not affect equivalent  
7 performance or the risk profile.

8           [Slide.]

9           A PMA is defined as a process where the  
10 FDA evaluates Class III medical devices. Class III  
11 devices are usually those that support or sustain  
12 human life, are of substantial importance in  
13 preventing impairment of human health or which  
14 present a potential unreasonable risk of illness or  
15 injury.

16          [Slide.]

17          Now, to put this submission into that  
18 context, the EMBOL-X aortic filter originally  
19 underwent a clinical study to demonstrate the  
20 safety and effectiveness of the device in support  
21 of a PMA application.

22          However, in June of 2001, the PercuSurge  
23 device, which is also an embolic protection device,  
24 was cleared through the 510(k) regulatory pathway  
25 opening the doors for the EMBOL-X aortic filter to

1 be reviewed under the 510(k) regulations.

2 The PercuSurge device, which has a similar  
3 intended use to EMBOL-X, in conjunction with  
4 cardiopulmonary bypass arterial line blood filters,  
5 which has similar technology to the EMBOL-X device,  
6 will be used as a combination predicate for the  
7 EMBOL-X device in the determination of substantial  
8 equivalence under the 510(k) regulations.

9 [Slide.]

10 To go over a little bit of the history of  
11 how the endpoints for this study were developed, at  
12 the beginning, the sponsor wanted a nonclinical  
13 effectiveness endpoint--that is, to capture  
14 debris--and an equivalence safety study. The FDA  
15 consistently expressed concerns regarding the  
16 interpretability of the proposed endpoints.

17 FDA, however, agreed to the proposed  
18 effectiveness endpoint assuming that the safety  
19 endpoint, which included some neurologic outcomes  
20 and other embolic-related events, would capture the  
21 clinical effectiveness of the device; and that the  
22 device labeling would be restricted to only the  
23 facts from the study. No clinical implications  
24 could be made from the capture of debris since none  
25 was evaluated.

1 [Slide.]

2 And briefly, just one slide of  
3 engineering--that is my background; I performed the  
4 engineering review of the submission.

5 Overall, on the bench studies, there were  
6 some design concerns and/or test method concerns  
7 that remain that may be related to the endothelial  
8 injuries observed with this device. These concerns  
9 are presently being addressed.

10 Both the biocompatibility and  
11 sterilization, packaging, and shelf-life had no  
12 further questions; they were all fine.

13 I think this is the point where I turn it  
14 over to Julie for her clinical review.

15 DR. SWAIN: Thank you for the opportunity  
16 to present.

17 Let me make a comment first, that I am at  
18 somewhat of a disadvantage in that we traditionally  
19 exchange presentations with the sponsors before the  
20 talk so we can mold our presentations, and we  
21 provided our slides to the sponsor, and the sponsor  
22 chose not to provide theirs, so I think that some  
23 of the comments that I will make are a little bit  
24 off-the-cuff in response to some of the  
25 presentations that I had no knowledge that these

1 items were going to be presented.

2 [Slide.]

3 The clinical review was done by both Wolf  
4 Sapirstein and myself, and we are both  
5 cardiothoracic surgeons. I am a consultant to the  
6 FDA.

7 [Slide.]

8 The study design, as you have seen, was  
9 randomized, which is important in I think one of  
10 the discussions that we will have about the  
11 neuropsych--it is a randomized, multicenter trial  
12 and one of the largest trials done--control arm,  
13 patients without filter; and an interim data  
14 analysis was planned at 50 percent of the patients.

15 [Slide.]

16 In the study plan, it was said that "If  
17 the hypothesis tests performed at the interim are  
18 statistically significant, indicating emboli  
19 capture and equivalent safety, the study will be  
20 terminated."

21 However, the study was continued with the  
22 attempt to show safety superiority, and that is  
23 some of the data that we will discuss.

24 [Slide.]

25 Inclusion/exclusion criteria are patients

1 with elective operations; isolated coronary bypass  
2 or valve; greater than age 60; and there were a  
3 total of 24 exclusion criteria. One was  
4 neurological deficit; one was a history of major  
5 stroke as defined by the clinical history of a  
6 fixed, focal neurological deficit attributable to  
7 stroke; redo operations; and renal failure on  
8 dialysis.

9 [Slide.]

10 Neurological evaluation was essentially  
11 gross neurologic testing--history, physical exam;  
12 NIH stroke score; and no neuropsychological  
13 testing. And I have to say that I disagree with  
14 some of the comments made. Dr. Kuntz was talking  
15 about reviewing this 3 years. I have spent  
16 probably a quarter of a century as my major  
17 interest in the neurological effects on cardiac  
18 surgery, and the consensus conference, the Key West  
19 Conference in 1995, published in the Annals of  
20 Thoracic Surgery when Dr. Ferguson was the editor,  
21 and then the updates published when Dr. Edmunds is  
22 now the editor, listed the problems and the  
23 suggestions of the tests that could be done. And  
24 the comment was made that cardiac surgeons now  
25 still don't agree.

1           There was a very nice conference about 5  
2 months ago, sponsored by the NIH, where the leaders  
3 in the field were invited to discuss this problem.  
4 I attended the conference, and we essentially had  
5 subgroups, and the subgroup that I was in was a  
6 neuropsych group--I don't believe anyone else in  
7 this room that I recognize was at that conference  
8 in that particular area--and I disagree that there  
9 were no neuropsychological tests 3 years ago and  
10 that there is none now.

11           And again, this is a randomized study, so  
12 we know that a lot of things cause changes after  
13 cardiac surgery or knee operations or whatever--but  
14 that's the beauty of a randomized study, that one  
15 can then look at the changes.

16           And when you look at a device that perhaps  
17 you have difficulty showing efficacy, it may be  
18 that it is not efficacious or that you didn't  
19 measure the most sensitive measures. And that may  
20 be relevant to the discussion here.

21           [Slide.]

22           What are the endpoints? Efficacy is that  
23 greater than 75 percent of the filters would  
24 capture at least one particle. And there was a  
25 composite primary safety endpoint composing several

1 items, essentially saying that it wouldn't be worse  
2 than normal cardiac surgery, more than 5 percent  
3 worse.

4 One of the secondary safeties was aortic  
5 injury, which was not really part of the composite  
6 safety primary endpoint.

7 [Slide.]

8 On patient demographics, there are really  
9 no statistical differences in the baseline. It was  
10 a well-randomized study. The treatment group was  
11 73 percent male, 91 percent Caucasian, an average  
12 age of 71, and 84 percent of patients had an  
13 isolated coronary bypass operation.

14 [Slide.]

15 The composite safety endpoint comprised  
16 the items that are seen here. Renal were an  
17 elevation of creatinine, and then, a new dialysis  
18 requirement. Neurological divided into stroke,  
19 TIAs, nonmetabolics. Cardiac is Q-wave MI and  
20 non-Q-wave MI.

21 [Slide.]

22 We looked at the number of particles  
23 trapped, and the average was, I believe, 5.6 mean  
24 particles per filter. The problem is--a  
25 denominator has been mentioned by the panel



1 members--you don't know how many particles are  
2 liberated--there was no middle cerebral doppler or  
3 carotid measurements made, and that probably  
4 wouldn't help a lot.

5 We also know that particulate matter and  
6 gaseous emboli are the two main causes of  
7 neurological dysfunction. This does nothing to  
8 gaseous; we are talking about particulate, as Dr.  
9 Kouchoukos said.

10 [Slide.]

11 When you look at the maximum number  
12 trapped in some filters, it was 25 in the regular  
13 study and 38 in the roll-in patients, so maybe that  
14 gives you an idea of what the denominator is,  
15 because you would love to know the percentage of  
16 particles trapped, but you really can't know that  
17 information. That may also have an influence on  
18 the efficacy of this device or the clinical  
19 utility.

20 [Slide.]

21 I picked out just selected events--in  
22 neurological, I picked out stroke; in renal  
23 failure, I picked out dialysis; in MIs, I picked  
24 out Q-wave MIs. And I think as the comment was made  
25 yesterday by Dr. Kuntz, you like to see a trend of

1 everything in the same direction, and you really  
2 don't see a lot of trends here in these events.

3 Now, the comment was made by two of the  
4 three speakers, and certainly not Dr. Kouchoukos,  
5 that MI is irrelevant here. Well, it is highly  
6 relevant--it is actually the first adverse event  
7 that I would think of.

8 When you have that--I don't know where  
9 that device is--but the filter is between the  
10 outflow cannula and the aortic cross-clamp. Well,  
11 we all know that you get retrograde flow in the  
12 aorta. You get it in many instances; that's how we  
13 close aortic valves and get aortic  
14 insufficiency--but in cardiac surgery, when you  
15 take an aortic cross-clamp off, you've got an  
16 essentially normally pressured aorta, or we drop  
17 the pressure transiently, and you've got very often  
18 a sucked-on aorta that has collapsed with very  
19 minimal pressure, so you always get retrograde  
20 flow.

21 So I would look at myocardial infarction  
22 as a physiologist as being one of the more  
23 interesting adverse events in this. And then, the  
24 first branch of the aorta is the coronary artery,  
25 the second branch is the neuro-feeding vessels. So

1 those are the areas that I tend to be most  
2 concerned with.

3 [Slide.]

4 We looked at--I use the term  
5 "manipulation-related aortic injury." It is used  
6 in the literature, and in fact it is used in the  
7 literature quoted by the sponsor. So when you look  
8 at manipulation-related aortic injury, we can see  
9 that there are changes, and we really don't know  
10 what this means, as the sponsor pointed out.

11 In this acute study, the patients were  
12 followed an average--a median followup was 7.0  
13 days--it was during their hospitalization--and what  
14 that implies to whether you find an injury.

15 This occurred in 9.2 percent of the filter  
16 patients--42--in the regular study. Three of the  
17 filter patients, as has been said, required aortic  
18 repair. I know that the dissection ones that the  
19 control group had were the only ones placed up  
20 there in the table. However, three of the  
21 patients, the surgeons did choose to do an aortic  
22 repair, which is a fairly major procedure to add  
23 onto an isolated coronary bypass. Whether that was  
24 a correct decision or not, it is the data.

25 [Slide.]

1           The study protocol showed 30-day followup  
2 or hospital discharge, whichever occurred first; so  
3 the median was 7.0 days to study these patients.

4           And then there is the post hoc study that  
5 was a telephone followup that followed up 43 of the  
6 49 aortic injury patients and 18 of them were  
7 followed for greater than one year by telephone  
8 followup.

9           So I think it was a good effort by the  
10 sponsor to see if you could find in this type of  
11 study whether there were adverse events, and again,  
12 none was found. Also, there was no apparent  
13 training effect in that you didn't get more of  
14 these at the beginning. They were pretty evenly  
15 distributed throughout the study.

16           And again, they were not associated with  
17 the adverse events that were measured in this  
18 study.

19           [Slide.]

20           When we looked at the post hoc data  
21 analysis--I was not present at the FDA when this  
22 was discussed, but I don't think that the exact  
23 type of analysis was specified--in fact, I know it  
24 wasn't. The problem of a post hoc data analysis is  
25 that they are not planned actively in the

1     investigational plan, so the statistical treatment  
2     that Dr. Gray will talk about is somewhat difficult  
3     in that a .05 P value is probably not the P value  
4     that would be interesting for you. And nominal P  
5     values do not account for multiplicity; you can do  
6     multiple analyses and find something that has a P  
7     value. And there is really no way to know how much  
8     adjustment should be applied when judging the  
9     significance of this.

10           [Slide.]

11           We looked at Higgins score, and Higgins  
12     score was chosen because it is generally used at  
13     one institution; it is a lot less used than some of  
14     the more common, SDS database or New York Heart,  
15     and although those two don't look at adverse events  
16     so much, they really do tell you how sick a patient  
17     you have, because it is an estimation of mortality.

18           I have spent 24 years with the STS  
19     database, and I think that that is probably an  
20     interesting way to look at the data, and we don't  
21     have the results of that analysis.

22           However, when we look at the Higgins score  
23     greater than 5, again you can see that when you  
24     look at death, slightly favored in the filter  
25     group, and stroke, slightly favored in the control

1 group, there is really no consistent trend in these  
2 data.

3 [Slide.]

4 So in summary, the filter trapped at least  
5 one particle in 97 percent of the cases. The  
6 composite safety endpoints were really not  
7 different between the two, as the hypothesis was.  
8 The individual safety events, there were no  
9 significant differences; and the only difference  
10 was in the manipulation-related aortic injury.

11 [Slide.]

12 Conclusions: The filter traps particles;  
13 a correlation with clinical improvement was not  
14 shown; there were additional concerns raised by the  
15 occurrence of aortic injuries.

16 Thank you.

17 DR. GRAY: Good morning. My name is Gerry  
18 Gray, and I was the statistical reviewer for this  
19 submission.

20 [Slide.]

21 I am going to just address a few issues.  
22 First, I am going to talk a little bit about  
23 judging the results of the trial, because I think  
24 that's really the crux of any kind of disagreement  
25 we might have with this trial.

1           Next, I'm going to talk a bit about  
2   particulate capture as a surrogate endpoint; and  
3   finally, I will finish up with the subgroup  
4   analyses and the Higgins risk scores.

5           [Slide.]

6           The first issue is how are we going to  
7   judge the results of the trial, and of course, the  
8   bottom line is always the tradeoff between a  
9   probable benefit and a potential risk. I would say  
10   that both of those things are kind of hazy in this  
11   case.

12          And the other thing you have to think  
13   about is what the appropriate set of endpoints to  
14   be using, and what is the appropriate  
15   control/comparator group.

16          [Slide.]

17          So, really, in this case, there are three  
18   main ways that you might judge the results. The  
19   first one being the most compelling is that the  
20   results are judged internally to the study, with  
21   clinical outcomes compared to a control group in a  
22   randomized trial, and at that level, you can really  
23   make pretty sound causal inferences.

24          The next level might be to judge the  
25   results in comparison to other similar devices

1 and/or studies that you might have.

2 And the third way that you might think  
3 about it is judging the results narrowly, based  
4 purely on proposed claims, thinking of the device  
5 as a "tool," in effect, and comparing that to some  
6 predefined criteria you might have for  
7 effectiveness of that tool. And in this scenario,  
8 really, there is no demonstration of clinical  
9 effectiveness.

10 [Slide.]

11 So, internally, using the adverse event  
12 composite as the endpoint, again, we had 1,289  
13 patients randomized to the EMBOL-X versus control,  
14 the outcome being the major composite adverse event  
15 rate. We had 17.1 percent of events in the EMBOL-X  
16 arm and 18.9 percent in the control.

17 The first statistical test is for  
18 noninferiority, and that was with an equivalence  
19 delta of 5 percent. So in other words, we are  
20 testing whether the EMBOL-X is no more than 5  
21 percent worse than the control. And that is  
22 strongly rejected, so that certainly we can say  
23 that the EMBOL-X device is equivalent to the  
24 control if you measure that as 5 percent.

25 And contrary to what I heard in the



1 sponsor's presentation, there as an amendment to  
2 the IDE that specified that we would do another  
3 test for superiority for this adverse event  
4 composite, and you can do that without any worries  
5 about alpha penalty, but that test is not at all  
6 significant; the P value is 0.38.

7 So the bottom line is, as you have  
8 probably already come to a similar conclusion, for  
9 the adverse event rates, the devices are equivalent  
10 but not superior in terms of this endpoint.

11 [Slide.]

12 Internally to the study, using the  
13 endpoint of particulate capture, the outcome is the  
14 proportion of the filters that capture at least one  
15 particle. And for this endpoint, there is no way  
16 to judge internally, because there was no  
17 comparison group that we had.

18 [Slide.]

19 For the other safety endpoint, to me, it  
20 looks like the results for safety are actually  
21 remarkably similar. Of all the types of serious  
22 adverse events, all 32 that we saw that the sponsor  
23 presented in one of their tables, there was none  
24 that came even close to being significant in either  
25 direction.

1           So for all the other safety endpoints,  
2 there is no evidence of any difference between the  
3 device arm and the control arm.

4           And finally, there was a secondary  
5 endpoint, as has been discussed, on aortic  
6 endothelial injury that was significantly higher in  
7 the EMBOL-X arm, but again, there is no detectable  
8 effect or outcome from those aortic injuries on any  
9 clinical adverse events. For short-term followup,  
10 there were 42 randomized patients.

11         [Slide.]

12           You might think you could compare the  
13 results of this study to some other device or a  
14 similar device. The predicate for the EMBOL-X  
15 device is the PercuSurge balloon aspiration  
16 catheter for SVG patients combined with the meshes  
17 in the CPB arterial filters.

18           Unfortunately, though, the two devices are  
19 really fairly different in terms of their mechanism  
20 action, and they are fairly different in terms of  
21 the patient populations that were studied.

22           So really, from my point of view, this is  
23 sort of a dead end; you can't really make a  
24 comparison with this device.

25         [Slide.]

1           Finally, based on the proposed claims, the  
2 evaluation of the device as a tool, the claim says  
3 "to contain and remove particulate emboli," and  
4 certainly the device is successful in that regard  
5 because it was successfully deployed in about 96  
6 percent of the patients, and depending on the  
7 denominator that you use--either the number of  
8 filters or the number of patients--it captured one  
9 or more particle in either 97 or 92.5 percent of  
10 the time.

11           So that easily meets the predefined  
12 criteria of particulate capture in 75 percent of  
13 the filtered patients.

14           [Slide.]

15           So to summarize, the internal evidence is  
16 for equivalent safety, and there is no evidence of  
17 any effect on adverse event rates.

18           Externally, it is very difficult to make  
19 any comparison, and as a tool, the device certainly  
20 captures particulate material, and from the  
21 sponsor's summary slide, they said "Clinical  
22 efficacy can be reasonably extrapolated." So if  
23 you choose that route, you may get extrapolation.

24           [Slide.]

25           Let me just talk a little bit about

1 particular capture as a surrogate, because one way  
2 you might justify the particulate capture is to  
3 think it is a surrogate for some clinically  
4 meaningful endpoint. And the question here in  
5 bullet number two--is particulate capture, however  
6 you measure that, a valid surrogate for clinical  
7 adverse events?

8 Really, to be a valid surrogate, the  
9 endpoint has to be somehow correlated with the  
10 outcome of interest, and somehow it has to capture  
11 the effect of the treatment on that outcome.

12 [Slide.]

13 So just going down that path a little bit,  
14 here is a two-by-two table that shows whether or  
15 not particles were captured and then whether or not  
16 a composite event was observed. And as you can  
17 see, the correlation coefficient is quite small  
18 there. There is really no obvious correlation that  
19 I can see between particulate capture and whether  
20 or not there was a composite event.

21 It is the same if you do number of  
22 particles captured, whether it captured any  
23 particle or not, or if you measured the total  
24 surface area of particles captured.

25 [Slide.]

1           Looking at that a little further, this is  
2   called an Q-Q [phonetic] plot, and it plots the  
3   quantiles of two distributions.

4           Here, on the X axis is the number of  
5   particles captured in patients who had no composite  
6   event, and the Y axis shows the number of particles  
7   captured in patients who did have a composite  
8   event. And if the two distributions are the same,  
9   we would expect that to be a straight 45-degree  
10   lines, and indeed, it is almost entirely a straight  
11   45-degree line. So there is no real evidence that  
12   there is any difference, except potentially out in  
13   the fair tail there, where you have more than 12 or  
14   15 particles captured.countries

15           [Slide.]

16           So I did a little bit of--well, actually,  
17   first of all, here is the same kind of plot that is  
18   using particle area instead of number of particles,  
19   and it looks the same, visually.

20           The only thing that would make you think  
21   there might be a relationship would e what is going  
22   on in the extreme tails here, where you have more  
23   than 10, 12, 15 particles captured. So I did a  
24   little bit of data-dredging of my own to try to  
25   figure out if there was anything going on out

1     there, and that represents about 46 patients, I  
2     think it was, out of the 1,200 in the trial. It is  
3     relatively low numbers, and from a statistical  
4     point of view, you can't really draw any  
5     conclusions from those.

6             [Slide.]

7             So in regard to particulate capture as a  
8     surrogate, it really doesn't meet the condition  
9     that it is correlated with the outcome of interest.

10            And similar results hold if you use other  
11    endpoints, or other kinds of adverse events not  
12    included in the composite.

13            It is possible, of course, as Rick Kuntz  
14    pointed out, that there could be effects that are  
15    so subtle that they were not measured in this  
16    trial, and therefore, we have no way of knowing  
17    whether there is any effect on them.

18            [Slide.]

19            The third and final topic is just covering  
20    the basis on the subgroup analyses. The sponsor  
21    acknowledged this, that on Table 6 of their panel  
22    package, they have 36 different subgroup analyses  
23    that they performed, and certainly, when you look  
24    at them through the statistical viewpoint, the P  
25    values are small, but given the number of subgroups

1 that we have gone through, that really is not  
2 surprising at all.

3 [Slide.]

4 And finally, for the preoperative Higgins  
5 risk score, really, that is in a sense another  
6 subgroup analysis, because we have used the Higgins  
7 risk--we tried several cuts on the Higgins risk  
8 score and found one where the P value was slightly  
9 less than .05. But in order to make any strong  
10 statistical conclusion from that, we would need to  
11 know how small the P value has to be to be  
12 significant, and .047 really isn't it if you do any  
13 reasonable adjustment for multiplicity.

14 [Slide.]

15 So to summarize the last subject, the  
16 subgroup analyses really don't provide any evidence  
17 of superiority in terms of adverse event rates.  
18 And again, I heard this pretty clearly from the  
19 sponsor's presentation as well.

20 That concludes my presentation.

21 MS. WENTZ: Thanks very much.

22 At this point, I just want to review the  
23 Questions to the Panel, and I believe there are  
24 six.

25 Question 1. The primary safety endpoint

1 for this study was a composite of seven clinical  
2 adverse events including death, neurologic deficit  
3 mild and severe, renal insufficiency, perioperative  
4 myocardial infarction, gastrointestinal  
5 complications, and limb-threatening peripheral  
6 embolism, evaluated at hospital discharge or 30  
7 days, whichever was shorter. The median followup  
8 time was seven days.

9 Some facts from the study are: The  
10 observed overall composite event rates were 17.1  
11 percent in the EMBOL-X arm and 18.9 percent in the  
12 control.

13 The composite event rate for the EMBOL-X  
14 arm was shown to be equivalent--not more than 5  
15 percent higher--than the control.

16 Also as specified in the protocol, a  
17 separate test for a lower event rate in the EMBOL-X  
18 arm was not statistically significant.

19 The EMBOL-X arm demonstrated a  
20 significantly higher incidence of aortic  
21 endothelial injury--9.2 percent versus 2.0 percent.  
22 Although these patients did not appear to have any  
23 short-term clinical sequelae resulting from the  
24 injuries, the long-term effects are unknown.

25 And the final question being: Do these



1 data support the safety of the EMBOL-X intra-aortic  
2 filter?

3 Question 2. The primary effectiveness  
4 endpoint in this trial was to demonstrate that 75  
5 percent of the devices would capture at least one  
6 particle during elective CABG or single-valve  
7 procedures. This was demonstrated in the study.

8 There was no demonstrated reduction in any  
9 category of clinical adverse event in this  
10 well-controlled 1,289-patient trial. Please  
11 address the following concerns:

12 1) Can this method of embolic entrapment,  
13 from this study or elsewhere, be extrapolated to  
14 clinical efficacy?

15 2) Do these data support the effectiveness  
16 of the EMBOL-X intra-aortic filter?

17 Question 3. Do the study data support an  
18 appropriate risk/benefit profile?

19 Question 4. One aspect of the 510(k)  
20 review of a new product is the review of its  
21 labeling. The labeling must indicate which  
22 patients are appropriate for treatment, identify  
23 potential adverse events with the use of the  
24 device, and explain how the product should be used  
25 to maximize benefits and minimize adverse effects.

1           Please address the following questions  
2   regarding product labeling:  
3       1) Do the Indications for Use adequately  
4   define the patient population studied? For  
5   example, should the patient population receiving  
6   this device be limited to the same patient  
7   population utilized in the study--for example,  
8   nonemergent, patients over age 60, and first-time  
9   isolated valve or CABG patients.  
10       2) Are there any other restrictions that  
11   should be placed on the patient population  
12   receiving this device?  
13       3) Based on the clinical experience,  
14   should there be additional Contraindications,  
15   Warnings, and Precautions for the use of the  
16   EMBOL-X intra-aortic filter?  
17       4) Should the labeling include specific  
18   study information such as: no reduction of  
19   clinical events were noted in a 1,289-patient  
20   clinical study; and the EMBOL-X device appears to  
21   increase the rate of endothelial injury?  
22       5) What should the labeling include  
23   regarding the use of ultrasound both before--for  
24   assessment of the aorta--and after--monitoring of  
25   injury--the use of the device?

1           Question 5. Please provide any other  
2 recommendations or comments regarding the labeling  
3 of this device.

4           Question 6. If the data provided are not  
5 adequate to support safety and/or effectiveness,  
6 what additional data, analyses, or study would you  
7 require?

8           Thank you.

9           Questions and Answers

10          DR. TRACY: Does that complete the FDA  
11 presentation?

12          MS. WENTZ: Yes.

13          DR. TRACY: Does the panel have any  
14 questions for the FDA before we move on?

15          DR. LASKEY: I have one question to the  
16 engineer. Maybe it is trivial, maybe not.

17          Nitinol and its thermal properties--there  
18 is a nitinol frame here. Patients are generally  
19 cooled when they are put on bypass, hearts are  
20 cooled, and so on and so forth. Is there anything  
21 happening with the--should we be concerned about  
22 any change in function or configuration of the  
23 frame here?

24          MS. WENTZ: That is a very good question,  
25 and actually, that was brought up in a few of the

1 stent studies a few years ago, when you place  
2 stents that have nitinol in them because patients  
3 are cooled down as well. And I believe from that,  
4 I did not review that from a material  
5 standpoint--our OST scientists did--but that was  
6 looked at, and the temperature that the patients  
7 are cooled down to does not affect the nitinol.

8 DR. TRACY: Dr. Marler?

9 DR. MARLER: In Dr. Swain's discussion of  
10 the myocardial infarction data, I may not have been  
11 able to see the complete slide, but I didn't have  
12 the impression of the same difference between the  
13 groups and the incidence of MI, and I was wondering  
14 where that information came from.

15 DR. SWAIN: I think in your pack, you can  
16 talk about total MIs or Q-wave versus non-Q-wave  
17 MIs. So that is the difference. I believe the  
18 sponsor's presentation was total MIs, and my  
19 presentation was to pick out Q-wave MIs and stroke  
20 versus the lesser injury.

21 Do you want that slide back up? I might  
22 be able to do that.

23 DR. MARLER: I might be looking at the  
24 wrong table; I am looking on page 35, Table 62.

25 [Slide.]

1 DR. MARLER: Okay. It was a problem of  
2 not being able to read the slide, because I  
3 couldn't read the X axis. So it is clear now. I'm  
4 sorry.

5 DR. TRACY: Dr. Krucoff?

6 DR. EDMUNDS: Julie, when they did these  
7 epi-aortic ecograms, did they look at the area of  
8 the cross-clamp with the core laboratory look,  
9 where they were able to see things that nobody else  
10 could see?

11 DR. TRACY: Can we ask that during the  
12 open committee discussion, please? We haven't  
13 quite gotten there yet. I just don't want to mix  
14 the FDA with the sponsor.

15 DR. EDMUNDS: I just asked Dr. Swain. I  
16 don't know why--

17 DR. TRACY: If she doesn't know the  
18 answer, then, let's just hold the question.

19 DR. SWAIN: Yes, it was looked at. And  
20 the injuries are not seen by the surgeon--I didn't  
21 look at my screens very much when I was busy  
22 closing up and getting up pumps. It was identified  
23 at the institution, I guess.

24 DR. TRACY: Dr. Krucoff?

25 DR. KRUCOFF: I have a question for Dr.

1 Gray.

2 Gerry, in your particles analysis--I  
3 probably just didn't connect when it started--but  
4 to me somehow, the elimination of particles by the  
5 filter would ostensibly be associated with a  
6 reduction in clinical events; if you got them out  
7 with the filter, presumably, you are saving the  
8 patient that avalanche effect.

9 And somehow as I looked at these--is that  
10 where the divergence in total number of particles  
11 captured--well, I'm confused.

12 DR. GRAY: Let's go back to Slide 2.

13 DR. KRUCOFF: Yes, because that's where it  
14 started.

15 [Slide.]

16 DR. KRUCOFF: So is your expectation--is  
17 what you are testing here that if you get more  
18 particles out with the filter, you are more likely  
19 to have a composite event?

20 DR. GRAY: We have a problem here with the  
21 endpoint, because we all would like to see the  
22 clinical endpoint, and what we have is the  
23 particulate capture. My line of reasoning here was  
24 let's see if somehow we can justify using  
25 particulate capture as a surrogate for some

1 clinical outcome that we are interested in, in this  
2 case being the composite adverse events.

3 So I am trying to see if there is some  
4 correlation between whether or not particles were  
5 captured in a patient and did that patient have an  
6 adverse event or not.

7 DR. KRUCOFF: Okay, but the assumption  
8 here is what I am trying to get at, Gary--

9 DR. GRAY: I guess the assumption is that  
10 if removing particles does anything, we would  
11 expect that where the particles were removed, we  
12 would be reducing the composite adverse event rate.  
13 That is my assumption.

14 DR. KRUCOFF: Well, I would suggest that  
15 the assumption is if particles are a surrogate for  
16 bad things happening, that that is what happens in  
17 a control population--if you have 1,000 patients  
18 with no protection, some of them are going to have  
19 very few particles, and they would have fewer  
20 events; others are going to have showers of  
21 particles or big particles, and those would have  
22 clinical events.

23 The trouble is that as you start removing  
24 particles, if you capture very few particles, those  
25 may be patients who have very few particles, and if

1 you capture a lot of particles, those might be  
2 patients who would be high risk whom you are  
3 protecting.  
4 I just don't see how this would even begin  
5 to test that, because the assumption is that  
6 somehow particles removed would correlate with  
7 badness in the same patient population, and--  
8 DR. GRAY: I agree that--we wish we knew  
9 the denominator here, right; we wish we knew for  
10 any patient how many particles really there were  
11 present, and then we would have some idea of the  
12 effectiveness of removing those particles, so we  
13 could say that somehow, the amount of particles  
14 released in that patient is some indication of  
15 their potential risk for an adverse event, and then  
16 we could try to figure out, okay, if we remove a  
17 certain proportion of them or if this device can  
18 remove some proportion of them, what effect would  
19 that have on the outcome. That's what I wish we  
20 knew. But we don't know that. We don't really know  
21 at all--and you are right, we don't know for any  
22 patient--if the device captured one particle, that  
23 could be the only one that was released, or it  
24 could be one out of thousands. There is no way to  
25 tell that with the data that we have.



1           So I admit I was on a bit of a stretch  
2 here to try to figure out is there any way to take  
3 the data that we do have, which is purely the  
4 number of particles captured, and relate that to  
5 whether or not there was an adverse event. That's  
6 why I went on to the next slide, which was just  
7 that.

8           What is missing here is--you are  
9 right--what we really would love to know is the  
10 underlying information for each patient as to how  
11 many particles were actually released and how much  
12 risk was that patient exposed to. But we have  
13 nothing to tell us that as far as I know.

14          DR. KRUCOFF: Okay. Can I ask  
15 you--because when I actually walked through these  
16 slides, what I ended up sitting here thinking,  
17 which I took as different from what you were  
18 suggesting, is that although I agree it is a  
19 stretch, this might be taking as an imputation that  
20 you can take patients who are at much higher risk,  
21 i.e., the higher -particle group, and pull them  
22 down to a line of identify with patients who are at  
23 much lower risk, i.e., patients who have fewer  
24 particles, in a population where you are not  
25 allowing these particles into the systemic

1     circulation--you are removing them.  
2             Is that wrong?  
3             DR. GRAY: I can't say whether that is  
4     wrong or right, because if I understand correctly,  
5     what you are thinking is that somehow this mesh  
6     puts a limit on the amount of particles that are  
7     released, that actually escape through into the  
8     circulation, and therefore, that would be nothing  
9     but a good thing.  
10            Is that a correct interpretation of what  
11   you said?  
12            DR. KRUCOFF: But that would be one  
13   notion, I think, of the whole generation of distal  
14   protection devices, that basically, the more you  
15   get out with the device represents some sort of  
16   surrogate incremental protection afforded the  
17   patient.  
18            DR. GRAY: Yes. And that sort of gets  
19   back to my first set of ways of judging the results  
20   of the trial. And you can infer in your own mind  
21   that removing particles is undoubtedly a good thing  
22   and that the device only needs to be shown  
23   effective as a tool that removes particles, and  
24   that's all they care about, therefore, I'm happy  
25   with.

1           On the other hand, we have the internal to  
2 the study comparison between the treatment and the  
3 control group, where it was remarkably uniform,  
4 remarkable similarity in the adverse event rates  
5 across the board.

6           So how do we make that judgment--that's  
7 why I started out with that, because I think that  
8 is really the whole crux of the--

9           DR. TRACY: The less kind interpretation,  
10 Mitch, would be that it doesn't matter if you  
11 remove particles--the risk is the same--

12          DR. KRUCOFF: Yes.

13          DR. GRAY: That's right.

14          DR. KRUCOFF: Understood.

15          DR. TRACY: Okay. Can I ask Dr. Wentz to  
16 clarify--you mentioned that bench study, some  
17 design concerns and/or test method concerns  
18 remained that may be related to the endothelial  
19 injuries. Could you expand on that just a little  
20 bit?

21          MS. WENTZ: First of all, it is not  
22 "Doctor" but thank you.

23          When this submission first came in, we  
24 looked at the test methods and the results and  
25 procedures and all that, and everything looked

1    okay, and we let the study go on.  It wasn't until  
2    we started focusing on these endothelial injuries  
3    that we backtracked and said, okay, what could some  
4    of the possibilities be for these injuries.

5            Dr. Sapirstein and myself re-reviewed all  
6    of those test methods and found that there were a  
7    number of them that could possibly be related to  
8    those endothelial injuries.  So we just sent those  
9    questions to the company in the form of a 510(k)  
10   Additional Information Letter--and did that come  
11   back already--no--they are still formulating the  
12   answers to those.

13           Does that answer your question?

14           DR. TRACY:  I guess so.  I'm not sure what  
15   the design questions are--whether it is that the  
16   thing is too stiff or is too--is there some  
17   fundamental problem with this thing that you are  
18   asking them to clarify?

19           MS. WENTZ:  Yes, that's basically it.  When  
20   we tried to repeat their test methods using the  
21   sample device that we had, some of the forces that  
22   we felt were not anywhere near some of the forces  
23   that were on the paper that they said they had  
24   recorded.  So we just asked them to clarify some of  
25   their test methods and procedures in light of the

1     endothelial injuries.

2             DR. TRACY: Thank you.

3             Are there any other questions?

4             Dr. Aziz?

5             DR. AZIZ: Let me ask Julie a question.

6             Julie, these endothelial disruptions--I'm  
7     sure the company will focus on that later--in your  
8     review of the data, where were they occurring? Was  
9     it at the tip of where the net is? Where is it,  
10    and what do you think is causing it?

11            DR. SWAIN: Right. You know what you can  
12    see on a TEE; essentially, you are blocked off  
13    because of the airways so that in the area examined  
14    of the ascending aorta, they occurred. Some of  
15    them occurred proximal, where the clamp was, but I  
16    didn't see it, or don't recall it being broken out  
17    that much; but they certainly did occur in the area  
18    where you have aortic manipulation.

19            And I think also, in answer to Dr.  
20    Krucoff's question, the patients who had the  
21    biggest amount of atherosclerosis were screened  
22    out--that group wasn't studied. So the kind of  
23    catch-22 is that maybe they would benefit more, but  
24    again, aortic manipulation in the presence of known  
25    atherosclerosis, from all the data and the work of

1 20 years from Dr. Kouchoukos, is what we all  
2 learned we shouldn't do.

3 DR. AZIZ: And then actually looking at  
4 the ECO that was presented earlier, it seems like  
5 there were two different types. You had this  
6 fibrinous strand sort of waving at you, and then  
7 you had like an intramural hematoma on the wall.  
8 Maybe we'll look at that later.

9 DR. SWAIN: Yes. You may ask the sponsor  
10 about the intramural. I didn't remember seeing  
11 that. It's kind of like in surfing, we use the  
12 term "dings"--it is an aortic "ding"--no clinical  
13 consequence as evaluated in this short-term study.

14 Open Committee Discussion

15 DR. TRACY: At this point, let's move on  
16 to the open committee discussion. I think there  
17 are lots of questions waiting to be asked.

18 I would just like to remind everyone that  
19 this is a premarket notification or a 510(k)  
20 submission that is being brought to the panel at  
21 this time. And at the end, the FDA is asking for  
22 recommendations and advice. There will not be a  
23 final vote. And the two lead reviewers were Dr.  
24 Marler and Dr. Edmunds.

25 Dr. Marler, if you would like to lead off

1 with questions for the sponsor.

2 DR. MARLER: Okay. I guess this very  
3 technical--I still have this question comparing  
4 table 7-18 and 7-19. I tried to read and  
5 understand, and I'm sure there is an explanation,  
6 but for adverse event under NIH Score greater than  
7 4 in the control group, there are 13 in table 7-18  
8 with 644 patients, and then, when the sample size  
9 is reduced to 620, there are more--16. Is that  
10 because you are including events that occurred  
11 after the first exam?

12 DR. ALLEN: I think it's a very good  
13 observation. There were very little times when the  
14 initial NIH score was applied. In the initial  
15 design of the study, we had hoped to have a 24-hour  
16 evaluation on every patient, but it became quite  
17 obvious as the study progressed that that wasn't  
18 practical. Patients were [inaudible] and so forth.

19 So discussions then allowed us to do our  
20 first evaluation at 3 plus or minus one day. So  
21 that initial evaluation is variable in time.

22 When you look at the 7-day evaluation  
23 which is applied evenly among both groups--call it  
24 the end evaluation of the New York Stroke  
25 Scale--the rates were essentially identical between

1 the groups--2.6 was [inaudible].

2 DR. MARLER: So, then, which one do you  
3 think most accurately reflects the strokes that  
4 were due to the cardiopulmonary bypass procedure  
5 and the surgery? Is it Table 7-18?

6 DR. ALLEN: That's a very interesting  
7 question, because when you look at actual frank  
8 stroke after cardiopulmonary bypass, Dr. Kuntz in  
9 his presentation outlined the multiplicity of  
10 reasons for why patients have strokes. The stroke  
11 rate overall in our study was about 2.5 percent.  
12 It is interesting that if you look at the time  
13 course as to when those strokes occurred, only  
14 about half of them actually occurred greater than  
15 24 hours--the patients woke up neurologically  
16 intact, and at day 2 or day 3 had an event.

17 So the device's potential to impact frank  
18 stroke is with the operation. When we did an  
19 analysis, for example, on the impact of atrial  
20 fibrillation, 60 percent of patients who had stroke  
21 also had atrial fibrillation, which we know is a  
22 potential indicator for stroke.

23 Dr. MARLER: The neurologist was looking  
24 throughout this for any description of the strokes  
25 or any further analysis, even breakdown as to



1 hemorrhagic or ischemic--did I miss it? Is it  
2 somewhere in the writeup, or was there data  
3 available to me on how the strokes were diagnosed  
4 as to their type?  
5 DR. ALLEN: I don't believe we broke the  
6 strokes down, and I don't--  
7 DR. MARLER: For severity?  
8 DR. ALLEN: --for severity as far as  
9 whether it is a hemorrhagic stroke or--  
10 MS. CHANG: These are all ischemic  
11 strokes.  
12 DR. ALLEN: --but they were adjudicated by  
13 a blinded events committee and felt to be related  
14 to surgery.  
15 MS. CHANG: We provided narratives in the  
16 510(k) filing on these.  
17 DR. MARLER: Okay. But those aren't in  
18 the packet here.  
19 MS. CHANG: No.  
20 DR. MARLER: Okay. They are not very  
21 interesting reading to the cardiovascular surgeon,  
22 I'm sure.  
23 [Laughter.]  
24 DR. MARLER: So, what I am looking for is  
25 an argument as to the logic of--I mean, we have

1 said that there is not a surrogate marker for  
2 safety. What is the clinical efficacy--what is the  
3 benefit to the patient--I mean, what's the talk on  
4 this? What are you expecting--why do this? It  
5 seems to me the study has shown that it is as good  
6 as doing nothing. But why is it better, and what  
7 are you thinking?

8 I was concerned that the discussion about  
9 the cognitive outcomes indicated--which I thought  
10 would be an obvious possible benefit--it was stated  
11 that that wasn't really thought of as a potential  
12 benefit. So I am a little unclear on the thinking  
13 of really what this study means to the patient.

14 DR. ALLEN: We grappled with that, and I  
15 thought Dr. Kuntz tried to outline that in his  
16 presentation with regard to study design.

17 You know, intuitively, reducing the  
18 particulate emboli load is a good thing, but we  
19 grappled in the design of the study with the very  
20 question that you are asking: How can we  
21 demonstrate an efficacy endpoint?

22 And the conclusion was that, for example,  
23 unlike the SAFER trial where you had a very  
24 specific marker--CPK isoenzymes that affected a  
25 very specific end organ--with the exception of

1 perhaps serum creatinine, we didn't have specific  
2 sensitive markers that might detect subtle clinical  
3 changes in patients' outcome. So we ended up with  
4 a design trial that essentially looked at safety  
5 equivalency to show that the device wasn't causing  
6 harm and that the particulate capture was the  
7 efficacy endpoint, and that capture of particles  
8 was a good thing.

9       You know, the difficulty of designing a  
10 trial comes down to what can be practically applied  
11 and logistically applied across multiple centers,  
12 and the one power calculation that Dr. Kuntz  
13 did--if you look at just, for example, frank stroke  
14 and assume you have a 3 percent incidence of frank  
15 stroke, not all of those strokes occurred in the  
16 operating room, so you wouldn't even expect that  
17 the device would prevent all of those strokes, but  
18 let's assume for argument that the event rate was 3  
19 percent. A 20 percent reduction in that 3 percent  
20 rate would require a sample size of a little over  
21 22,000 patients to demonstrate that.

22       So you weigh what seems clinically  
23 intuitive with the practical aspects of designing a  
24 trial that demonstrates that clinical efficacy.  
25 And I think the additional analysis that we

1 provided--and I make full disclaimer--I take my mea  
2 culpa in that I don't make claims of superiority  
3 when we look at that additional analysis, but it  
4 does provide some element of risk-benefit as to  
5 what population this may truly benefit, and it is  
6 the higher-risk group that we know as surgeons have  
7 an increased risk for morbidity and mortality  
8 postoperatively that, intuitively, reducing that  
9 embolic load in those patients seems very  
10 reasonable, and that additional analysis, when you  
11 looked at an endpoint that has a specific marker,  
12 i.e., serum creatinine, you began to see clinical  
13 ticks that, yes, there is something maybe going on  
14 there.

15 But I agree with the FDA, and I don't want  
16 to make claims of superiority. We simply rely on  
17 the clinician's intuitiveness that a reduction of  
18 this embolic load is a good thing.

19 DR. MARLER: What can you say that would  
20 reassure me or the committee--I mean, we have  
21 knowledge that there are instances where something  
22 that is really intuitively obvious--blood pressure,  
23 arrhythmias, I hesitate to mention, but ECIC bypass  
24 I am pretty confident of, in which there is a real  
25 obvious case in which the intervention did what it

1 was supposed to do but it wasn't clear at all, and  
2 some people still think it may have actually been  
3 harmful but hiding underneath the obvious clinical  
4 benefit.

5 Is there something different about this  
6 that would--is there any reassurance you can offer?

7 DR. ALLEN: No, and I think the panel  
8 members are grappling with the same things that I  
9 grapple with when I think about this data.

10 You are absolutely right. There are many  
11 instances where your intuition tells you something  
12 is good, and a well-designed trial tells you that  
13 now your intuition wasn't as good as you thought it  
14 was.

15 In this particular trial, it was designed  
16 as a safety study to demonstrate that the device,  
17 compared to current cannulation techniques that we  
18 use every day in open heart operations, isn't  
19 worse, and that particulate capture was the  
20 clinical efficacy endpoint.

21 DR. MARLER: Thank you.

22 DR. TRACY: Dr. Edmunds?

23 DR. EDMUNDS: Keith, are you going to be  
24 the one who responds, or someone else? On this  
25 injury, if we discount the three, one of which was

1 the scalpel and the other one, a surgeon took a  
2 stitch or two, and just concentrate on the 42, most  
3 of which were not surgeon-noticed at the time, or  
4 anesthesiologist reading the ecocardiogram, whoever  
5 it was, how many of those were more than just the  
6 endothelium? How many actually got into the media?

7 DR. ALLEN: I think Dr. Kouchoukos is  
8 experienced in this field, so I'll let him answer  
9 that question.

10 DR. KOUCHOUKOS: These were all basically  
11 endothelial disruptions. In other words, they are  
12 just small fragments of intima.

13 DR. EDMUNDS: And do you think they were  
14 scratches from the nitinol wire?

15 DR. KOUCHOUKOS: Well, the question was  
16 raised earlier about in the filter group, where  
17 these intimal disruptions were located, and they  
18 were distributed throughout the ascending aorta.  
19 Some of them were clearly related to the filter,  
20 but others occurred in the mid-aorta or perhaps in  
21 the more proximal part. So they would be expected  
22 to have resulted from other manipulations of the  
23 aorta, and that's basically why they occurred in  
24 the control group.

25 DR. EDMUNDS: And that's why they occurred

1 in the control group, too.

2 DR. KOUCHOUKOS: That's what I said;  
3 that's why they were present in the--they were  
4 present in 2 percent of the control group. And I  
5 think--

6 DR. EDMUNDS: right. So how long do you  
7 think it took for that to heal?

8 DR. KOUCHOUKOS: Well, we don't have  
9 followup ecocardiograms or epi-aortic images to  
10 know the answer to that. It is also important to  
11 note that only one of these was detected by  
12 transesophageal ecocardiography. They were all  
13 detected for the most part by epi-aortic scanning.  
14 So you wouldn't see one of these with a  
15 transesophageal, and certainly not with a  
16 two-dimensional surface ECO.

17 DR. EDMUNDS: Isn't it a stretch to call  
18 this an "injury"?

19 DR. KOUCHOUKOS: Well, it's a good  
20 question. We termed it a "disruption," but I think  
21 others would term it an "injury."

22 DR. EDMUNDS: Gosh, I would consider it an  
23 overinterpretation of the ecocardiogram myself, if  
24 it were one of my cases.

25 Does anyone really think that this will

1 progress to any problem downstream for the patient?

2 DR. KOUCHOUKOS: I think, again, it is  
3 important to put it in historical context. These  
4 endothelial disruptions have been occurring since  
5 we started doing cardiac surgery.

6 DR. EDMUNDS: Exactly.

7 DR. KOUCHOUKOS: They have been there  
8 forever. And from what we know about the outcomes  
9 of patients who have cardiac surgery, they are  
10 probably of no significance. We know that an  
11 intra-aortic dissection is a catastrophic event,  
12 and we know how frequently that occurs, and it is  
13 very rare. And I think to extrapolate to what  
14 happens to the endothelial disruptions is hard,  
15 because they are not as significant as the others,  
16 and we really have no way of following what happens  
17 to them. We would surmise that they probably heal  
18 eventually, but we have no hard data to support  
19 that.

20 DR. EDMUNDS: In your experience as a very  
21 busy cardiac surgeon over the long period, which is  
22 the greater injury--the cross-clamp injury to the  
23 endothelium, or produced by this filter--in your  
24 opinion? I know you don't have data, but you have  
25 a lot of clinical experience.



1 DR. KOUCHOUKOS: Well, I think that what,  
2 if anything, we have learned from this study, and  
3 as we have learned from our own clinical  
4 experience, is that we want to manipulate the  
5 ascending aorta as little as possible. And  
6 certainly a cross-clamp is a major insult, if you  
7 will, to the ascending aorta. It is exposed to a  
8 lot of surface of the aorta, with the potential for  
9 dislodgment of atheromatous debris, and a  
10 side-butting clamp is the same.

11 DR. EDMUNDS: So you think the clamps are  
12 a bigger injury?

13 DR. KOUCHOUKOS: I do.

14 DR. EDMUNDS: Now, as I understand it, the  
15 company does not intend to make any statement on  
16 the labeling about--am I out of order already--

17 DR. TRACY: No--not yet.

18 [Laughter.]

19 DR. EDMUNDS: --okay--about clinical  
20 benefit; is that correct?

21 MS. CHANG: To both of them, yes.

22 DR. EDMUNDS: Okay. There is no evidence  
23 that particulate emboli to the brain is good, so it  
24 is logical to assume that reducing it is at least  
25 not bad and is probably good.

1           What percentage of the atherosclerotic  
2 emboli to the brain would you guess this filter  
3 when it is deployed decreases from some unknown  
4 hole? What would be your clinical estimate? I  
5 have one in my head.

6           DR. KOUCHOUKOS: I'm not sure I understand  
7 your question, Dr. Edmunds.

8           DR. EDMUNDS: What percentage of all the  
9 emboli that go to the head from surgical  
10 manipulation doing a case do you think this filter  
11 catches?

12          DR. KOUCHOUKOS: Well, it catches a  
13 different amount of material from each patient, and  
14 I think you saw that. There are patients who  
15 release small numbers and small sizes of  
16 particulate matter and others who release large  
17 amounts.

18          From what we know from Dr. Barbut's  
19 studies and from our own experience, I think about  
20 20 percent of those have the potential to go to the  
21 cerebral circulation, and a percentage of those  
22 would probably be dispersed to the brachial  
23 arteries and not enter the brain, but in her study,  
24 I think on overage, about 9 percent of the emboli  
25 in one small study that were released went to the

1 brain.

2 DR. EDMUNDS: Well, I would answer the  
3 question a little bit differently. Fourteen  
4 percent of the cardiac output goes to the brain.  
5 That is physiology. So we would presume that the  
6 amount of emboli would be the same, unless there is  
7 some streaming.

8 We also know that the injury to the brain  
9 is due to a whole lot more than atherosclerotic  
10 emboli--complement activation, cytokines, regional  
11 perfusion differences, temperature movement, and  
12 all that sort of thing go into cumulative brain  
13 injury--but from what I know, I think you are only  
14 catching a fraction of the total exposure of the  
15 brain to atherosclerotic emboli. I have no idea  
16 what the fraction is exactly, but I suspect it is  
17 less than 50 percent. Would you disagree strongly  
18 with that, any of you?

19 DR. KOUCHOUKOS: No.

20 MS. CHANG: No.

21 DR. EDMUNDS: My case rests.

22 DR. TRACY: Thank you.

23 We'll go around the table with panel  
24 members to allow them to ask any questions they  
25 have for the sponsor, and we'll begin with Dr.

1 Pina.

2 DR. PINA: I have a question about the  
3 renal dysfunction. I see your definition of renal  
4 insufficiency being an increase of greater than 2  
5 or 50 percent increase, and I may have missed it  
6 here, but do you actually have the values of the  
7 creatinines? Do you have the mean values--because  
8 so many things happen around surgery with drugs  
9 that we give that can alter renal function back and  
10 forth, and yet to those of us who take care of  
11 these patients afterward, that is a very  
12 significant point, the renal function.

13 DR. ALLEN: I think the important aspect  
14 of that is that it is 50 percent above baseline. I  
15 think one of the things that the investigators  
16 wanted to put into this study is that if you start  
17 out at a creatinine of 1.8 and go to 2.0, it's not  
18 fair to count that as a patient who has renal  
19 insufficiency, but it is a 50 percent increase from  
20 baseline or any increase above 2 that is important.

21 DR. PINA: But I can also do that if I  
22 give a lot of diuretics to a patient in a  
23 perioperative period.

24 I would like to know what happens to those  
25 patients later. Do you have any followup after

1 those 7 days about the renal function?

2 DR. ALLEN: Actually, what you see is that  
3 renal function, as you well know, worsens after  
4 cardiac surgery. There were some patients who  
5 required dialysis, but renal function returns to  
6 normal.

7 The beauty of this is that this was in a  
8 randomized trial, so the same variables like did  
9 you start Altase postoperatively, at a time when  
10 you are diuresing a patient. You make the  
11 assumption--and it is the reason you do the  
12 randomized trial, to allow for those variables to  
13 be adjudicated.

14 DR. PINA: My point about followup with  
15 the renal function is it may help to differentiate  
16 the things that are strictly just the drugs that we  
17 do, or is it really emboli events to the kidneys,  
18 which may not result.

19 DR. ALLEN: I guess I don't know that our  
20 data can help you answer that. All I know is that  
21 in a randomized trial, when you look at the safety  
22 endpoint, one of the endpoints in the composite was  
23 renal insufficiency, and you can't make a statement  
24 that renal insufficiency was significantly reduced.  
25 Only when you look at the higher-risk patients do

1 you begin to see trends or ticks in favor of a  
2 reduction of renal insufficiency. But you have to  
3 assume that if you did everything identical in the  
4 two groups except one got a filter and one didn't,  
5 and you see the impact that renal insufficiency has  
6 on length of stay--and you know that from your  
7 clinical experience--in our situation, if you had  
8 renal insufficiency in our study, your length of  
9 stay was almost 15 days compared to 7 if you didn't  
10 have renal insufficiency. So it does have a  
11 dramatic impact.

12 DR. PINA: I have no further questions.

13 DR. TRACY: Dr. Ferguson?

14 DR. FERGUSON: First, I want to  
15 congratulate the presenters, both your group and  
16 the FDA, for very lucid presentations.

17 I have a couple of questions that relate  
18 to the particulate matter. The difficult is, as  
19 has been mentioned many times before, that we don't  
20 know what the denominator is, whether in the total  
21 spectrum of open heart surgery on a person who has  
22 atherosclerosis or some clot in the ventricle,  
23 whether a total screen would capture 1,000  
24 particles, 500,000, whatever. So that is of a  
25 little concern, and I will get back to it in a

1 second, relative to the time that you deploy.

2 But first I want to ask about in the  
3 Higgins above-5 group, did you note that there were  
4 more numbers of particles in that than you would  
5 expect? I missed that; I'm sorry.

6 DR. ALLEN: I think that's a very good  
7 question. to be honest, I'm not sure we did that  
8 analysis. We looked at the Cleveland Clinic score  
9 that was prespecified in those papers, picked their  
10 number of 5 and used that number.

11 DR. FERGUSON: You would expect that would  
12 be the case. And that gets to my second question,  
13 which is that the instrument was not stressed to  
14 the max, if you will, because it was not purposely  
15 p ut in the kinds of aortas that everybody is  
16 seeing today. I think that's a fair statement--or  
17 is it?

18 DR. ALLEN: Yes, sir, and part of the  
19 exclusion--although we didn't specifically exclude  
20 patients with, for example, Grade 4 aortas--the  
21 exclusion criterion is that if you did your  
22 stronotomy and opened the patient, and it was an  
23 aorta that the surgeon did not feel that he could  
24 clamp or wanted to clamp, then, patients were  
25 excluded.

1           So you are absolutely right, it did not  
2 necessarily even apply to the worst patients.

3           DR. FERGUSON: The issue, then, for us is  
4 that--I don't know how we would approach this, and  
5 FDA will tell us--but this obviously is going to be  
6 used in the very severe aortas all the way up to  
7 the porcelain aorta.

8           Nick, do you want to respond?

9           DR. KOUCHOUKOS: Well, there is certainly  
10 the potential to use it in those patients. I think  
11 it would depend on the comfort level of the  
12 individual surgeon. But I see no reason why it  
13 would not be used in severely atherosclerotic  
14 aortas. We didn't encounter many patients who fell  
15 into that category by virtue of the patient group  
16 that we were elected to study.

17          DR. FERGUSON: It gets to the disruption  
18 issue and whether there are going to be more  
19 disruptions in that group. I suspect there will  
20 be, because you say they--that gets to my next  
21 question, if I could go on to that, and that is you  
22 have had more experience with epi-aortic ECO than  
23 anybody in the world, I suspect. Have you seen  
24 these disruptions in the series that you did prior  
25 to this study?



1 DR. KOUCHOUKOS: One of the interesting  
2 things about epi-aortic scanning is that we have  
3 been using it for a long time, but we never for the  
4 most part until we began this study or until we  
5 became aware of some other publications did another  
6 scan after the completion of the procedure. You  
7 see, that's the difference with this study and how  
8 we have come to identify these endothelial  
9 disruptions.

10 The point I made earlier is that it is  
11 quite likely that if we did epi-aortic scans on  
12 patients after the procedure, we would have found  
13 these endothelial disruptions a long time ago.

14 DR. FERGUSON: So my next extension of  
15 that would be in your opinion, the group's opinion,  
16 should epi-aortic scanning before and after use of  
17 this device be recommended in the use of the  
18 device.

19 I know that a lot of people don't use  
20 epi-aortic scanning, and I understand the  
21 ramifications there, but from a safety standpoint,  
22 I am just bringing that up.

23 DR. KOUCHOUKOS: Currently, epi-aortic  
24 scanning is not standard of care, and it is my  
25 impression that it probably won't be for the

1 foreseeable future. And based on what we know  
2 about the outcomes of the patients who develop  
3 endothelial disruptions, I would say it would not  
4 be necessary.

5 DR. ALLEN: I think as a side that Dr.  
6 Ferguson did not use epi-aortic scanning--we do  
7 1,800 pumps a year at our hospital, and epi-aortic  
8 scanning certainly is not standard of care by any  
9 means. My personal belief and how I would use this  
10 device, epi-aortic scanning is not additive.

11 The thing that I would--

12 DR. FERGUSON: If the disruptions are only  
13 seen with epi-aortic scanning, and if most people  
14 in the study didn't use it, you don't know what the  
15 real incidence of that is. That's what I'm getting  
16 to.

17 DR. ALLEN: I think, though, that the  
18 corollary to that is that we did see endothelial  
19 disruptions, and--

20 DR. FERGUSON: You did with--

21 DR. ALLEN: --with epi-aortic  
22 scanning--and you are right, those did occur--but I  
23 think you have to put that in the clinical context  
24 of what those mean, and the two endothelial  
25 disruptions that were repaired at the single center

1 very early on in the series were--as Dr. Kouchoukos  
2 said, there was really no historical background as  
3 to what those meant, and they acted upon their--I  
4 won't say clinical inexperience--but their lack of  
5 historical background about them.

6 The 10 endothelial disruptions that were  
7 subsequently identified by surgeons, none of those  
8 surgeons acted on those, because they had been  
9 educated and kind of knew now what they were  
10 seeing, and they didn't overreact to--I won't use  
11 Dr. Edmunds' term--but as an overinterpretation of  
12 a very sophisticated imaging technique.

13 DR. FERGUSON: I understand the data very  
14 well. The next question is what percent of the  
15 sites used epiaortic scanning before and after,  
16 because that to me would be the gold standard to  
17 really define whether this is going to turn out to  
18 be significant.

19 DR. KOUCHOUKOS: Someone can provide me  
20 with the exact number of patients who had scanning  
21 in this study, but I think it is over--is it 500 or  
22 thereabouts--

23 MS. CHANG: It's over 500.

24 DR. KOUCHOUKOS: So about 500 of the  
25 patients had epiaortic scanning. And again, I

1 think it is important to emphasize that 78 percent  
2 of these endothelial disruptions were not seen  
3 either by the operating surgeon or by the  
4 anesthesiologist, who is taking perhaps a little  
5 closer look at these ecocardiograms  
6 intraoperatively. They were not recognized. They  
7 were only recognized by the core laboratory.  
8 DR. FERGUSON: I see. Thank you.  
9 That's all I have.  
10 DR. KOUCHOUKOS: Four hundred and nineteen  
11 patients had epiaortic scanning.  
12 DR. TRACY: It's exactly 12 o'clock now,  
13 so at this point, let's take an intermission for  
14 lunch and resume at 1 o'clock.  
15 [Whereupon, at 12 o'clock p.m., the  
16 proceedings were recessed, to reconvene at 1:06  
17 p.m. this same day.]

AFTERNOON SESSION

[1:06 p.m.]

DR. TRACY: If everybody is ready, I'd like to resume the open committee discussion.  
Dr. Ferguson, were you finished with your questions?

DR. FERGUSON: Yes, thank you.

DR. TRACY: Okay. I'll pass it on, then, to Dr. DeMets.

Open Committee Discussion - Continued

DR. DeMETS: Thank you.

Some of the questions I had have either been addressed earlier or addressed in the questions, but I still have a couple more.

Could you tell me a bit more about the rationale for the particular delta, the 5 percent that was decided? Obviously, that's very critical in the size of the study you came up with and the goal that you were after. So could you comment on how that rationale went?

DR. ALLEN: Dr. DeMets, I apologize, but if I could have Dr. Kuntz answer your questions, I would appreciate it.

Thanks.

MS. CHANG: Dr. Kuntz can provide more

1 detail, but that delta was decided on after several  
2 meetings with the FDA, and it was mutually agreed  
3 upon.

4 DR. DeMETS: Okay. So this was not  
5 something that was based on clinical considerations  
6 or what would be important to rule out as a safety  
7 issue, or--

8 MS. CHANG: Not being a statistician, I--

9 DR. DeMETS: Well, it's not a statistical  
10 question. It's actually what clinical difference  
11 do you want to rule out, and I'm just trying to  
12 understand how the 5 percent was arrived at. There  
13 are a lot of statistical implications about that,  
14 but how you got to that decision is what I'm trying  
15 to understand.

16 MS. CHANG: The delta of 5 percent.

17 DR. DeMETS: Why 5 percent.

18 DR. KUNTZ: The deltas are always  
19 inexactly determined in general, and I think that  
20 in our decision, with the baseline rate of 15  
21 percent as established, 5 percent is already 33  
22 percent delta, which is kind of on the high end of  
23 deltas to begin with, but has been in the range for  
24 devices in the cardiovascular arena. But I think  
25 that overall, the final arbitrator was that the

1 clinicians felt that if they could remove emboli  
2 and have a plus or minus 5 percent overall event  
3 rate, they would accept that the device would  
4 remove the emboli, and that was the thing that we  
5 passed around, and that seemed to be the logical  
6 background to the 5 percent decision.

7 DR. DeMETS: Okay. I asked you about the  
8 [inaudible] issue earlier; perhaps I jumped the  
9 gun. I am trying to understand as a non-surgeon  
10 what percent of patients who have this surgery in  
11 fact release particles. Is it 100 percent of them,  
12 or is it half of them?

13 DR. KOUCHOUKOS: We don't know the answer  
14 to that question because we have had until now no  
15 way to assess that. Dr. Barbet, who is here, has  
16 done some studies with ecocardiography and  
17 transcranial doppler suggesting that there is a  
18 large number of particles. The issue there is that  
19 some of these particles are gaseous and some of  
20 them are particulate, so it is difficult to tell.

21 DR. DeMETS: Well, my question actually  
22 has two parts. One, when it is released, how much  
23 is there in a patient, but how many patients is  
24 it--is it almost always? Is it rarely? I don't  
25 know as a non-surgeon.

1 DR. KOUCHOUKOS: Well, this study would  
2 suggest that almost all the patients release--or,  
3 at least over the age of 60--release some  
4 particulate matter, and it is a spectrum,  
5 obviously, depending to a great extent on the  
6 severity of atheromatous disease in the ascending  
7 aorta would determinate how many particles and  
8 their size are released.

9 DR. DeMETS: The second part of the  
10 question which is important is given that they are  
11 released, what percent does this device capture.  
12 And at least in the FDA review, there was a  
13 suggestion that if you average 5 to 5.6, whatever  
14 it was, and it was 25 to 30 particles, that  
15 suggests a 20 percent or so capture rate. Is it  
16 higher or lower than that, because if there is a  
17 lot of it, and you aren't getting much of it, then,  
18 what are you really accomplishing, I guess is what  
19 I am trying to understand.

20 DR. KOUCHOUKOS: Well, again, it is hard  
21 to say. This filter is occlusive, so  
22 theoretically, at least, it should capture all of  
23 the particles that are released proximal to where  
24 the filter is located. So we would think that the  
25 capture rate should be high.



1 DR. DeMETS: Okay, thank you.

2 Another question is that this trial  
3 obviously is not blinded, and when you are trying  
4 to establish equivalence, or at least safety  
5 equivalency, one of the challenges is always that  
6 you have to do a high-quality study. If you don't  
7 do a high-quality study, then it is easy to show or  
8 easier to show two things being equivalent,  
9 whatever you define as equivalent.

10 So my question is given that this is  
11 clearly an ongoing study, what comments can you  
12 make about that there wasn't some bias between the  
13 two procedures, if you will. I'm not sure it is  
14 possible to introduce bias, but at least the  
15 potential seems to me to be to do that. So can you  
16 help me on that a little bit?

17 DR. KUNTZ: Yes, sir. It's an excellent  
18 question.

19 In any study where we are using a device,  
20 especially a surgical study, it is impossible to  
21 blind because the ethics would make it a sham and  
22 impossible.

23 So most of the time--and this goes to the  
24 question about our endpoint per se--the endpoint  
25 had always been focused on safety. And we talked

1 earlier about the fact that we had included the  
2 myocardial infarction part because that is a safety  
3 endpoint that you would be interested in--not  
4 necessarily an efficacy endpoint, because the  
5 filter is north of the heart there.

6 But the bottom line is that the components  
7 of the endpoints were all hard endpoints, that is,  
8 they could be determined by an external  
9 adjudication committee that would hardly be  
10 malleable by someone who had a conflict of  
11 interest. That is, myocardial infarction is a new  
12 enzyme elevation or change in the EKG, death is  
13 death, stroke is stroke. These are very hard,  
14 nonsubjective endpoints, and they tend to help  
15 minimize the influence of bias of unblindedness.

16 So we tried to make sure that the  
17 constellation of [inaudible] endpoints were in fact  
18 hard endpoints, none of which would be too  
19 subjective or that would lend itself to too much  
20 bias.

21 DR. KOUCHOUKOS: Can I amplify on that for  
22 just a moment?

23 DR. DeMETS: Sure.

24 DR. KOUCHOUKOS: The examiners for the  
25 neurologic events were blinded as were the

1 patients. So at least the neurological assessment  
2 as blinded.

3 DR. DeMETS: Well, I think that to your  
4 credit, you worked very hard on that end of the  
5 process. Again, I ask the question as a  
6 non-surgeon, but one can imagine that if you had a  
7 bias about a device, you could be more careful or  
8 more careless, if you will, in actually doing the  
9 surgery and therefore artificially introducing one  
10 group looking better than the other. And like I  
11 said, maybe that's an ignorant question for a  
12 non-surgeon, but the issue is how you deliver it  
13 also affects--even if you have everything blinded,  
14 and the ascertainment bias is minimized--how you  
15 deliver the therapies can also introduce bias, and  
16 I'm just trying to understand that process.

17 DR. ALLEN: Let me give you a real example  
18 where what you are saying could be totally true.  
19 Let's say, for example, as a surgeon, I randomize  
20 my patient to not a filter. I could change my  
21 technique, for example, of the operation and do  
22 less proximal anastomosis or, for example, not use  
23 a side-biting clamp, or do things that might do  
24 what you are saying; but you will recall that in  
25 the demographics, we really specifically looked at

1 those things, so that the things that the operator  
2 could vary that might impact the outcome, such as  
3 doing T-graphs off the mammary instead of putting  
4 them on the aorta, were not occurring, or not using  
5 a side-biting clamp, were not occurring.

6 It is a very valid question, and I think  
7 the size and scope of the study, we did the very  
8 best job we could, not only from a design  
9 standpoint but then from an analysis standpoint of  
10 the operative data, to ensure that that wasn't  
11 going on.

12 DR. DeMETTS: Okay. I think that answers  
13 most of the questions I had--and I still struggle  
14 with the issue of the clinical relevancy in the  
15 surrogate, but I'm not sure what you can say that  
16 you haven't already said.

17 Thanks.

18 DR. TRACY: Dr. Aziz?

19 DR. AZIZ: I just have a few questions,  
20 some sort of technically related.

21 When I look at the cannula, you have the  
22 side arm, and when you cannulate, do you get air in  
23 that side arm, and how do you de-air that?

24 DR. ALLEN: The concern about air comes  
25 from two things, primarily in that the filter is

1 set in a heparin solution and then retracted up  
2 into the device. That still doesn't ensure that  
3 there couldn't be air within this cannula--

4 DR. AZIZ: But when you initially  
5 cannulate, when you put that cannula in right away,  
6 the first time around--

7 DR. ALLEN: You de-air by taking--there is  
8 an operator that goes in the sideport, and it is  
9 de-aired through that, through that one-way valve,  
10 so when you pull that out, the operator allows the  
11 air to flush out, and then, the cannula itself is  
12 de-aired as I described earlier, with venting  
13 through this air release plug.

14 The operator also has the little--I call  
15 it an air release plug--it is the little wet plug  
16 that allows the air to go through it, so it  
17 actually vents through that plug when you put it  
18 in.

19 DR. AZIZ: And that filter comes up at the  
20 time that you have taken the aorta cross-clamp off,  
21 so whatever is there in the ascending aorta, you  
22 capture.

23 DR. ALLEN: Right. The filter is inserted  
24 and deployed right before you take the cross-clamp  
25 off.

1 DR. AZIZ: So at the time that the  
2 protamine is being given, this filter is down, or  
3 is the filter mesh still up?  
4 DR. ALLEN: No. The filter has been  
5 withdrawn, and most of the time, patients have been  
6 decannulated. In my center, we decannulate before  
7 we give protamine. But if you leave your cannula in  
8 when you are giving protamine, the filter has been  
9 withdrawn.  
10 DR. AZIZ: It has been withdrawn. All  
11 right.  
12 Let me just ask a few more questions,  
13 then. When you were giving your talk, you showed,  
14 obviously, two extremes--one with the flap of the  
15 aortic dissection, which I think anybody could see.  
16 The other side-by-side sort of TE, the surface ECO  
17 that you had, there were two--one that had these  
18 strands sort of waving at you--it could be  
19 fibrin--but the other one--and maybe we could have  
20 one of the ECO guys look at that with us--it seemed  
21 to me like there wasn't a disruption in the intima,  
22 but that there was a gap or a gray zone in the  
23 actual media itself.  
24 Could we look at that?  
25 DR. ALLEN: Dr. Weissman, who was

1 our--while we are teeing that up, if you would like  
2 us to show that, we can have Dr. Weissman go over  
3 that specifically, and Dr. Kouchoukos might comment  
4 on that since he does an extensive amount of  
5 epiaortic imaging.

6 DR. KOUCHOUKOS: All of these ECOs were  
7 reviewed by Dr. Weissman, and he indicated to me  
8 and I think will indicate to you if there is any  
9 question about it that there are no medial injuries  
10 at all that were identified. These were all  
11 endothelial or intimal injury.

12 DR. WEISSMAN: That is correct.

13 I am Neil Weissman, and I was the director  
14 of the ECO core lab for this study. I am a  
15 cardiologist at Washington Hospital Center.

16 I have no financial conflict of interest;  
17 they gave a grant to the hospital for my work on  
18 this.

19 There have been a number of different  
20 points brought up, and I think they have been  
21 answered very well, methodological issues and the  
22 extent of these endothelial disruptions and what  
23 they look like. So as we boot this up, let me just  
24 go through a couple of those things.

25 DR. AZIZ: Sure.

1 DR. WEISSMAN: People are asking about the  
2 methodology, and the methodology--and I wrote the  
3 protocol--was to do the scanning in a transverse  
4 manner, starting right proximally and capturing  
5 proximally at least 5 beats, and then moving one  
6 centimeter at a time, capturing 5 beats one  
7 centimeter, and so forth. So it was pretty  
8 methodological. And then, you do transverse  
9 imaging across the ascending aorta.

10 As that was done, it had to be annotated  
11 on the screen or verbally to let me know where they  
12 are, and that's how we got location information.

13 What you saw--were you referring to this  
14 picture or earlier on--

15 [Slide.]

16 DR. AZIZ: There was another one.

17 DR. WEISSMAN: Yes. I think one of the  
18 things--I don't know if it came through  
19 completely--was that this is the more typical  
20 thing, which I have to admit I have trouble seeing  
21 here. There is a little wiggly right over  
22 there--and you have got to turn the lights down.

23 These images in the core lab were reviewed  
24 three times--once by a technician, who would write  
25 down their preliminary results; then, independently



1 by me; then, after I reviewed it, I looked at what  
2 the technician said to see if I missed anything and  
3 went to the spot where the technician thought they  
4 might have seen something to see if I missed  
5 anything.

6 That is why I think 78 percent of these  
7 things were not seen by the anesthesiologists or  
8 the surgeons at the time.

9 So I think the terms used here--"strands"  
10 and "dings" and "footprints"--are all pretty  
11 accurate.

12 [Slide.]

13 DR. WEISSMAN: So this is worst case  
14 scenario here. Where you see that thing sort of  
15 flipping is definitely among the worst case  
16 scenario. These, you could see from across the  
17 room on a projection with the lights on, okay?  
18 This was not the typical thing.

19 So, show me what area you are concerned  
20 about?

21 [Dr. Aziz indicating.]

22 DR. WEISSMAN: Actually, you can tell that  
23 that is coming away from the wall. Right there,  
24 you see it is coming away from the wall. That is  
25 not even part of the wall. The intima is that very

1 light ECO within it. The intima there is probably  
2 on the order of 2 or 3 millimeters thick, and it is  
3 the ECO-density of the flap that is being lifted  
4 up.

5 That is essentially a monolayer, because  
6 the intima is so thin that you aren't seeing it.  
7 That is extra-aortic that you are seeing that  
8 little lifting.

9 DR. AZIZ: Okay, good. Again, now that in  
10 a sense we have identified that you do have these  
11 endothelial disruptions, in my own mind, apart from  
12 the sites where you have an aortic cross-clamp on  
13 them, it could happen either in the control group  
14 or in the other group, but in the group where the  
15 filter is in place, how do you think that is  
16 causing that? Do you think it is the tip of the  
17 sheath that you are putting in, and is it occurring  
18 at the posterior wall of the aorta?

19 For me, that is an important issue.

20 DR. WEISSMAN: And I'm going to  
21 defer--since I was not in the operating room, I'm  
22 going to defer that to the surgeons to comment on.  
23 Again, the results show that there were these  
24 little disruptions distributed along the whole  
25 ascending aorta. To conjecture how they arose, I'm

1 not going to do that; I just read the images.

2 DR. KOUCHOUKOS: One has to assume that  
3 some of these were caused by the filter itself,  
4 although clearly not all of them were, even in the  
5 filter groups, because we saw them distributed in  
6 areas where the filter was not located.

7 The metal rim of the filter, it is  
8 conceivable, could create a small intimal  
9 disruption, and that is probably the explanation  
10 for why they occurred.

11 DR. AZIZ: But clearly, the goal of the  
12 filter is to prevent the emboli going upstream  
13 when you take the cross-clamp off. But I think as  
14 was mentioned in the FDA presentation, you could  
15 envision where particulate matter is caught in that  
16 mesh. You take the aortic cross-clamp off, and  
17 blood from upstream obviously hits the mesh on the  
18 other side and dislodges particles going downward,  
19 and maybe that is what is responsible for the  
20 Q-wave MI.

21 What do you think about that?

22 DR. KOUCHOUKOS: It is theoretically  
23 possible, and I think that is one--it is fortuitous  
24 that we didn't look at myocardial infarction, but  
25 again, we found no difference in the prevalence of

1 myocardial infarction--

2 DR. AZIZ: But theoretically, it doesn't  
3 protect against that, and it could predispose to  
4 something going in the reverse direction.

5 DR. ALLEN: Actually, I spoke with Dr.  
6 Ferguson, and the design of the filter is a  
7 windsock design. Myocardial infarction was  
8 specifically put in as a safety endpoint for  
9 concerns for that very point. But in designing  
10 that windsock which drapes down over--it is like my  
11 son does when fishing for tadpoles--you are  
12 catching them in the windsock that falls down  
13 below, and you wouldn't expect when you have that  
14 pressure change for that to blow it out of that  
15 windsock as you would, for example, if it were a  
16 flat filter, like a seine. And that's why it was  
17 designed like that.

18 DR. AZIZ: Okay. I have just a couple of  
19 other questions.

20 Clearly, the majority of patients--we all  
21 do bypass cases, and again, the case mix here was  
22 mainly patients who came for bypass  
23 surgery--patients who come in for valve operations,  
24 particularly aortic valve operations, obviously,  
25 you had more calcium and bits of material that

1 could come off there. When you analyzed your  
2 particulate trapping, did you find that it was  
3 higher in patients having valve operations, aortic  
4 valve versus mitral versus--did you look at that  
5 subset?

6 DR. ALLEN: We did, and what we primarily  
7 found was that the vast majority of the histologic  
8 material that we have treated was atheromatous.

9 There is a wide range of material such as  
10 calcific material, organized clot that looked like  
11 it came from LV or left atrial appendage. We  
12 didn't specifically see a correlation between if  
13 you had a valve and you had more, for example,  
14 calcium versus atheromatous.

15 DR. AZIZ: Particularly in aortic valve.

16 The other thing--and I know you can't do  
17 it now--but in the study design, the reason you  
18 decided not to use TCD monitoring was because--

19 DR. KUNTZ: I am not an expert in TCD  
20 monitoring, but we have discussed this with other  
21 trials. It is not clear even in carotid  
22 interventions that TCD monitoring can be very  
23 helpful, because the high-intensity transience that  
24 occurs with that occurs, for example, in every  
25 operation for carotid enterectomy, so it is clearly

1 possibly an overly sensitive measure of  
2 high-intensity transience, whatever that is, with  
3 respect to ultrasound as a measure of emboli.

4       Clearly, they do measure emboli, but they  
5 may be measuring other things as well, because it  
6 is so frequent. Now, there is a lot of interest in  
7 looking at transcranial dopplers, and you have to  
8 do bilateral transcranial dopplers during that, and  
9 I think that might have also been logistically a  
10 little bit difficult during the operation. But I  
11 think because of the lack of a good sensitivity  
12 specificity profile for that test per se, it wasn't  
13 used.

14       DR. AZIZ: Looking at it, you could have  
15 seen, for example, compared with the control  
16 group--just take the CABs, where you are operating  
17 up the aorta--you might have seen less numbers.  
18 Clearly, it has been shown that there is a  
19 correlation between the number of hits you get on  
20 TCD and cognitive dysfunction.

21       DR. ALLEN: I think the difference is it  
22 is hard to know--we can actually see that certainly  
23 the filter doesn't capture air or gaseous emboli;  
24 it captures particulate matter. And I think that's  
25 the hard thing with transcutaneous dopplers, that

1 it lumps everything as to specks on a spectrum, and  
2 whether it is a particulate material or a gaseous  
3 emboli, they all look the same.

4 DR. AZIZ: One other thing--did you see  
5 any correlation in people who did get these  
6 endothelial let's say injuries--was there a  
7 correlation between the thinness of the aortic wall  
8 and the size of the aorta--in other words, big,  
9 dilated aortas were more prone to getting it?

10 DR. ALLEN: That's a great question, and  
11 actually, what we looked at was whether the size of  
12 the filter, which obviously corresponds to the size  
13 of the aorta, correlated to an increase or decrease  
14 in endothelial disruptions, and it didn't. It just  
15 wasn't correlated. So was a larger filter size  
16 more prone to causing endothelial disruptions--no.  
17 Was a smaller filter size less prone--no.

18 DR. AZIZ: Thank you.

19 DR. TRACY: Dr. Krucoff?

20 DR. KRUCOFF: Let me just ask a couple of  
21 quick questions. First, you mentioned at the very  
22 beginning of your presentation, I believe, a number  
23 for the percentage of patients who were screened  
24 relative to those actually enrolled.

25 DR. ALLEN: Yes. About 15 percent of

1 patients who were screened were eventually enrolled  
2 in the study.

3 DR. KRUCOFF: Do you have any sense of who  
4 the other 85 percent were or why they were--were  
5 they eligible but just didn't want to be in a  
6 research protocol, or are we really talking about a  
7 patient population that comes from 15 percent of  
8 the open heart surgery universe?

9 DR. ALLEN: I think there are a lot of  
10 reasons, and you touched on both of them. I think  
11 Dr. Kouchoukos' slide showing the types of patients  
12 that cardiac surgeons are operating on any more,  
13 the cardiologists just don't send us patients who  
14 are low-risk, and this was a safety study looking  
15 specifically at low-risk patients, and although we  
16 enrolled a lot of patients, it took us 20 months to  
17 do that. Even in the low-risk group, there  
18 certainly would be some patients who opted not to  
19 do it, but quite honestly, the enrollment--as  
20 surgeons became familiar with this device and saw  
21 what they were capturing, enrollment in the study  
22 was pretty accelerated, and surgeons wanted to  
23 participate in the study.

24 DR. KRUCOFF: I also just wanted to  
25 ask--and thank you for passing the model around,



1 because that helped me compared to the pictures--as  
2 I look at the actual retraction process--and I want  
3 to ask you a little bit about the windsock  
4 design--when that pulls back in, it seems to me  
5 that there is a point when it is partially  
6 retracted where not the tip of the windsock but the  
7 upper part is actually just kind of flattened. And  
8 I would worry about whether that was capable of  
9 dumping debris that was not down in the windsock  
10 but that was higher up. Have you all--and I just  
11 wasn't aware, at least in our panel pack, of any  
12 sort of bench-testing or preclinical modeling that  
13 has been done to see at what point or what size  
14 particles would be dumped rather than captured.

15 DR. ALLEN: It is a preclinical test, and  
16 I'll let Jean speak to that.

17 MS. CHANG: Yes. We did extensive  
18 preclinical tests, and our [inaudible] with the  
19 panel package includes the clinical information  
20 there. Our preclinical test was with little  
21 polyester beads, polystyrene beads. They are like  
22 little pinballs, so that when you do this, when you  
23 capture, you measure percent capture. And our  
24 capture rate was well above 80 percent.

25 DR. KRUCOFF: Eighty percent. Are these

1 sticky beads?

2 MS. CHANG: No. So it is worst case.

3 DR. KRUCOFF: Because one thing that I  
4 take--and this is from looking at your own  
5 pictures--is that a lot of the particles that you  
6 photographed are not down in the windsock; they are  
7 up in the sort of billowing part of the material.  
8 When I look at these pictures, and just thinking  
9 about what sticky particles, lipid particles or  
10 thrombus particles--there is no question that when  
11 you get the big one in here, that's down at the tip  
12 of the windsock, but a lot of these others are not.

13 DR. KOUCHOUKOS: When these were removed,  
14 the technician who was in the operating room was  
15 responsible to collect it, to flatten it out, and  
16 to display it; and I suspect that part of that is a  
17 flattening effect that was done so that we could  
18 get a photograph of the material in the filter. So  
19 it is partially related to that.

20 To the question about possible loss of  
21 material, the filter is removed after 20 or 30  
22 minutes, so we would surmise that most of the  
23 embolization would have occurred, so it's possible  
24 that we missed some of the material that might have  
25 passed through the filter as you are removing it,

1 but I think the probability of that resulting in  
2 the loss of a large number of particles would be  
3 very, very strong.

4 DR. KRUCOFF: Okay, I take your point, Dr.  
5 Kouchoukos. I was just sitting here looking at  
6 this and was thinking about, for example, some of  
7 the pressure shifts that Dr. Swain described, and  
8 if that was actually not trapped down in the tail  
9 but was sitting up in the higher, whether a sudden  
10 shift of pressure would dislodge it and do  
11 something else.

12 I think most of my comments have been  
13 mentioned. I think the real issue here is a  
14 denominator one. In a very complex array of end  
15 organ problems, and even the precedents that Dr.  
16 Kuntz mentioned in the SAFER study, where the  
17 PercuSurge device was used, actually, there was no  
18 actual or even attempt to measure particulate  
19 capture in that study; that was driven entirely by  
20 a clinical measure of an end organ whose effect  
21 could be imputed to probably particle capture, but  
22 actually, it was purely a clinical measure--and the  
23 IIb/IIIa is the same thing. We all sort of suspect  
24 in the angioplasty environment that we may create  
25 particulate matter that is responsible for end

1 organ wounds to the heart, but the fact that a  
2 IIb/IIIa inhibitor corrects the wounds to the heart  
3 to some degree, or that a distal protection system  
4 in a vein graft protects the wounds to the heart,  
5 we have never really directly measured the role of  
6 particles. And what you all deal with in the brain  
7 and renal failure and all of the end organ effects,  
8 again, I think has been clearly recognized as  
9 multifactorial. Some of that is probably  
10 particulate embolization, and some of it is  
11 probably noncirculatory arrest and predisposing  
12 factors and transient hypertension and everything  
13 else in the world that comes with that.

14 So I think that even if we start with the  
15 end organ denominator where, obviously, from a  
16 patient misery point of view, you would love to  
17 find a better way of bringing patients through open  
18 heart, the particulate component may only be a  
19 subset of that.

20 And then, my difficulty with the other  
21 denominator is that, based on your data, it sounds  
22 to me like at least 96.8 percent of patients who  
23 undergo open heart surgery have something  
24 capturable in a filter, and certainly, 96.8 percent  
25 don't have measurable end organ effects. And

1 again, there is no question when you look at a  
2 thrombus as huge as the one you have pictured in  
3 one of these pictures, you've got to believe  
4 intuitively that you have benefitted the patient by  
5 pulling that out of his body rather than letting it  
6 go wherever the heck it was going to go.

7 Ultimately, sorting this out is tough, and  
8 I think that point has been made.

9 And then, the safety issues become more  
10 preeminent, because understand what is efficacy is  
11 hard.

12 One suggestion that I would like to  
13 amplify as my last comment is if you could in fact  
14 correlate--and this is just speaking from my own  
15 seat--but if you could correlate that descriptors  
16 that you would think, prospectively, for a patient  
17 preoperatively, would identify a higher likelihood  
18 of having an embolic untoward event, whether they  
19 are morphologic descriptors of the aorta or low EFs  
20 or whatever you would think would be the ones that  
21 would say this population is likely to have more  
22 particles, or more frequency of particles, or  
23 whatever, and correlate that in your own dataset to  
24 a higher capture rate of particles with the filter,  
25 to me, that would be at least a first step toward

1 saying maybe what this filter is able to do is to  
2 take higher-risk patients--and again, this is my  
3 question to the FDA statistician--one of the things  
4 that I think is at least a possible interpretation  
5 of what Gerry actually did was to say that what  
6 you're doing is taking the highest-risk patients  
7 and creating a more linear risk by capturing more  
8 particles. It's just that we don't have the  
9 descriptors match or really know what the  
10 descriptors model would be to say are those really  
11 the higher-risk patients or is that just sort of an  
12 [inaudible] finding.

13 I think if you could build that model  
14 somewhat, that higher risk of whatever you think  
15 would predict embolic events and a higher capture  
16 of particles is a correlation, then you could start  
17 to think about what would you look for in a smaller  
18 venue as a way of really showing a benefit to  
19 putting a filter in.

20 I just want to acknowledge that what is  
21 very clear, particularly from our two surgeons  
22 presenting, is a clear desire to try to prevent  
23 these types of untoward events with an intuitively  
24 obvious kind of mechanical approach, but in the  
25 face of a real difficult trial planning environment

1 to the point where I'm not sure that this one study  
2 achieves everything that you would want out of it.

3 DR. KUNTZ: Mitch, can I make a few points  
4 about your comments? I think I would actually like  
5 to amplify some of your points.

6 If we look at the SAFER trial per se,  
7 there were actually two studies that demonstrated a  
8 relationship between the emboli and the outcome.  
9 There is the initial John Webb study done, where he  
10 actually counted the number of particles, and that  
11 was associated with the amount of CK-MB that was  
12 elevated. That was the pilot study before the  
13 SAFER study.

14 The second one was done in my institution  
15 by Campbell Rogers and compared the amount of  
16 emboli removed by the EPI device compared to that  
17 by the SAFER device, showing that the devices were  
18 equivalent once you controlled for the amount of  
19 particles and their enzyme elevations.

20 So we do start to see some connection  
21 between emboli and what is causing damage to the  
22 heart. The other thing about the SAFER trial was  
23 that the device that was used to remove particles,  
24 the PercuSurge device, compared to doing nothing at  
25 all, the main and only difference was that one

1 device removed particles and one didn't, and there  
2 was a 50 percent reduction in MI.

3 So I think that that is actually pretty  
4 solid proof that the emboli were measurable in  
5 their impact on the heart. It's about as solid as  
6 you can get, I think, from an 800-patient  
7 randomized trial.

8 So we do know that when using the heart as  
9 a surrogate, as an organ that has small vessels,  
10 like all organs do, that can be clogged up by  
11 emboli that have necrosis and damage as  
12 demonstrated by IIb/IIIa inhibitor trials and so  
13 on, that means something to the organ with a  
14 readily available measurable outcome, that we  
15 actually do prevent cell death by removing emboli,  
16 at least in that [inaudible].

17 So the next transition to say that  
18 actually capturing emboli in the body in general is  
19 not so much of a high-falluting notion or theory as  
20 something, as Dr. Marler mentioned earlier, like  
21 ECIC bypass or other things that are intuitive but  
22 maybe don't have as much connection.

23 So I think the evidence is growing, and  
24 there are lots of different venues now for  
25 investigations in which the notion of putting a



1 filter on a device is actually already intuitively  
2 being planned. For example, many carotid stents  
3 are now being packaged with filters where there is  
4 no attempt to look at the filter component part of  
5 safety; it just makes sense to put a filter on  
6 there, because emboli don't make any sense if they  
7 go to the brain, for example.

8 The same is starting to be done with  
9 studies in the renal area as well.

10 So I think there is a growing body of  
11 evidence that small emboli are not good for organs;  
12 in the heart, I think it is established; and as we  
13 start to look at other organs, that notion is  
14 starting to grow. And this trial was caught in the  
15 middle of that in having technology available for  
16 surgery using background information, like Dr.  
17 Kouchoukos had shown for years that this was a  
18 problem, and we were caught in the crosshairs of  
19 being able to demonstrate the emboli being removed  
20 with the growing idea that emboli can be measured  
21 in some organs but not in all organs, and  
22 potentially, if this trial were to be repeated in a  
23 year or two, maybe we would have much more  
24 sensitive measures. But that's kind of how we put  
25 that in perspective.

1 DR. KRUCOFF: I understand, and I take  
2 many of your points, Rick. On the coronary distal  
3 protection, as you well know, there is still  
4 ongoing dialogue about occlusive and nonocclusive  
5 and whether there are other elements besides  
6 particles. I don't want to intimate for a second  
7 that there is any proof that particles are good for  
8 anything.

9 I think the reverse side here, and  
10 particularly in a large vessel that feeds virtually  
11 every other vessel in the body, is what size of  
12 particles, how many of them, and at what cost is  
13 unfortunately where I think "caught in the middle"  
14 is probably a good phrase.

15 DR. TRACY: Dr. Laskey?

16 DR. LASKEY: When one gets to this end of  
17 the table, one had better be brief--or insightful.  
18 I'll be both.

19 [Laughter.]

20 Rick, you can't compare PercuSurge and  
21 IIb/IIIa inhibitors in the same breath. I mean,  
22 IIb/IIIa inhibitors don't do anything for large  
23 embolic goobers, yet they decrease the rate of  
24 necrosis.

25 So this is, I think, without being

1 dismissive, a rare outcome of a very prevalent  
2 disease. Atherosclerosis is diffuse, and without  
3 simplifying this any further, if you do diagnostic  
4 catheterization, you take a wire out of the body  
5 that has been in the body for 30 seconds, you are  
6 wiping of thrombi, you are wiping off platelets,  
7 you are wiping off clots--and yet the rate of  
8 stroke or embolization or other horrible things  
9 during cath is so acceptable that we don't even  
10 think about it. So we are not thinking about  
11 putting filters on our diagnostic caths, but those  
12 clots are there--let there be no debate about  
13 that--and it is a question of how sensitive the  
14 test is to look for them. So if you do scanning EM  
15 on your guidewires, you are going to find it.

16 And if we parlay that to where we are  
17 today, atherosclerosis in the ascending aorta is  
18 virtually present in 100 percent of the patients  
19 that you all operate on, yet the event rate--thank  
20 goodness--is acceptably low--1 to 2 percent, maybe  
21 3 percent adverse embolic-type event rate. Now, it  
22 would be great if that were zero, but I don't think  
23 that's why we are here today. But I think we do  
24 need to be careful about signal and noise and  
25 reducing an event rate which agreeably is low but

1     could be lower.

2             So I am not entirely sure what we are all

3     about here, and we are putting an instrument to the

4     ascending aorta purportedly with the aim of

5     collecting debris, but it is in there transiently,

6     it is in there at a moment in time that you just

7     sort of arbitrarily said is the moment of risk, and

8     yet introducing the trocar [phonetic] into the

9     ascending aorta could just as well release debris.

10    Case-in-point--virtually every, single brachial

11    arteriotomy I have ever done in a patient over 50

12    [inaudible] catheterization, you open the artery,

13    and there is plaque right there. I mean, it is a

14    universally present disease in these patients, and

15    yet the dread sequelae are fairly infrequent, and

16    developing sensitive tools is critical--you have

17    heard that; I don't need to repeat that. Assays

18    for efficacy are sorely needed in this. And

19    certainly, capturing the universe of the period of

20    risk is critical. I think you need to be in there

21    for the whole period at which the patient is at

22    risk, and that includes from the moment you

23    instrument or manipulate the aorta to the time that

24    you go on or off bypass and give protamine.

25             Those are just more editorial-type

1    comments, but that is what we are all grappling  
2    with here, and I'm sure it has a lot to do with why  
3    we are not and may not ever, until we have  
4    developed an incredibly sensitive test, be able to  
5    demonstrate the efficacy of these tools.

6            But my one question to you is why 75  
7    percent. Why did you pick that? Why didn't you  
8    pick 95 percent? Seventy-five percent is so little  
9    to my mind, given the prevalence, the universal  
10   prevalence, of this stuff in these aortas. Why not  
11   go for a higher figure?

12           DR. ALLEN: We certainly could have gone  
13   for a higher figure, and if it had been anything  
14   under 96.8, we would have met that higher figure.  
15   There isn't a historical background. We have never  
16   had the ability to place a filter in the ascending  
17   aorta before, so we don't know what those numbers  
18   are.

19           I think Dr. Edmunds early on was quite  
20   astute when he said the device isn't designed to  
21   capture all emboli, but is it better to have a  
22   device that captures some emboli or just not have  
23   the device at all, and we let those emboli go.

24           I guess that's the crux of the  
25   philosophical debate. I guess as a surgeon, having

1 put this device in over 100 patients and seen what  
2 it pulls out, I'd rather have something that I can  
3 use to pull out some of those emboli. You're  
4 right--I'd love to have a device that captured  
5 everything, but that's not what I have. I've got a  
6 device that captures a lot of emboli, and that has  
7 intuitively got to be good for the patient, and as  
8 you pointed out, the event rates that we are  
9 measuring are so small that to power a study to  
10 show pertinent and important reductions in those  
11 event rates would require such a huge study that  
12 from a practicality standpoint, it is not  
13 reasonable. So you design the trial to capture  
14 particles and show the device is safe. And I think  
15 that we have accomplished that by meeting both of  
16 our primary endpoints.

17 DR. KUNTZ: Just one statistical thing,  
18 too--you have to have a little bit of room for your  
19 [inaudible] about the number that you're trying to  
20 show so that you can demonstrate that the number  
21 you have has a lower battery that is above 75  
22 percent.

23 DR. LASKEY: I understand that, Rick, but  
24 really, a device which is 75 percent efficient is  
25 nothing that I would want to fiddle around with.

1 I think that what we're talking about is  
2 efficiency on two levels, but certainly, the  
3 efficiency of retrieval--we have no idea what the  
4 efficiency of retrieval is. Is it 2 percent? Is  
5 it 20 percent? Is it close to something on the  
6 order of 90 percent?

7 Don't equate the measure of efficacy that  
8 you have here with the measure of efficiency of the  
9 device. We don't know how much of that stuff at  
10 risk of embolization is actually retrieved. Yes,  
11 96 percent of your device have done so, but that's  
12 not the same as how much of the stuff which is at  
13 risk of embolizing is actually retrievable ergo how  
14 much do you lower the risk of embolization.

15 Dr. White, it's all yours.

16 DR. WHITE: I can't be insightful, so I  
17 will be brief.

18 [Laughter.]

19 I am intrigued by emboli protection. I am  
20 involved with emboli protection in multiple organs,  
21 as Rick knows, and I like the intuitive argument  
22 that I never saw an emboli that I liked. I think  
23 the problem is--and you guys probably know this  
24 better than I do, but for the rest of our panel  
25 members--as Dr. Laskey just said, taking five

1 emboli out of circulation is a great idea, but not  
2 if 500 get by. And I think there is a threshold at  
3 which we would decide that there is efficiency or  
4 efficacy in taking those--even partial prevention  
5 is better than no prevention--but there is a  
6 threshold where the partial prevention meets the  
7 road, and that is, I think--we have heard from  
8 multiple people who keep trying to get to this  
9 denominator.

10 The only thing I found in the whole  
11 pack--and tell me if I am wrong about this--it is  
12 under agency review for us; I don't know where you  
13 guys have it--it is under "Summary of FDA Methods"  
14 on page 3. At the very bottom, it says that you  
15 did an in vitro study with these 120-micron beads,  
16 and that your acceptance criterion was to capture  
17 50 percent of those beads. Am I--I don't want to  
18 go faster than you can go--it is Number 5 in the  
19 Agency's summary. I don't know where it came from  
20 in the primary pack. Do they have the Agency  
21 summary?

22 DR. ZUCKERMAN: Yes. It is in Dr. Wentz'  
23 initial review.

24 DR. WHITE: On page 3--do you see what I  
25 am referring to? It is Number 5 at the bottom of



1 that page. And what it shows is that in vitro, I  
2 guess this is, if you said 50 percent was captured,  
3 and that was your goal, that was your acceptance  
4 rate. Is that what you set up?

5 MS. CHANG: That was the lower threshold.

6 DR. WHITE: So you would have been happy  
7 with a 50 percent capture rate?

8 MS. CHANG: These were polystyrene beads.  
9 We were testing the worst case. So they were  
10 literally like pinballs.

11 DR. WHITE: Yes.

12 MS. CHANG: Now, in the body, the  
13 particulates would be more sticky. So yes, we  
14 chose that lower threshold--this is based on  
15 initial bench-testing. Our values ranged from the  
16 large filters to small filters, so it was sometimes  
17 as high as 80, 90 percent.

18 DR. WHITE: Right. The design itself  
19 looks to me like it would be better than 50  
20 percent. If I put that in a plastic tube and blew  
21 balls to it, just as you mentioned with the  
22 windsock, we ought to capture those balls. So when  
23 you start to talk about the failure mechanism, the  
24 reason that you would fail to capture them, is that  
25 because you are not getting uniform deployment of

1 the ring in the aorta? Are the balls sneaking  
2 around the ring? Why do you fail to capture a  
3 ball?

4 MS. CHANG: There is a little bit of a  
5 teeny gap right here, so again, on the smaller  
6 filter, this gap percentage is smaller. In a  
7 larger filter, it is a lower percentage. And it is  
8 just the way the filter is deployed, because again,  
9 there is a little bit of--

10 DR. WHITE: Because the next question I  
11 have for you that again goes to the clinical arena  
12 is could the surgeons judge the adequacy of the  
13 deployment. For example, the only experience I  
14 have like this is the EPI filter, which is a  
15 nitinol ring. And we image that radiologically and  
16 actually find a reasonable number of times that we  
17 have to adjust that filter to get it to actually  
18 oppose the wall; otherwise, it cants and tilts, and  
19 we don't have apposition.

20 Do you guys have any direct control over  
21 the apposition of this filter? Do you know if it  
22 is cocked, or do you know if it is--do you know  
23 what I mean--canted in one way?

24 DR. ALLEN: I understand what you are  
25 asking, and I think the answer to that is that

1    there is a tactile sense to the device that you can  
2    tell when it is certainly deployed, and when you  
3    are having problems with deployment, if you had a  
4    problem with deployment, that tactile sensation  
5    gives you that feedback, and you need to make  
6    adjustments.

7           DR. WHITE: But you don't actually look at  
8    the ring. You are not seeing it; it is going  
9    through the wall of the aorta.

10          DR. ALLEN: No.

11          DR. WHITE: And your ultrasound--Dr.  
12    Kouchoukos, when you image these with your imaging,  
13    can you actually see the ring on ultrasound?

14          DR. KOUCHOUKOS: It would be very  
15    difficult, because the probe is a 7-megahertz probe  
16    that sits on top of the aorta, and it would be  
17    impossible, really, to effectively image that area.  
18    And TEE doesn't give you a good image of that  
19    particular part of the aorta, as you well know.

20          DR. WHITE: And then--do you want to say  
21    something else?

22          DR. ALLEN: I just wanted to come back to  
23    your comment about do things go around. I think  
24    Jean talked about the small area at the top, but  
25    you all know when you look at femoral arteries or,

1 as surgeons, you look at the ascending aorta, that  
2 these aren't perfect manufactured tubes,  
3 particularly when you have disease in them--they  
4 have nooks and crannies and stalactites and  
5 stalagmites--and placing the device down in here,  
6 you wouldn't expect to get perfect apposition of  
7 the nitinol ring to the inside diameter of the  
8 vessel wall.

9 DR. WHITE: I agree, but that's the whole  
10 purpose of reaching the threshold where--your  
11 argument which you are trying to make it, which is  
12 that without measurable efficacy, less emboli are  
13 better than more emboli. And I would like to get  
14 comfortable that we are taking out most of the  
15 emboli, because taking out one out of 300  
16 million--I wouldn't agree with your argument if you  
17 were taking out a vast minority of the emboli.

18 The other thing is in sizing--I notice  
19 that your device comes in 3 mm and 4 mm increments,  
20 and you size the aorta with a device that measures  
21 the outside diameter of the aorta. Have you looked  
22 at or measured any internal consistency among  
23 operators at being able to fit the aorta and get it  
24 right? Are you able to judge the right size of  
25 that aorta? How much do you miss? What is the

1 variability there?

2 DR. ALLEN: Actually, we make estimates of  
3 the thickness of the ascending aorta, and when you  
4 are sizing the device, you step down so that you  
5 take into account the internal diameter--

6 DR. WHITE: My question is is there any  
7 measurement of how accurate any given surgeon is,  
8 or between surgeons, at making the right choice for  
9 the filter in order to fit that aorta, because if  
10 you are going to have high-efficiency capture, you  
11 really want to be measuring very carefully, or at  
12 least be very on-the-money about the right size.

13 Do you have any measurement of that  
14 consistency or accuracy or the ability to correctly  
15 deploy the filter?

16 DR. ALLEN: I don't have specific measures  
17 where we measured how effective the surgeons were  
18 at measuring the device. The devices that you use  
19 to measure the size of the ascending aorta are  
20 graft-sizers that vascular and cardiac surgeons use  
21 every day, and to measure the size of the ascending  
22 aorta, it's not rocket science.

23 DR. WHITE: But with a 3 mm sensitivity, a  
24 little bit of mistakes make for incomplete loops,  
25 and again, if we are talking about that less emboli

1 are better than more emboli, I would want to know  
2 that you are accurately doing the best you can to  
3 screen all of those out.

4 DR. ALLEN: I can just tell you that as a  
5 surgeon, I think I do a pretty accurate job of  
6 using the vascular sizers to tell me what the size  
7 of the ascending aorta is and in choosing the  
8 appropriate filter size.

9 DR. WHITE: It would be an interesting  
10 experiment to actually do it in a model, even an  
11 animal model, and measure your emboli three  
12 different times or two different surgeons or take a  
13 couple of your fellows--you could tell us that, "Do  
14 you know what--we have five guys do it, and it's an  
15 easy thing to do; this thing fits no problem," or  
16 you could tell us that there is a tremendous amount  
17 of variability between surgeons and your ability to  
18 capture these or fit this device appropriately.

19 MS. CHANG: Dr. White, actually, I'm  
20 sorry--in Europe when we did our first cases in the  
21 early 1990s, we did that correlation, and that's  
22 how we came up with the aortic sizer. So we would  
23 have the surgeons--these were the first 20 or so  
24 cases--do the aortic sizers and also do MJ  
25 [phonetic], and they correlated.

1 DR. WHITE: How did you measure that?  
2 MS. CHANG: I think they did imaging.  
3 DR. WHITE: They imagined the loop?  
4 MS. CHANG: The cross-section.  
5 DR. WHITE: With what?  
6 MS. CHANG: Epi-aortic.  
7 DR. WHITE: With what--ultrasound?  
8 MS. CHANG: Yes.  
9 DR. WHITE: Okay.  
10 The only other issue is that, again going  
11 back to the efficiency, the filters only filter  
12 everything bigger than they are, and that's a  
13 debate that we all have about distal protection  
14 devices. If you choose a filter which is nice to  
15 use, and you can have flow and all those things,  
16 you have to give up everything smaller than the 120  
17 microns if that is the size you pick, or if you  
18 pick a smaller pore size, you get problems with  
19 that as well.  
20 But that goes to the efficiency of fewer  
21 emboli, and that is that we really don't know that  
22 the big emboli are the problem. As you have shown  
23 in your graph, the smaller emboli block the smaller  
24 brain arteries. So you may be picking out the big  
25 chunks, and the little stuff is still causing a

1 problem, which is why I think we get back to if  
2 there is a threshold where taking emboli out is  
3 good, then we would like to get some measure of  
4 that efficacy. That is why I think just the  
5 rationale for me that taking out some emboli is  
6 better than no emboli is difficult to get  
7 enthusiastic about.

8 That's all I have.

9 DR. TRACY: I have just a couple of quick  
10 questions.

11 This study was specifically not done in  
12 people with very severe aortic disease, and yet it  
13 did show, whether these things are clinically  
14 relevant or not, more evidence for aortic  
15 disruption than not having a device deployed.

16 What is there that tells us that if we  
17 move into people with much more severe aortic  
18 disease that we won't have greater consequences of  
19 an increased number of aortic disruptions? Why  
20 would it be safe?

21 DR. ALLEN: If you recall the odds ratio  
22 table that I showed you at the very end of the  
23 slides, for example, in the higher-risk patients,  
24 there was a specific component there that looked at  
25 imaged plaque what grade the aorta was, and there



1 was a correlation between grades of plaque versus  
2 whether or not an EDS was occurring.

3 I can give you that assurance, that there  
4 didn't appear to be an increased incidence of EDS  
5 in patients who had worse aortas. What I can't  
6 come back to is that it is an assumption that more  
7 atheroembolism is generated in patients who have  
8 worse aortas--but I can't give you the denominator,  
9 as Dr. Laskey and Dr. White have both asked.

10 DR. TRACY: And the other question or  
11 comment I have is that there are other things that  
12 cause neurologic events. How do we know that the  
13 things that aren't getting by aren't the things  
14 that would be causing problems? I am struck by  
15 this lack of any kind of endpoint to look at that  
16 with.

17 DR. ALLEN: I think the whole flavor of  
18 particularly the last several questions illustrates  
19 the struggle that the panel has with the intuitive  
20 notion that particle removal is bad, but the study,  
21 because it is a safety equivalence study, doesn't  
22 show this dramatic reduction in events. And it  
23 relates, as we tried to go through, and it is--you  
24 want this device to be able to demonstrate a  
25 reduction of events, but as Dr. Laskey pointed out,

1 the events that fortunately occur with excellent  
2 cardiac surgical care today are not that high. So  
3 to power and design studies that can appropriately  
4 measure those events is very difficult.

5 I think the issue is that the composite  
6 events, or the events that comprise that composite,  
7 were chosen to look at the device as far as the  
8 safety standpoint, and that's how the study was  
9 powered. That's how the study was done, and as an  
10 investigator, I am pretty proud of how that study  
11 was done.

12 DR. KOUCHOUKOS: I think it's also  
13 important to recognize that the stroke rate in  
14 80-year-old patients is not 2 percent or 3 percent.  
15 It is more like 8 or 10 percent. And the  
16 prevalence of significant renal dysfunction is also  
17 higher. We did not have a large percentage of our  
18 patients in this study, for obvious reasons, who  
19 were in that category, but one can assume that it  
20 might be possible to demonstrate efficacy in this  
21 high-risk group because of the higher prevalence of  
22 both of these major complications.

23 DR. TRACY: Dr. Aziz?

24 DR. AZIZ: The size of the filter in the  
25 heart-lung machine is usually about 20 microns.

1 The size of the pores here is about 120 microns, if  
2 I am right. In your testing and design before you  
3 came up with this, did you try filters with smaller  
4 pores?

5 MS. CHANG: Yes, we did. We looked at  
6 85-micron pores, and the issue is the  
7 back-pressure, which then starts to create arterial  
8 resistance. So the 120 allows for basically almost  
9 virtually no pressure drop between the filter, and  
10 also to catch particulates of a size that-- again,  
11 Dr. Barbut and Yao [phonetic] did a study where  
12 they looked at embolic size with regard to  
13 neurologic outcomes, and there seemed to be a  
14 collection of larger sizes at about the 120 mark.

15 DR. AZIZ: This is has obviously been  
16 available in Europe since 1998 or so. Outside the  
17 group of patients that we have discussed today, is  
18 there a general feeling that in the older patients  
19 that are being done there--over 80--that there has  
20 been a clinical benefit?

21 DR. ALLEN: The short answer to that is  
22 yes, the European data would suggest that there is  
23 a risk reduction particularly in high-risk  
24 patients, but you'll note we haven't shown any of  
25 that data, because it is not randomized data. I

1 think it is not appropriate data, and I think you  
2 have to stand on the randomized control data, but  
3 you asked the question, and--

4 DR. AZIZ: And that's published, or was  
5 that just an impression?

6 DR. ALLEN: That's published  
7 data--European Journal of Cardiovascular and  
8 Thoracic Surgery.

9 DR. AZIZ: This is actually a question for  
10 future trials that involve neuro-protective sorts  
11 of mechanisms. Can the S-100 protein be used as a  
12 marker? Some people have done that for brain  
13 injury on cardiopulmonary bypass.

14 DR. ALLEN: I can't--Dr. Kouchoukos and  
15 Dr. Kuntz are both eager to answer that question.

16 DR. KUNTZ: And I'd love to hear Dr.  
17 Marler on that, because we would love to use  
18 something for carotid studies as well. I don't  
19 know if he knows about that.

20 DR. MARLER: No.

21 DR. KOUCHOUKOS: There is data in the  
22 cardiopulmonary literature, and Dr. Edmunds, who is  
23 the current editor of one of our journals, can  
24 probably address it. But it is a very insensitive  
25 marker in patients who are undergoing

1 cardiopulmonary bypass with neurologic injury.

2 That's the short answer.

3 DR. EDMUNDS: The S-100 protein is a  
4 marker of neurologic injury, but it is also a  
5 marker of macrophage [inaudible] activation, and  
6 since these operations are all contaminated with  
7 field suction, reclaiming field blood, the marker  
8 is not a reliable index of neurologic injury.

9 DR. TRACY: Are there any other questions  
10 from the panel for the sponsor?

11 DR. MARLER: Could I ask one more  
12 question?

13 DR. TRACY: Yes.

14 DR. MARLER: I haven't heard much about  
15 the indications and the precautions, and we are  
16 asked about that. Could you walk me through your  
17 thinking on going from the selection criteria in  
18 the trial to what you are recommending as  
19 indications for use?

20 It seems to me that the trial selected  
21 patients at, at least neurologically, a lower risk  
22 of events, and yet, it seems you are actually  
23 intending to use this for a much broader range of  
24 patients. Is that correct, and could you walk me  
25 through at least some of the exclusion criteria to

1 explain why you wouldn't continue to apply them  
2 when it is actually used?

3 DR. ALLEN: I think from a practical  
4 standpoint, the device, except in aortas that  
5 couldn't be cross-clamped or that the surgeon chose  
6 not to cross-clamp, you are right, the device would  
7 probably be more broadly applied.

8 I think you do have an opportunity when  
9 you look at the subset of patients that were  
10 considered higher-risk that you saw some mitigating  
11 effect in those patients with the filter, so those  
12 patients were deriving a benefit.

13 It is an inference, and I don't draw  
14 superiority in those patients--I don't make that  
15 claim at all--but it does allow us to show that  
16 there is a subset of higher-risk patients whom we  
17 certainly didn't harm, and actually, some of the  
18 data suggests that we saw some benefit.

19 DR. MARLER: But I think you made  
20 reference that those patients at higher risk were  
21 excluded for obvious reasons from the trial, and  
22 they are not so obvious to me if you then intend to  
23 use the device in them.

24 DR. ALLEN: The patients--we specifically  
25 wanted to look at--when you design a trial for

1 safety, to demonstrate safety of a device, I don't  
2 choose patients who are going to have tons and tons  
3 of complications. If I want to demonstrate  
4 superiority, I choose a population that is going to  
5 have a lot of complications so I don't have to  
6 enroll as many patients, and I can demonstrate  
7 that.

8 So the trial, for all the reasons that we  
9 have discussed--and it is still a very large trial  
10 just to demonstrate that equivalency--was designed  
11 in that fashion.

12 There were patients, though, who were at  
13 moderate or high risk as measured by the Cleveland  
14 Clinic Score, 20 percent of our population, and in  
15 those patients, certainly the device was safe, and  
16 in that subset analysis, there may have been some  
17 benefit shown.

18 So I think surgeons are going to have to  
19 use their judgment as to whom they are going to use  
20 this in. The only patients from an aortic  
21 standpoint who were excluded were if you couldn't  
22 cross-clamp the ascending aorta. If they met  
23 inclusion and exclusion criteria, and you got to  
24 the operating room, and the patient unfortunately  
25 had a porcelain ascending aorta, those patients

1 weren't included.

2       So I think some of the precautions you are  
3 asking me to come to, clearly, if you can't put a  
4 clamp on the aorta, you are probably going to have  
5 to try to figure out some other method to  
6 revascularize these patients, and the filter isn't  
7 going to be appropriate in that population, because  
8 it is obviously attached to a cardiopulmonary  
9 bypass cannula.

10       DR. EDMUNDS: If I might just comment, I  
11 think the Higgins Score is almost irrelevant to  
12 this problem, because there are lots of ways to  
13 die, and the Higgins Score will be influenced by  
14 whether or not somebody has emphysematous  
15 [phonetic] lungs, and I can't see how this device  
16 is going to affect that--and a lot of other  
17 factors. I think that local factors are the ones  
18 that are relevant here--things that directly  
19 produce particulate emboli.

20       DR. ALLEN: I don't disagree with you, Dr.  
21 Edmunds.

22       DR. EDMUNDS: I hope not.

23       [Laughter.]

24       DR. TRACY: Dr. Krucoff?

25       DR. KRUCOFF: I just have one question,



1 really, to Dr. Allen and Dr. Kouchoukos, and  
2 accepting that this is not going to be a data-based  
3 answer. But if this device were approved and came  
4 on the market commercially, as two of the  
5 individuals who have obviously had their hands on  
6 it in human application more than anybody else, can  
7 you tell me just what ball park of your total  
8 clinical practice of open heart surgery you think  
9 you would use this thing in? In what percentage of  
10 patient would you actually pull this off the shelf  
11 and use it?

12 DR. KOUCHOUKOS: Based on what we know  
13 from epi-aortic scanning, as I showed you, we know  
14 that the prevalence of atherosclerosis begins to go  
15 up at age 60 or 65. The cut-off in this study was  
16 60. I would consider using it in every patient  
17 over the age of 60.

18 In this study, we retrieved emboli in 96  
19 percent of the patients. There is every  
20 expectation that as you apply this to older and  
21 older patients, we would retrieve more debris. So  
22 I would use it in any adult patient over the age of  
23 60 undergoing a cardiac procedure.

24 DR. KRUCOFF: And again, Dr. Kouchoukos,  
25 in a very broad sense, is that 30 percent, 50

1 percent, 70 percent of your practice?

2 DR. KOUCHOUKOS: Well, it depends on an  
3 individual's practice, but in an average adult  
4 practice, that would probably encompass probably 85  
5 to 90 percent of patients who have cardiac surgical  
6 procedures.

7 DR. KRUCOFF: Dr. Allen?

8 DR. ALLEN: I'm sure it's the same in your  
9 institutions. I don't have cardiologists referring  
10 me too many young, healthy patients anymore. Most  
11 of my patients are over the age of 64. The median  
12 age in my practice is 72.

13 I would agree with Dr. Kouchoukos. About  
14 20 percent of my patient are done off-pump. I am  
15 not a huge advocate of off-pump, but I use it  
16 selectively in appropriate patients, so you can  
17 already take my number down to about 80 percent. I  
18 am going to use it in a lot of patients.

19 DR. KRUCOFF: So most of the patients in  
20 whom you would cannulate the aorta, you would use  
21 this device.

22 DR. ALLEN: Patients that I would put on  
23 cardiopulmonary bypass and cannulate the ascending  
24 aorta, I think the device is very safe, and the  
25 stuff you see on the filter, it's hard to say that

1 it's not a good thing to take these things out; so  
2 I would be using it pretty frequently.  
3 DR. TRACY: It's hard to say it's not a  
4 good thing to take it out, but it's not easy to say  
5 that it is a good thing to take it out.  
6 DR. ALLEN: Yes.  
7 DR. TRACY: I mean, we have very little  
8 that says it is a good thing to take it out. I'd  
9 just like to make that point.  
10 DR. ALLEN: Dr. Tracy, you are absolutely  
11 right, and I think that gets back to the whole  
12 issue that it is a safety study and not an efficacy  
13 study.  
14 DR. TRACY: Right.  
15 Are there any other questions from the  
16 panel at this point?  
17 [No response.]  
18 DR. TRACY: If not, we'll take a 15-minute  
19 break and then reconvene.  
20 [Break.]  
21 DR. TRACY: If everybody would take their  
22 seats, we can reconvene.  
23 There have been a number of questions  
24 regarding the efficiency of this device at  
25 collecting whatever "goobers" it is collecting, and

1 I believe the sponsor may have some additional  
2 information that might help us understand the  
3 efficiency of this device.

4 MS. CHANG: Yes. We just got our bench  
5 test results from California, and the average in  
6 the bench test is 80 to 90 percent.

7 DR. WHITE: Is that the experiment with  
8 the 120-micron beads?

9 MS. CHANG: Yes--the pinballs flying all  
10 over.

11 DR. TRACY: And there were no studies done  
12 with something that was more similar to  
13 atheromatous or to blood clots; is that correct?

14 MS. CHANG: Yes.

15 DR. TRACY: All right. At this point,  
16 we'll start going through the questions that were  
17 posed to us by the FDA.

18 Questions for the Panel

19 DR. TRACY: The first question: The  
20 primary safety endpoint for this study was a  
21 composite of seven adverse clinical events detailed  
22 on this slide. The median followup was 7 days.  
23 Some facts from the study are: The observed  
24 overall composite event rates were 17.1 percent in  
25 the EMBOL-X arm and 18.9 percent in the control;

1 the composite event rate for the EMBOL-X arm was  
2 shown to be equivalent or not different from that  
3 in the control; also as specified in the protocol,  
4 a separate test for a lower event rate in the  
5 EMBOL-X arm was not statistically significant; the  
6 EMBOL-X arm demonstrated a significantly higher  
7 incidence of aortic endothelial injury--9.2 percent  
8 versus 2.0 percent. Although these patients did  
9 not appear to have any short-term clinical sequelae  
10 resulting from these injuries, long-term effects  
11 are unknown.

12 So the first question posed to us is: "Do  
13 these data support the safety of the EMBOL-X  
14 intra-aortic filter?"

15 I am supposed to summarize the discussion,  
16 and I really have no idea. I think in terms of  
17 being equivalent, if that is equivalent to doing  
18 nothing, then, I suppose it is equivalent to doing  
19 nothing. I really don't know how to answer that  
20 question. I'll have to ask the other panel members  
21 if they can be more articulate than I on this  
22 question.

23 DR. KRUCOFF: Isn't safety and issue of  
24 doing harm? I think from what I heard discussed,  
25 other than the scanning finding of the "ding" or

1    whatever we are calling the little flap of tissue,  
2    which it seemed pretty unclear is related to any  
3    kind of clinical sequelae, it sounded like there  
4    was reasonable information to support that no harm  
5    was being done. The non-inferiority statistic  
6    actually also, to my understanding--and David,  
7    maybe you can chime in--implies that at least no  
8    harm is being done.

9           DR. DeMETS: Yes. My assessment of that  
10   would be that it certainly met the criteria that  
11   were established for the clinical delta that you  
12   were after, and moreover, the rates are actually  
13   lower in the treatment arm; overall, the composite  
14   is lower; and when you examine the individual  
15   components, most of them are at least in the  
16   direction--a few are in the wrong direction, but  
17   just by a little bit. So we would like to have  
18   more data on that, of course, but this is what you  
19   have.

20           So for the issue of safety within the  
21   criteria that were set up, I think they have met  
22   those goals.

23           DR. TRACY: Yes. I guess the thing that  
24   I'm struggling with is that I question whether this  
25   was an appropriate safety endpoint, but it was what

1 was predetermined to be the endpoint of the study.  
2 Bram?  
3 DR. ZUCKERMAN: Right. Dr. Tracy, I think  
4 multiple people have commented that in the year  
5 2002, they might design a different-type trial, but  
6 we have to appreciate how FDA and the sponsor  
7 designed the safety primary endpoint when the trial  
8 was first designed. And I believe safety was  
9 really designed the way Dr. DeMets just summarized.  
10 They met the delta. The trends were in the right  
11 direction. And our third concern, which the panel  
12 has commented on, was that the aortic disruptions  
13 did not have significant clinical sequelae, and if  
14 that is the agreement of this panel, then, for our  
15 purposes, it has met a safety definition.  
16 DR. TRACY: I think it has met it in the  
17 patients that it was tested in. I don't think you  
18 can extrapolate beyond the patients who were  
19 tested. I don't think there are data that would  
20 support expanding into a different group of people,  
21 for example, people with greater degrees of aortic  
22 disease. There are no data that support that.  
23 DR. ZUCKERMAN: And I am glad that you  
24 mention that, because multiple panelists asked  
25 about that this morning, and that question of

1 labeling is a critical one that we will get to in  
2 Question 4.

3 DR. EDMUNDS: Yes, but this was a random  
4 sample of a population of patients that has been  
5 well-characterized and defined. The mean age was  
6 71. So I do think that if statistics apply to  
7 anything, they apply to a set of patients that fit  
8 these descriptors.

9 DR. TRACY: It applies to the patients  
10 that were included in this study. It doesn't apply  
11 to another type of patient. So I would agree with  
12 you. But that has--

13 DR. EDMUNDS: Other type of patients.

14 DR. TRACY: Right--that were not studied.

15 Are there any other comments regarding  
16 this first question?

17 [No response.]

18 DR. TRACY: If not, we'll move to Question  
19 2.

20 "The primary effectiveness endpoint in  
21 this trial was to demonstrate that 75 percent of  
22 the devices would capture at least one particle  
23 during elective CABG or single-valve procedures.  
24 This was demonstrated in the study. There was no  
25 demonstrated reduction in any category of clinical



1 adverse event in this well-controlled patient  
2 trial."

3 "Please address the following concerns:

4 1) Can this method of embolic entrapment, from  
5 this study or elsewhere, be extrapolated to  
6 clinical efficacy?"

7 I think that the answer to that is that we  
8 probably cannot extrapolate beyond what data we  
9 have. We do see that it is efficiency in  
10 retrieving material from patients undergoing  
11 surgery. We don't have data that would support  
12 clinical advantage to that. We don't have data  
13 that it is harmful to retrieve this material. But  
14 I would be cautious about extrapolating on the  
15 basis of intuition.

16 DR. LASKEY: I'm not sure that we know it  
17 is efficient. We know that it does. We don't know  
18 the efficiency. So perhaps it's best just not to  
19 use that term.

20 DR. WHITE: The confusion comes from the  
21 fact that 90 percent of the patients had some; so I  
22 think that's where you say that a high rate of  
23 recovery. But the number of total emboli that was  
24 recovered is unknown, and I tend to think of that  
25 as the efficiency of the filter, so it is an

1 unknown efficiency rate.  
2 DR. TRACY: Okay. Good point.  
3 Are there any other comments on that  
4 particular bullet?  
5 DR. KRUCOFF: Well, there are sort of two  
6 parts--clinical efficacy and efficiency. On  
7 clinical efficacy, it would seem like we were  
8 pretty consistently clear there is no demonstration  
9 of an obvious relationship other than thinking that  
10 emboli are bad. There is no relationship in the  
11 data to clinical efficacy.  
12 And the second part is really about  
13 effectiveness or efficiency of thrombi, and there,  
14 I think we have, as Chris said, the denominator  
15 issues.  
16 DR. TRACY: And the second part: "Do  
17 these data support the effectiveness of the EMBOL-X  
18 intra-aortic filter?"  
19 I think you are hearing comments on that  
20 that we don't know how much was missed that was not  
21 captured. There is no way to know what was not  
22 captured by the device.  
23 DR. FERGUSON: Does the fact that it  
24 captures a known quantity that we know about give  
25 them a plus? I think it does.

1 DR. TRACY: I think the problem again is  
2 that we don't know what it doesn't capture. There  
3 is no bench data that tells us exactly what it  
4 doesn't capture with something that would be  
5 bioequivalent to human atheromatous material or  
6 clot.

7 We know that it captures something, but we  
8 don't know if there is 1,000 times or one time or  
9 10 times as much that is getting past the filter.

10 DR. FERGUSON: But you are posing that as  
11 a negative, and what I am looking at is the  
12 positive, which is that it captures a known  
13 quantity that we know it captures over the non-use  
14 of the filter.

15 DR. EDMUNDS: Yes--you have distorted the  
16 question.

17 DR. FERGUSON: No.

18 DR. MARLER: I don't think he has. My  
19 concern is that my common sense and intuition tell  
20 me that using the word "effectiveness" and not  
21 meaning "clinical benefit" is distorting the  
22 question. Whether the filter pulls back objects or  
23 not, I don't think there is any question.

24 DR. TRACY: Maybe the FDA can clarify what  
25 they are actually asking us here.

1 DR. ZUCKERMAN: In an ideal world,  
2 certainly, we would have liked to have been able to  
3 say yes to clinical efficacy; it makes it a much  
4 easier decision for everyone concerned. But I  
5 think everyone agrees that those aren't the data  
6 that we have in front of us, so what we are looking  
7 for is expert clinical opinion. For example, the  
8 fact that we have this device, this tool, that can  
9 be used in cardiac surgery to take out a certain  
10 number of particles, for the cardiac surgeons on  
11 the panel, is that an effective device? We won't  
12 be able to give you the denominator, but we would  
13 like your clinical impression.

14 Certainly for those who don't do this  
15 procedure every day, we may have a different  
16 impression, but we are interested especially in the  
17 cardiac surgical perspective.

18 DR. EDMUNDS: Well, the term "clinical  
19 efficacy" is a very ambiguous term. It's a bad  
20 term. Does this method of embolic entrapment  
21 remove particulate emboli from the circulation?  
22 That is unambiguous. That's what the question  
23 should be.

24 The second question is does it have  
25 clinical effectiveness, if that's what you want to

1 ask. But you have got to have an unambiguous  
2 question.  
3 DR. ZUCKERMAN: Okay, then, how would you  
4 answer those two questions?  
5 DR. EDMUNDS: You would like me to answer?  
6 DR. ZUCKERMAN: Yes.  
7 DR. EDMUNDS: I'd be delighted to answer.  
8 Do you want it in half an hour or less?  
9 DR. ZUCKERMAN: Yes.  
10 DR. EDMUNDS: It does remove emboli, and  
11 they have not shown any clinical effectiveness in  
12 this study.  
13 DR. FERGUSON: I would agree with that. I  
14 don't think we have data to support the clinical  
15 effectiveness. They have said that, and we have  
16 said it here, so that's easy to answer.  
17 But I do think that we still have to say  
18 that in the sum total of taking care of patients in  
19 the operating room, anything that you can do where  
20 you can prove that you are taking out this  
21 material, but anything you can do like putting a  
22 filter in the arterial pressure line of the  
23 heart-lung machine, and we see the material that is  
24 trapped in that filter, we know that that's a good  
25 thing for the patient. And I view this as the same

1 sort of thing.  
2 DR. ZUCKERMAN: And, Dr. Aziz, do we have  
3 consensus there?  
4 DR. AZIZ: Yes, I think I would agree with  
5 the other two surgical members of the team.  
6 DR. TRACY: Okay. So I think you have  
7 heard the spectrum of answers on that. We don't  
8 know how much is left behind, but it's a good thing  
9 to take away something.  
10 DR. WHITE: Do you want to hear the other  
11 side of that argument?  
12 DR. TRACY: Absolutely. Dr. White?  
13 DR. WHITE: I think that that is not  
14 clear. I could put anything into this patient's  
15 body and have a few bits stick to it and say,  
16 "Look, I got one." Is it a tool that takes out  
17 emboli? Yes--I got one out of 150,000 emboli.  
18 My problem is--I won't argue that 97  
19 percent of these patients had emboli removed--what  
20 I am concerned about is that I have no idea whether  
21 that is doing the patient any good. I think that I  
22 would feel much more happy with partial capture and  
23 a clinical benefit than to support the efficacy of  
24 a tool that may not have any benefit.  
25 So the question then is are we denying the

1 surgeons a tool that they might want to use. I  
2 think that's a much more difficult question. But I  
3 don't want anybody to get confused that I have some  
4 understanding that the tool that removes some  
5 benefits is something that I could say is a good  
6 thing to use in the next patient that has bypass.  
7 I think that's the difference between that and the  
8 bypass filter. The bypass filter is very small; it  
9 captures lots of bits. The question of efficacy  
10 there is much more easily satisfied, because you  
11 can look at the other side of the filter and see  
12 what it misses. What these guys can't tell us is  
13 what is on the other side of the filter, what is  
14 being missed--plus the filter is partially  
15 filterable, and there is disease distal to the  
16 filter that causes these events, both renal and  
17 neurologic, multifactorial disease. So I think we  
18 don't have a very good handle on this.

19 DR. TRACY: I think these are very  
20 difficult questions to answer, and the problems are  
21 the difference between clinical efficacy, and I  
22 think the answer is pretty clear on that part that  
23 we don't have demonstrate of clinical efficacy, and  
24 the problem with effectiveness of the device--it  
25 removes something, and the surgical feeling is that

1 it is a good thing to remove something, and the  
2 other side of the coin is that we don't know if it  
3 is removing one one-thousandth or 90 percent of  
4 what is there available to be retrieved. That is  
5 the ambiguity that I think--  
6 DR. LASKEY: Can't we do better than that?  
7 It is so disingenuous to let us go on the record as  
8 saying that taking something out is a good thing.  
9 It's just so disingenuous, it makes me very  
10 uncomfortable.  
11 Suppressing PVCs is a good thing. It also  
12 kills people.  
13 I just can't accept that. Can we change  
14 the language? Taking something out is a good  
15 thing--no.  
16 DR. EDMUNDS: Warren, we have got to put  
17 this problem in the context that it really is.  
18 Embolization to the brain has been known since  
19 Lee's paper in 1960. It has been a huge problem in  
20 cardiothoracic surgery with bypass since that time.  
21 The improvements that we made that allow  
22 us to have the cognitive deficits that we have  
23 today are small, incremental improvements--heparin  
24 dosage, antifibrinolytics, protenine--all those  
25 sorts of things. It is going to be incremental.



1 This is just a little baby step, perhaps, but it is  
2 a step in the right direction because you are  
3 dragging out some garbage.

4 I think you have to look at it in that  
5 context. That means that there is a lot of garbage  
6 still in there. We know that. But we are going to  
7 be taking it out spoonful-by-spoonful. That's what  
8 it has been for the last 42 years.

9 DR. TRACY: I have the feeling there is  
10 not going to be consensus on this.

11 DR. MARLER: No. I think that the  
12 question about whether it is good or not to remove  
13 the emboli in a way is kind of independent of this  
14 discussion. I mean, that could be answered by  
15 different studies or looking at the literature and  
16 forming an opinion that way, which we really have  
17 not done.

18 I thought we were looking at the results  
19 of a particular trial.

20 DR. TRACY: If we remove the question of  
21 whether it is good or not, then I think the answer  
22 becomes even more difficult, because if we don't  
23 assume that there may be value to removing it, we  
24 have an unknown percentage of something that is  
25 being removed, and it becomes even more difficult

1 to answer the second bullet.

2 So I think I respect the surgical opinion  
3 that removal of material that would otherwise have  
4 gone somewhere probably is a good thing, although  
5 we do not have that clinical efficacy answer.

6 It is unfortunate that there is not a  
7 better endpoint in this study to look at. I'm not  
8 sure we're going to get much farther than where we  
9 are with this question.

10 DR. KRUCOFF: I'll just make one comment  
11 and probably end up having to change specialties,  
12 because I actually lean toward the surgical group  
13 on this one.

14 I think it would concern me if we felt  
15 that these little intimal "dings" or that some  
16 other significant safety issue or that some  
17 technical element that really made you redo or do  
18 differently the basic procedure of cannulating the  
19 aorta were a part of this device. Then I would  
20 feel very conservative about all the issues that  
21 have been so much discussed, about whether pulling  
22 grunge out of the [inaudible] meant anything.

23 I guess, based on our consensus on the  
24 first point, if there is really no significant  
25 safety issue--and my understanding of the

1 discussion this morning, including all of us  
2 surgeons on both sides of the table, is that you  
3 put this thing in, you cannulate the aorta the same  
4 way you would for another purpose, and you  
5 basically slide the filter in with a pretty  
6 straightforward--essentially, with no real  
7 technical change from how you would do a routine  
8 open heart procedure--then, if debris is bad, and  
9 you are pulling it out, is that the first  
10 spoonful--it sounds to me like the surgeons from  
11 both sides of the table feel like it probably is.  
12 And it is very clear that we aren't going to see  
13 any answer to that question in the dataset. But if  
14 there is not a significant safety issue, then, to  
15 me, a lot of the judgment about whether or not to  
16 use this probably ought to come not from the panel  
17 but from a community of surgeons who do this in  
18 live patients, and the only way they can do that is  
19 if it is available.

20 DR. TRACY: Other comments on this  
21 question?

22 [No response.]

23 DR. TRACY: If not, we'll move to Question  
24 3.

25 "Do the study data support an appropriate

1 risk/benefit profile?"

2 I think you have--I'm not sure what the  
3 question is, really.

4 DR. ZUCKERMAN: Okay. Dr. Krucoff just  
5 gave an answer to the question. He is looking at  
6 the data in this study as well as general  
7 experience, literature, et cetera, to try to come  
8 up with a risk/benefit profile which he judges to  
9 be positive.

10 We are asking for other comments. At the  
11 end of the day, one needs to cut to the chase. Is  
12 there enough data within this study and external  
13 data to agree with Dr. Krucoff's comments?

14 DR. TRACY: I think you're going to hear  
15 the same kind of debate back and forth among the  
16 panel members. I think within the confines of the  
17 type of patient who was involved in this particular  
18 study, Dr. Krucoff's answer probably is the correct  
19 answer. But I would caution again not to  
20 extrapolate too widely, and I am having the sense  
21 that there would be wide extrapolation if the  
22 device were clinically available.

23 DR. MARLER: I guess I can say that it's  
24 hard to compute a risk-benefit when you have no  
25 indication of any benefit. I am extremely cautious

1 about, quote, "approving," whatever that means,  
2 something for which there is no evidence of any  
3 benefit. And there is some risk--we have discussed  
4 it--and there were a lot of questions about whom it  
5 is going to be applied in.

6 So I would say that you can't say that the  
7 study data support a risk-benefit profile because  
8 they don't support any benefit.

9 DR. EDMUNDS: I think I can disagree with  
10 that. It does not support any clinically  
11 demonstrable benefit except that it removes a  
12 filter full of garbage, and that's a benefit.

13 DR. TRACY: Accepting that on,  
14 basically--at this point, we have been asked to  
15 accept that on intuition, and I have a hard time  
16 accepting that on the basis of intuition, because  
17 we still struggle with not knowing what got by. So  
18 it is very difficult, I agree.

19 DR. WHITE: The problem is it only takes  
20 one piece to cause a stroke. You can get an  
21 endpoint with one piece. That is why it is hard to  
22 know that if you got 60 percent of them, there is a  
23 clinical benefit associated with that, when it only  
24 takes one to cause an endpoint.

25 DR. EDMUNDS: I am just really frustrated.

1 For 42 years, we have been doing open heart  
2 surgery, and we have known that we have been  
3 circulating emboli. We're going to do it tomorrow,  
4 with or without the filter, circulating emboli. So  
5 there is no argument that there are no more emboli  
6 going to circulate if you use this device. There  
7 are going to be lots, but there are going to be  
8 less. The amount less is what is on the filter.  
9 That's the benefit. It is not a clinically  
10 demonstrable benefit.

11 They didn't measure clinical benefit in  
12 this study. They didn't say they measured it, they  
13 didn't intend to measure it. We can't hold them to  
14 that standard.

15 DR. TRACY: I think it's a little more  
16 than that. I think it's the whole issue of what is  
17 not known about what gets by. That's a very, I  
18 think, unknowable thing. There is something  
19 appealing about removing debris, but we don't know  
20 how much is being removed.

21 Dr. Laskey?

22 DR. LASKEY: I think what Hank says is  
23 absolutely on the money. I would agree with it 100  
24 percent. It is just unfortunate that most people  
25 when they hear "risk-benefit" think about clinical

1 benefit. They think about the risk of harm, and  
2 they think about the risk of clinical benefit.

3 So perhaps again, the language could be  
4 softened here somewhat. But "benefit" here needs  
5 to be strictly qualified that this is not an  
6 artifactual benefit, but it is a benefit in terms  
7 of the study, which is strictly defined as catching  
8 stuff on the filter. But it is not a clinical  
9 benefit. When people see "risk-benefit," that's  
10 what they think of.

11 DR. WHITE: Do you need to consider the  
12 potential risks because you didn't measure a  
13 clinical risk in these patients; are there other  
14 theoretical or potential risks? We had a lot of  
15 discussion about air. We had discussion about  
16 disruptions, use in people with more sick aortas.

17 If we are going to hypothesize about a  
18 possible benefit, should we hypothesize about how  
19 potentially dangerous this could be if the device  
20 were misused?

21 DR. LASKEY: Thank you, Chris.

22 One point I wanted to mention is that we  
23 don't know that these disruptions are the triponine  
24 [phonetic] of cardiac surgery. We didn't know  
25 about triponine until we started looking at these

1 very, very subtle markers of injury. I'm not sure  
2 we know--none of us wants to believe that they are  
3 bad, but they may be, and it may be a very subtle  
4 marker of injury that we just don't have a handle  
5 on, just like triponine was in the early days.

6 So I wouldn't dismiss it, and I would  
7 certainly keep it in the mix. I agree with you.

8 DR. TRACY: Dr. Ferguson?

9 DR. FERGUSON: I agree totally with the  
10 concept that this question is bad in the sense that  
11 "risk-benefit" does convey something other than  
12 what we did.

13 I think we have all agreed that the risk  
14 is no greater than the control in this, and the  
15 benefit is that it takes out some clot, which we  
16 know is bad. And that's a better way to put it, I  
17 think, than the way it is written.

18 DR. TRACY: Are there any other comments  
19 on this question?

20 [No response.]

21 DR. TRACY: Does that satisfy your  
22 question, Dr. Zuckerman?

23 DR. ZUCKERMAN: Well, I think it is  
24 important for the record to know if the other two  
25 cardiac surgeons agree with Dr. Ferguson's last



1 statement, that if we more precisely define what we  
2 are trying to get at here, if they would agree.

3 DR. EDMUNDS: Well, we are splitting  
4 hairs, and I have on my loop, so that first of all,  
5 it is not clots that tare being taken out. It is  
6 atherosclerotic debris, principally. And I don't  
7 think that we have demonstrated clinical benefit,  
8 but we have demonstrated a benefit, and we have  
9 discussed that.

10 There is a risk. The risk is exceedingly  
11 low. It is the EDS risk. And while it isn't zero,  
12 it is the number next to zero.

13 DR. AZIZ: I think--how can I put  
14 it--clearly, I think that those two things are both  
15 in a sense true, and I hope they are connected in  
16 the sense that we do believe that removing clot  
17 will give us reduced brain injury, and I agree with  
18 what Chris is saying that we don't know how much is  
19 getting through on the other hand, and you really  
20 don't need to have a lot going through.

21 But to sort of crystallize it, I think the  
22 wording of this statement is a little fuzzy. I  
23 think that what one should say is that clearly,  
24 this device has some risks, but at the same time,  
25 it reduces some other risks--namely, the risks

1 associated with the lot.

2 So I think I would agree with what Dr.  
3 Ferguson is saying.

4 DR. TRACY: Okay. We'll move on, then, to  
5 the fourth question.

6 We are being asked in this question to  
7 review the labeling. "The labeling must indicate  
8 which patients are appropriate for treatment,  
9 identify potential adverse events with the use of  
10 the device, and explain how the product should be  
11 used to maximize benefits and minimize adverse  
12 effects. Please address the following questions  
13 regarding product labeling: Do the Indications for  
14 Use adequately define the patient population  
15 studied? For example, should the patient  
16 population receiving this device be limited to the  
17 same patient population utilized in the study. For  
18 example, non-emergent; patients over the age of 60;  
19 first-time isolated valve or CABG patients."

20 We'll take that first piece first. I  
21 think that there may be some sense that it would be  
22 applicable in other patient populations, but we  
23 have pretty scanty information as it stands, and I  
24 would be very cautious about expanding beyond the  
25 population that was studied in this protocol.

1 DR. AZIZ: Well, I don't know quite how to  
2 take that, because we know that certain patients  
3 are at higher risk of getting this sort of problem,  
4 particularly the guys who are over 85, and although  
5 that didn't form a large percentage of the cohort,  
6 I think it comes back to what happened in a sense  
7 yesterday. Although you may say that, I think that  
8 in clinical practice, if I wanted to use this  
9 device, I would probably want to use it in a  
10 patient group that I know is at increased risk of  
11 having a neurological event.

12 So no matter what we say, I think that in  
13 clinical reality, you probably would target the  
14 higher-risk patient anyway.

15 DR. TRACY: Yes. My only point is not to  
16 extend it to the porcelain aortas or beyond the  
17 scope of this particular study. There did seem to  
18 be a group in whom there might be greater benefit,  
19 and that was the higher-risk patient population  
20 within this study. I think we have no data beyond  
21 this patient population, which did include some  
22 higher-risk but not the extraordinarily high-risk  
23 patients.

24 I think you have what you have.

25 DR. EDMUNDS: Dr. Tracy, I would like to

1 take two words or two phrases out of that. I would  
2 like to take "non-emergent" and "first-time" out of  
3 that statement. Otherwise I can live with it. But  
4 there is no sense in handicapping the surgeons when  
5 they are doing more difficult cases who are  
6 otherwise over 60 and at risk of this problem.

7 DR. TRACY: Do either of the other  
8 surgeons have any comments?

9 DR. FERGUSON: We're talking about  
10 labeling here, and the question to me would be--and  
11 you have a good point, I think, Cynthia--or what I  
12 am wrestling with is should the labeling say that  
13 this device was tested under these conditions, and  
14 put those conditions in, which would work for the  
15 FDA and work for our consciences and so forth.  
16 Now, the way the device is going to be used is of  
17 concern to us, but it is of no concern, because it  
18 is going to be used--Dr. Kouchoukos already said  
19 the ones he is going to define and use--

20 DR. ZUCKERMAN: Maybe we can have a  
21 time-out here and talk briefly about what we are  
22 getting at in this question.

23 Certainly the agency doesn't regulate the  
24 practice of medicine, and if the device is  
25 approved, there will be surgeons who will use it as

1 they want to. But what we are talking about today  
2 is truthful and accurate labeling in an indications  
3 statement, and then, part and parcel, we  
4 traditionally describe the clinical trial that was  
5 performed in the Clinical Trials section.

6 I guess the main question that FDA has is  
7 when we look at the indications and intended use,  
8 it says that "The EMBOL-X aortic filter is  
9 indicated for use with the EMBOL-X aortic cannula in  
10 cardiac surgery procedures to contain and remove  
11 particulate emboli."

12 Based on the data that we discussed this  
13 morning, should there be additional qualifiers that  
14 better describe who was actually studied?

15 DR. TRACY: You can refer to that in the  
16 "Proposed Labeling" section; page 2 of 10 at the  
17 top has the proposed indications for intended use.  
18 And then, on page 4 of 10, it begins the  
19 description of the patient population.

20 DR. FERGUSON: What was the first  
21 reference?

22 DR. TRACY: It is in the section titled,  
23 "Proposed Labeling," page 2 of 10, down at the  
24 bottom. And it is at the very top, Number 2, and  
25 it is exactly what Dr. Zuckerman said.

1           "The EMBOL-X aortic filter is indicated  
2 for use with the EMBOL-X aortic cannula in cardiac  
3 surgery procedures to contain or remove particulate  
4 emboli."

5           And the question is do you then add the  
6 phraseology "in patients over the age of 60 who are  
7 undergoing first-time surgery for isolated valve or  
8 CABG"--I think that's the question, or is it  
9 adequate on page 4 to state the description of the  
10 patient population that was studied.

11          My instinct would be that in general, we  
12 state the indication, unless there is a particular  
13 reason to put a qualification on it, and usually  
14 put those caveats in in the description of the  
15 patient population. I'm not sure that it is  
16 critical to put that up front in the labeling, but  
17 I think it has to be somehow there that that is the  
18 patient population that was studied.

19          DR. MARLER: So we are talking about  
20 indications now?

21          DR. TRACY: Right.

22          DR. MARLER: I guess it depends on whether  
23 you think your recommendations need to be driven by  
24 data, or not. I would personally be concerned,  
25 because my knowledge is that if you exclude

1 patients who have had prior stroke in the study,  
2 you are very unlikely to get any estimate of the  
3 risk of stroke from a procedure, because that's one  
4 of the highest-risk groups of patients.

5 And we heard from the sponsor that the  
6 reason this population was chosen--at least, my  
7 interpretation of the response--was because they  
8 wanted to select a group of patients in whom it was  
9 safe, or most likely to be safe.

10 So I think to take it beyond that without  
11 any data to support it would be against my  
12 understanding of what we are doing here today,  
13 which should be driven by data from the study that  
14 we are presented.

15 DR. TRACY: So are you supporting adding  
16 the phrases in "Indications and Intended Use,"  
17 "non-emergent; patients over age 60; first-time  
18 isolated valve or CABG patients"? Would you  
19 propose putting that in the indications statement?

20 DR. MARLER: Well, because on some of  
21 those exclusion criteria, I don't have the  
22 expertise to interpret them, I would certainly  
23 think that you are changing the game if you include  
24 patients who had prior stroke or carotid stenosis.  
25 I don't know how the filter could possibly relate

1 to that.

2 What I would say is that I think each of  
3 the exclusions need to be discussed carefully in  
4 terms of whether there is any indication that it is  
5 also safe in those patients.

6 DR. TRACY: Maybe a more detailed  
7 description of the exclusion criteria would be  
8 helpful in there. That doesn't seem to be  
9 particularly well-detailed in the proposed  
10 labeling. And I think perhaps in the section  
11 entitled "Indications and Intended Use," if there  
12 were some reference specifically to "Please see  
13 below for specific patient population inclusion and  
14 exclusion criteria" and a statement that the device  
15 was tested in these patients only, would be  
16 appropriate.

17 DR. FERGUSON: Yes. As data?

18 DR. MARLER: Right--whatever. I just  
19 think that if you make the decision that intuition  
20 is going to drive this whole process, there is  
21 almost no purpose to even do the trial beyond the  
22 first few number of patients where you show it  
23 catches some emboli.

24 So I think that the basis of approval for  
25 anything should have to do with data that show that



1 it is safe and effective.

2 DR. TRACY: Okay?

3 DR. ZUCKERMAN: Well, there are still a  
4 lot of question marks for these other potential  
5 patient populations. Sometimes, what we do in our  
6 labeling is after the "Clinical Trials" section, in  
7 what would be before Section 8, we have a  
8 discussion of individual patient considerations  
9 where we could talk about some of these  
10 sub-populations and the lack of known data right  
11 now.

12 Do you think that it would be appropriate  
13 to add that section?

14 DR. TRACY: I think that that would  
15 perhaps be helpful to somehow, if you can, capture  
16 some of the questions about the other patient  
17 populations. I think that would be appropriate,  
18 and reference to that in the original Number 2,  
19 "Indications and Intended Use."

20 DR. EDMUNDS: I think it would be more  
21 concise just to say what the study was done on, and  
22 that is what is up there, and then put that in as  
23 data, and then put a disclaimer that the  
24 manufacturer does not extrapolate these data to  
25 anybody. They probably wouldn't say it just that

1 way.

2 [Laughter.]

3 DR. TRACY: Let's move on to the other  
4 bullets, and maybe it will help us clarify.

5 The second bullet is: " Are there any  
6 other restrictions that should be placed on the  
7 patient population receiving this device?" And I  
8 guess that means in terms of contraindications. I  
9 think rather than contraindications, simply stating  
10 that other patient populations were not studied  
11 would be the more appropriate way of stating it.

12 DR. FERGUSON: Contraindications are in  
13 the next bullet.

14 DR. TRACY: Right.

15 DR. KRUCOFF: Yes, and one that has come  
16 up that I didn't see in their labeling is the  
17 porcelain aorta as at least one morphologic  
18 descriptor that has come up a couple of times today  
19 that at least in my version is not specifically  
20 listed as someone who would probably not want to--I  
21 wonder, is there a broader range, or is it worth  
22 saying the obvious, which is that in patients in  
23 whom you would simply not want to cannulate the  
24 aorta, those are patients who are obviously not  
25 candidates for this?

1 DR. TRACY: The next section is, "Based on  
2 the clinical experience, should there be additional  
3 Contraindications, Warnings, and Precautions for  
4 the use of the EMBOL-X intra-aortic filter?"

5 DR. MARLER: I want to back up a little  
6 bit. On these indications and contraindications,  
7 aren't we more deciding not what an individual  
8 physician can do in his or her practice, but what  
9 the company can advertise as what the FDA has  
10 looked at and approved as an indication,  
11 presumably, following the meaning of FDA approval,  
12 that it is driven by data and good evidence?

13 MR. MORTON: Madam Chair, a couple of  
14 points, not to supersede Dr. Zuckerman.

15 DR. TRACY: Yes.

16 MR. MORTON: There is a difference between  
17 approving and clearing a 510(k), and that's what  
18 the FDA action will be; it will be a clearance. So  
19 the sponsor or the manufacturer will represent the  
20 device as "cleared," not as "approved," and that is  
21 significantly different.

22 Additionally, there is a tremendous  
23 difference in the regulatory burden on indications  
24 and changing indications once those are locked into  
25 a clearance. It requires data and could even

1     require another 510(k).

2             DR. MARLER: But the impact of the  
3     indications is primarily the way the product can be  
4     advertised and sold, not the way it is used by the  
5     physician; is that correct?

6             MR. MORTON: That is a practice of  
7     medicine issue; correct.

8             DR. ZUCKERMAN: Okay, but the point that  
9     we are trying to get at is what is a truthful way  
10    to describe the dataset that has been discussed  
11    here and other external data that is potentially of  
12    importance, Dr. Marler. So that the indications  
13    statement doesn't need to necessarily follow, dot  
14    by dot, the clinical trial guidelines if it is  
15    reasonable to extrapolate farther. That is where  
16    we need your help; it may not be.

17            DR. TRACY: Yes. I think it is very hard  
18    to extrapolate the way this particular study is  
19    constructed. It is hard to extrapolate, to expand  
20    on the indications. I think that the indications  
21    as stated in the--I'm not sure I would want to be  
22    more restrictive than what is stated here, other  
23    than to say that somewhere in the subsequent body,  
24    there has to be a statement of exactly who was  
25    studied and exactly who was not included. I think

1 the exclusion criteria have to be more clearly  
2 stated than they are in the current proposed  
3 labeling.

4 As far as placing additional restrictions,  
5 I don't think we need to place additional  
6 restrictions beyond clearly reiterating what was  
7 studied and what was not studied.

8 Dr. Pina?

9 DR. PINA: Let me go beyond the inclusion  
10 criteria. I would like to see a table with the  
11 baseline values of the patients. They were over  
12 the age of 60, but the age was up there, so I think  
13 that whomever is going to use this needs to see the  
14 mean values of the population that was actually  
15 studied. And it can be very simple--age, gender,  
16 type of surgery, even number of vessels involved.  
17 That information should be available. But that's a  
18 descriptor that I think a surgeon needs to look at  
19 to make a decision about whether they want to use  
20 this or not.

21 It's not just the inclusion criteria; it  
22 is what the data actually are.

23 DR. TRACY: So that is there in paragraph  
24 4 on page 4 of 10 and perhaps would benefit from  
25 expanding that into a more inclusive table.

1 DR. MARLER: And I guess I would want to  
2 add that I thought I pretty clear asked for an  
3 explanation of how to extrapolate, and I did not  
4 hear anything except that emboli are bad.

5 DR. EDMUNDS: I don't think we can  
6 extrapolate, and I don't think we can say the  
7 negatives, either, because the negative list will  
8 be long. You start out with dissecting aneurysms,  
9 marfands [phonetic], airlos-dolos [phonetic],  
10 Siamese twins--you can keep going in an endless  
11 list. I think we have got to just stick to what  
12 this trial was. The use of this product is based  
13 on the demonstration that this filter captures  
14 embolic material when used in these patients, and  
15 that's where I would recommend that you stop,  
16 because that's all the data that we have.

17 DR. LASKEY: And it is exceedingly  
18 difficult to stretch things beyond the equivalence.  
19 If we had evidence of benefit, it could conceivably  
20 stretch this to another sample, but with  
21 equivalence, that's a long way.

22 DR. KRUCOFF: I would even amplify that to  
23 emphasize that since part of what I think makes  
24 this whole consideration reasonable is the absence  
25 of safety concerns, pushing the envelope in

1 directions to new, untested populations, I think  
2 whether you encounter safety issues would be an  
3 important question that would be outside of this  
4 discussion.

5 DR. ZUCKERMAN: Good. So you know what  
6 the present indications statement reads. How do you  
7 make it in a concise fashion more applicable to the  
8 data that have been presented?

9 DR. TRACY: I think you have to simply  
10 state, "See patient selection criteria for patients  
11 involved in this protocol." I think you have to  
12 refer to some other section. Otherwise, you end up  
13 with a 15-paragraph--you have to refer to other  
14 areas in the labeling and then perhaps a statement  
15 that this device simply was not tested in other  
16 patient populations.

17 DR. ZUCKERMAN: Okay. Another option that  
18 we sometimes use--the one that you have suggested  
19 is just to put, "(See Clinical Trials section)" in  
20 parentheses. Another option that we sometimes use  
21 is to indicate some key clinical parameters right  
22 in the indications statement of which there aren't  
23 any right now.

24 Are there any that are real show-stoppers  
25 that should be up front?

1 DR. TRACY: I think they are all relevant,  
2 and I think you just have to refer to who was  
3 involved in this study and specifically who was not  
4 involved in this study.

5 Dr. Ferguson?

6 DR. FERGUSON: Are we through with that?  
7 I want to bring up another point.

8 DR. TRACY: Yes.

9 DR. FERGUSON: I think it belongs here;  
10 maybe not. That is that as I read through  
11 these--and correct me if I am wrong--I see nothing  
12 at all in the deployment and use of the device that  
13 indicates that you should use imaging as a guide.  
14 I know that you don't need it to put the instrument  
15 in, but is there any reason--I am just bringing it  
16 up as a question, because half of the patients had  
17 either TEE or the epiaortic.

18 I would like to have some discussion about  
19 whether that should be somehow included here.

20 MS. WENTZ: That's actually the last part  
21 of the question.

22 DR. FERGUSON: Oh, okay. I jumped the  
23 gun. It's not in there now; right?

24 DR. KRUCOFF: My understanding from the  
25 comments from the investigators was that it just



1 would not be conceivably or technically feasible to  
2 do that on a routine basis.

3 DR. TRACY: I think their point on the TEE  
4 was that you wouldn't pick up these small  
5 disruptions.

6 DR. FERGUSON: TEE isn't available for  
7 everybody.

8 DR. TRACY: It is available, but it may  
9 not pick up these small disruptions that we don't  
10 know, as Dr. Laskey said, what the clinical  
11 relevance of these are. We know that they are  
12 occurring at a higher rate, and short-term, we  
13 don't see any increased adverse outcome related to  
14 it.

15 DR. FERGUSON: Excuse me. That jumps the  
16 point. I'm not talking about that. I was  
17 referring to the fact that it is very useful when  
18 you are getting ready to put an aortic cannula in,  
19 as Dr. Kouchoukos has done in monumental studies,  
20 it is very useful to have some sort of  
21 visualization of the arteriosclerotic aorta before  
22 you put the cannula in. I would like to have some  
23 discussion of that.

24 DR. EDMUNDS: Well, that adds a little bit  
25 of a burden and delays the operation by about 10

1 minutes, because first of all, you have to get the  
2 probe sterilized, and then you have to find an  
3 anesthetist who knows where it is, and it take  
4 about 10 minutes to do the study.

5 DR. FERGUSON: And probably if you do  
6 that, and you put the cannula in a spot that is not  
7 arteriosclerotic, you are probably doing the  
8 patient more good maybe even than the filter.  
9 That's my point.

10 DR. PINA: Dr. Tracy, let me go back to  
11 Dr. Zuckerman's point, because I don't think we  
12 answered that first statement for you.

13 In your statement of indications, I think  
14 you can very easily say, just like it says now,  
15 for cardiac surgical procedures that are  
16 non-emergent and in patients over the age of 60,  
17 for either bypass or valve surgery--you can say all  
18 of those indications in one sentence, and that will  
19 describe generally the population, and then say  
20 "Refer to Table such-and-such."

21 DR. TRACY: I think you are going to have  
22 some argument from your surgical colleagues that  
23 that is appropriate to limit the surgeons' ability  
24 to do this in emergency cases. So that was why I  
25 was holding back from making that particular

1 statement. I think it is important to say that  
2 that is the patient population that was involved,  
3 but I don't know that I would put that pu front.

4 DR. PINA: Yes, but the FDA will not  
5 regulate medical practice. They are going to do  
6 whatever they want anyway. The surgeons are going  
7 to use it any way they want, but I think the  
8 truth-in-labeling has to be the population that was  
9 studied. And the surgeons will make their  
10 decisions clinically, as we have always done with  
11 everything we do.

12 DR. TRACY: I think that is going beyond  
13 the scope of--I think to put that in there in the  
14 indications creates some liability issues that I  
15 would rather not open up to. I think we need to  
16 have it very clear that this is the patient  
17 population that was studied, and these are the  
18 people who were excluded. I'm not sure that it  
19 belongs as an additional sentence in the  
20 indications for usage.

21 DR. EDMUNDS: I agree with you. I think  
22 that's the data-supported course.

23 DR. TRACY: But the reference needs to be  
24 there.

25 Just to try to get through these other

1 issues, "Are there any other restrictions that  
2 should be placed on the patient population  
3 receiving this device?"

4 No, there are no other--I cannot think of  
5 any other restrictions that we have discussed here  
6 that need to be placed. Again, it has got to be  
7 clear who was excluded from the patient population.

8 And finally, on the third bullet of this  
9 question, "Should there be additional  
10 Contraindications, Warnings, and Precautions for  
11 the use of the EMBOL-X intra-aortic filter?"

12 I think, with the exception of the  
13 porcelain aorta, which has come up a couple of  
14 times, the contraindications as stated are fair,  
15 and I think somewhere, the idea of the porcelain  
16 aorta has to come up, either as a warning or as a  
17 precaution.

18 Dr. Krucoff?

19 DR. KRUCOFF: I have a question, and there  
20 may not be a precedent to make this helpful, but I  
21 wonder if it would be worth, given the whole  
22 spirit of this discussion, separating out technical  
23 effectiveness from clinical effectiveness and to  
24 label this as a device that, with reasonable  
25 safety, has been shown to be technically effective

1 at retrieving particulate matter in the setting of  
2 cannulation of the aorta, but in a precaution sort  
3 of environment, or make a statement that while this  
4 has been shown to be technically effective at  
5 removing particular matter, it has not been shown  
6 to be clinically beneficial, and therefore, caution  
7 in the use of this device would be warranted.

8 DR. TRACY: I think it is difficult to go  
9 in that direction since we don't even know that it  
10 has been shown to be effective at removing  
11 particulate material. We don't know what  
12 percentage it is removing. I think that to  
13 introduce that into the labeling would be  
14 confusing, at the very least.

15 DR. KRUCOFF: I was actually just wondering  
16 from a precedent, because I can think back as far  
17 as a day where an indication for reducing ischemia  
18 during angioplasty was sort of a technical  
19 achievement that, actually, devices and biologics  
20 were both ultimately approved for, even though the  
21 full clinical ramifications of reducing ischemia as  
22 a technical feat were never implied or  
23 demonstrated. And I just wondered whether there  
24 was a precedent that we could get clearer language  
25 for docs who are considering using this in their

1 patients to make it very clear why this is on the  
2 market, but also the limitations of what we can  
3 understand about its clinical utility.

4 DR. TRACY: The study results are pretty  
5 clear as they are stated here. I think it doesn't  
6 have to be expanded beyond the study results. I  
7 think it is explained in here what the study  
8 results are, and I think it is fair to leave that  
9 open to interpretation for the operator whether  
10 that study results supports the use for them in  
11 their individual patients.

12 DR. AZIZ: If I could just make one point  
13 about the porcelain aorta, I think that shouldn't  
14 be put as a contraindication. I think I would just  
15 use it as a precaution, because you could probably  
16 see a case where you do a beating heart on-pump  
17 case where you won't cross-clamp the aorta, so you  
18 could still put the mesh inside.

19 So I would say that that should be a  
20 precaution but shouldn't be a contraindication.

21 DR. LASKEY: How about if you're just not  
22 willing to cross-clamp the aorta if it's a  
23 contraindication?

24 DR. FERGUSON: And I'll get back to the  
25 point again about the ECO--nobody wants to talk

1 about it--many times, you can't tell if you've got  
2 a porcelain aorta or not until you actually do the  
3 TEE. You don't want to do this, you don't want to  
4 do that, because that's going to be bad.

5 DR. TRACY: Shall we flip the page?

6 DR. MARLER: So, at least at a minimum, I  
7 hear that as a precaution, we could specifically  
8 notify the clinician that certain patients were  
9 excluded from the study, and there is no evidence  
10 of the safety in patients that--and then list the  
11 clinical exclusions.

12 DR. TRACY: I think it's very clear that  
13 more detailed description of clinical exclusion  
14 needs to be included in the labeling.

15 DR. MARLER: Rather than just referring to  
16 the protocol.

17 DR. TRACY: Right.

18 The next bullet, then, is: "Should the  
19 labeling include specific study information such as  
20 no reduction of clinical events were noted in a  
21 1,289-patient clinical study; and the EMBOL-X  
22 device appears to increase the rate of endothelial  
23 injury?"

24 The study results on page 4 of 10 do  
25 indicate that none of the surgical procedure

1 differences between the randomized group achieved  
2 statistical significance. So that is there, I  
3 believe. And yes, I think it is important that  
4 those pieces of information be conveyed very  
5 clearly and in a fashion that can be readily picked  
6 up in the labeling.

7 And then, to grapple a little more with  
8 "What should the labeling include regarding the use  
9 of ultrasound both before--for assessment of the  
10 aorta--and after--monitoring of injury--the use of  
11 the device?" there are data provided on the use of  
12 TEE or epiaortic imaging, and that is presented on  
13 page 5. Is there something peculiar about this  
14 device that would make it necessary to mandate the  
15 use of TEE?

16 DR. FERGUSON: My suggestion would be  
17 something in the instructions for use of how  
18 helpful TEE can be in assessing the aorta both at  
19 the time the cannula is put in and also assessing  
20 the aorta before and after, but not to make it a  
21 mandatory part of the situation.

22 DR. TRACY: I think that's reasonable. It  
23 is stated here, but perhaps that could be clarified  
24 a little bit as to exactly what was seen with the  
25 two modalities of assessment, and then, certainly



1 in the instructions for use, that should be  
2 emphasized.

3 DR. LASKEY: What is the standard of care?  
4 Are all patients generally getting TEEs now?

5 DR. EDMUNDS: No. Could we go back to a  
6 previous slide? Number 2 of the first paragraph,  
7 you kind of went over, but I take objection to that  
8 statement, because "endothelial injury" is not  
9 defined. If you are going to use the term, I think  
10 you have to define it. We have discussed this, and  
11 we have been unable to demonstrate that this is a  
12 harmful finding, and we have shown, or the study  
13 showed, that 78 percent of OR personnel were unable  
14 to detect it at all.

15 So I think we have got to back off a  
16 little bit about that unless you start to raise a  
17 whole bunch of thorns that really don't need to be  
18 raised.

19 DR. TRACY: The thorns are there, though,  
20 unfortunately, and that was found, and it is  
21 defined on page 5 of 10, the presence of  
22 ecocardiographically-evident endothelial disruption  
23 is noted" and the statement is made "did not put  
24 the patient at a statistically greater risk for  
25 composite endpoint event." I think that's fairly

1     stated. It was there. There was no evidence that  
2     it increased the composite risk.

3             So it has to be there, and perhaps  
4     something that needs to be observed over time,  
5     because we certainly have very limited information  
6     about what the prognostic significance of this is,  
7     and in particular as more diseased aortas are  
8     approached with this device, I wouldn't be  
9     surprised if there were a greater risk of  
10    disruption in those patients. So we need to be  
11    tracking something like this.

12            Are there any other comments on this?

13            [No response.]

14            DR. TRACY: I think we are on to Number 5.

15            "Please provide any other recommendations  
16    or comments regarding the labeling of this device."

17            I think we touched on--Dr. Pina?

18            DR. PINA: I just want to go back and ask a  
19    question of the sponsor. Somewhere along the way,  
20    you stated that the endothelial injury was almost  
21    center-specific and operator-specific, or you saw  
22    it several times in the same operator, so that  
23    perhaps experience may have a lot to do with lack  
24    thereof.

25            Am I correct? Did I hear that right?

1 DR. TRACY: I don't think that was my  
2 impression.

3 DR. ALLEN: Actually, a good recollection.  
4 The two that were repaired were from the same  
5 center very early on in the experience, and then,  
6 after an historical basis for what we were seeing,  
7 and realizing they weren't causing clinical events,  
8 the additional 10 that were actually observed by  
9 surgeons weren't intervened on.

10 DR. PINA: Again, I don't know if it would  
11 pay to say something in there that, first of all,  
12 they are not that common. We don't know how many  
13 times this happens, as Dr. Edmunds said, and we  
14 don't even know about them. So something to take  
15 a little bit away from the fear, even though I know  
16 it is there, but not to cause undo alarm may have  
17 something to do with the experience of the surgeon  
18 or the surgeon's ability to see this, just to kind  
19 of temper a little bit the fear of the endothelial  
20 injury.

21 DR. TRACY: Wasn't it the repairs were  
22 done early on because people didn't understand the  
23 lack of clinical importance. It is not that the  
24 number of disruptions decreased over time.

25 DR. MARLER: Wait a minute, now.

1 Intuition is ruling today; right?  
2 DR. TRACY: I'm sorry?  
3 DR. MARLER: Emboli are bad. Endothelial  
4 injury is likewise bad. I don't see why there is  
5 any discussion of this.  
6 DR. KRUCOFF: I think there is a reason,  
7 because I think you have to be fundamentally  
8 consistent about our assessment of the safety of  
9 the safety of this thing. To me, that's actually  
10 the much more rigorous part of this than the  
11 efficacy issue. And if we have all reached the  
12 consensus that this is safe, part of that is  
13 clearly based on our assessment that for the data  
14 presented, the significantly increased incidence of  
15 this finding in fact doesn't translate into  
16 significant clinical sequelae once, at least, you  
17 get enough competence in the surgical group to stop  
18 putting stitches into the darn things.  
19 DR. TRACY: But I certainly wouldn't  
20 advocate removing page 5 of 10.  
21 DR. PINA: No, no, I'm not advocating  
22 that, either.  
23 DR. KRUCOFF: No, no. I think this is  
24 just consistent with what was brought up before  
25 about the statement if we are going to say

1 "endothelial injury," that implies something that  
2 is unsafe or bad. I think if we say this is a  
3 finding that has not appeared to translate into a  
4 clinically significant finding, that to me is  
5 consistent with our saying this is safe. And I  
6 think that's one where we can and need to be  
7 rigorous, particularly if we are unanimous.

8 DR. EDMUNDS: If you use the term  
9 "endothelial injury," detectable only on  
10 post-repair epiaortic ecocardiography one-quarter  
11 of the time that it is there by the OR personnel.  
12 In other words, I think this is totally  
13 impractical.

14 DR. TRACY: But it is there, and it is  
15 part of the description of the patient study. The  
16 phrase, I believe, is fair the way it is described,  
17 and I believe it should be left in the labeling.  
18 Okay.

19 Question 5. "Please provide any other  
20 recommendations or comments regarding the labeling  
21 of this device."

22 I think we along the course made other  
23 comments. Unless anybody else has additional  
24 comments to make regarding the labeling, I think we  
25 have covered this one.

1           Okay. And Question 6: "If the data  
2 provided are not adequate to support safety and/or  
3 effectiveness, what additional data analysis or  
4 study would you require?"

5           I think we are going to get back into the  
6 same discussion about what we need by safety and  
7 effectiveness. I think everybody would be happy if  
8 there were some other cognitive endpoint that could  
9 be analyzed at some point, obviously not in this  
10 dataset since it was not collected, but I think  
11 that would be something that we would be looking  
12 for in future studies, other measures that might be  
13 more appropriate than the composite endpoint that a  
14 priori is going to miss the thing that you are  
15 looking for, or you are hoping to reduce.

16          DR. ZUCKERMAN: Dr. Marler, can you be  
17 more specific regarding what measurements of  
18 neurocognitive dysfunction you would be looking for  
19 in future studies?

20          DR. MARLER: Given a menu, can I pick my  
21 favorite neuropsychological tests? Yes, I can. I  
22 like several tests because they are easy to  
23 administer, take little time, and cause minimal  
24 irritation to the patient and presumably the  
25 surgeon.

1           But I think that we have found it possible  
2 in the research that I have sponsored in this from  
3 Wake Forest and now Johns Hopkins that if you put  
4 them in a room and don't let them out, it gives  
5 them specified limits that a neuropsychologist can  
6 reach some agreement on how to do  
7 neuropsychological evaluation in an efficient and  
8 cost-effective way.

9           It is not answering your question, but to  
10 say Trailmaking B or Trailmaking A or this or that  
11 test I don't think is helpful in this situation. I  
12 think there is a way to come up with cognitive  
13 evaluation, and I think it has been done before,  
14 and I think it would move the field forward.

15           You read different things in the newspaper  
16 from year to year. Sometimes it has cognitive  
17 effects, and I think most recently in the  
18 newspaper, it doesn't, at least, long-term effects.

19           Sorry to be so unhelpful. I think it is a  
20 question that can be answered; I'm not going to  
21 answer it now.

22           DR. TRACY: The other piece of  
23 effectiveness is what percentage of material is the  
24 device capturing. And it seems like there might be  
25 some other bench test that could be better designed

1 other than non-sticky things being passed through  
2 the system. There must be some other biologic way  
3 of testing to get a better idea of what percentage  
4 of material is being missed or picked up by the  
5 device.

6 I would suggest that being part of the  
7 mechanical effectiveness assessment that should be  
8 done. And I am still a little troubled by your  
9 original question about is there some design  
10 problem here with the device that is resulting in  
11 whatever these disruptions are. We need to  
12 somewhere along the line satisfy the FDA on that  
13 with this device and certainly with any future  
14 device.

15 Dr. Krucoff?

16 DR. KRUCOFF: I would also suggest at  
17 least for future work, because I think one of the  
18 things that Chris mentioned that we have seen in  
19 other filters, since these are circular devices,  
20 and whether they are aligned, whether they are  
21 rotationally aligned and actually transverse across  
22 the aorta or whether they are cockeyed, that at  
23 that level, if there were a way--imaging or  
24 otherwise--to get a sense of how frequently these  
25 things simply are or are not aligned the way you



1 would ideally envision them, to me, that would be  
2 useful information somewhere along the line in the  
3 evolution of these things.

4 DR. MARLER: I wanted to add to my comment  
5 that my intuition--since intuition is important  
6 today--is that the way to find effects on the brain  
7 is not to look at low-risk patients but to include  
8 patients who have had prior stroke, have atrial  
9 fibrillation, have high risk of stroke, and they  
10 are usually the ones who show the effects of  
11 interventions. It is easier to see an effect.

12 I would say that what you have is the  
13 baseline stroke rate that goes with the whole  
14 procedure here, and there was no increase or  
15 decrease, but in particularly risky patients, you  
16 might be more likely to see the cognitive effects  
17 and the stroke effects. That is based on a number  
18 of trials that I am saying that.

19 DR. AZIZ: I think the problem with that  
20 would be to read lighting up of the stroke in those  
21 patients may be unrelated to emboli. It may be  
22 hypertensive episodes. So the protective effects--

23 DR. MARLER: We are having increasing  
24 evidence--well, okay--stroke is a systemic problem,  
25 and what triggers it varies from individual and

1 from time to time, and what we are seeing here is  
2 these emboli may be more of a trigger of a whole  
3 cascade of events rather than the entire event  
4 itself.

5 DR. EDMUNDS: I think this will come out  
6 if this is approved and used over time, but we  
7 don't have the data to say anything about it now,  
8 and that's why I think we have to wait. But I'm  
9 sure that clinicians will want to use it in the  
10 high-risk patients; it's just horse sense.

11 DR. MARLER: I'm just trying to provide  
12 advice where I would go if I had to find clinical  
13 benefit. I am saying that the higher-risk  
14 patients, certainly with cognitive measures, but go  
15 where the things are happening.

16 DR. TRACY: Are there any other comments  
17 on this question?

18 DR. FERGUSON: I have one question. With  
19 the PMAs, we talk a lot about post-market approval  
20 studies. That is not an issue with the 510(k) as I  
21 understand; right?

22 DR. ZUCKERMAN: That's correct.

23 DR. FERGUSON: Thank you very much.

24 DR. TRACY: Are there additional questions  
25 or comments that the FDA would like to make?

1 DR. ZUCKERMAN: No.  
2 DR. TRACY: Does the sponsor have any  
3 additional comments or questions at this time?  
4 MS. CHANG: No, thank you.  
5 DR. TRACY: Mr. Dacey, any comments or  
6 questions?  
7 MR. DACEY: After over 4 years, this is my  
8 first 510(k) experience, and I certainly can't  
9 bring any clinical experience or intuition to this  
10 process, but speaking for the consumer, I am  
11 awfully glad that this process is taking place.  
12 That's all that I have to say.  
13 DR. TRACY: Thank you.  
14 I guess Mr. Morton, by the fact that he is  
15 not here, apparently has no additional comments.  
16 So at this point, we will have another  
17 open public hearing.  
18 Is there anyone in the audience who wishes  
19 to address the panel on today's topic?  
20 [No response.]  
21 DR. TRACY: If not, we'll close the open  
22 public hearing.  
23 Are there any final recommendations from  
24 the panel?  
25 [No response.]

1 DR. TRACY: Dr. Zuckerman?  
2 DR. ZUCKERMAN: Geretta, do you have  
3 something to read about finding out about each  
4 panel participant's view on this topic?  
5 MS. WOOD: No.  
6 DR. ZUCKERMAN: I stand corrected. Thank  
7 you.  
8 DR. TRACY: I think you heard them.  
9 The meeting is adjourned.  
10 Thank you all very much.  
11 [Whereupon, at 3:42 p.m., the proceedings  
12 were concluded.]  
13 - - -