

# TSE Decontamination: Studies Relevant to Facility and Equipment Cleaning

Prof. Dr. Christoph Kempf

ZLB

for PPTA TSEWG

- ✓ NaOH and NaClO most commonly used in industry
- ✓ Inconsistent results in literature


Compound	Time/ Temperature	Agent tested	Log ID <sub>50</sub> change	Comment(s)
1 N sodium hydroxide	1 hr/RT	Hamster Scrapie	≥5.5	Effective
2 N sodium hydroxide	2 hr/RT	Hamster Scrapie	≥5.1	Effective
0.09 N sodium hydroxide	2 hr/RT/121 °C/1 hr	Hamster Scrapie	≥7.4	Effective
0.9 N sodium hydroxide	2 hr/RT/121 °C/1 hr	Hamster Scrapie	≥7.4	Effective
1 N sodium hydroxide	1 hr/RT	Hamster Scrapie	6.0	Effective
0.25 N sodium hydroxide	2 hr/RT	Mouse CJD	~1	Ineffective
1 N sodium hydroxide	2 hr/RT	Mouse CJD	2.9	Moderately effective

Reviewed by: Douglas Lee\*, Wytold Lebing, Christopher Stenland, Stephen Petteway, Jr.

Downstream 27


- ....The only methods that appear to be completely effective under worst case conditions are strong sodium hypochlorite solutions or hot solutions of sodium hydroxide ... Taylor, Vet J. 159, 10 -17, 2000
- ...Residual infectivity has been detected following treatment of 263K with 1M sodium hydroxide...(Diringer & Braig, Lancet 1, 439-40, 1989; Ernst & Race, J Virol Meth 41,193-202, 1993.
- etc...

## Denaturation/Inactivation of Prions by NaOH-Treatment



263K hamster brain homogenate was incubated for various times at different temperatures and NaOH concentrations. PrP<sup>res</sup> was determined by WB using the 3F4 antibody


## Denaturation/Inactivation of Prions by NaOH-Treatment



Microsomes purified from scrapie-infected hamster brain were incubated in 0.1 N NaOH at 22°C for 15 or 60 min


# Denaturation/Inactivation of Prions by NaOH-Treatment

Kinetics of PrP<sup>res</sup> disappearance after exposure of scrapie-positive hamster brain homogenate to NaOH at RT



# Denaturation/Inactivation of Prions by NaOH-Treatment

Kinetics of PrP<sup>res</sup> disappearance after exposure of scrapie-positive hamster brain homogenate to NaOH at RT






ZLB



## Denaturation/Inactivation of Prions by NaClO-Treatment

- PrP<sup>c</sup> and PrP<sup>res</sup> are both not detectable by the standard WB after short NaClO (100ppm) treatment even if no Proteinase K is used.
  - The same was observed with recombinant PrP<sup>c</sup>.
  - Antibodies other than the 3F4 also do not result in a signal.
  - Model proteins are not detectable after NaClO treatment either by Coomassie, silver or auro dye staining.
- 
- Protein fragmentation is known to occur upon hypochlorite treatment by chloramine formation and nitrogen-centered radicals (Hawkins & Davies Biochem J 332, 617-625, 1998)

Inactivation of 263K by NaClO



Data:  
Rohwer R.

# Denaturation/Inactivation of Prions by NaOH/NaClO-Treatment

Treatment	Inactivation temperature or final concentration of chemical (wt/vol)	Decrease in log <sub>10</sub> LD <sub>50</sub> after exposure for			
		15 min	60 min	15 min	60 min
		CJD		Scrapie	
NaOH	0.1 <i>N</i>	4.8	4.8	5.0	6.0
	1.0 <i>N</i>	4.5	≥5.0†	6.0	≥6.8†
NaOH	0.25 <i>N</i>	4.7	4.8		
NaOH + 1% SDS	0.25 <i>N</i>	≥4.0	4.0		
NaOCl	0.5%			4.5	4.5
	1.0%			6.0	6.0
	2.5%		3.3	≥6.0	≥7.0

PAUL BROWN, ROBERT G. ROHWER,  
D. CARLETON GAJDUSEK

THE JOURNAL OF INFECTIOUS DISEASES • VOL. 153, NO. 6 • JUNE 1986



Sanitization of stainless steel surfaces was evaluated

- Hamster brain homogenate dried on stainless steel coupons
- Shaken in 0.1 N NaOH for 30 min, then purified water
- Recovery of PrP<sup>Sc</sup> by swabbing
- Experiment performed in triplicate

### Results:

- Quantitative recovery of PrP<sup>Sc</sup> from untreated coupon
- No residual PrP<sup>Sc</sup> recovered from 0.1 N NaOH treated coupons

	Input Load [log <sub>10</sub> ]	Output Load [log <sub>10</sub> ]	Reduction [log <sub>10</sub> ]
Experiment 1	5.6	< 1.6	> 4.0
Experiment 2	5.1	< 1.6	> 3.5
Experiment 3	6.1	< 1.6	> 4.5




ZLB



## Evaluation of Sanitization

Iron powder previously loaded with scrapie-positive hamster brain homogenate was exposed to NaOH. Where indicated, iron powder samples were subsequently treated with PK (+). Prion proteins were dissolved in loading buffer and detected by Western blot analysis (equivalent of approx. 50 cm<sup>2</sup> / lane)




## Infectivity of surface-bound mouse prions after various treatments

Inoculation	Sick/Total
<b>1. Uninfected wires</b>	
Untreated	0/3
<b>2. Infectious wires</b>	
Untreated	6/6
Sodium hydroxide (1M, 1 h, 25°C)	0/6
Formaldehyde (10%, 1 h, 25°C)	6/6
Guanidinium thiocyanate (4M, 16 h, 25°C)	0/6

**E. Flechsig et al.**

Molecular Medicine 7(10): 679–684, 2001

## Evaluation of sanitization of DEAE-sepharose loaded with 263K microsomes by 0.1N NaOH or WFI (mock) for 5.5 hrs



- 
- ✓ NaOH and NaClO treatment both destroy PrP<sup>res</sup>
  - ✓ Kinetics showed a rapid destruction of PrP<sup>res</sup> even at low NaOH and NaClO concentrations
  - ✓ Reduction of PrP<sup>res</sup> in the range of up to >4.5 Log could be demonstrated
  - ✓ These results are in agreement with many single point measurements based on infectivity data